


NOVA



User Manual

2.1.4 / November, 2018



Metrohm Autolab B.V.
Kanaalweg 29/G
3526 KM, Utrecht
The Netherlands
+31302893154
autolab@metrohm.com
<https://www.metrohm.com/en/products/electrochemistry/>

NOVA

2.1.4

User Manual

Metrohm Autolab Teachware
Metrohm Autolab B.V.
3526 KM, Utrecht

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above or at autolab@metrohm.com.

This documentation is protected by copyright. All rights reserved.

Table of contents

1	NOVA installation	1
1.1	Software compatibility	1
1.2	Hardware compatibility	2
1.3	Software installation	2
1.4	External devices	5
1.4.1	Metrohm Devices support	6
1.4.2	Metrohm Devices installation	7
1.4.3	Spectrophotometer support	7
1.4.4	Spectrophotometer installation	8
1.4.5	Autolab RHD Microcell HC support	8
1.4.6	Autolab RHD Microcell HC installation	9
1.5	Powering the instrument	9
1.6	Autolab hardware installation	10
1.7	Software license	12
1.8	Intended use	13
1.9	Options	13
2	Conventions	16
2.1	Scientific conventions	16
2.2	Software conventions	16
2.3	Numbering conventions	17
2.4	Warning label conventions	18
2.5	NOVA information, warnings and errors	18
2.6	NOVA menus and controls	20
3	Release notes	23
3.1	Version 2.1.4 release	23
3.1.1	Chrono charge discharge galvanostatic	24
3.1.2	Editable overlay plot properties	25
3.1.3	Disable event logging	26
3.1.4	Import table in the repeat loop during a running procedure ...	27
3.2	Version 2.1.3 release	28
3.2.1	Procedure information	28
3.2.2	Export to RelaxIS	29
3.2.3	OCP command	30
3.2.4	Import table in Repeat command	31
3.2.5	Light source shutter control	31
3.2.6	Demo database file	32



3.2.7	Colored data file command tiles	33
3.3	Version 2.1.2 release	34
3.4	Version 2.1.1 release	34
3.4.1	Signal names, identity and locations	34
3.4.2	Current range logging	37
3.4.3	Event logging	38
3.4.4	Export options for Spectrophotometer control panel	40
3.4.5	Export options for Spectrophotometer control panel	41
3.4.6	Spectroelectrochemistry procedure	42
3.5	Version 2.1 release	44
3.5.1	Search function	44
3.5.2	Check cell	46
3.5.3	Current interrupt	46
3.5.4	Spectrophotometer manual control	47
3.5.5	Spectroelectrochemical measurements	48
3.5.6	Repeat number in Repeat command	50
3.5.7	Custom command name	51
3.5.8	Zoom function	52
3.5.9	Electrochemical Frequency Modulation	53
3.5.10	Corrosion rate analysis	54
3.5.11	New Plots frame controls	55
3.5.12	Device drivers installation	57
3.6	Version 2.0.2 release	58
3.6.1	Managed schedules	58
3.6.2	New color picker	60
3.6.3	Data handling command shortcut button	61
3.6.4	Library filters	62
3.6.5	Extended Sampler information	63
3.6.6	Zoom function for data analysis commands	64
3.6.7	Custom name for Build signal command	65
3.7	Version 2.0.1 release	66
3.7.1	Procedure and data tags	67
3.7.2	New plot options	68
3.7.3	Import data	70
3.7.4	Number of recent items	70
3.7.5	Plot preview	71
3.7.6	Print plot	72
3.7.7	Region insensitivity	74
3.7.8	Electrochemical interface toggle	74
3.7.9	ECl10M measurements	75
3.7.10	Library column display	77
3.7.11	Data grid column display	77
3.7.12	Estimated duration	78
3.7.13	Interpolate command	79
3.7.14	Hydrodynamic analysis	79
3.8	Version 2.0 release	80
3.8.1	Dynamic data buffers	80

3.8.2	Value of Alpha	81
3.8.3	Autolab RHD Microcell HC support	81
3.8.4	PGSTAT204 and M204 combination with Booster10A support	82
3.8.5	ECI10M module support	82
3.8.6	AC voltammetry	83
4	Dashboard	84
4.1	Actions	85
4.2	Recent items	86
4.3	What's going on	88
4.4	Instruments panel	90
5	Instruments panel	92
5.1	Change the default instrument	95
5.2	Autolab control panel	96
5.2.1	Instrument information panel	97
5.2.2	Tools panel	98
5.2.3	Autolab display panel	127
5.3	Autolab RHD Microcell HC control panel	134
5.3.1	Autolab RHD Microcell HC hardware setup	134
5.3.2	Autolab RHD Microcell HC manual control panel	135
5.4	Autolab Spectrophotometer control panel	136
5.4.1	Autolab Spectrophotometer hardware setup	137
5.4.2	Autolab Spectrophotometer manual control panel	140
5.5	Metrohm devices control panel	149
5.5.1	Metrohm Dosino control panel	149
5.5.2	Metrohm Sample Processor control panel	155
5.5.3	Metrohm Stirrer control panel	164
5.5.4	Metrohm Remote box control panel	169
6	Library	173
6.1	Default procedures	175
6.2	Add location	177
6.3	Default save Location	179
6.4	Moving files to a new location	180
6.5	Remove location	180
6.6	Load from Library	182
6.7	Edit name and remarks	183
6.8	Rating and tagging	184
6.9	Preview plot	186
6.10	Column visibility	187



6.11	Filtering the Library	188
6.12	Sorting the Library	191
6.13	Rearranging Library columns order	192
6.14	Locating files	193
6.15	Delete files from Library	194
6.16	The data repository	195
6.17	Merge data	197
6.18	Search function	200
7	NOVA commands	203
7.1	Control commands	204
7.1.1	Message	204
7.1.2	Send email	205
7.1.3	Repeat	206
7.1.4	Increment	224
7.1.5	Play sound	227
7.1.6	Build text	228
7.1.7	.NET	229
7.1.8	Procedure Information	233
7.2	Measurement - general	235
7.2.1	Autolab control	236
7.2.2	Apply	239
7.2.3	Cell	240
7.2.4	Wait	240
7.2.5	OCP	246
7.2.6	Set pH measurement temperature	247
7.2.7	Reset EQCM delta frequency	248
7.2.8	Autolab R(R)DE control	250
7.2.9	MDE control	251
7.2.10	Synchronization	254
7.3	Measurement - cyclic and linear sweep voltammetry commands	257
7.3.1	CV staircase	257
7.3.2	CV linear scan	260
7.3.3	LSV staircase	262
7.4	Measurement - voltammetric analysis commands	264
7.4.1	Sampled DC voltammetry	265
7.4.2	Normal pulse voltammetry	267
7.4.3	Differential pulse voltammetry	270
7.4.4	Differential normal pulse voltammetry	273
7.4.5	Square wave voltammetry	276
7.4.6	PSA (Potentiometric stripping analysis)	279
7.4.7	AC voltammetry	283

7.5	Measurement - chrono methods commands	286
7.5.1	Record signals	286
7.5.2	Chrono methods	290
7.6	Measurement - impedance commands	302
7.6.1	FRA measurement	303
7.6.2	FRA single frequency	305
7.6.3	Additional properties	307
7.6.4	Electrochemical Frequency Modulation	327
7.7	Data handling commands	332
7.7.1	Windower	332
7.7.2	Build signal	337
7.7.3	Calculate signal	345
7.7.4	Get item	358
7.7.5	Import data	358
7.7.6	Export data	362
7.7.7	Generate index	364
7.7.8	Shrink data	365
7.8	Analysis - general commands	367
7.8.1	Smooth	368
7.8.2	Peak search	372
7.8.3	Regression	373
7.8.4	Derivative	376
7.8.5	Integrate	378
7.8.6	Interpolate	379
7.8.7	FFT analysis	379
7.8.8	Convolution	381
7.8.9	Calculate charge	385
7.8.10	Hydrodynamic analysis	386
7.8.11	ECN spectral noise analysis	387
7.8.12	iR drop correction	391
7.8.13	Baseline correction	392
7.8.14	Corrosion rate analysis	398
7.9	Analysis - impedance	403
7.9.1	Electrochemical circle fit	403
7.9.2	Fit and simulation	405
7.9.3	Kronig-Kramers test	449
7.9.4	Include all FRA data	452
7.9.5	Potential scan FRA data	453
7.10	Metrohm devices commands	454
7.10.1	Dosino	455
7.10.2	Sample Processor	460
7.10.3	Stirrer	468
7.10.4	Remote	469
7.11	External devices commands	472
7.11.1	Spectroscopy	472
7.11.2	External device control	479
7.11.3	RHD control	483



8	Default procedures	486
8.1	Cyclic voltammetry	487
8.1.1	Cyclic voltammetry potentiostatic	487
8.1.2	Cyclic voltammetry galvanostatic	490
8.1.3	Cyclic voltammetry potentiostatic current integration	493
8.1.4	Cyclic voltammetry potentiostatic linear scan	496
8.1.5	Cyclic voltammetry potentiostatic linear scan high speed	499
8.2	Linear sweep voltammetry	503
8.2.1	Linear sweep voltammetry potentiostatic	503
8.2.2	Linear sweep voltammetry galvanostatic	506
8.2.3	Linear polarization	509
8.2.4	Hydrodynamic linear sweep	512
8.2.5	Hydrodynamic linear sweep with RRDE	518
8.2.6	Spectroelectrochemical linear sweep	525
8.3	Voltammetric analysis	529
8.3.1	Sampled DC polarography	530
8.3.2	Normal pulse voltammetry	535
8.3.3	Differential pulse voltammetry	539
8.3.4	Differential normal pulse voltammetry	543
8.3.5	Square wave voltammetry	547
8.3.6	AC voltammetry	551
8.4	Chrono methods	555
8.4.1	Chrono amperometry ($\Delta t > 1$ ms)	556
8.4.2	Chrono coulometry ($\Delta t > 1$ ms)	559
8.4.3	Chrono potentiometry ($\Delta t > 1$ ms)	561
8.4.4	Chrono amperometry fast	564
8.4.5	Chrono coulometry fast	566
8.4.6	Chrono potentiometry fast	570
8.4.7	Chrono amperometry high speed	573
8.4.8	Chrono potentiometry high speed	576
8.4.9	Chrono charge discharge potentiostatic	580
8.4.10	Chrono charge discharge galvanostatic	583
8.5	Potentiometric stripping analysis	587
8.5.1	Potentiometric stripping analysis	587
8.5.2	Potentiometric stripping analysis constant current	588
8.6	Impedance spectroscopy	590
8.6.1	FRA impedance potentiostatic	591
8.6.2	FRA impedance galvanostatic	594
8.6.3	FRA potential scan	597
8.6.4	FRA current scan	601
8.6.5	FRA time scan potentiostatic	605
8.6.6	FRA time scan galvanostatic	608
8.6.7	Electrochemical Frequency Modulation	612
9	Additional measurement command properties	615
9.1	Sampler	616

9.2	Automatic current ranging	618
9.3	Cutoffs	620
9.3.1	Cutoff configuration	621
9.3.2	Combining cutoffs	623
9.4	Counters	624
9.4.1	Counter configuration	625
9.4.2	Counter action - Pulse	627
9.4.3	Counter action - Autolab control	629
9.4.4	Counter action - Shutter control	630
9.4.5	Counter action - Get spectrum	631
9.4.6	Combining counters	632
9.5	Plots	633
9.5.1	Default plots	634
9.5.2	Custom plots	634
9.5.3	Plot options	636
9.6	Automatic integration time	641
9.7	Value of Alpha	642
10	Procedure editor	644
10.1	Creating a new procedure	645
10.2	Global options and global sampler	647
10.3	End status Autolab	651
10.4	Procedure tracks	652
10.5	Procedure wrapping	653
10.6	Procedure zooming	654
10.7	Command groups	655
10.7.1	Grouping commands	655
10.7.2	Ungrouping commands	656
10.7.3	Renaming groups	657
10.8	Enabling and disabling commands	658
10.8.1	Disabling commands	658
10.8.2	Enabling commands	659
10.9	Adding and removing commands	660
10.9.1	Adding commands	660
10.9.2	Removing commands	667
10.10	Moving commands	668
10.10.1	Moving commands using the drag and drop method	669
10.10.2	Using the drag and drop method to move commands to a command group or a sub-track	670
10.11	Moving multiple commands	672
10.12	Stacking commands	674
10.12.1	Creating command stacks	675



10.12.2	Remove commands from stacks	677
10.13	Links	678
10.13.1	Viewing links	679
10.13.2	Creating links	681
10.13.3	Editing links	689
10.14	My commands	692
10.14.1	Saving a My command	693
10.14.2	Editing My commands	695
11	Running measurements	698
11.1	Starting procedure	698
11.2	Procedure validation	701
11.3	Procedure cloning	702
11.4	Plots frame	704
11.4.1	Displaying multiple plots	706
11.5	Real time modifications	708
11.5.1	Real-time properties modification	708
11.5.2	Procedure control	711
11.5.3	Reverse scan direction	712
11.5.4	Display the Manual control panel	713
11.5.5	Enable and disable plots	715
11.5.6	Q+ and Q- determination	718
11.6	End of measurement	719
11.6.1	Procedure time stamp	719
11.6.2	Post validation	720
11.7	Specify plot preview	721
11.8	Detailed plot view	722
11.8.1	Plot properties	723
11.8.2	Toggle the 3D view	725
11.8.3	Toggle the step through data mode	725
11.8.4	Add an analysis command	727
11.8.5	Zooming options	728
11.8.6	Print plot	729
11.8.7	Export plot to image file	731
11.8.8	Relocate plots	733
11.9	Viewing the data grid	737
11.9.1	Current range logged in the data grid	739
11.9.2	Events logged in the data grid	739
11.9.3	Formatting the data grid	740
11.9.4	Sorting the data grid	741
11.9.5	Changing the order of the columns in the data grid	742
11.9.6	Exporting the data from the data grid	743
11.10	Convert data to procedure	745

12 Data analysis	747
12.1 Smooth analysis	748
12.1.1 SG mode	750
12.1.2 FFT mode	752
12.2 Peak search	754
12.2.1 Automatic search mode	757
12.2.2 Manual peak search	757
12.2.3 Manual adjustments	768
12.2.4 Results	770
12.3 Regression analysis	771
12.4 Integrate	775
12.5 Interpolate	779
12.6 Hydrodynamic analysis	784
12.7 Baseline correction	788
12.7.1 Zooming in/out	793
12.7.2 Fine tuning the baseline correction	794
12.8 Corrosion rate analysis	796
12.8.1 Tafel Analysis	799
12.8.2 Polarization Resistance	803
12.9 Electrochemical circle fit	805
12.9.1 Zooming in/out	809
12.9.2 Fine tuning the baseline correction	810
12.9.3 Copy as equivalent circuit	812
12.10 Fit and simulation	813
12.10.1 Direct fitting or simulation	814
12.10.2 Fitting or simulation using the dedicated editor	816
12.10.3 Viewing the result	823
13 Data handling	827
13.1 Get item	827
13.2 Shrink data	830
14 Data overlays	834
14.1 Create an overlay	834
14.2 Adding data to an overlay	836
14.3 Changing overlay plot properties	838
14.4 Changing overlay plot settings	840
14.5 Hiding and showing plots	842
14.6 Remove data from overlay	845
14.7 Additional Overlay controls	847



15 Procedure scheduler	850
15.1 Remove instrument from schedule	852
15.2 Creating a procedure schedule	853
15.2.1 Open procedures	853
15.2.2 Recent procedures	855
15.2.3 Search Library	856
15.2.4 Remove procedure	857
15.3 Using synchronization points	858
15.4 Naming and saving the schedule	861
15.5 Running the schedule	864
15.5.1 Starting the complete procedure schedule	865
15.5.2 Starting the schedule sequentially	866
15.5.3 Procedure schedule control	867
15.6 Inspecting procedures or data	869
15.7 Schedule zooming	871
16 Hardware description	873
16.1 General considerations on the use of the Autolab potentiostat/galvanostat systems	874
16.1.1 Electrode connections	874
16.1.2 Operating principles of the Autolab PGSTAT	877
16.1.3 Environmental conditions	897
16.1.4 Noise considerations	897
16.1.5 Cleaning and inspection	899
16.2 Instrument description	900
16.2.1 Autolab N Series (AUT8) instruments	900
16.2.2 Autolab F Series (AUT8) instrument	912
16.2.3 Autolab MBA N Series (AUT8) instruments	925
16.2.4 Autolab Compact Series (AUT4/AUT5) instruments	927
16.2.5 Multi Autolab Series (MAC8/MAC9) instruments	941
16.2.6 Autolab 7 Series (AUT7) instruments	951
16.2.7 μ Autolab Series instruments	963
16.3 Module description	972
16.3.1 Common modules	973
16.3.2 Optional modules	998
17 Diagnostics	1181
17.1 Connecting the instrument	1181
17.2 Running the Diagnostics	1183
17.3 Integrator calibration	1186
17.4 Diagnostics options	1186
17.5 Firmware update	1187



18 Warranty and conformity	1189
18.1 Warranty	1189
18.2 Spare part availability	1190
18.3 Safety information	1190
18.4 EU Declaration of conformity	1192
18.4.1 EU Declaration of Conformity	1192
18.4.2 EU Declaration of Conformity	1195
18.4.3 Declaration of Conformity	1197
18.4.4 Declaration of Conformity	1199
18.5 Environmental protection	1201
Index	1202

1 NOVA installation

This chapter describes how to install the NOVA software on the host computer and to how connect Autolab and external devices to the host computer.

The NOVA installation package is supplied on CD-ROM or USB support provided with the Autolab instrument. It can also be downloaded from the Metrohm Autolab webpage.



NOTICE

Leave the Autolab disconnected from the computer when installing NOVA for the first time.

1.1 Software compatibility

NOVA requires Windows 7 or later as operating systems in order to run properly. NOVA can be installed on 32 bit and 64 bit versions of Windows.



NOTICE

Previous versions of Windows are not supported.

The minimum and recommended specifications are reported in *Table 1* and *Table 2*, respectively.

Table 1 Overview of the minimum specifications for NOVA

CPU	1 GHz or faster 32-bit (x86) or 64-bit (x64) processor
RAM	2 GB RAM
HD	20 GB available hard disc space
GPU	DirectX 9.0c compliant display adapter with 64 MB RAM

Table 2 Overview of the recommended specifications for NOVA

CPU	Intel Core i5 or equivalent AMD processor
RAM	8 GB RAM




Figure 1 The installation wizard

Click the **INSTALL** button to start the installation. The files will be copied on the computer. If needed, the installation folder can be changed using the installation wizard.

When prompted to do so, please click the **Install** button provided in the Metrohm Autolab Driver installation window (see figure 2, page 4).




Figure 2 Install the Metrohm Autolab device drivers



NOTICE

Make sure that the Always trust software from Metrohm Autolab B.V. check box is ticked.

When prompted to do so, please click the **Install** button provided in the Avantes Driver installation window (see figure 3, page 4).




Figure 3 Install the Avantes device drivers



NOTICE

Make sure that the Always trust software from Avantes BV check box is ticked.

The installer will indicate when the installation is completed, as shown in Figure 4.




Figure 4 The installation is complete

1.4 External devices

The following additional external devices can be connected to the host computer:

- **Metrohm liquid handling devices:** these devices can be used to handling liquid samples and to automate the handling thereof.
- **Spectrophotometers:** the Autolab or the supported Avantes spectrophotometers can be used to perform spectroelectrochemical measurements in combination with the Autolab potentiostat/galvanostat.
- **Autolab RHD Microcell HC:** this device can be used to perform temperature-controlled measurements.

1.4.2 Metrohm Devices installation

Connecting a **USB** controlled Metrohm device to the host computer will trigger the installation of the instrument (see figure 5, page 7).




Figure 5 The Metrohm Device Driver installer

The installation will complete automatically.



NOTICE

In order to control the supported Metrohm Sample Processors, an additional Windows component must be present on the computer. The controls for the Metrohm Sample processors use the **Microsoft msxml6.0.dll** library for the configuration files (XML file format). This component may not be preinstalled on every Microsoft operating system. Please ensure the availability of this dll on the operating. If this package is missing, please download the installation package from the Microsoft website.

1.4.3 Spectrophotometer support

NOVA provides support for Autolab spectrophotometers.

The following Autolab spectrophotometer are supported through a **USB** connection to the host computer:

- Autolab Spectrophotometer UA
- Autolab Spectrophotometer UB

Additionally, compatible Avantes spectrophotometers are also supported when connected to the host computer through a **USB** connection. The following devices are supported:

- AvaSpec ULS2048-USB2

1.4.6 Autolab RHD Microcell HC installation

No driver is required to control the Autolab RHD Microcell HC. When this type of device is connected to the host computer, it is immediately recognized by NOVA and listed in the **Instruments** panel.



NOTICE

The Autolab RHD Microcell HC is only detected by NOVA if a stage is connected to the controller with a cell mounted on the stage.



NOTICE

The Autolab RHD Microcell HC is connected to the host computer through a serial port. If no serial port is present on the computer, a USB to Serial port adapter can be installed. The drivers required for this adapter are not included in the installation package of NOVA and need to be installed separately.

1.5 Powering the instrument

In order to use the instrument, it must be connected to the mains using the mains connection socket, located on the back plane of the instrument. Before connecting the instrument to the mains make sure that the mains output voltage matches the value indicated on the main voltage indicator, located above the connector (*see figure 7, page 9*).



Figure 7 The required mains voltage is indicated above the connector




Figure 9 The Driver Manager application

The Driver Manager can be used at any time to select the driver to use to control the Autolab. Two drivers are available:

- **Nova only (recommended setup):** this is the latest driver for the Autolab, allowing up to 127 instruments to be connected to the host computer. This driver is compatible with 32 bit and 64 bit versions of Windows.
- **Legacy driver:** this is an older driver version which can be used in combination with the GPES or FRA software. No further developments are planned for this driver. The maximum number of devices connected to the host computer is 8. Data transfer may be slower than with the NOVA only driver. This driver is **only** compatible with 32 bit versions of Windows.



NOTICE

The **Nova only** driver will not work with previous versions of NOVA (version 1.10 and older). In order to use previous versions of NOVA, it is necessary to start the Driver Manager application provided with the previous version and select one of the available drivers provided with this previous version (please refer to the **Getting Started** manual of the previous version of NOVA for more information).

Click the installation button for the required driver to change the device driver and follow the instructions on screen. The selected driver will be installed for all connected Autolab instruments. New instruments connected to the host computer will be configured using the selected driver.



NOTICE

The **Driver Manager** application can be used to change the device driver at any time.

1.7 Software license

The Autolab NOVA software, and all its components, provided in conjunction with the Metrohm Autolab potentiostat/galvanostat instruments is copyrighted and owned by Metrohm Autolab.

The software is provided as a **Free Licensed Closed-Source** product with limited warranty. The software can be installed on any computer without specific authorization from Metrohm Autolab.

Metrohm Autolab retains the copyright to the software. You may neither modify nor remove references to confidentiality, proprietary notices or copyright notices. Modifications of the software in part or as a whole is not permitted.

Metrohm Autolab warrants that the software, when operated properly, is suitable for the specified use with the electrochemical instrumentation from Metrohm Autolab or compatible external instrumentation.

Metrohm Autolab is exempt from further warranty or liability. Metrohm Autolab is neither liable for third-party damages or consequential damage not for loss of data, loss of profits or operating interruptions, etc.

1.8 Intended use

All Metrohm Autolab products are designed for electrochemical research and development within the normal environment of a laboratory. The instrumentation shall therefore only be used for this purpose and within the specified environmental conditions. All other uses fall out of the scope of the instrumentation and may lead to voiding of any warranty.

1.9 Options

The application options can be defined by selecting the *Options* from the **Edit** menu. A window will be displayed, showing two different sections (see figure 10, page 13).




Figure 10 The application Options window

The following properties are available in the **General** section (see figure 10, page 13):

- **Auto save measured data:** specifies if measured data should be saved automatically at the end of each measurement, using the provided toggle. This option is on by default.
- **Embedded executable IF030:** specifies the path to the embedded application for instrument fitted with the **IF030** controller, using the provided button. This is a system property, do **not** change this unless instructed by Metrohm Autolab.
- **Embedded executable IF040:** specifies the path to the embedded application for instrument fitted with the **IF040** controller, using the provided button. This is a system property, do **not** change this unless instructed by Metrohm Autolab.



- **Logging:** specifies if error logging should be used and which level of logging should be used, if applicable, using the provided drop-down list. This is a system property, do **not** change this unless instructed by Metrohm Autolab.
- **Check for updates at:** specifies the URL for version checks of NOVA. This is a system property, do **not** change this unless instructed by Metrohm Autolab.
- **Number of recent items:** defines the number of recent items shown in the **Recent items** panel of the dashboard. The default value is 5. Please refer to *Chapter 4.2*, for more information on the recent items.
- **Defined tags:** provides a list of tags used in NOVA. This list is empty by default and is automatically populated by user-defined tags through the tagging feature of NOVA. If needed, tags can be removed or added to this list directly. More information on the use of tags can be found in *Chapter 6.8*.
- **Event logging:** specifies if event logging should occur during the measurement, using the provided toggle. This option is on by default. For more information about **Event logging** please refer to *Chapter 11.9.2*



CAUTION

Modifying the system properties shown in the **General** section can interfere with the operation of the instrument. Do **not** change these properties unless instructed by Metrohm Autolab.

The **Plots** section displays the default plot options used in NOVA for all plots (see *figure 11, page 15*).




Figure 11 The plot options

In this section all the default options for plots can be specified. Clicking the **Reset** button will reset all the options to the factory default values.



NOTICE

Please refer to *Chapter 9.5.3* for more information on the Plot options.

- To select any item on screen, click the item.
- To select consecutive items on screen, click the first item, press and hold down **[SHIFT]**, and then click the last item.
- To select nonconsecutive items on screen, press and hold down **[CTRL]**, and then click each item.

2.3 Numbering conventions

All **numeric** values are defined in NOVA according to the local culture defined for the Windows operating system. Depending on these settings, the decimal separator symbol can either be **.** or **,**.

Improper use of the local culture settings defined in Windows may lead to wrong values. For example, typing 0,3 in NOVA on a computer which uses the **.** as decimal separator will be validated as 3.



NOTICE

It is recommended to consult the local culture settings defined in Windows before using the NOVA software.

Scientific (exponential) numbering is done using the **e** or **E** symbol. A value of **1e2** or **1E2** is converted to **100**.

The following prefixes are using in NOVA for engineering notation:

- **T**, for *Tera* (1000000000000).
- **G**, for *Giga* (1000000000).
- **M**, for *Mega* (1000000).
- **k**, for *Kilo* (1000).
- **m**, for *Milli* (0.001).
- **μ**, for *Micro* (0.000001).
- **n**, for *Nano* (0.000000001).
- **p**, for *Pico* (0.000000000001).




Figure 12 Information is highlighted in blue

- **Warning:** any item highlighted in yellow indicates that an issue has been identified and that user intervention is recommended in order to resolve the issue (see figure 13, page 19). Whenever possible, the cause and a possible solution will be offered. It is possible to ignore the warning and continue working with the software however this may lead to invalid data.




Figure 13 Warnings are highlighted in yellow

- **Error:** any item highlighted in red indicates that a problem has been identified and that user intervention is required in order to resolve the error (see figure 14, page 20). Whenever possible, the cause and a possible solution will be offered. It is **not** possible to ignore the error. No measurements are possible until the error is resolved.

Redo 'Action name'	Redoes the specified action. Keyboard shortcut: [CTRL] + [Y]
Cut	Cuts the selected item(s) to the clipboard. Keyboard shortcut: [CTRL] + [X]
Copy	Copies the selected item(s) to the clipboard. Keyboard shortcut: [CTRL] + [C]
Paste	Pastes the items in the clipboard at the specified location. Keyboard shortcut: [CTRL] + [V]
Select All	Selects all visible items. Keyboard shortcut: [CTRL] + [A]
Options	Specifies the default options used in the application.
<hr/>	
View	
Zoom in	Zooms in on a plot. Keyboard shortcut: [CTRL] + [=]
Zoom out	Zooms out on a plot. Keyboard shortcut: [CTRL] + [-]
Fit all	Adjusts the plot area to the best possible scale. Keyboard shortcut: [F4]
Manual control	Displays the Manual control panel for the default instrument. Keyboard shortcut: [F10]
<hr/>	
Measurement	
Run	Starts the procedure defined in the selected tab on the default instrument. Keyboard shortcut: [F5]
Run on ►	Starts the procedure defined in the selected tab on the specified instrument.
Instrument #1	
Instrument #2	
...	



Pause	Pauses the running command in the selected tab.
Skip	Skips the running command in the selected tab.
Stop	Stops the measurement running in the selected tab.
<hr/>	
Help	
User manual	Displays the NOVA User Manual. Shortcut key: [F1]
About	Displays the About dialog.

3 Release notes

This chapter describes the release notes of the current and previous versions of NOVA. The release notes are provided in reverse chronology. The following version have been released:

- **Version 2.1.4:** minor update of NOVA 2.1.
- **Version 2.1.3:** minor update of NOVA 2.1. *Version 2.1.3 release (see chapter 3.2, page 28)*. This version was released on March 9th, 2018.
- **Version 2.1.2:** minor update of NOVA 2.1. *Version 2.1.2 release (see chapter 3.3, page 34)*. This version was released on May 5th, 2017.
- **Version 2.1.1:** minor update of NOVA 2.1 *Version 2.1.1 release (see chapter 3.4, page 34)*. This version was released on March 24th, 2017.
- **Version 2.1:** the current major release of NOVA *Version 2.1 release (see chapter 3.5, page 44)*. This version was released on November 15th, 2016.
- **Version 2.0.2:** minor update of NOVA 2.0 *Version 2.0.2 release (see chapter 3.6, page 58)*. This version was released on July, 6th, 2016.
- **Version 2.0.1:** minor update of NOVA 2.0 *Version 2.0.1 release (see chapter 3.7, page 66)*. This version was released on April, 1st, 2016.
- **Version 2.0:** the original major release of NOVA 2 *Version 2.0 release (see chapter 3.8, page 80)*. This version was release on October, 7th, 2015.

3.1 Version 2.1.4 release

Version 2.1.4 adds the following new functionality:

1. A new **Default Procedure** "Chrono charge discharge galvanostatic" was added, along with a new **Demo Database** file that demonstrates this procedure. *Chrono charge discharge galvanostatic (see chapter 3.1.1, page 24)*.
2. Possibility to edit plot properties within the **Overlay** tab. *Editable overlay plot properties (see chapter 3.1.2, page 25)*.
3. Option to disable event logging. *Disable event logging (see chapter 3.1.3, page 26)*
4. Possibility to import a table in the Repeat for multiple values command while the procedure is running. *Import table in the repeat loop during a running procedure (see chapter 3.1.4, page 27)*




Figure 16 Demo data file 22 shows galvanostatic charge discharge of a cylinder cell battery

3.1.2 Editable overlay plot properties

The possibility to edit plot properties has been added to the **Overlay** tab. This new functionality allows the plotted data series to be configured, including the point or line styles and colors. Each data series in the overlay plot can now be given a title that is propagated automatically to the legend. The overlay plot can be given a title as well. For more detailed information, please proceed to *Changing overlay plot properties* (see chapter 14.3, page 838).




Figure 17 Overlay tab with a new plot properties panel

3.1.3 Disable event logging

The possibility to disable event logging has been added to the **Application options**. When event logging is disabled, events such as cutoffs, user intervention, and the current range in which data points were collected will not be logged. Disabling event logging reduces slightly the size of the measurement data .nox file and reduces the RAM used when creating the .nox data file. Event logging is enabled by default. It is not recommended to disable event logging unless the measurement is expected to create a very large data file and RAM availability during data file creation is expected to be a problem. For more detailed information, please proceed to *Options* (see chapter 1.9, page 13).




Figure 18 The application options window



NOTICE

Modifying the options shown in the **General** section of the application options can interfere with the operation of the instrument. Do **not** change these properties unless instructed by Metrohm Autolab.

3.1.4 Import table in the repeat loop during a running procedure

The possibility to import a table while the procedure is running has been added to the **Repeat** command. This new functionality allows the repeat values to be imported at the moment the Repeat command is encountered in the procedure sequence. Therefore, changes can be made to the values in the table during the measurement sequence, until the moment the command is executed. Use of this feature is recommended for advanced NOVA users only because validation cannot be performed on the values within the table when it is imported during run time. For more detailed information, please proceed to: *Import a table using the Import table option (see chapter 7.1.3.2.4, page 214)*.

Figure 19 The repeat for multiple values command properties with a toggle to repeat a table at runtime



NOTICE

Use of the **Procedure information** command tile is optional.



NOTICE

More information on the **Procedure information** command is available in *Chapter 7.1.8*.

3.2.2 Export to RelaxIS

NOVA 2.1.3 offers the possibility to export electrochemical impedance spectroscopy data in the RelaxIS file format. This functionality has been added to the existing **Export data** command *Chapter 7.7.6*.



NOTICE

RelaxIS is a comprehensive software solution for the analysis of impedance spectra. This software is a product of **rhd instruments**. For more information, please visit <https://www.rhd-instruments.de/index.php/en/relaxis-en.html>.

The RelaxIS file format option is located in the File format drop down menu of the **Export data** command, as shown in (*see figure 21, page 29*).





Figure 21 The export data command with RelaxIS file format

When the RelaxIS file format is selected, the list of link-able properties is pre-configured to those supported by the RelaxIS file type, as shown in (*see figure 22, page 30*). Data from a **FRA measurement** can be linked

3.2.4 Import table in Repeat command

The **Repeat** for multiple values command now provides the possibility to import a table. The table can be imported from a CSV file, as shown in (see figure 23, page 31).



The screenshot shows a control panel titled 'Values' with a table and two configuration sections.

	Scan Rate	Concentration	Temperature
1	0,05	0,1	25
2	0,01	0,1	25
3	0,005	0,1	25
4	0,001	0,1	25

Add range

Begin value

End value

Number of values

Distribution

Import Table

Column delimiter

Decimal separator

Rows to skip

Figure 23 Import a table in the Repeat for multiple values command



NOTICE

More information on how to import a table in the Repeat command can be found in .

3.2.5 Light source shutter control

NOVA 2.1.3 adds the possibility to open and close the shutter of an Auto-lab light source from the Spectrophotometer control panel. A toggle is provided that allows the light source shutter to be opened and closed using a TTL pulse from the default connected Autolab instrument.




Figure 24 The light source shutter open toggle is located in the Spectrophotometer control panel



NOTICE

More information on the Spectrophotometer manual control panel of NOVA 2.1.3 can be found in *Chapter 5.4.2*.

3.2.6 Demo database file

A new file has been added to the **Demo database**. Demo data file 21 demonstrates linear sweep voltammetry with in-situ UV/Vis spectroscopy, as shown in *Demo database file* (see chapter 3.2.6, page 32).




Figure 25 Demo data file 21 shows LSV of ferrocyanide oxidation with in-situ UV/Vis spectroscopy

3.2.7 Colored data file command tiles

The appearance of the command tiles within a data file has been altered to provide a visual contrast to the command tiles within a procedure file. This style change aims to provide a visual cue that a data file is open in the **Procedure editor**. The tiles in the data file are shaded green. In addition, the word **DATA** now appears at the top left corner of the **Procedure editor**, as shown in (see figure 26, page 33).




Figure 26 Command tiles within a data file are shaded green



3.3 Version 2.1.2 release

Version 2.1.2 several bugs are corrected; no new functionality is added.

3.4 Version 2.1.1 release

Version 2.1.1 adds the following functionality:

1. New signal names and locations when using selected data analysis commands *Signal names, identity and locations* (see chapter 3.4.1, page 34).
2. Current range logged for all measurement commands *Current range logging* (see chapter 3.4.2, page 37)
3. Event logging for all measurement commands *Event logging* (see chapter 3.4.3, page 38).
4. Export possibility for Spectrophotometer manual control panel *Export options for Spectrophotometer control panel* (see chapter 3.4.4, page 40).
5. Step through data option added to the Spectrophotometer manual control panel *Export options for Spectrophotometer control panel* (see chapter 3.4.5, page 41).
6. A procedure for spectroelectrochemical measurements has been added to the **Default** procedures *Spectroelectrochemistry procedure* (see chapter 3.4.6, page 42).

3.4.1 Signal names, identity and locations

NOVA 2.1.1 changes the way the following analysis commands work:

- Smooth
- Derivative
- Integrate
- Baseline correction

In previous versions of NOVA, these commands created two result signals regardless of the number of signals in the source data. For example, when using the **Integrate** command in previous versions of NOVA, the calculated signals were called *Integration result X* and *Integration result Y*.

These calculated signals, only available in the analysis command itself, no longer had units or identity.

In NOVA 2.1.1, when either one of these four commands is used on data provided by a parent measurement command, the nature of the signals involved in the analysis command and their units are retained and the resulting signals are duplicated in the original parent measurement command.

For example, applying a **Smooth** command on i vs E data (WE(1).Current vs Potential applied), as shown in *Figure 27*, now produces two new signals called *Smoothed WE(1).Current* and *Potential applied*, as shown in *Figure 28*.




Figure 27 Adding a Smooth command to the i vs E plot

The calculated signals are automatically plotted (see *figure 28*, page 36).




Figure 28 Smoothed WE(1).Current vs Potential applied plot is created. Additionally, the calculated data is copied to the parent measurement command (CV staircase), as shown in Figure 29.




Figure 29 The calculated signals are also available in the parent command

The calculated *Smoothed WE(1).Current* signal is known by the NOVA procedure as a valid current signal, obtained by applying a **Smooth** command on the WE(1).Current signal of the **CV staircase** command.

3.4.2 Current range logging

NOVA 2.1.1 now logs the active current range for each data point recorded in all measurement commands. This information is stored in the data file and is reported in the data grid (see figure 30, page 37).



Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Index	Current range
0,288086	10,4533	2,86011E-7	0,287781	118	100 nA
0,290527	10,4777	2,88239E-7	0,290253	119	100 nA
0,292969	10,5021	2,97638E-7	0,292664	120	100 nA
0,29541	10,5265	2,98279E-7	0,295074	121	100 nA
0,297852	10,5509	2,96173E-7	0,297607	122	100 nA
0,300293	10,5753	2,99164E-7	0,299957	123	100 nA
0,302734	10,5997	3,02643E-7	0,302338	124	100 nA
0,305176	10,6241	3,08868E-7	0,305298	125	1 µA
0,307617	10,6485	3,0954E-7	0,30777	126	1 µA
0,310059	10,6729	3,07281E-7	0,310394	127	1 µA
0,3125	10,6973	3,1076E-7	0,312897	128	1 µA
0,314941	10,7217	3,21289E-7	0,315277	129	1 µA
0,317383	10,7461	3,21442E-7	0,317657	130	1 µA
0,319824	10,7705	3,18298E-7	0,32016	131	1 µA
0,322266	10,7949	3,19824E-7	0,322601	132	1 µA

Figure 30 The active current range is now reported in the data grid



NOTICE

This new feature only applies to measurements carried out with NOVA 2.1.1 or later. Measurements carried out with earlier versions of NOVA may not display the active current range properly.



NOTICE

More information on current range logging is available in *Chapter 11.9.1*.

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range	Overload	Cutoffs
0.478516	12.0811	4.81323E-7	0.479523	1	196	1,51173E-6	0	100 nA	Current overload	
0.480957	12.1055	4.83459E-7	0.481964	1	197	1,51173E-6	0	100 nA	Current overload	
0.483398	12.1299	4.85382E-7	0.484436	1	198	1,51173E-6	0	100 nA	Current overload	
0.48584	12.1543	4.89227E-7	0.486816	1	199	1,51173E-6	0	100 nA	Current overload	
0.488281	12.1787	4.91089E-7	0.489258	1	200	1,51173E-6	0	100 nA	Current overload	
0.490723	12.2031	4.92767E-7	0.49176	1	201	1,51173E-6	0	100 nA	Current overload	
0.493164	12.2275	4.93378E-7	0.49408	1	202	1,51173E-6	0	100 nA	Current overload	
0.495605	12.2519	4.9762E-7	0.496552	1	203	1,51173E-6	0	100 nA	Current overload	
0.498047	12.2763	4.99329E-7	0.499115	1	204	1,51173E-6	0	100 nA	Current overload	
0.500488	12.3007	5.03754E-7	0.501495	1	205	1,51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A
0.50293	12.3251	5.05707E-7	0.503937	1	206	1,51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A
0.505371	12.3495	5.05615E-7	0.506256	1	207	1,51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A
0.507813	12.3739	5.10681E-7	0.508759	1	208	1,51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A
0.510254	12.3983	5.12939E-7	0.511047	1	209	1,51173E-6	0	100 nA	Current overload	WE(1).Current > SE-07 A

Figure 32 The details of the cutoff condition are reported in the data grid

The same applies to user intervention, where the action performed by the user is reported in the grid (see figure 33, page 39).

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range	User events
0.710449	16.7508	7.02209E-7	0.71106	1	529	1,72456E-5	-1,44562E-5	1 µA	
0.708008	16.763	7.06177E-7	0.708679	1	530	1,72456E-5	-1,44562E-5	1 µA	
0.705566	16.7752	6.9458E-7	0.706207	1	531	1,72456E-5	-1,44562E-5	1 µA	
0.703125	16.7874	6.99463E-7	0.703766	1	532	1,72456E-5	-1,44562E-5	1 µA	
0.700684	16.7997	7.01599E-7	0.701324	1	533	1,72456E-5	-1,44562E-5	1 µA	
0.698242	16.8119	6.95801E-7	0.698853	1	534	1,72456E-5	-1,44562E-5	1 µA	
0.695801	16.8241	6.9519E-7	0.696411	1	535	1,72456E-5	-1,44562E-5	1 µA	
0.693359	16.8363	6.93359E-7	0.694061	1	536	1,72456E-5	-1,44562E-5	1 µA	
0.690918	16.8485	6.86035E-7	0.691528	1	537	1,72456E-5	-1,44562E-5	1 µA	Reverse scan direction toggled
0.688477	16.8607	6.91833E-7	0.689056	1	538	1,72456E-5	-1,44562E-5	1 µA	
0.686035	16.8729	6.84204E-7	0.686707	1	539	1,72456E-5	-1,44562E-5	1 µA	
0.683594	16.8851	6.88782E-7	0.684204	1	540	1,72456E-5	-1,44562E-5	1 µA	
0.681152	16.8973	6.81763E-7	0.681793	1	541	1,72456E-5	-1,44562E-5	1 µA	
0.678711	16.9095	6.72302E-7	0.679382	1	542	1,72456E-5	-1,44562E-5	1 µA	

Figure 33 User events are reported in the data grid



NOTICE

This new feature only applies to measurements carried out with NOVA 2.1.1 or later. Measurements carried out with earlier versions of NOVA may not display the recorded events properly.



NOTICE

More information on event logging is available in *Chapter 11.9.2*.

Furthermore, NOVA now provides indications whenever the measurement conditions can be improved, by highlighting the affected command in blue in the procedure editor (see figure 34, page 40).




Figure 34 Commands are highlighted in blue when the measurement conditions can be improved

3.4.4 Export options for Spectrophotometer control panel

NOVA 2.1.1 adds the possibility to export data measured in the **Spectrophotometer** control panel. The data can be exported to ASCII or Excel using the  button located in the top right corner of the control panel (see figure 35, page 40).




Figure 35 The measured data can be exported to ASCII or Excel

It is also possible to export the chart displayed in the **Spectrophotometer** control panel using the  button, as shown in Figure 36.




Figure 36 The chart can also be exported



NOTICE

More information on the manual control of the Autolab and Avantes spectrophotometers can be found in *Chapter 5.4*.

3.4.5 Export options for Spectrophotometer control panel

NOVA 2.1.1 adds the possibility to toggle the *Step through data* option on or off in the Spectrophotometer control panel using the  button in the top right corner of the control panel (see *figure 37, page 42*).




Figure 37 The Step through data option can be used in the Spectrophotometer control panel

When the *Step through data* mode is on, an additional indicator is added to the plot, showing the X and Y coordinates of the point indicated by the arrow. The indicator can be relocated anywhere in the plot area.



NOTICE

More information on the manual control of the Autolab and Avantes spectrophotometers can be found in *Chapter 5.4*.

3.4.6 Spectroelectrochemistry procedure

A procedure for spectroelectrochemical measurements in combination with support spectrophotometers has been added to the **Default procedures**. This procedure is based on a synchronized measurement, using a **LSV staircase** command. The procedure is shown in *Figure 38*.




Figure 38 The default procedure for spectroelectrochemical measurements

This procedure includes a counter that is used to trigger a spectroscopy measurement every points. The **Spectroscopy** command stacked on the **LSV staircase** command uses the data from the two preceding **Spectroscopy** commands to calculate the absorbance and transmittance automatically.



CAUTION

This procedure requires an Autolab spectrophotometer or a supported Avantes spectrophotometer.



NOTICE

More information on this procedure can be found in *Chapter 8*.



3.5 Version 2.1 release

Version 2.1 adds the following functionality:

1. Application-wide search function for Procedures, Data and Schedules *Search function (see chapter 3.5.1, page 44)*.
2. Check cell tool *Check cell (see chapter 3.5.2, page 46)*.
3. Cell off after current interrupt *Current interrupt (see chapter 3.5.3, page 46)*.
4. Manual control for Autolab and Avantes spectrophotometers *Spectrophotometer manual control (see chapter 3.5.4, page 47)*.
5. New command and command options for spectroelectrochemical applications *Spectroelectrochemical measurements (see chapter 3.5.5, page 48)*.
6. Repeat number added to **Repeat** command *Repeat number in Repeat command (see chapter 3.5.6, page 50)*.
7. Custom names for commands *Custom command name (see chapter 3.5.7, page 51)*.
8. Zoom function for procedure editor and schedule editor *Zoom function (see chapter 3.5.8, page 52)*.
9. New Electrochemical Frequency Modulation (EFM) measurement command available *Electrochemical Frequency Modulation (see chapter 3.5.9, page 53)*.
10. Corrosion rate analysis command expanded with Linear polarization analysis mode *Corrosion rate analysis (see chapter 3.5.10, page 54)*.
11. Improved **Plot** frame controls *New Plots frame controls (see chapter 3.5.11, page 55)*.
12. All device drivers are now included in the NOVA installer *Device drivers installation (see chapter 3.5.12, page 57)*

3.5.1 Search function

NOVA now provides the possibility to search for Procedures, Data or Schedules. A dedicated input field is located in the top right corner of the application (see figure 39, page 45).




Figure 39 A search box is provided in the top right corner

The search function can be used to specify a string, with or without wild-cards. When triggered, the search function will look for all Procedures, Data and Schedule items in all the *Locations* specified in the **Library**, except the Default procedures.

The results of the search will be reported in a dedicated tab, grouped by item type. The table controls used to display the results are the same as those used by the **Library**. The results can therefore be sorted or filtered as required.



NOTICE

The search function will look for all items that match the specified search string in the Name or Remarks.



NOTICE

More information on the search function can be found in *Chapter 6.18*.



3.5.2 Check cell

The Check cell tool is now available from the instrument control panel. This tool can be used to check the electrode connections and the noise level by performing five consecutive current or potential measurements and determining the average value and standard deviation of each measurement (see figure 40, page 46).




Figure 40 The Check cell tool can now be used to check the noise level

The tool can therefore be used to assess the instrument noise pickup and optimize the measurement conditions.



NOTICE

More information on the **Check cell** tool can be found in *Chapter 5.2.2.4*.

3.5.3 Current interrupt

The current interrupt tool has been modified and now allows the possibility to switch the cell off at the end of the measurement (see figure 41, page 47).






Figure 41 The current interrupt tool now provides the possibility to set the cell end state

The **Cell state after measurement** toggle, located in the **Settings** panel, can be used to specify the state of the cell at the end of the measurement. This toggle is off by default.



NOTICE

More information on the current interrupt tool can be found in *Chapter 5.2.2.2*.

3.5.4 Spectrophotometer manual control

NOVA now provides a complete manual control interface for Autolab and Avantes spectrometers. This interface can be used to setup the hardware configuration of the connected spectrophotometer and manually control the spectrophotometer (see figure 42, page 48).




Figure 42 Autolab and Avantes spectrophotometers can be manually controlled

Using this interface it is possible to acquire spectra using the specified properties. It is also possible to save measured spectra as dark and reference (blank) spectra and convert the measured data to absorbance, transmittance or reflectance.



NOTICE


More information on the manual control of the Autolab and Avantes spectrophotometers can be found in *Chapter 5.4*.

3.5.5 Spectroelectrochemical measurements

New measurement command and command options have been added to NOVA in order to facilitate spectroelectrochemical measurements. The **Avantes** command is now replaced with the **Spectroscopy** command, which can be used to control Autolab (and Avantes) spectrophotometers.

It is no longer necessary to initialize and close this type of device in a procedure and the new **Spectroscopy** command now supports a software acquisition mode, which can be used at any time without triggers (*see figure 43, page 49*).





The screenshot shows a 'Properties' dialog box with a 'Software trigger' section. The settings are as follows:

Property	Value	Unit
Command name	Software trigger	
Action	Software trigger	
Device name	ASM80001	
Start wavelength	175	nm
Stop wavelength	1333	nm
Integration time	500	ms
Number of averages	1	
Enable light source shutter control	<input checked="" type="checkbox"/>	
DIO connector	P1	
Shutter open	<input type="checkbox"/>	

A 'More' button is located at the bottom right of the dialog.

Figure 43 Software and hardware control is now possible to Autolab and Avantes spectrophotometers

New measurement options are available for all measurement commands that support them. These options can be used to control the light source shutter position or to acquire a spectrum on the spectrophotometer connected to the DIO port (see figure 44, page 50).




Figure 44 New control options are available for measurement commands



NOTICE

More information on the **Spectroscopy** command and the spectrophotometer control options are available in *Chapter 7.11.1* and *Chapter 9*, respectively.

3.5.6 Repeat number in Repeat command

The **Repeat** command now provides a new signal, Repetition number, that can be used in combination with other commands. This new signal is a single value that is incremented at the beginning of each repetition (see *figure 45, page 51*).




Figure 45 The Repetition number is now available in the Repeat command



NOTICE

More information on the Repeat command can be found in *Chapter 7.1.3*.

3.5.7 Custom command name

For improved readability in the procedure editor, it is now possible to specify a name for all commands in a procedure. Providing a custom name will overrule the default name of the command (see figure 46, page 52).




Figure 46 Custom names can now be given to all commands



NOTICE

When a custom name is provided for a command, the content of this command is no longer updated during a measurement, if applicable.

3.5.8 Zoom function

The procedure editor and the schedule editor now offer the possibility to zoom in or out at any time to increase or decrease the size of the items shown on screen. The controls for this new zoom function are located in the top right corner of the editor frame (see figure 47, page 52).




Figure 47 Zoom controls are now available

Using this function will either scale the size of the items and the text up or down (between 200 % and 50 % of the original size), as shown in Figure 48.




Figure 48 Zooming in on the procedure editor

The following zooming controls are available:

- **Zoom out:** decreases the scaling of the items and text shown on screen. The  button or **[CTRL] + [-]** keyboard shortcut can be used to do this.
- **Zoom to 100%:** resets the scaling of the items and text shown on screen to the default size. The  button or **[F4]** keyboard shortcut can be used to do this.
- **Zoom in:** increases the scaling of the items and text shown on screen. The  button or **[CTRL] + [=]** keyboard shortcut can be used to do this.



NOTICE

More information on the zoom controls of the procedure editor and the schedule editor can be found in *Chapter 10.6* and *Chapter 15.7*, respectively.

3.5.9 Electrochemical Frequency Modulation

This version of NOVA provides support for Electrochemical Frequency Modulation (EFM) measurements. These measurements are based on the application of a small amplitude voltage perturbation and recording of the electrochemical response of the cell. Using the measured data, corrosion rate information can be determined.

EFM measurements use a special two component sinewave modulation. During this type of measurements, the response from the cell at the applied frequency, higher harmonics of these frequencies and intermodulated frequencies are recorded. *Figure 49* shows a typical measurement.




Figure 49 Example of an EFM measurement



CAUTION

Electrochemical Frequency Modulation measurements require a **FRA32M** module.



NOTICE

More information on Electrochemical Frequency Modulation command can be found in *Chapter 7.6.4*.

3.5.10 Corrosion rate analysis

The **Corrosion rate analysis** command has been complemented with a new mode: **Polarization Resistance**. This analysis method is based on the **ASTM G59** standard and it uses the *Stern-Geary* equation to determine the corrosion current and the corrosion rate (see figure 50, page 55).




Figure 50 The Linear polarization method has been added to the Corrosion rate analysis command

Provided that the analysis is carried out in a low overpotential range with respect to the corrosion potential, the Linear polarization analysis method can provide a direct estimation of the corrosion current and corrosion rate, using user-defined Tafel slopes.



NOTICE

More information on the **Corrosion rate analysis** command can be found in *Chapter 7.8.14*.

3.5.11 New Plots frame controls

The **Plots** frame now provides new controls that can be used to disable plots (see figure 51, page 56).




Figure 51 Disabling plots in the Plots frame



NOTICE

Disabling plots can be done at any time.

It is also now possible to relocate the plot order or overlay plots by dragging the plots in the frame (see figure 52, page 57).




Figure 52 Rearranging the plot order



NOTICE

It is not possible to relocate plot during a measurement.



NOTICE

More information on the disabling of plots and the relocation of plots can be found in *Chapter 11.5.5* and *Chapter 11.8.8*, respectively.

3.5.12 Device drivers installation

The installation package of NOVA 2.1 now installs all required device drivers during the installation process, as described in *Chapter 1.3*.

The following drivers are installed:

- **Autolab device drivers:** required for using the Autolab potentiostat/galvanostat.
- **Metrohm device driver:** required for using any supported Metrohm liquid handling instrument.
- **Spectrophotometer device driver:** required for using any supported Autolab (or Avantes) spectrophotometer.




Figure 53 Schedules are now managed through the Library

The most recent **Schedules** are now also listed in the **Recent items** panel on the **Dashboard**, as shown in Figure 54.




Figure 54 The most recent Schedules are now listed in the Recent items panel



NOTICE

More information on the **Library** can be found in *Chapter 6*.




Figure 56 The Custom tab provides additional controls for specifying the color

On the Custom tab, colors can be specified using RGB values or by changing the hue of the selected color or by selecting any available color in the provided RGB color matrix.

3.6.3 Data handling command shortcut button

A shortcut button has been added to this version of NOVA allowing data handling commands to be added to a procedure or data. The  shortcut button, located in the top right corner of the procedure editor, works in the same way as the data analysis shortcut button already available in the NOVA (see figure 57, page 62).




Figure 57 The data handling shortcut button can be used to add data handling commands easily



NOTICE

The data handling commands shown in the popout menu depend on the selected command.



NOTICE

More information on the use of the data handling shortcut button can be found in *Chapter 13*.

3.6.4 Library filters

To facility data management, the **Library** now provides filtering options that can be used to force the **Library** to display items that fit within the specified filter conditions. *Figure 58* shows an example using two filter conditions, one on the instrument serial number and one on the rating.

Demo Database

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 04 - Hydrodynamic linear sweep	Fe ²⁺ /Fe ³⁺ , NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	14-6-2016 12:47:46	★★★★	<input type="checkbox"/> ★☆☆☆☆
Demo 05 - Fe(II) - Fe(III) on pcPt	Fe ²⁺ /Fe ³⁺ Reversibility Test - LSV with i	AUT71848	31-8-2015 14:40:14	14-6-2016 12:47:47	★★★★	<input type="checkbox"/> ★☆☆☆☆
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvanostatic	AUT71848	31-8-2015 11:27:11	14-6-2016 12:47:48	★★★★	<input checked="" type="checkbox"/> ★★★★★

Figure 58 The Library now provides filtering options for better data handling



NOTICE

More information on the **Library** filters can be found in *Chapter 6.11*.

3.6.5 Extended Sampler information

To provide more information on how and when the **Sampler** records the signals during any electrochemical measurement, the **Sampler** editor has been extended with a table that provides more details on the sampling conditions. In the **Sampler** editor, shown in *Figure 59*, an additional [More](#) button is now available.

CV staircase

Sampler

Signal	Sample	Average	d/dt
WE(1).Current	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
WE(1).Potential	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Power	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Resistance	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Charge	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
External(1).External 1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Integrator(1).Charge	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Integrator(1).Integrated Current	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Time	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Sample alternating

[More](#)

Figure 59 The Sampler editor now provides more information




Figure 61 Zooming in and out is now possible for relevant analysis commands



NOTICE

More information on this new option can be found in *Chapter 12.7.1* and in *Chapter 12.9.1* for the **Baseline correction** and **Electrochemical circle fit** commands, respectively.

3.6.7 Custom name for Build signal command

It is now possible to define a custom name to the signals generated by the **Build signal** command, as shown in *Figure 62*.




Figure 62 The Build signal command now offers the possibility to provide custom names

3.7.1 Procedure and data tags

It is now possible to assign tags to procedures and data files. This is a convenient tool for bookkeeping purposes. The controls for tagging data or procedures are provided in **Tags** sub-panel of the **Properties** panel of the procedure editor (see figure 63, page 67).

The screenshot shows the 'Properties' panel for a procedure named 'Demo 20 - Iron screw in seawater'. It features several input fields and buttons:

- Procedure name:** A text box containing 'Demo 20 - Iron screw i...'.
- Remarks:** A text area containing 'Linear sweep', 'voltammetry', and 'potentiostatic'.
- Estimated duration:** A numeric input field with '0' and a unit 's'.
- Buttons:** 'Edit End status' and 'More'.
- Tags section:** A 'Rating' field with five stars and an 'Add' button.

Figure 63 The Tags sub-panel provides controls for tagging data and procedure files

Two types of tags can be assigned to data or procedures:

- **Rating:** a rating based on a stars system can be assigned to each data or procedure file. By default, no stars are assigned to a file, but it is possible to change this at any time.
- **Tags:** text tags can be added to each data or procedure file. By default, no tags are assigned to a file, but it is possible to change this at any time.



NOTICE

The rating and tags are updated when the file is saved.

It is also possible to provide a rating and tags directly from the Library (see figure 64, page 68).

Demo Database							+ -	
Name ▲	Remarks	Instrument	Rating	Tags	Measurement date	Last modified		
Demo 01 - Copper	CuSO4 0.01 M	AUT71848	☆☆☆☆☆	Add	31-8-2015 11:07:50	11-2-2016 14:3		
Demo 02 - Lead de	Pb(ClO4)2 0.0		★★★★☆	Add	4-2-2009 11:04:15	11-2-2016 14:29		
Demo 03 - Bipoten	RRDE measurr	MAC80064#3	☆☆☆☆☆	Add	15-7-2013 13:45:21	11-2-2016 14:30		
Demo 04 - Hydrod	Fe2+/Fe3+, N	AUT71848	☆☆☆☆☆	Add	31-8-2015 13:53:57	11-2-2016 14:3		
Demo 05 - Fe(II) - F	Fe2+/Fe3+ Re	AUT71848	☆☆☆☆☆	Add	31-8-2015 14:40:14	11-2-2016 14:3		
Demo 06 - Galvanc	Lead depositic	AUT71848	☆☆☆☆☆	Add	31-8-2015 11:27:11	11-2-2016 14:3		
Demo 07 - Chrono	Example of fa:	AUT71848	☆☆☆☆☆	Add	1-9-2015 13:20:24	11-2-2016 14:3		
Demo 08 - Superca	Supercapacito	AUT71848	☆☆☆☆☆	Add	1-9-2015 13:29:23	11-2-2016 14:3		
Demo 09 - Superca	Supercapacito	AUT50229	☆☆☆☆☆	Add	1-9-2015 13:50:29	11-2-2016 14:3		
Demo 10 - Differen	Differential pu	AUT50477	☆☆☆☆☆	Add	18-8-2015 15:11:45	11-2-2016 14:3		

Figure 64 Tags can also be defined in the Library



NOTICE

More information on the rating and tagging of procedures and data can be found in *Chapter 6.8*.

3.7.2 New plot options

The plot options have been expanded in order to allow for additional control of the plotting of data (see *figure 65, page 69*).




Figure 65 New plot options have been added to the Axes sub-panel

The following options have been added to the **Axes** sub-panel:

- **Fixed scale (on/off)**: defines if the axis should be automatically scaled or if a fixed scale should be used. When this property is switched on, it is possible to define a minimum and maximum value for the axis.
- **Custom ticks (on/off)**: defines if custom major and minor ticks should be used for the axis. When this property is switched on, it is possible to define the distribution of major and minor ticks.




Figure 67 The number of recent items can be edited in the NOVA options

The default number of items is 5 and can be edited at any time.

3.7.5 Plot preview

NOVA now offers the possibility to assign one of the plots of a data file as a preview plot to display in the Library. This provides a quick preview of the data contained in each data file. The preview plot is shown in a tooltip (see figure 68, page 72).




Figure 68 A plot preview is now added to each data item in the Library



NOTICE

Data measured with previous versions of NOVA will create a preview plot when changes to the file are saved in the current version.



NOTICE

More information on the plot previews can be found in *Chapter 6.9*.

3.7.6 Print plot

It is now possible to print plots using the provided  button (see figure 69, page 73).




Figure 69 Plots can now be printed

A print preview dialog will be displayed, allowing finetuning of the print output (see figure 70, page 73).




Figure 70 A print preview dialog is shown




Figure 72 A tooltip shows the active electrochemical interface in bold



NOTICE

More information on the optional ECI10M module can be found in *Chapter 16.3.2.8*.

3.7.9 ECI10M measurements

In order to more easily identify measurements carried out with the **ECI10M** as the active electrochemical interface, the **(ECI10M)** suffix will be shown in the procedure editor, next to the serial number of the active instrument, below the procedure title (*see figure 73, page 76*).




Figure 73 The ECI10M suffix is shown in the procedure editor

For measurements that have been carried out with the **ECI10M**, the same suffix will be added to the instrument serial number in the **Library**(see figure 74, page 76).




Figure 74 The ECI10M suffix is shown in the Library



NOTICE

This suffix is not shown for measurements carried out with previous version of NOVA.

3.7.10 Library column display

The display settings of the grids used in the **Library** are now non-volatile. This change affects the following settings of the Library:

- **Column order:** the order in which the columns appear.
- **Column visibility:** the visibility of the available columns.
- **Sorting options:** the sorting options used in the Library.



NOTICE

The display settings used in the Library are stored on the local computer and can be defined for each type of Library location.



NOTICE

More information on the display settings used in the Library can be found in *Chapter 6.10* and *Chapter 6.13*.

3.7.11 Data grid column display

The display settings of the grids used in the **Data grid** are now non-volatile. This change affects the following settings of the data grid:

- **Column order:** the order in which the columns appear.
- **Column formatting:** the data formatting used in each column.
- **Sorting options:** the sorting options used in the data grid.



NOTICE

These settings are stored for each command in the data file.



NOTICE

More information on the display settings used by the data grid can be found in *Chapter 11.9*.




Figure 76 When commands are stacked, the Estimated duration takes underlying commands into account



NOTICE

The Estimated duration is determined based on the interval time and the estimated number of points as well as the duration of underlying commands, if applicable.

3.7.13 Interpolate command

The new Interpolate command, available in the Analysis - general group of commands, is now available. This command can be used to determine Y or X value based on a user-defined X or Y value, by linear interpolation.



NOTICE

More information on the Interpolate command can be found in *Chapter 7.8.6*.

3.7.14 Hydrodynamic analysis

The Hydrodynamic $i \propto \sqrt{\omega}$ command has been renamed to **Hydrodynamic analysis** and it has been extended with the Koutecký-Levich analysis method.

3.8.2 Value of Alpha

For the **CV staircase** command and the **LSV staircase** command, the *Alpha value* advanced property is now available (see figure 77, page 81).




Figure 77 The Alpha value property is available for the CV staircase and LSV staircase command



NOTICE

More information on the use of the value of Alpha can be found in Chapter 9.7.

3.8.3 Autolab RHD Microcell HC support

NOVA 2.0 introduces support for the Autolab RHD Microcell HC support (see figure 78, page 81).



Figure 78 The Autolab RHD Microcell system



Figure 80 NOVA introduces the support of the ECI10M module



NOTICE

More information on the ECI10M module can be found in *Chapter 16.3.2.8*.

3.8.6 AC voltammetry

This new version of NOVA provides a new command for AC voltammetry measurements *AC voltammetry* (see *chapter 7.4.7, page 283*). A default procedure for this electrochemical method is also available *AC voltammetry* (see *chapter 8.3.6, page 551*).

4 Dashboard

The **Dashboard** is the home screen of NOVA. Whenever NOVA starts, the **Dashboard** is always shown to the user (see figure 81, page 84).




Figure 81 The Dashboard

At any time, when NOVA is used, it is possible to show the display by clicking the home tab (🏠). The Dashboard provides four different panels:

- **Actions:** this panel provides a list a shortcut buttons to trigger a common task in NOVA *Actions* (see chapter 4.1, page 85).
- **Recent items:** this panel provides a list of the last procedures, data and schedule items *Recent items* (see chapter 4.2, page 86).
- **What's going on:** this panel provides messages to the user about ongoing or finished events in NOVA *What's going on* (see chapter 4.3, page 88).
- **Instruments:** this panel provides a list of connected instruments *Instruments panel* (see chapter 4.4, page 90).

4.1 Actions

The **Actions** panel provides a series of buttons that can be used to quickly trigger a common action or control of the NOVA software (see figure 82, page 85).




Figure 82 The Actions panel provides shortcut to the most common actions in NOVA

The following shortcut buttons are provided:

- **Open library:** this button open the **Library**. More information on the **Library** can be found in *Chapter 6*.
- **New procedure:** creates a new blank procedure. More information on the procedure editor can be found in *Chapter 10.1*.
- **Import procedure:** imports a procedure from a *.nox* file in the **Library**. More information on the **Library** can be found in *Chapter 6*.
- **Import data:** imports NOVA data from a *.nox* file in the **Library**. More information on the **Library** can be found in *Chapter 6*.
- **Import command:** imports a command from a *.noi* file in the **My Commands** group of command. More information on the **My Commands** can be found in *Chapter 10.14*.
- **New schedule:** creates a new procedure schedule. More information on the **Procedure scheduler** can be found in *Chapter 15*.
- **Import schedule:** imports a procedure schedule from a *.nos* file in the **Library**. More information on the **Procedure scheduler** can be found in *Chapter 15*.



NOTICE

The **Import data** action can also be used to directly import data from the **GPES** and **FRA** software into NOVA.

4.2 Recent items

The **Recent items** panel lists the most recent procedures, data and schedules items (*see figure 83, page 86*).




Figure 83 The Recent items panel shows the last procedures, data and schedules items



NOTICE

By default, the **five** most recent items are shown in the **Recent items** panel. This number can be adjusted in the NOVA Options *Options* (*see chapter 1.9, page 13*).

The last items that are saved are automatically updated in the **Recent items** panel each time an item is saved (*see figure 84, page 87*).




Figure 84 The Recent items panel is automatically updated when data or procedures are saved

It is possible to remove items from the **Recent items** panel by right-clicking an item and selecting the *Remove from recent items* option from the context menu (see figure 85, page 88).




Figure 86 The What's going on panel is used to provide messages to the user

At any time, it is possible to clear the **What's going on** panel using the **X** button (see figure 87, page 89).




Figure 87 Clearing the panel content



Figure 90 Clicking the refresh button will update the content of the Instruments panel

The content of the Instrument panel is updated when the  button is clicked (see figure 91, page 91).




Figure 91 The Instruments panel is refreshed

 <p>VIRT00001</p>	<p>This symbol is used to identify all Autolab PGSTAT204 instruments. These instruments have a serial number starting with AUT5.</p>
 <p>VIRT00001</p>	<p>This symbol is used to identify all Autolab PGSTAT101 instruments. These instruments have a serial number starting with AUT4.</p>
 <p>VIRT00001</p>	<p>This symbol is used to identify all Multi Autolab Series instruments (M101 and M204). These instruments have a serial number starting with MAC8 (for the M101 Multi Autolab systems) and MAC9 (for the M204 Multi Autolab systems).</p>
 <p>VIRT00001</p>	<p>This symbol is used to identify all µAutolab type II and µAutolab type III instruments. The µAutolab type II instruments are identified by a serial number starting with µ2AUT7 and the µAutolab type III instruments are identified by a serial number starting with µ3AUT7.</p>
 <p>VIRT00001</p>	<p>This symbol is used to identify all Autolab 7 Series instruments (PGSTAT302, PGSTAT30, PGSTAT12, PGSTAT100) as well as the older Autolab 9 Series instruments (PGSTAT30, PGSTAT20, PGSTAT10 and PGSTAT100). These instruments are identified by a serial number starting with AUT7 or USB7.</p>

The following tiles are used to identify the connected **Autolab RHD Microcell HC** controllers:

 <p>12</p>	<p>This symbol is used to identify all Autolab RHD Microcell HC controllers connected to the computer through a RS232 connection. These instruments are identified by their serial number (or device name).</p>
---	--

The following tiles are used to identify the connected **Autolab or Avantes** spectrophotometers:



 ASM80002	This symbol is used to identify all Autolab or Avantes spectrophotometers connected to the computer through a USB connection. These instruments are identified by their serial number (or device name).
---	--

The following tiles are used to identify the connected **Metrohm** devices:

 2429	This symbol is used to identify all Metrohm 800 Dosino devices connected to a USB controlled Metrohm device. These instruments are identified by their serial number (or device name).
 2358	This symbol is used to identify all Metrohm 801 Magnetic Stirrers or 804 Titration Stands with a stirrer connected to it (either Metrohm 802 Rod Stirrer or Metrohm 741 Magnetic Stirrer). These instruments are identified by their serial number (or device name).
 4589	This symbol is used to identify all Metrohm 814, 815 or 858 Sample Processor devices connected by USB to the host computer. These instruments are identified by their serial number (or device name).
 41774	This symbol is used to identify all Metrohm 6.2148.010 Remote Box devices connected to a USB controlled Metrohm Device. These instruments are identified by their serial number (or device name).

In *Figure 92*, two instruments are connected (a Multi Autolab system with Serial Number MAC91234 and an Autolab N Series instrument with serial number AUT81234). A virtual instrument, with serial number VIRT00001 is also connected. This instrument is identified as a PGSTAT204.

The serial number of the N series instrument is shown in bold underlined font (**AUT81234**) indicating that this is the *default* instrument.



NOTICE

When a measurement is started, it will always be executed on the *default* instrument, unless otherwise specified.

The following actions can be performed in the **Instruments** panel:

- Change the *default* **Autolab** instrument.
- Open the instrument control panel.

5.1 Change the default instrument

The default instrument, displayed in the bold underline in the **Instruments** panel, is the instrument used in any measurement, unless otherwise specified. This is also the instrument used for procedure validation purposes.

To change the default instrument, right-click any instrument tile in the **Instruments** panel and select the *Make [Instrument serial number] the default instrument* from the context menu (see figure 93, page 95).




Figure 93 Defining the default instrument



NOTICE

Only one instrument can be set as default instrument.

- **Instrument information panel:** this panel displays information about the instrument.
- **Tools panel:** this panel provides quick access to the hardware setup and a number of direct measurement tools like current interrupt and positive feedback.
- **Autolab display panel:** this panel provides a number of manual controls of the instrument.



NOTICE

The available channels in a multi channel Autolab are listed in the **Channels** sub-panel. Each channel is identified by a letter or a number. More information is provided in *Chapter 16.2.5*.

5.2.1 Instrument information panel

The **Instrument information** panel shown in the instrument control panel provides information on the selected instrument (*see figure 96, page 97*).




Figure 96 The Instrument information panel

This information is updated in real time and is provided for information only. The following items are listed:

- **State:** indicates the state of the instrument (idle or measuring).
- **Instrument type:** indicates the type of instrument.
- **Modules:** shows the extension modules of the instrument.
- **Embedded processor:** shows the type of embedded processor installed in the instrument (IF030 or IF040).
- **Embedded software:** shows the embedded application name and version number.



NOTICE

The i-Interrupt, Positive feedback, pH calibration and Reset integrator drift tools are only shown on if the instrument provides the functionality used by these tools.

5.2.2.1 Hardware setup

The  button can be used to edit the **Hardware setup**. The hardware setup screen shows three panels (*see figure 98, page 99*).




Figure 98 The hardware setup

The following panels are provided:

- **Autolab module panel:** used to specify the type of Autolab and additional properties of this instrument.
- **Additional modules panel:** provides a list of compatible extension modules which can be installed in the instrument or connected to the instrument.
- **Properties panel:** provides additional parameters for the extension modules.

5.2.2.1.1 Autolab module panel

The **Autolab module** panel can be used to specify the following properties (*see figure 99, page 100*):

- **Main module:** specifies the type of Autolab using the provided drop-down list.
- **Power supply frequency:** specifies if the mains frequency is 50 or 60 Hz.




Figure 100 Optional module can be selected in the optional modules panel



NOTICE

Only the FRA2, BA, FI20 - Integrator, Booster10A, Booster20A, External Devices, IME303, IME663 and MUX module have additional properties to display in the **Properties** panel.

5.2.2.2 Current interrupt

The button can be used to perform a current interrupt (i-Interrupt) measurement. This tool can be used to determine the uncompensated resistance, R_U .



NOTICE

This tool is not available for μ Autolab type II and type III instrument as well as the Autolab PGSTAT10.

During a current interrupt measurement, a constant potential is applied on the cell before the current interrupt circuit is triggered. This circuit interrupts the current flow in the cell and measures the potential decay. From the measured potential decay, the uncompensated resistance (R_u) value is determined, using a linear and an exponential regression.

Two values of the uncompensated resistance, R_u , are determined automatically at the end of the measurement:

- **Ru linear:** this value is obtained from a linear regression performed on the initial segment of the voltage decay.
- **Ru exponential:** this value is obtained from an exponential regression performed on the initial segment of the voltage decay.

Proper determination of this value requires an accurate measurement of the current. The measurements must therefore be carried out at a potential value where the current is high enough to be measured properly and the current range must be adjusted in accordance.



NOTICE

For accurate measurements, the current should be at least in the order of 1 mA.

When the i-Interrupt tool is used, the control screen for this tool will be displayed (*see figure 101, page 103*). The control screen provides two panels and one plot area.




Figure 101 The *i-Interrupt* tool

The **Settings** panel shows the properties used in the current interrupt measurement (see figure 102, page 103).




Figure 102 The *i-Interrupt* Settings panel

The following properties and controls are available:

- **High speed:** a toggle that can be used to switch the high speed ADC module (ADC10M or ADC750) off or on (default off).




Figure 103 The measured and fitted data

The measured data points are shown as a point plot. The linear regression is shown using a green line and the exponential regression is shown as a red line. The start and end value of the two regression methods are shown using vertical lines with matching colors.

It is possible to hide or show the measured data or the regression data by checking or unchecking the check boxes shown in the legend (see figure 104, page 105).




Figure 104 Using the check boxes to show or hide the data measured during the current interrupt



NOTICE

The copied value can be pasted in a suitable property field in NOVA.

5.2.2.3 Positive feedback

The Positive feedback button can be used to perform a positive feedback measurement. This tool can be used to determine the uncompensated resistance, R_u .



NOTICE

This tool is not available for μ Autolab type II and type III instrument as well as the Autolab PGSTAT10.

During a positive feedback measurement, a potential pulse is applied on the cell and the potential is recorded. The iR compensation value can be adjusted upwards manually until its value is close to the actual value of the uncompensated resistance, R_u . When the compensated resistance reaches a value close to the actual value of R_u , potentiostatic loop will start to ring. When the compensated resistance exceeds the R_u value, the potentiostatic loop is no longer stable and the instrument will oscillate.

When the positive feedback tool is used, the control screen for this tool will be displayed (see figure 108, page 107). The control screen provides two panels and one plot area.




Figure 108 The positive feedback tool

The **Settings** panel shows the properties used in the current interrupt measurement (see figure 109, page 108).



Settings

Current range

iR compensation value

DC potential V

Pulse potential mV

Step duration ms

Figure 109 The positive feedback Settings panel

The following properties and controls are available:

- **Current range:** the current range in which the positive feedback measurement is performed.
- **iR compensation value:** the value of the compensated resistance, Ω .
- **DC potential:** the start and stop potential applied during the positive feedback measurement, in V.
- **Pulse potential:** the potential value applied in the pulse during the positive feedback measurement, in V.
- **Step duration:** the duration of the pulse applied during the positive feedback measurement, in s.

Clicking the button initiates the positive feedback measurement, using the specified properties.



CAUTION

The positive feedback tool switches the cell on and applies a potential pulse. It is highly recommended to specify the measurement properties carefully before starting the measurement.

While the positive feedback measurement is running, the spinning symbol  is shown. The instrument cannot be used until the measurement is finished. When the measurement is finished, the measured data is displayed next to the **Settings** panel (see figure 110, page 109).




Figure 110 The measured data

The measured data shows the potential profile applied on the cell. Since the positive feedback tool uses an iterative approach, it is possible to adjust the value of the iR compensation value property and repeat the measurement (see figure 111, page 109).




Figure 111 The measured data and previous results shown in overlay

The Show previous results toggle provided in the **Results** panel can be used to show or hide the data from the previous measurement (see figure 112, page 109).



NOTICE

The tool only stores the current measured data and the data from the previous measurement.

Results

Show previous results

Figure 112 The previous data can be enabled and disabled

When the check cell tool is used, the control screen for this tool will be displayed *Figure 114*. The control screen provides two panels and one plot area.




Figure 114 The check cell tool

The **Settings** panel shows the properties using the cell check measurement (see *figure 115, page 111*).

Figure 115 The check cell Settings panel

The following properties and controls are available:

- **Current range:** the current range in which the cell check is performed.
- **Bandwidth:** a drop-down control that can be used to specify the bandwidth of the instrument (high stability, high speed or ultra-high speed).
- **Potential/Current:** a numeric field that can be used to specify the applied potential (in potentiostatic mode) or the applied current (in galvanostatic mode).

$$\frac{\bar{\sigma}}{[CR]}$$

If the ratio exceeds a value of 0.25, a message is shown, indicating that the measured noise is too high (see figure 117, page 113).

Cell noise level

Noise level may be too high for this current range



Figure 117 A message is shown if the noise levels are too high

The check cell tool also evaluates the stability of the feedback loop by testing if the measured potential, E_{measured} , is within the acceptable limits of the applied potential, E_{applied} , according to:

$$|E_{\text{measured}} - E_{\text{applied}}| \geq 2 \frac{20}{2^{16}} + 0.005$$

If this inequality is true, then the measured potential does not correspond to the applied potential and a message is displayed (see figure 118, page 113).

Check electrodes

The measured data is invalid. Measured 9,699V does not match the applied 0,000V. Please check the electrodes and the cables.



Figure 118 A message is shown if the data is invalid

In this case, the measurement is invalid and no data will be displayed. This problem usually occurs when the connections to the cell are not correct or when the reference electrode is not functional.



NOTICE

This test is only carried out when the check cell tool is used in **Potentiostatic** mode.

If the current range can (in case of a current underload) or must be (in case of an overload) optimized, a message will also be displayed at the end of the measurement (see figure 119, page 114).





Figure 120 The pH calibration tool

The **Settings** panel shows the properties used for the pH measurement (see figure 121, page 115).

Figure 121 The pH calibration Settings panel

The following properties can be specified using the drop-down lists:

- **Input mode:** defines how the pH electrode is connected (single, differential). This setting depends on the specifications of the pH sensor. For pH sensors fitted with an internal reference electrode, the single input mode is used. For pH sensors using an external reference electrode, the differential mode is used. For all supported Metrohm sensors, the input mode is single.
- **Temperature mode:** defines if the temperature is measured through the pH sensor, if possible, or if the temperature is specified manually. It is only possible to measure the temperature if the pH sensor is fitted with an internal temperature sensor.
- **Temperature:** defines the temperature at which the pH sensor is calibrated, in °C. This value can only be specified if the **Temperature mode** property is set to manual control.



Calibration

Calibration temperature 25,00 °C

pH	E (V)
4,000	0,177
7,000	0,000
9,000	-0,118

Print report

Figure 123 The Calibration panel

Predefined calibration data points are available. These data points are stored on the module (for the pX1000 module) or in a file locally stored on the computer (for the pX module).

Finally, a plot is shown on the right-hand side of the control panels. This plot shows the three stored calibration points and the regression line. Below the plot, the equation of the regression line and the slope is displayed, in V/pH units, as well as the offset and the correlation coefficient.

The **Print report** button located below the plot. This button can be used to create a printable calibration report *Printing calibration report* (see chapter 5.2.2.5.5, page 123).

The pH calibration tool can be used to perform the following tasks:

- **Clear all calibration points** *Remove all previous calibration points (see chapter 5.2.2.5.1, page 117)*
- **Add a new calibration point** *Adding calibration points (see chapter 5.2.2.5.2, page 119)*
- **Edit a calibration point** *Editing a calibration point (see chapter 5.2.2.5.3, page 122)*
- **Remove a calibration point** *Removing calibration points (see chapter 5.2.2.5.4, page 122)*
- **Save the calibration data** *Saving calibration data (see chapter 5.2.2.5.6, page 125)*

5.2.2.5.1 Remove all previous calibration points

Clicking the **X** button in the **Calibration** panel removes all calibration points from the table (see figure 124, page 118).




Figure 124 Removing all the calibration data



NOTICE

The calibration points are still stored in the on-board memory of the pX1000 module or on the computer for the pX module. The change to the calibration data points is only finalized when the pH calibration tool is closed.

After clearing all the calibration data points, the plot area is also cleared (see figure 125, page 118).





Figure 125 The cleared plot

5.2.2.5.2 Adding calibration points

The calibration of a pH sensor requires at least two data points. It is possible to add data points to the calibration data in two different ways:

- By manually adding data points to the calibration data.
- By measuring the pH of a buffer of known pH value with the connected pH sensor.


To manually add data points to the calibration data, click the first available cell in the **Calibration** panel and directly type the pH and corresponding potential value (see figure 126, page 119).



Calibration	
Calibration temperature 25,00 °C	
pH	E (V)
9,000	-0,118

Figure 126 Manually adding calibration data

To add a measured data point to the calibration data, select the required pre-defined pH buffer value from the drop-down list provided in the **Settings** panel (see figure 127, page 119).



Settings	
Input mode	Single
Temperature mode	Measured
Temperature	25 °C
pH calibration buffer	4,000
	4,000
	7,000
	9,000

Figure 127 Selecting the pH buffer value

It is also possible to directly type the pH value of the buffer in the **Settings** panel (see figure 128, page 120).



Settings

Input mode

Temperature mode

Temperature °C

pH calibration buffer

Figure 128 Manually specifying the pH buffer value

With the buffer value specified in the **Settings** panel, click the button in the **Measurement** panel when a stable and correct voltage is displayed (see figure 131, page 121).

Settings

Input mode

Temperature mode

Temperature °C

pH calibration buffer

Measurement

Measured pH -
dpH/dt -
Voltage E -0,119 V
dE/dt 0,000 V/s
Measured temperature 21,301 °C
dT/dt 0,000 °C/s

Figure 129 Accepting a calibration point

As soon as the button is clicked, the measured value and the specified buffer value are added to the calibration data (see figure 130, page 121).




Figure 130 The updated calibration data



NOTICE

The calibration data is automatically plotted on the right-hand side of the **Settings** panel when two or more calibration data points are specified.

It is possible to add more calibration points, using the method described above.

If the temperature is measured using the built-in temperature sensor, a validation message may be shown when accepting a new value if the temperature at which this new data point is measured differs by more than 0.5 °C from the existing calibration data (see figure 131, page 121).

Confirm temperature


The current calibration temperature differs more than 0.5 °C from the stored temperature. Are you sure you want to replace the calibration temperature?

Yes

No

Figure 131 A warning is shown when the temperature deviates by more than 0.5 °C

Clicking the **Yes** button validates the new data point despite the temperature difference. Clicking the **No** button cancels the validation of the new calibration point.




Calibration

Calibration temperature 25,00 °C

pH	E (V)
9,000	-0,118
7,000	0,000
4,000	0,177

Figure 134 Removing a calibration point

Press the **[Delete]** key on the keyboard to remove the complete row from the calibration data (see figure 135, page 123).



Calibration

Calibration temperature 25,00 °C

pH	E (V)
9,000	-0,118
4,000	0,177

Figure 135 The calibration point is removed

5.2.2.5.5 Printing calibration report

When the calibration is complete, it is possible to generate a printable report for bookkeeping purposes. To do this, click the **Print report** button located below the plot area (see figure 136, page 124).




Figure 136 Generating a calibration report

A report will be generated (see figure 137, page 125).




Figure 137 The calibration report

This report contains all the calibration data and regression data, as well as the date of the report.

5.2.2.5.6 Saving calibration data

When the pH calibration is finished, close the pH calibration tool. If the calibration data was modified, you will be prompted to save or discard the data (see figure 138, page 125).

Write calibration data

Calibration data has changed. Write calibration data to module?

Yes No

Figure 138 The data can be saved when closing the pH calibration tool

Saving the data will overwrite the existing calibration data stored in the on-board memory of the **pX1000** module or in the calibration file stored on the computer. The previous data can no longer be used after it is overwritten.



NOTICE

It is recommended to reset the integrator drift each time a current range is changed.

5.2.3 Autolab display panel

The **Autolab display** panel, shown in *Figure 140*, provides basic manual control of the Autolab and the extension modules that can be manually controlled.



Figure 140 The Autolab display panel



NOTICE

The actual content of the **Autolab display** panel depends on the hardware setup. In *Figure 140* manual control of the IME663 module and the R(R)DE electrode is available.

The **Autolab display** panel always shows the **Instrument** panel at the left-most position. This panel provides manual control of the Autolab potentiostat/galvanostat as well as an overview of the real-time values measured by the instrument.

The information is provided in three panels:



NOTICE

The iR compensation control is only available in potentiostatic mode.

5.2.3.2 Instrument Signals sub-panel

The **Instrument Signals** sub-panel provides the real time values measured by the Autolab (see figure 142, page 129).




Figure 142 The Instrument Signals sub-panel

The following signals are shown in the **Instrument Signals** sub-panel:

- **Potential:** the potential difference, in V, measured between the reference electrode (RE) and the sense electrode (S) or between the reference electrode (RE) and the working electrode (WE), depending on the type of instrument. This is a measured value and the noise level for this signal is reported using the bars located below the value.
- **Current:** the current, in A, flowing between the counter electrode (CE) and the working electrode (WE). This is a measured value and the noise level for this signal is reported using the bars located below the value.
- **Resistance:** the resistance value of the cell, calculated from the Potential and Current values, in Ω .
- **Power:** the power value of the cell, calculated from the Potential and Current values, in W.

The measured values reported in this sub-panel are shown with a noise level, represented by a number of small bars (between 0 and 8 bars). The bars are determined from the standard deviation, σ . Table 4 shows the correspondence between noise bars and the standard potential of potential and current (where [CR] is the current range of the instrument).

Table 4 Overview of the noise bars

Number of bars	Potential noise	Current noise
0	$\sigma = 0$	$\sigma = 0$
1	$\sigma < \frac{0.5}{2^7}$	$\sigma < \frac{0.5[\text{CR}]}{2^7}$



Number of bars	Potential noise	Current noise
2	$\sigma < \frac{0.5}{2^6}$	$\sigma < \frac{0.5[CR]}{2^6}$
3	$\sigma < \frac{0.5}{2^5}$	$\sigma < \frac{0.5[CR]}{2^5}$
4	$\sigma < \frac{0.5}{2^4}$	$\sigma < \frac{0.5[CR]}{2^4}$
5	$\sigma < \frac{0.5}{2^3}$	$\sigma < \frac{0.5[CR]}{2^3}$
6	$\sigma < \frac{0.5}{2^2}$	$\sigma < \frac{0.5[CR]}{2^2}$
7	$\sigma < \frac{0.5}{2^1}$	$\sigma < \frac{0.5[CR]}{2^1}$
8	$\sigma \geq \frac{0.5}{2^1}$	$\sigma \geq \frac{0.5[CR]}{2^1}$

5.2.3.3 Instrument Warnings sub-panel

The **Instrument Warnings** sub-panel provides the real time instrument warnings (see figure 143, page 130).




Figure 143 The Instrument Warning sub-panel

The following warnings can be displayed are displayed in the **Instrument Warnings** sub-panel:

- **Current:** this indicator will be lit when a current overload is detected. The current overload warning will be triggered whenever the measured current exceeds the measurable range of the active current range.
- **Potential:** this indicator will be lit when a potential overload is detected. The potential overload warning will be triggered whenever the output potential of the instrument reaches the compliance voltage limit.
- **Temperature:** this will be lit when a temperature overload is detected. The temperature overload warning will be triggered whenever the operating temperature of the instrument exceeds the maximum allowed value.



NOTICE

The **PGSTAT204** and **M204** module are fitted with an unrecoverable temperature overload circuit. When this warning is triggered, the instrument needs to be switched off completely in order to recover from the temperature overload.

- **Oscillation:** this indicator will be lit when the feedback loop cannot be regulated properly and oscillation is detected.



NOTICE

When the oscillation warning is triggered, the cell is automatically switched off for safety reasons. On instruments fitted with a **Cell enable** button on the front panel, this button must be engaged in order to recover from the oscillation warning.

5.2.3.4 Docking and undocking Autolab display panels

For convenience, it is possible to undock the **Autolab display** panel and display its content in a separate window. To do this, click the  button in the top right corner of the **Autolab display** panel (see figure 144, page 131).




Figure 144 Undocking the Autolab display panel




Figure 146 Closing parts of the undocked Autolab display panel window

The closed sub-panel will be removed from the window and the window will be resized, if applicable (see figure 147, page 133).




Figure 147 The updated Autolab display panel window

- **Temperature range:** this slider can be used to specify the minimum and maximum allowed temperature for the Autolab RHD Microcell HC controller. The default values are -50 °C and 100 °C by default. By clicking and dragging the black ends of the slider, the minimum and maximum temperature can be adjusted, as shown in *Figure 149*.
- **Hold time:** specifies a holding time to use during a measurement after the temperature of the Autolab RHD Microcell HC controller has stabilized. The default value is 120 s.
- **Equilibration condition:** specifies the minimum value of the first derivative of the temperature versus time to reach a stable temperature. The default value is 0.5 °C/min.
- **Equilibration time:** specifies the minimum time during which the equilibration condition must be valid in order to consider the temperature of the Autolab RHD Microcell HC controller stable, in s.
- **Temperature timeout:** specifies a maximum time which is allowed to pass for the Autolab RHD Microcell HC controller to stop adjusting the temperature, in min.

The temperature regulation of the Autolab RHD Microcell HC controller works in the following way:

1. After setting the new temperature, the actual temperature is measured.
2. The temperature is considered stable when the derivative of the temperature versus time is smaller than the **Equilibration condition** for a duration equal or longer than the **Equilibration time**.
3. If no stable temperature can be reached, the controller will stop regulating the temperature after the specified **Temperature timeout**.

5.3.2 Autolab RHD Microcell HC manual control panel

The **Manual control** panel of the Autolab RHD Microcell HC can be used to read the current temperature of the controller and can be used to set a new temperature of the controller (see *figure 150*, page 135).




Figure 150 The Autolab RHD Microcell HC Manual control panel

The following properties are available in the **Manual control** panel:




Figure 153 The Spectrophotometer control panel

The **Autolab Spectrophotometer control** panel has two sub-panels:

- **Hardware setup panel:** this panel displays instrument settings for the Autolab Spectrophotometer.
- **Spectrophotometer display panel:** this panel provides manual control of the Autolab Spectrophotometer.

5.4.1 Autolab Spectrophotometer hardware setup

The configuration of the connected **Spectrophotometer** can be adjusted in the **Hardware setup** panel (see figure 154, page 137).

Figure 154 The Spectrophotometer Hardware setup panel

The **Hardware setup** panel displays information or properties of the connected Spectrophotometer. The following properties are available:

- **Name:** an input field which can be used to give a dedicated name to the instrument. By default, the name of the instrument corresponds to the instrument serial number.



- **Serial number:** a read-only field that provides the serial number of the instrument.
- **Number of pixels:** a read-only field that provides the number of pixels of the detector of the instrument.
- **Minimum wavelength:** a read-only field that provides the lowest measurable wavelength of the detector the instrument.
- **Maximum wavelength:** a read-only field that provides the highest measurable wavelength of the detector the instrument.
- **Firmware version:** a read-only field that provides the firmware version of the instrument.
- **FPGA version:** a read-only field that provides the FPGA version of the instrument.
- **Integration delay:** an input field which can be used to specify the integration delay in ms.
- **Enable dark correction:** a toggle which can be used to enable or disables the dark correction (default OFF).
- **Dark correction percentage:** an input field which can be used to specifies the percentile value of dark correction (0-100 %).
- **Number of smoothing pixels:** an input field which can be used to specify the number of pixels used in the smoothing algorithm. When this value is set to 0, no smoothing is used. The optimal value depends on the fiber diameter, pixel size and type of spectrophotometer.
- **Use high resolution ADC:** a toggle which can be used to enables or disables the high resolution ADC of the spectrometer. When enabled, the measured values are resolved using a 16 Bit ADC, when disabled a 14 Bit ADC is used instead (default ON).

NOVA supports all Autolab Spectrophotometers and all Avantes USB 2.0 AvaSpec Spectrophotometers with a suitable firmware installed. The following firmware versions are supported:

- **000.031.000.000** or **009.031.000.000:** these two versions of the firmware support all options provided in NOVA.
- **009.028.000.000:** this firmware version supports all options provided in NOVA except spectrum averaging. When this firmware is detected, a warning symbol is shown in the **Hardware setup** panel, with an indication that an outdated firmware is detected, as shown in *Figure 155*.




Figure 155 A warning is shown when an outdated firmware is detected

- Other versions:** all other firmware versions are not supported in NOVA. When an unsupported firmware is detected, a warning error symbol is shown in the **Hardware setup** panel, with an indication that an unsupported firmware is detected, as shown in Figure 156. In this case, the spectrophotometer cannot be used.




Figure 156 An error is shown when an unsupported firmware is detected



NOTICE

Please contact Metrohm Autolab for information on the update process of the installed firmware.

Table 5 provides an overview of the optimal number of **Smoothing pixels** for the Autolab spectrophotometers in function of the fiber diameter.

Table 5 Optimal smooth pixel settings for the different optical fiber diameters

Fiber diameter (µm)	Optimal smoothing pixels
10	0



Fiber diameter (μm)	Optimal smoothing pixels
25	1
50	2
100	3
200	7
400	14
500	17
600	21



NOTICE

For information on the optimal number of **Smoothing pixels** for compatible Avantes spectrophotometer, please refer to the Avantes user manual

5.4.2 Autolab Spectrophotometer manual control panel

The **Spectrophotometer display** panel provides manual control of the connected Autolab (or Avantes) spectrophotometer (see figure 157, page 140).




Figure 157 The Spectrophotometer display panel

The following properties are available (see figure 158, page 141):

- **Start wavelength:** an input field which can be used to specify the start wavelength used in the measurement, in nm.
- **Stop wavelength:** an input field which can be used to specify the stop wavelength used in the measurement, in nm.

- **Integration time:** an input field which can be used to specify the integration time, in ms. The smallest integration time depends on the type of spectrophotometer used. For the Autolab Spectrophotometer instruments, the smallest possible value is 1.05 ms.
- **Number of averages:** an input field which can be used to specify the number of averages, as an integer.
- **Measurement mode:** a drop-down control that can be used to specify the measurement mode (continuous or single). In continuous mode, the spectrophotometer will acquire spectra until stopped by the user. In single mode, the spectrophotometer will acquire a single spectrum.
- **Light source shutter open:** a toggle that allows the light source shutter to be opened and closed with a TTL pulse from the default Autolab instrument.

The image shows a control panel for measurement properties. It includes the following fields and controls:

- Start wavelength:** 175 nm
- Stop wavelength:** 1333 nm
- Integration time:** 500 ms
- Number of averages:** 1
- Measurement mode:** A dropdown menu currently set to "Continuous".
- Light source shutter open:** A dropdown menu currently set to "Continuous".
- Start:** A button located at the bottom right of the panel.

Figure 158 The measurement properties

To start the acquisition of a spectrum, the **Start** button can be pressed. Depending on the Measurement mode property, the spectrophotometer will acquire one or more spectra and display the measured data in the plot on the right hand side (see figure 159, page 142).




Figure 159 Measured spectra are displayed in the plot on the right-hand side



NOTICE

The measured data is displayed in arbitrary units.

If needed, the measurement properties can be adjusted while spectra are acquired (see figure 160, page 142).




Figure 160 Measurement properties can be adjusted



NOTICE

While spectra are being acquired, the **Hardware setup** of the spectrophotometer cannot be adjusted.

In continuous measurement mode, the acquisition of data can be stopped by pressing the **Stop** button.

After stopping the acquisition, it is possible to save the last measured spectrum as a *Dark spectrum* or as a *Reference spectrum*, by clicking the  button or  button, respectively (see figure 161, page 143).




Figure 161 Saving a measured spectrum as Reference spectrum

When a *Dark* or *Reference* spectrum is saved, a check mark ( or ) will be visible in the top right corner of the **Spectrophotometer display** window (see figure 162, page 143).




Figure 162 Saved spectra are indicated by a check mark



NOTICE

It is possible to overwrite the saved *Dark* or *Reference* spectrum by clicking the associated buttons again.



NOTICE

Changing the acquisition properties will discard the saved *Dark* and *Reference* spectrum.

5.4.2.1 Display modes

The **Spectrophotometer display** panel provides the possibility to toggle between different display modes, using the buttons (**S**, **S_d**, **A**, **T**) located in the top right corner (see figure 163, page 144).




Figure 163 Controlling the display mode of the measured data

The following display modes are available:

- **Scope mode** (**S**): this mode shows the raw data from the spectrophotometer in arbitrary units. This display mode is always available.
- **Dark corrected scope mode** (**S_d**): this mode shows the raw data (S_{Measured}) from the spectrophotometer corrected with the stored *Dark* spectrum (S_{Dark}), in arbitrary units. This display mode is only available if a *Dark* spectrum is saved. The dark corrected scope data is calculated according to:

$$S_{-d} = S_{\text{Measured}} - S_{\text{Dark}}$$

- **Absorbance mode (A):** this mode shows the absorbance values calculated from the measured data (S_{Measured}), the stored *Dark* spectrum (S_{Dark}) and the stored *Reference* spectrum ($S_{\text{Reference}}$). This display mode is only available if a *Dark* and a *Reference* spectrum are saved.

$$A = -\log\left(\frac{S_{\text{Measured}} - S_{\text{Dark}}}{S_{\text{Reference}} - S_{\text{Dark}}}\right)$$

- **Transmittance mode (T):** this mode shows the transmittance values calculated from the measured data (S_{Measured}), the stored *Dark* spectrum (S_{Dark}) and the stored *Reference* spectrum ($S_{\text{Reference}}$). This display mode is only available if a *Dark* and a *Reference* spectrum are saved.

$$T = 100 \cdot \left(\frac{S_{\text{Measured}} - S_{\text{Dark}}}{S_{\text{Reference}} - S_{\text{Dark}}}\right)$$



NOTICE

The modes can be toggled while spectra are acquired.

5.4.2.2 Step through data

The **Spectrophotometer display** panel provides the possibility to toggle the *Step through data* mode on or off using the  button in the top right corner (see figure 164, page 145).




Figure 164 The *Step through data* option can be used in the Spectrophotometer control panel

When the *Step through data* mode is on, an additional indicator is added to the plot, showing the X and Y coordinates of the point indicated by the arrow, as shown in *Figure 164*.



NOTICE

The indicator is always shown for the first data point of the plot.

It is possible to relocate the indicator in the following ways (see *figure 165, page 146*):

- By clicking anywhere in the plot area: the indicator is relocated to the closest data point of the plot.
- Using the [←]/[→]: the indicator can be moved by 1 point at a time.
- Using the [←]/[→] and [CTRL]: the indicator can be moved by 10 points at a time.
- [←]/[→] and [CTRL] and [SHIFT]: the indicator can be moved by 100 points at a time.




Figure 165 It is possible to relocate the indicator using the mouse or keyboard

5.4.2.3 Export data and plot

The **Spectrophotometer display** panel provides the possibility to export the measured data. Measured value can either be exported to ASCII or Excel format or as an image, using the provided  and  buttons in the top right corner of the panel (see *figure 166, page 147*).




Figure 166 The measured data can be exported

Clicking the  button displays a pop-out menu providing controls of the format of the exported file (see figure 167, page 147).




Figure 167 The data can be exported to ASCII or Excel

The data can be exported as ASCII or to Excel. The following properties can be specified:

- **File format:** specifies the format of the output file (ASCII or Excel), using the provided drop-down list.
- **Write column headers:** a toggle that can be used to indicate if the names of the signals need to be added to the output file.
- **Column delimiter:** specifies the symbol used as a column separator, using the provided drop-down list. This property is only available for ASCII output.
- **Decimal separator:** specifies the decimal separator symbol used in the output file, using the provided drop-down list. This property is only available for ASCII output.

Clicking the  button displays a save dialog window which can be used to specify the filename and location (see figure 168, page 148).




Figure 168 Specifying the filename and location



NOTICE

All of the available data is exported to the file.

Clicking the  button displays a pop-out menu providing controls of the format the size of the exported image file (see figure 169, page 148).




Figure 169 Exporting the data as image

Two types of image types can be used when exporting plots:

- **Pixel based output:** the data is exported to a pixel based file format, with or without compression (*.bmp, *.png, *.jpg, *.tiff, *.gif). When this type is used, the size of the image is specified in pixels.
- **Vector based output:** the data is exported to a vector file format (*.emf, *.svg, *.wmf). When this type is used, the size of the image is specified in arbitrary units.

Clicking the  button displays a Windows explorer dialog which can be used to specify the path, name and file type used to create the output image file (see figure 170, page 148).




Figure 170 Specifying the name, location and type of output file

5.5 Metrohm devices control panel

Double clicking a Metrohm device tile in the **Instruments** panel opens the **Metrohm device control** panel in a new **tab**. Depending on the type of device, the content of the Metrohm device control panel will be different. Four categories of Metrohm devices are supported in NOVA:

- Metrohm 800 Dosino with 807 Dosing Cylinder
- Metrohm 814, 815 and 858 Sample Processor
- Metrohm 801 Magnetic Stirrer and Metrohm 804 Titration Stand with 741 Magnetic Stirrer or 802 Rod Stirrer
- Metrohm 6.2148.010 Remote Box

5.5.1 Metrohm Dosino control panel

The **Metrohm Dosino control** panel opens in a new tab when a **Dosino** tile, shown in the **Instruments** panel, is double clicked (*see figure 171, page 149*).




Figure 171 The Metrohm Dosino control panel

The **Metrohm Dosino control** panel shows three different sub-panels:

- **Dosino information panel:** this panel displays information about the Dosino.
- **Tools:** this panel provides quick access to the hardware setup.
- **Metrohm display:** this panel provides manual control of the Dosino.



NOTICE

The **Name** must be unique.

5.5.1.2 Dosino hardware setup

The configuration of the connected **Dosino** can be adjusted in the Hardware setup. To open the Hardware setup, click the dedicated button in the **Tools** panel (see figure 174, page 151).




Figure 174 Click the Hardware setup button in the Tools panel to open the Dosino Hardware setup

A new window will be displayed, showing the settings for the selected **Dosino** (see figure 175, page 152).

- **Port 4**
 - **Active:** specifies if the port is active or not, using the provided  toggle. Port 4 is active by default. The active state is set to off, the port will be skipped by **Prepare** and **Empty** actions executed on the Dosino.
 - **Rate:** the rate used by the port, in ml/minute. The maximum value is 150 ml/minute.

5.5.1.3 Dosino manual control

The **Metrohm display** panel provides controls which can be used to manually operate the selected **Dosino**. These controls can be used at any time (see figure 176, page 153).




Figure 176 The manual controls of the Dosino

The following controls are provided:

- **Prepare:** starts a single prepare cycle on the Dosino.
- **Fill:** fill the Dosing cylinder completely, using the specified fill port.
- **Empty:** starts an empty cycle on the Dosino.
- **Port:** selects the active port, using the provided drop-down list.
- **Volume:** an input field which can be used to specify a volume to manually dose using the Dosino, in ml.
- **Dosed:** a read-only field which shows the dosed volume.
- **Dose:** starts a dosing action, using the specified *Volume* and using the selected *Port*.




Figure 178 The undocked Metrohm display panel window

5.5.2 Metrohm Sample Processor control panel

The **Metrohm Sample Processor control** panel opens in a new tab when a **Sample Processor** tile, shown in the **Instruments** panel, is double clicked (see figure 179, page 156).




Figure 179 The Metrohm Sample Processor control panel

The **Metrohm Sample Processor control** panel shows three different sub-panels:

- **Sample Processor information panel:** this panel displays information about the Sample Processor.
- **Tools:** this panel provides quick access to the hardware setup.
- **Metrohm display:** this panel provides manual control of the Sample Processor.



CAUTION

The controls provided for the **Metrohm Sample Processor** depend on the configuration and the options installed on the device. The hardware setup needs to be adjusted in order to match the instrument configuration. For more information, please refer to the user manual provided with the instrument.

5.5.2.1 Sample Processor information panel

The **Sample Processor information** panel shown in the **Instrument control** panel provides information on the selected instrument (*see figure 180, page 157*).

4589

Name

Type 858

State Idle

Parallel execution

Figure 180 The Instrument information panel for the Metrohm Sample Processor

The following items and controls are listed:

- **Name:** an input field which can be used to give a dedicated name to the instrument. By default, the name of the instrument corresponds to the last four or five digits of the instrument serial number.
- **Type:** indicates the type of instrument. For the **Sample Processor**, the type can be 814, 815 or 858.
- **State:** indicates the state of the instrument.
- **Parallel execution:** a toggle that can be used to specify if the parallel execution is allowed for this device (off by default).

The **Name** can be edited if required. This is convenient for identifying the **Sample Processor** in NOVA. If a specific name is provided, this name will be used throughout the whole NOVA application to identify the **Sample Processor** (see figure 181, page 157).

4589

Name

Type 858

State Idle

Parallel execution

Figure 181 The Sample Processor name can be modified if required



NOTICE

The **Name** must be unique.



5.5.2.2 Sample Processor hardware setup

The configuration of the connected **Sample Processor** can be adjusted in the Hardware setup. To open the Hardware setup, click the dedicated button in the **Tools** panel (see figure 182, page 158).




Figure 182 Click the Hardware setup button in the Tools panel to open the Sample Processor Hardware setup

A new window will be displayed, showing the settings for the selected **Sample Processor** (see figure 183, page 158).




Figure 183 The Metrohm Sample Processor hardware setup

The following settings can be specified for each **Sample Processor**:

- **Rack specification:** defines the sample rack mounted on the **Sample Processor**, using the provided drop-down list. The available sample racks are identified by their Metrohm part numbers.


■ **Tower 1**

- **Active:** specifies if the tower of the **Sample processor** is active or not, using the provided  toggle. Towers are active by default. If the active state is set to off, the tower will not be available for use.
- **Lift rate:** the rate used by the lift, in mm/s. The maximum value is 25 mm/s.
- **Shift rate:** the rotation rate used by the sample rack, in °/s. The maximum value is 20 °/s.
- **Swing rate:** the swing rate used by the swing arm, if available, in °/s. The maximum is 55 °/s.
- **Work position:** the work position used by the lift of the tower, in mm with respect to the top of the tower. The maximum value is 235 mm.
- **Position limit:** the maximum position used by the lift of the tower, in mm with respect to the top of the tower. The maximum value is 235 mm. The position limit must always be larger or equal to the work position.
- **Pumps:** specifies if pumps are installed on the back of the tower or connected to the back of the tower.
- **Valves:** specifies if valves are installed on the back of the tower.
- **Stirrer:** specifies if a stirrer is connected to the back of the tower.
- **Peristaltic pump:** specifies if a peristaltic pump is installed on the side of the tower. This option is only available with the **Metrohm 858 Professional Sample Processor**.
- **Injection valve:** specifies if an injection valve is installed on the side of the tower. This option is only available with the **Metrohm 858 Professional Sample Processor**.



■ Tower 2

- **Active:** specifies if the tower of the Sample processor is active or not, using the provided  toggle. Towers are active by default. If the active state is set to off, the tower will not be available for use.
- **Lift rate:** the rate used by the lift, in mm/s. The maximum value is 25 mm/s.
- **Shift rate:** the rotation rate used by the sample rack, in °/s. The maximum value is 20 °/s.
- **Swing rate:** the swing rate used by the swing arm, if available, in °/s. The maximum is 55 °/s.
- **Work position:** the work position used by the lift of the tower, in mm with respect to the top of the tower. The maximum value is 235 mm.
- **Position limit:** the maximum position used by the lift of the tower, in mm with respect to the top of the tower. The maximum value is 235 mm. The position limit must always be larger or equal to the work position.
- **Pumps:** specifies if pumps are installed on the back of the tower or connected to the back of the tower.
- **Valves:** specifies if valves are installed on the back of the tower.
- **Stirrer:** specifies if a stirrer is connected to the back of the tower.



NOTICE

Tower 2 is only available with the **Metrohm 814** and **815 Sample Processors**.



NOTICE

The settings defined in the **Sample Processor** hardware setup affect the controls provided in the **Metrohm display**.

5.5.2.3 Sample Processor manual control

The **Metrohm display** panel provides controls which can be used to manually operate the selected **Sample Processor**. These controls can be used at any time (*see figure 184, page 161*).




Figure 184 The manual controls of the Sample Processor

The following controls are provided:



- **Tower 1**
 - **Rack position:** sets the position of the sample rack with respect to tower 1.
 - **Lift position:** sets the position of the lift on tower 1, in mm with respect to the top of the tower. Shortcut buttons are provided for the *Work* position (**W**), *Shift* position (**S**) and *Home* position (**H**).
 - **Pump 1:** switches pump 1 on or off using the provided  toggle.
 - **Pump 2:** switches pump 2 on or off using the provided  toggle.
 - **Valve 1:** switches valve 1 on or off using the provided  toggle.
 - **Valve 2:** switches valve 2 on or off using the provided  toggle.
 - **Stirrer:** sets the rotation rate of the stirrer, from -15 to 15. When the rotation rate is set to 0, the stirrer will stop.
 - **Peristaltic pump:** sets the rotation rate of the peristaltic pump, from -15 to 15. When the rotation rate is set to 0, the pump will stop.
 - **Injection valve:** sets the state of the injection valve, using the provided drop-down list (*Fill* or *Inject*).
- **Tower 2**
 - **Rack position:** sets the position of the sample rack with respect to tower 2.
 - **Lift position:** sets the position of the lift on tower 2, in mm with respect to the top of the tower. Shortcut buttons are provided for the *Work* position (**W**), *Shift* position (**S**) and *Home* position (**H**).
 - **Pump 1:** switches pump 1 on or off using the provided  toggle.
 - **Pump 2:** switches pump 2 on or off using the provided  toggle.
 - **Valve 1:** switches valve 1 on or off using the provided  toggle.
 - **Valve 2:** switches valve 2 on or off using the provided  toggle.
 - **Stirrer:** sets the rotation rate of the stirrer, from -15 to 15. When the rotation rate is set to 0, the stirrer will stop.
- **Hold/Continue:** holds the current action, if possible. This button is only enabled when the Sample processor is not idle and when the action carried out by the Sample processor can be held. When the Sample processor is held, the **Hold** button switches to a **Continue** button which can be clicked again to resume the action.

- **Stop:** stops the current action, if possible. This button is only enabled when the Sample processor is not idle and when the action carried out by the Sample processor can be stopped.



NOTICE

The controls shown in the **Metrohm display** panel depend on the hardware setup of the selected **Sample Processor**.

For convenience, it is possible to undock the **Metrohm display** panel and display its content in a separate window. To do this, click the  button in the top right corner of the **Metrohm display** panel (see figure 185, page 163).




Figure 185 Undocking the Metrohm display panel for Sample Processor control

The content of the **Metrohm display** panel will be duplicated in a new window on top of the main NOVA software window. This new window can be moved next to the main NOVA window, or to another computer display if available (see figure 186, page 164).




Figure 186 The undocked Metrohm display panel window

5.5.3 Metrohm Stirrer control panel

The **Metrohm Stirrer control** panel opens in a new tab when a **Stirrer** tile, shown in the **Instruments** panel, is double clicked (see figure 187, page 165).





Figure 187 The Metrohm Stirrer control panel

The **Metrohm Stirrer control** panel shows three different sub-panels:

- **Stirrer information panel:** this panel displays information about the stirrer.
- **Metrohm display:** this panel provides manual control of the stirrer.



NOTICE

The  button, located in the **Tools** panel, is disabled for the Metrohm Stirrer.

5.5.3.1 Stirrer information panel

The **Stirrer information** panel shown in the **Instrument control** panel provides information on the selected instrument (see figure 188, page 165).




Figure 188 The Instrument information panel for the Metrohm Stirrer

The following items and controls are listed:

- **Name:** an input field which can be used to give a dedicated name to the instrument. By default, the name of the instrument corresponds to the last four or five digits of the instrument serial number.



- **Type:** indicates the type of instrument. For the **Stirrer**, the type can be 801, 802 or 741, depending on the type of device.
- **State:** indicates the state of the instrument.

The **Name** can be edited if required. This is convenient for identifying the **Stirrer** in NOVA. If a specific name is provided, this name will be used throughout the whole NOVA application to identify the **Stirrer** (see figure 189, page 166).




Figure 189 The Stirrer name can be modified if required



NOTICE

The **Name** must be unique.

5.5.3.2 Stirrer manual control

The **Metrohm display** panel provides a control which can be used to manually operate the selected **Stirrer**. This control can be used at any time (see figure 190, page 166).




Figure 190 The manual controls of the Stirrer

The following control is provided:

- **Speed:** specifies the rotation rate of the Stirrer.

The **Speed** value can be adjusted between -15 and 15, with integral increments. Negative values will force the stirrer to rotate in the anti-clockwise direction while positive values will translate into a clockwise rotation. The higher the value, the higher the rotation rate. The speed value can be specified in two ways:

- **Numerically:** by typing the value directly in the **Metrohm display** panel (see figure 191, page 167).




Figure 191 Adjusting the Speed numerically

- **Slider:** by clicking and dragging the slider control provided in the **Metrohm display** panel (see figure 192, page 167).




Figure 192 Adjusting the Speed with the slider

In both cases, the **Metrohm display** panel will be updated after the **Speed** value is adjusted (see figure 193, page 167).




Figure 193 The Speed value is adjusted




Figure 195 The undocked Metrohm display panel window

5.5.4 Metrohm Remote box control panel

The **Metrohm Remote box control** panel opens in a new tab when a **Remote box** tile, shown in the **Instruments** panel, is double clicked (see figure 196, page 169).




Figure 196 The Metrohm Remote box control panel

The **Metrohm Remote box control** panel shows three different sub-panels:

- **Control box information panel:** this panel displays information about the control box.
- **Metrohm display:** this panel provides manual control of the control box.



NOTICE

The Hardware setup button, located in the **Tools** panel, is disabled for the Metrohm Control box.

5.5.4.1 Remote box information panel

The **Remote box information** panel shown in the **Instrument control** panel provides information on the selected instrument (see figure 197, page 170).




Figure 197 The Instrument information panel for the Metrohm Remote box

The following items and controls are listed:

- **Name:** an input field which can be used to give a dedicated name to the instrument. By default, the name of the instrument corresponds to the last four or five digits of the instrument serial number.
- **Type:** indicates the type of instrument. For the **Remote box**, the type is 770.
- **State:** indicates the state of the instrument.

The **Name** can be edited if required. This is convenient for identifying the **Remote box** in NOVA. If a specific name is provided, this name will be used throughout the whole NOVA application to identify the **Remote box** (see figure 198, page 170).




Figure 198 The Remote box name can be modified if required



NOTICE

The **Name** must be unique.

5.5.4.2 Remote box manual control

The **Metrohm display** panel provides controls which can be used to manually operate the selected **Remote box**. These controls can be used at any time (see figure 199, page 171).




Figure 199 The manual controls of the Remote box

The following controls are provided:

- **Output:** specifies the state of the 14 output lines (numbered OUT13 to OUT0) of the Remote box. The state of each output line can be set to either *low* or *high* state, represented by a **0** or a **1**, respectively. The state of the 14 output lines is specified as a 14 character string, consisting of **0** and **1**, representing the state of the output lines, from **OUT13 to OUT0**.
- **Input:** specifies the state of the 8 input lines (numbered IN7 to IN0). The state of each input line can be either *low* or *high*, represented by a **0** or a **1**, respectively. This is a read-only control.

For convenience, it is possible to undock the **Metrohm display** panel and display its content in a separate window. To do this, click the  button in the top right corner of the **Metrohm display** panel (see figure 177, page 154).




Figure 200 Undocking the Metrohm display panel for Remote box control

The content of the **Metrohm display** panel will be duplicated in a new window on top of the main NOVA software window. This new window can be moved next to the main NOVA window, or to another computer display if available (see figure 178, page 155).




Figure 201 The undocked Metrohm display panel window

6 Library

The **Library** provides an interface to NOVA procedures, data and schedules. The Library can be accessed at any time from the **Dashboard**, by clicking the  button (see figure 202, page 173).




Figure 202 Opening the Library



NOTICE

The **Recent items** panel in the **Dashboard** shows the most recent procedures (Recent procedures) and data files (Recent data) when NOVA is used *Recent items* (see chapter 4.2, page 86). Both lists are updated whenever a new procedure or data file is saved or updated.

The **Library** opens in a new tab, represented by the  symbol. This tab is always located to the immediate right of the **Dashboard** tab (Home tab, ) , as shown in Figure 203.




Figure 203 The Library tab is opened to the right of the Dashboard tab

The **Library** tab contains two panels. The panel on the left-hand side is a navigation panel that provides the possibility to select a **location** in which procedure, data or schedule files are located. The panel on the right-hand side lists the available procedure, data or schedule files for the selected location.

Four locations are always visible in the navigation panel:

- **Default procedures:** this location provides all the factory default procedures installed with NOVA. These procedures cannot be deleted or modified. They can be loaded and modified in the procedure editor and saved as new procedures.
- **My procedures:** this location contains all the user-defined procedures. This location maps all the procedure files located in the \My Documents\NOVA 2.1\Procedures folder.
- **My data:** this location contains all the user generated data files. This location maps all the procedure files located in the \My Documents\NOVA 2.1\Data folder.
- **My schedules:** this location contains all the user generated schedules. This location maps all the schedule files located in the \My Documents\NOVA 2.1\Schedules folder.

6.1 Default procedures

The **Default procedures** location is always visible in the **Library** panel. This location contains a series of factory default procedures. These procedures are intended to perform simple measurements and can be used for routine experiments or as templates for more elaborate procedures. These procedures are provided as read-only examples. They cannot be deleted or modified but it is possible to open these procedures and to save them as a modified version in one of the user-accessible locations.



NOTICE

The procedures located in the **Default procedures** location are generated by the NOVA software. None of these procedures is available as an individual file on the computer.

The **Default procedures** are listed in the panel on the right-hand side (see figure 204, page 175).

Default procedures	
Name	Remarks
Cyclic voltammetry potentiostatic	Cyclic voltammetry potentiostatic
Cyclic voltammetry galvanostatic	Cyclic voltammetry galvanostatic
Cyclic voltammetry current integration	Cyclic voltammetry current integration: requires Integrator module
Cyclic voltammetry linear scan	Cyclic voltammetry linear scan: requires SCAN250 or SCANGEN module
Cyclic voltammetry linear scan high speed	Cyclic voltammetry linear scan high speed: requires SCAN250 or SCANGEN and ADC10M or ADC750 module
Linear sweep voltammetry potentiostatic	Linear sweep voltammetry potentiostatic
Linear sweep voltammetry galvanostatic	Linear sweep voltammetry galvanostatic
Linear polarization	Linear polarization
Hydrodynamic linear sweep	Hydrodynamic linear sweep: requires an R(R)DE connected
Hydrodynamic linear sweep with RRDE	Hydrodynamic linear sweep with RRDE: requires an RRDE and a BA module connected
Sampled DC polarography	Sampled DC polarography: requires IME module
Normal pulse voltammetry	Normal pulse voltammetry: requires IME module
Differential pulse voltammetry	Differential pulse voltammetry: requires IME module
Differential normal pulse voltammetry	Differential normal pulse voltammetry: requires IME module

Figure 204 The Default procedures

It is possible to expand the **Default procedures** location in the panel on the left-hand side to show groups of procedures that use the same type of experimental conditions (see figure 205, page 176).




Figure 205 Expanding the Default procedures

Selecting one of the groups in the **Default procedures** will reduce the number of procedures shown in the panel on the right-hand side (see figure 206, page 176).




Figure 206 Selecting a group in the Default procedures reduces the number of procedures displayed

6.2 Add location

It is possible to add one or more locations for either procedures, data or schedules by clicking the **+** button and selecting the required type of location (procedure or data) from the popout menu (see figure 207, page 177).




Figure 207 Adding a location to the Library

A Windows folder selection window will be displayed. Using this control, it is possible to navigate to the folder to be added to the list of locations (see figure 208, page 178).




Figure 208 Any folder can be added to the list of locations

The new location will be added to the **Library** and the content of folder will be displayed in the frame on the right-hand side (see figure 209, page 178).




Figure 209 The Module test folder is added as a location to the Library



NOTICE

It is possible to add as many new locations as needed.





NOTICE

Any sub-folder containing NOVA procedures, data or schedules located in the folder added as location in NOVA will be displayed in the **Library**.

6.3 Default save Location

When more than one location is specified in the **Library**, one of these locations will be used as the default save location. This will be used to save procedures or data unless otherwise specified. By default, the **My procedures**, **My data** and **My schedules** locations are used. It is possible to assign another location as the default save location by selecting it in the **Library** panel and right-clicking it. The context menu can be used to set the selected location as the new default location (*see figure 210, page 179*).




Figure 210 Defining the default save location



The new location will be used as default save location.



NOTICE

The default save locations are indicated with an underline font in the **Library** panel.

6.4 Moving files to a new location

When two or more locations are specified in the **Library** panel, it is possible to move files from one location to another by selecting the files and drag and dropping them from the source location to the destination location (see figure 211, page 180).




Figure 211 Moving files to a new location using the drag and drop method

The moved files will be removed from the source location and copied to the destination location.

6.5 Remove location

It is possible to remove a location from the **Library**, by clicking the  button, located in the top right corner of the **Library** panel. This will remove the highlighted location from the **Library** (see figure 212, page 181).




Figure 212 Click the remove button to remove a location

A confirmation message is displayed before the location is removed. Clicking the button removes the location. Clicking the button cancels the remove action (see figure 213, page 181).

Remove library collection

Are you sure you want to remove "Module test" from the library?



Figure 213 A confirmation message is shown when a Location is removed from the Library



NOTICE

Only the location is removed from the **Library**. The content of the folder associated with this location is not deleted from the computer.



NOTICE

It is possible to open more than one item at the same time by using the multi selection method and clicking the button.

6.7 Edit name and remarks

Using the **Library**, it is possible to change the name and remarks of an item. To change the name or remarks, click the cell to be edited in the table shown in the right-hand side panel of the **Library** to select it and then click it again to go in edit mode (see figure 216, page 183).

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper	CuSO4 0.01 M, H2SO4 0.1 M, A	AUT71848	31-8-2015 11:07:50	11-2-2016 14:31:16	☆☆☆☆☆	Add
Demo 02 - Lead de	Pb(ClO4)2 0.01 M / HClO4 0.1 l		4-2-2009 11:04:15	11-2-2016 14:29:57	☆☆☆☆☆	Add
Demo 03 - Bipoten	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	11-2-2016 14:30:28	☆☆☆☆☆	Add
Demo 04 - Hydrod	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	11-2-2016 14:31:29	☆☆☆☆☆	Add
Demo 05 - Fe(II) - f	Fe2+/Fe3+ Reversibility Test - l	AUT71848	31-8-2015 14:40:14	11-2-2016 14:31:34	☆☆☆☆☆	Add
Demo 06 - Galvanc	Lead deposition on gold, galva	AUT71848	31-8-2015 11:27:11	11-2-2016 14:31:20	☆☆☆☆☆	Add
Demo 07 - Chrono	Example of fast options messu	AUT71848	1-9-2015 13:20:24	11-2-2016 14:31:42	☆☆☆☆☆	Add
Demo 08 - Supercz	Supercapacitor, 3.3 F, CV, differ	AUT71848	1-9-2015 13:29:23	11-2-2016 14:31:49	☆☆☆☆☆	Add
Demo 09 - Supercz	Supercapacitor, 3.3 f	AUT50229	1-9-2015 13:50:29	11-2-2016 14:31:55	☆☆☆☆☆	Add
Demo 10 - Differer	Differential pulse voltammetry:	AUT50477	18-8-2015 15:11:45	11-2-2016 14:31:11	☆☆☆☆☆	Add

Figure 216 Editing the name or remarks

Edit the name or remarks and click away from the cell or press the **[Enter]** key to validate the change.



CAUTION

Changing the name of an item in the **Library** only changes the display name. The name of the file on the computer remains unchanged.



6.8 Rating and tagging

It is possible to directly edit the rating and tags for items in the **Library**. To specify the rating of an entry in the **Library**, click the highest star in the rating field (see figure 217, page 184).

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M, Ag/	AUT71848	31-8-2015 11:07:50	11-2-2016 14:31:16	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQCM	Pb(ClO4)2 0.01 M / HClO4 0.1 M		4-2-2009 11:04:15	11-2-2016 14:29:57	★★★★★	Add
Demo 03 - Bipotentiostat measurement	RRDE measurement	MAC60064#3	15-7-2013 13:45:21	11-2-2016 14:30:28	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear sweep	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	11-2-2016 14:31:29	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility Test - LSV	AUT71848	31-8-2015 14:40:14	11-2-2016 14:31:34	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvanoc	AUT71848	31-8-2015 11:27:11	11-2-2016 14:31:20	☆☆☆☆☆	Add
Demo 07 - Chrono measurement with f	Example of fast options measurem	AUT71848	1-9-2015 13:20:24	11-2-2016 14:31:42	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic voltan	Supercapacitor, 3.3 F, CV, differen	AUT71848	1-9-2015 13:29:23	11-2-2016 14:31:49	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impedance :	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	11-2-2016 14:31:55	☆☆☆☆☆	Add
Demo 10 - Differential pulse measurem	Differential pulse voltammetry: ret	AUT50477	18-8-2015 15:11:45	11-2-2016 14:31:11	☆☆☆☆☆	Add

Figure 217 Rating data or procedure items in the Library

It is also possible to edit the tags for a **Library** item. To add a tag, click the **Add** button and specify the tag to add to the item (see figure 218, page 184).

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M, Ag/	AUT71848	31-8-2015 11:07:50	11-2-2016 14:31:16	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQCM	Pb(ClO4)2 0.01 M / HClO4 0.1 M		4-2-2009 11:04:15	11-2-2016 14:29:57	☆☆☆☆☆	Add
Demo 03 - Bipotentiostat measurement	RRDE measurement	MAC60064#3	15-7-2013 13:45:21	11-2-2016 14:30:28	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear sweep	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	11-2-2016 14:31:29	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility Test - LSV	AUT71848	31-8-2015 14:40:14	11-2-2016 14:31:34	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvanoc	AUT71848	31-8-2015 11:27:11	11-2-2016 14:31:20	☆☆☆☆☆	Add
Demo 07 - Chrono measurement with f	Example of fast options measurem	AUT71848	1-9-2015 13:20:24	11-2-2016 14:31:42	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic voltan	Supercapacitor, 3.3 F, CV, differen	AUT71848	1-9-2015 13:29:23	11-2-2016 14:31:49	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impedance :	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	11-2-2016 14:31:55	☆☆☆☆☆	Add
Demo 10 - Differential pulse measurem	Differential pulse voltammetry: ret	AUT50477	18-8-2015 15:11:45	11-2-2016 14:31:11	☆☆☆☆☆	Add

Figure 218 Adding tags to data or procedure items in the Library

A popout field will be displayed, allowing specification of a text used for tagging the data or procedure item in the Library (see figure 219, page 185).

Name ▲	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M, Ag/AgCl	AUT71848	31-8-2015 11:07:50	11-2-2016 14:31:16	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQCM	Pb(ClO4)2 0.01 M / HClO4 0.1 M		4-2-2009 11:04:15	11-2-2016 14:29:57	☆☆☆☆☆	Add
Demo 03 - Bipotentiostat measurement	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	11-2-2016 14:30:28	☆☆☆☆☆	Demo X
Demo 04 - Hydrodynamic linear sweep	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	11-2-2016 14:31:29	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility Test - LSV	AUT71848	31-8-2015 14:40:14	11-2-2016 14:31:34	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvanostatic	AUT71848	31-8-2015 11:27:11	11-2-2016 14:31:20	☆☆☆☆☆	Add
Demo 07 - Chrono measurement with fast options	Example of fast options measurement	AUT71848	1-9-2015 13:20:24	11-2-2016 14:31:42	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic voltammetry	Supercapacitor, 3.3 F, CV, differential	AUT71848	1-9-2015 13:29:23	11-2-2016 14:31:49	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impedance	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	11-2-2016 14:31:55	☆☆☆☆☆	Add
Demo 10 - Differential pulse measurement	Differential pulse voltammetry; reference	AUT50477	18-8-2015 15:11:45	11-2-2016 14:31:11	☆☆☆☆☆	Add

Figure 219 The tag can be specified in the popout field

The tag will be added to the Tags column in the Library (see figure 220, page 185).

Name ▲	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M, Ag/AgCl	AUT71848	31-8-2015 11:07:50	11-2-2016 14:31:16	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQCM	Pb(ClO4)2 0.01 M / HClO4 0.1 M		4-2-2009 11:04:15	11-2-2016 14:29:57	☆☆☆☆☆	Add Demo X
Demo 03 - Bipotentiostat measurement	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	11-2-2016 14:30:28	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear sweep	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	11-2-2016 14:31:29	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility Test - LSV	AUT71848	31-8-2015 14:40:14	11-2-2016 14:31:34	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvanostatic	AUT71848	31-8-2015 11:27:11	11-2-2016 14:31:20	☆☆☆☆☆	Add
Demo 07 - Chrono measurement with fast options	Example of fast options measurement	AUT71848	1-9-2015 13:20:24	11-2-2016 14:31:42	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic voltammetry	Supercapacitor, 3.3 F, CV, differential	AUT71848	1-9-2015 13:29:23	11-2-2016 14:31:49	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impedance	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	11-2-2016 14:31:55	☆☆☆☆☆	Add
Demo 10 - Differential pulse measurement	Differential pulse voltammetry; reference	AUT50477	18-8-2015 15:11:45	11-2-2016 14:31:11	☆☆☆☆☆	Add

Figure 220 The tag is added to the item in the Library



NOTICE

It is possible to remove a tag by clicking the small X next to tag

Demo X



6.9 Preview plot

All data items in the **Library** provide a plot preview in the tooltip (see figure 221, page 186).




Figure 221 A plot preview is displayed in a tooltip

The plot preview is automatically generated when the data set is saved. By default, the first plot of the data set is used to create the plot preview, however if needed the preview plot can be edited.



NOTICE

Measurements performed with older versions of NOVA do not have a preview plot. This plot can be generated when changes to older files are saved in the current version of NOVA.




NOTICE

More information on specifying the preview plot can be found in *Chapter 11.7*.

6.10 Column visibility


For each type of location, the visibility of the columns shown in the **Library** can be edited. To hide a visible column in the **Library**, right-click the column header and select the column to hide from the context menu (see figure 222, page 187).



Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper dep	CuSO4 0.01 M, H2SO4 0.1 M, Ag/AgCl	AUT71848	31-8-2015 11:07:50	11-2-2016 14:31:16	☆☆☆☆☆	Add
Demo 02 - Lead depos	Pb(ClO4)2 0.01 M / HClO4 0.1 M	AUT71848	4-2-2009 11:04:15	11-2-2016 14:29:57	☆☆☆☆☆	Add
Demo 03 - Bipotentiost	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	11-2-2016 14:30:28	☆☆☆☆☆	Add
Demo 04 - Hydrodynam	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	11-2-2016 14:31:29	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility Test - LSV	AUT71848	31-8-2015 14:40:14	11-2-2016 14:31:20	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvanostatic	AUT71848	31-8-2015 11:27:11	11-2-2016 14:31:20	☆☆☆☆☆	Add
Demo 07 - Chrono measurement with f	Example of fast options measurement	AUT71848	1-9-2015 13:20:24	11-2-2016 14:31:42	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic voltammetry	Supercapacitor, 3.3 F, CV, differential	AUT71848	1-9-2015 13:29:23	11-2-2016 14:31:49	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impedance	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	11-2-2016 14:31:55	☆☆☆☆☆	Add
Demo 10 - Differential pulse measurement	Differential pulse voltammetry; reverse	AUT50477	18-8-2015 15:11:45	11-2-2016 14:31:11	☆☆☆☆☆	Add

Figure 222 Right-click the column header to hide a visible column

The column will be hidden (see figure 223, page 187).



Name	Remarks	Instrument	Measurement date	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M, Ag/AgCl	AUT71848	31-8-2015 11:07:50	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQCM	Pb(ClO4)2 0.01 M / HClO4 0.1 M	AUT71848	4-2-2009 11:04:15	☆☆☆☆☆	Add
Demo 03 - Bipotentiostat measurement	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear sweep	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility Test - LSV	AUT71848	31-8-2015 14:40:14	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvanostatic	AUT71848	31-8-2015 11:27:11	☆☆☆☆☆	Add
Demo 07 - Chrono measurement with f	Example of fast options measurement	AUT71848	1-9-2015 13:20:24	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic voltammetry	Supercapacitor, 3.3 F, CV, differential	AUT71848	1-9-2015 13:29:23	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impedance	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	☆☆☆☆☆	Add
Demo 10 - Differential pulse measurement	Differential pulse voltammetry; reverse	AUT50477	18-8-2015 15:11:45	☆☆☆☆☆	Add

Figure 223 The column is hidden

To make a hidden column visible again, right-click the column header and select the hidden column from the context menu (see figure 224, page 188).

Demo Database

Name	Remarks	Instrument	Measurement date	Rating	Tags
Demo 01 - Copper	0.4 0.01 M, H2SO4 0.1 M, Ag/AgCl	AUT71848	31-8-2015 11:07:50	☆☆☆☆☆	Add
Demo 02 - Lead de	Pb(ClO4)2 0.01 M / HClO4 0.1 M		4-2-2009 11:04:15	☆☆☆☆☆	Add
Demo 03 - Bipoten	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	☆☆☆☆☆	Add
Demo 04 - Hydrod	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility Test - LSV	AUT71848	31-8-2015 14:40:14	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvan	AUT71848	31-8-2015 11:27:11	☆☆☆☆☆	Add
Demo 07 - Chrono measurement with f	Example of fast options measur	AUT71848	1-9-2015 13:20:24	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic voltan	Supercapacitor, 3.3 F, CV, differen	AUT71848	1-9-2015 13:29:23	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impedance	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	☆☆☆☆☆	Add
Demo 10 - Differential pulse measur	Differential pulse voltammetry: ret	AUT50477	18-8-2015 15:11:45	☆☆☆☆☆	Add

Figure 224 Hidden columns can be displayed again

NOTICE

It is not possible to hide the *Name* column.

6.11 Filtering the Library

The columns used to display the items in the **Library** can be used for filtering. To filter content of a column, click the button located in the right corner of the column header (see figure 225, page 188).

Demo Database

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M	AUT71848	31-8-2015 11:07:50	13-6-2016 14:59:51	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQC	Pb(ClO4)2 0.01 M / HClO4 0.1 M		4-2-2009 11:04:15	13-6-2016 14:59:54	☆☆☆☆☆	Add
Demo 03 - Bipotentiostat measu	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	13-6-2016 14:59:57	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	13-6-2016 14:59:59	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility	AUT71848	31-8-2015 14:40:14	13-6-2016 14:59:43	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galvan	AUT71848	31-8-2015 11:27:11	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 07 - Chrono measurement	Example of fast options	AUT71848	1-9-2015 13:20:24	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic	Supercapacitor, 3.3 F, CV, differen	AUT71848	1-9-2015 13:29:23	13-6-2016 14:59:37	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impec	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	13-6-2016 14:59:38	☆☆☆☆☆	Add
Demo 10 - Differential pulse mee	Differential pulse voltan	AUT50477	18-8-2015 15:11:45	13-6-2016 14:59:38	☆☆☆☆☆	Add

Figure 225 The columns displayed in the Library can be filtered

When the button is clicked, a menu will appear below the button, providing a list of filters options which can be selected to filter the content of the column based on the specified argument. Four type of filters are available:

- **Alphanumeric filter:** this filter provides the possibility to filter the content of the column based on items that start with a letter or number in the selected bracket(s). This filter is available for the **Name** and **Remarks** columns.

- **Enumeration filter:** this filter provides the possibility to filter the content of the column based on the list of available arguments. *Figure 226* shows an example of an enumeration filter, which displays all the available instrument serial numbers. This filter is available for the **Instrument** and **Tags** columns.
- **Date filter:** this filter provides the possibility to filter the content of the column based on a specific date or timeframe. This type of filter is available for the **Measured date** and **Last modified** columns.
- **Rating filter:** this filter provides the possibility to filter the content of the column based on the assigned rating. This type of filter is available for the **Rating** column.

In the example shown in *Figure 226* a list of instrument serial number is provided.

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 C	AUT71848	7:50	13-6-2016 14:59:51	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQCP	Pb(ClO4)2 0.01 M / HCl	AUT71848	:15	13-6-2016 14:59:54	☆☆☆☆☆	Add
Demo 03 - Bipotentostat measu	RRDE measurement	MAC80064#3	5:21	13-6-2016 14:59:57	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear	Fe2+/Fe3+, NaOH 0.2 N	AUT71848	3:57	13-6-2016 14:59:59	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility	AUT71848	0:14	13-6-2016 14:59:43	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold	AUT71848	31-8-2015 11:27:11	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 07 - Chrono measurement	Example of fast options	AUT71848	1-9-2015 13:20:24	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic	Supercapacitor, 3.3 F, C'	AUT71848	1-9-2015 13:29:23	13-6-2016 14:59:37	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impec	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	13-6-2016 14:59:38	☆☆☆☆☆	Add
Demo 10 - Differential pulse mee	Differential pulse voltan	AUT50477	18-8-2015 15:11:45	13-6-2016 14:59:38	☆☆☆☆☆	Add

Figure 226 It is possible to filter on the instrument serial number



NOTICE

The *Unspecified* filter check box can be used to filter the entries that have no associated value.

Selecting one or more of the available check boxes immediately removes all the entries that do not match the specified filter argument from view, as shown in *Figure 227*.



Demo Database								+
Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags		
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 C	AUT71848	<input checked="" type="checkbox"/> μ3AUT70530	7:50	13-6-2016 14:59:51	☆☆☆☆☆	Add	
Demo 04 - Hydrodynamic linear	Fe2+/Fe3+, NaOH 0.2 N	AUT71848	<input type="checkbox"/> AUT50229	3:57	13-6-2016 14:59:59	☆☆☆☆☆	Add	
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility	AUT71848	<input checked="" type="checkbox"/> AUT71848	0:14	13-6-2016 14:59:43	☆☆☆☆☆	Add	
Demo 06 - Galvanostatic CV	Lead deposition on golk	AUT71848	<input type="checkbox"/> AUT64146	7:11	13-6-2016 14:59:32	☆☆☆☆☆	Add	
Demo 07 - Chrono measurement	Example of fast options	AUT71848	<input type="checkbox"/> AUT85396	24	13-6-2016 14:59:32	☆☆☆☆☆	Add	
Demo 08 - Supercapacitor cyclic	Supercapacitor, 3.3 F, C'	AUT71848	<input type="checkbox"/> MAC80064#3	1-9-2015 13:29:23	13-6-2016 14:59:37	☆☆☆☆☆	Add	
Demo 11 - Hydrodynamic FRA w	With Autolab RDE at 10	μ3AUT70530	<input type="checkbox"/> Unspecified	8-4-2010 12:23:00	13-6-2016 14:59:42	☆☆☆☆☆	Add	
Demo 18 - Galvanic coupling Al/	Galvanic coupling, Al, Al	AUT71848		1-9-2015 15:02:51	13-6-2016 15:00:11	☆☆☆☆☆	Add	
Demo 19 - Cyclic voltammetry Fe	Cyclic voltammetry pote	AUT71848		25-11-2015 16:14:06	13-6-2016 15:00:13	☆☆☆☆☆	Add	

Figure 227 Applying the filter

When a column has a filter active, the symbol will be shown in the right-hand corner of the column header on which the filter is applied (see figure 228, page 190).

Demo Database								+
Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags		
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 C	AUT71848	31-8-2015 11:07:50	13-6-2016 14:59:51	☆☆☆☆☆	Add		
Demo 04 - Hydrodynamic linear	Fe2+/Fe3+, NaOH 0.2 N	AUT71848	31-8-2015 13:53:57	13-6-2016 14:59:59	☆☆☆☆☆	Add		
Demo 05 - Fe(II) - Fe (III) on pcPt	Fe2+/Fe3+ Reversibility	AUT71848	31-8-2015 14:40:14	13-6-2016 14:59:43	☆☆☆☆☆	Add		
Demo 06 - Galvanostatic CV	Lead deposition on golk	AUT71848	31-8-2015 11:27:11	13-6-2016 14:59:32	☆☆☆☆☆	Add		
Demo 07 - Chrono measurement	Example of fast options	AUT71848	1-9-2015 13:20:24	13-6-2016 14:59:32	☆☆☆☆☆	Add		
Demo 08 - Supercapacitor cyclic	Supercapacitor, 3.3 F, C'	AUT71848	1-9-2015 13:29:23	13-6-2016 14:59:37	☆☆☆☆☆	Add		
Demo 11 - Hydrodynamic FRA w	With Autolab RDE at 10	μ3AUT70530	8-4-2010 12:23:00	13-6-2016 14:59:42	☆☆☆☆☆	Add		
Demo 18 - Galvanic coupling Al/	Galvanic coupling, Al, Al	AUT71848	1-9-2015 15:02:51	13-6-2016 15:00:11	☆☆☆☆☆	Add		
Demo 19 - Cyclic voltammetry Fe	Cyclic voltammetry pote	AUT71848	25-11-2015 16:14:06	13-6-2016 15:00:13	☆☆☆☆☆	Add		

Figure 228 A filtered view of the location

It is possible to adjust the filter at any time by repeating the process described above. Each time a check box is either ticked or unticked, the information displayed in the **Library** will be automatically updated (see figure 229, page 190).

Demo Database								+
Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags		
Demo 10 - Differential pulse me	Differential pulse voltan	AUT50477	<input type="checkbox"/> μ3AUT70530	1:45	13-6-2016 14:59:36	☆☆☆☆☆	Add	
Demo 16 - Non linear FRA meas	Example of non-linear n	AUT50477	<input checked="" type="checkbox"/> AUT50477	17	13-6-2016 15:00:26	☆☆☆☆☆	Add	

Figure 229 Adjusting the filter

If needed, additional filters can be applied. In that case, the content of the **Library** is adjusted in order to only display the items that match all the filter conditions, as shown in *Figure 230*.




Figure 230 Adding additional filters



NOTICE

The specified filter(s) only apply to the active **Location** in the **Library**. For each **Location**, unique filter can be specified.



NOTICE

The specified filter(s) remain active until they are cleared or until **NOVA** is closed. To clear an active filter, uncheck all the check boxes in use by this filter.

6.12 Sorting the Library

The columns used to display the items in the **Library** can be used for sorting. To sort the data, click the column header. Clicking the header again toggles from ascending sorting to descending sorting (*see figure 231, page 192*).



Chrono methods

Name ▲	Remarks
Chrono amperometry ($\Delta t > 1$ ms)	Chrono amperometry ($\Delta t > 1$ ms)
Chrono amperometry fast	Chrono amperometry fast
Chrono amperometry high speed	Chrono amperometry high speed: requires ADC10M or ADC750 module
Chrono charge discharge	Chrono charge discharge
Chrono coulometry ($\Delta t > 1$ ms)	Chrono coulometry ($\Delta t > 1$ ms): requires Integrator module
Chrono coulometry fast	Chrono coulometry fast: requires Integrator module
Chrono potentiometry ($\Delta t > 1$ ms)	Chrono potentiometry ($\Delta t > 1$ ms)
Chrono potentiometry fast	Chrono potentiometry fast
Chrono potentiometry high speed	Chrono potentiometry high speed: requires ADC10M or ADC750 module

Figure 231 Sorting the columns in the Library



NOTICE

It is possible to sort the Library content using multiple columns by holding the **[SHIFT]** key and clicking the column headers.

6.13 Rearranging Library columns order

If necessary, it is possible to arrange the columns shown in the **Library** in whichever order necessary. To move a column in the **Library**, click the column header and while holding the mouse button, slide the column left or right in the **Library** panel (see figure 232, page 192).

Demo Database

Name ▲	Remark	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper	CuSO4 0.01 M, H2SO4 0.1 M, A	AUT71848	31-8-2015 11:07:50	11-2-2016 14:31:16	☆☆☆☆☆	Add
Demo 02 - Lead de	Pb(ClO4)2 0.01 M / HClO4 0.1 I		4-2-2009 11:04:15	11-2-2016 14:29:57	☆☆☆☆☆	Add
Demo 03 - Bipoten	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	11-2-2016 14:30:28	☆☆☆☆☆	Add
Demo 04 - Hydrod	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	11-2-2016 14:31:29	☆☆☆☆☆	Add
Demo 05 - Fe(II) - I	Fe2+/Fe3+ Reversibility Test - I	AUT71848	31-8-2015 14:40:14	11-2-2016 14:31:34	☆☆☆☆☆	Add
Demo 06 - Galvanc	Lead deposition on gold, galva	AUT71848	31-8-2015 11:27:11	11-2-2016 14:31:20	☆☆☆☆☆	Add
Demo 07 - Chrono	Example of fast options measu	AUT71848	1-9-2015 13:20:24	11-2-2016 14:31:42	☆☆☆☆☆	Add
Demo 08 - Superc	Supercapacitor, 3.3 F, CV, differ	AUT71848	1-9-2015 13:29:23	11-2-2016 14:31:49	☆☆☆☆☆	Add
Demo 09 - Superc	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	11-2-2016 14:31:55	☆☆☆☆☆	Add
Demo 10 - Differer	Differential pulse voltammetry:	AUT50477	18-8-2015 15:11:45	11-2-2016 14:31:11	☆☆☆☆☆	Add

Figure 232 Arranging the column order in the Library panel

Release the mouse button when the column is relocated.



NOTICE

The column order can be defined for the Default procedures, Procedures, Data and Schedules locations independently. The order will be used by all locations of the same type.

6.14 Locating files

The **Library** provides the option to quickly locating a file on the computer. Right-clicking an item in the **Library** displays a context menu that provides the choice to *Show in Windows Explorer*, as shown in *Figure 233*.

Demo Database							+ -	
Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags		
Demo 01 - Copper	CuSO4 0.01 M, H2SO4 0.1 M, A	AUT71848	31-8-2015 11:07:50	11-2-2016 14:31:16	☆☆☆☆☆	Add		
Demo 02 - Lead de	Pb(ClO4)2 0.01 M / HClO4 0.1 I		4-2-2009 11:04:15	11-2-2016 14:29:57	☆☆☆☆☆	Add		
Demo 03 - Bipo		MAC80064#3	15-7-2013 13:45:21	11-2-2016 14:30:28	☆☆☆☆☆	Add		
Demo 04 - Hydrod	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	11-2-2016 14:31:29	☆☆☆☆☆	Add		
Demo 05 - Fe(II) - f	Fe2+/Fe3+ Reversibility Test - l	AUT71848	31-8-2015 14:40:14	11-2-2016 14:31:34	☆☆☆☆☆	Add		
Demo 06 - Galvanc	Lead deposition on gold, galva	AUT71848	31-8-2015 11:27:11	11-2-2016 14:31:20	☆☆☆☆☆	Add		
Demo 07 - Chrono	Example of fast options messu	AUT71848	1-9-2015 13:20:24	11-2-2016 14:31:42	☆☆☆☆☆	Add		
Demo 08 - Supercz	Supercapacitor, 3.3 F, CV, differ	AUT71848	1-9-2015 13:29:23	11-2-2016 14:31:49	☆☆☆☆☆	Add		
Demo 09 - Supercz	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	11-2-2016 14:31:55	☆☆☆☆☆	Add		
Demo 10 - Differer	Differential pulse voltammetry:	AUT50477	18-8-2015 15:11:45	11-2-2016 14:31:11	☆☆☆☆☆	Add		

Figure 233 The Show in Windows Explorer option can be used to find a file on the computer

Using this option, a **Windows Explorer** window will be opened, showing the location of the file matching the selected item (see *figure 234*, page 194).




Figure 234 The selected file is shown in Windows Explorer

6.15 Delete files from Library

Through the **Library** interface, it is possible to delete one or more files from a location.

To delete a file from the active location, click the  button located in the top right corner of the right-hand side panel (see figure 235, page 194).




Figure 235 Deleting a file from the Library

A confirmation message is displayed before the file is deleted. Clicking the **Yes** button deletes the file. Clicking the **No** button cancels the delete action (see figure 236, page 195).

Remove library item

Are you sure you want to remove "Demo 01 - Copper deposition" from the library?



Figure 236 A confirmation is required to delete the file from the Library



CAUTION

Deleting a file from the **Library** also deletes the source file from the computer. The file is moved to the **Recycle Bin** and if needed, it can be restored (if possible).

6.16 The data repository

The **Library** provides access to a data repository. With the repository, it is possible to create one or more internal backups of a data item in the **Library**. This makes it possible to make one or more backups of the original data and recover the original data from one of backups, if required.



NOTICE

The repository can only be used for data files.

To store data in the repository, right-click the corresponding item in the **Library** and choose the *Store in repository* option from the context menu (see figure 237, page 195).

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO ₄ 0.01 M, H ₂ SO ₄ 0.1 M, A	AUT71848	31-8-2015 11:07:50	13-6-2016 14:59:51	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQC	Pb(ClO ₄) ₂ 0.01 M / HClO ₄ 0.1 l		4-2-2009 11:04:15	13-6-2016 14:59:54	☆☆☆☆☆	Add
Demo 03 - Bipotentiostat		MAC80064#3	15-7-2013 13:45:21	13-6-2016 14:59:57	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear	Fe ²⁺ /Fe ³⁺ , NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	13-6-2016 14:59:59	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcP	Fe ²⁺ /Fe ³⁺ Reversibility Test - l	AUT71848	31-8-2015 14:40:14	13-6-2016 14:59:43	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galva	AUT71848	31-8-2015 11:27:11	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 07 - Chrono measuremer	Example of fast options messu	AUT71848	1-9-2015 13:20:24	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic	Supercapacitor, 3.3 F, CV, differ	AUT71848	1-9-2015 13:29:23	13-6-2016 14:59:37	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impe	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	13-6-2016 14:59:38	☆☆☆☆☆	Add
Demo 10 - Differential pulse me	Differential pulse voltammetry:	AUT50477	18-8-2015 15:11:45	13-6-2016 14:59:38	☆☆☆☆☆	Add

Figure 237 Storing data in the repository

Demo Database + -

Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M, A	AUT71848	31-8-2015 11:07:50	13-6-2016 14:59:51	☆☆☆☆☆	Add
Demo 02 - Lead dep	...	HClO4 0.1 I	4-2-2009 11:04:15	13-6-2016 14:59:54	☆☆☆☆☆	Add
Demo 03 - Bipotent	Store in repository	MAC80064#3	15-7-2013 13:45:21	13-6-2016 14:59:57	☆☆☆☆☆	Add
Demo 04 - Hydrody	Revert from repository	#8	31-8-2015 13:53:57	13-6-2016 14:59:59	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcP	Fe2+/Fe3+ Reversibility Test - I	AUT71848	31-8-2015 14:40:14	13-6-2016 14:59:43	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galva	AUT71848	31-8-2015 11:27:11	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 07 - Chrono measuremer	Example of fast options messu	AUT71848	1-9-2015 13:20:24	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic	Supercapacitor, 3.3 F, CV, differ	AUT71848	1-9-2015 13:29:23	13-6-2016 14:59:37	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impe	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	13-6-2016 14:59:38	☆☆☆☆☆	Add
Demo 10 - Differential pulse me	Differential pulse voltammetry:	AUT50477	18-8-2015 15:11:45	13-6-2016 14:59:38	☆☆☆☆☆	Add

Figure 239 Deleting a repository backup



NOTICE

Deleting a repository backup does not remove the source data from **Library**.

6.17 Merge data

An advanced feature of the **Library** provides the means of merging items. When **Library** items are merged, a new item containing the procedures and the data from the merged items will be copied to the new **Library** item. This can be used to involve the data from two or more different measurements in a calculation or other data handling steps. This option also provides the means to merge different procedures into a single one.



NOTICE

It is only possible to merge items located in the same **Library** location.

To merge two or more items, select them by holding the **[CTRL]** key and clicking the items to merge (see figure 240, page 198).



Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M, A	AUT71848	31-8-2015 11:07:50	13-6-2016 14:59:51	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQC	Pb(ClO4)2 0.01 M / HClO4 0.1 I		4-2-2009 11:04:15	13-6-2016 14:59:54	☆☆☆☆☆	Add
Demo 03 - Bipotentiostat measu	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	13-6-2016 14:59:57	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	13-6-2016 14:59:59	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcP	Fe2+/Fe3+ Reversibility Test - I	AUT71848	31-8-2015 14:40:14	13-6-2016 14:59:43	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galva	AUT71848	31-8-2015 11:27:11	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 07 - Chrono measurer	Example of fast options measu	AUT71848	1-9-2015 13:20:24	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic	Supercapacitor, 3.3 F, CV, differ	AUT71848	1-9-2015 13:29:23	13-6-2016 14:59:37	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impe	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	13-6-2016 14:59:38	☆☆☆☆☆	Add
Demo 10 - Differential pulse me	Differential pulse voltammetry:	AUT50477	18-8-2015 15:11:45	13-6-2016 14:59:38	☆☆☆☆☆	Add

Figure 240 Select two or more items

With two or more items selected, click the  button in the top right corner (see figure 241, page 198).



Name	Remarks	Instrument	Measurement date	Last modified	Rating	Tags
Demo 01 - Copper deposition	CuSO4 0.01 M, H2SO4 0.1 M, A	AUT71848	31-8-2015 11:07:50	13-6-2016 14:59:51	☆☆☆☆☆	Add
Demo 02 - Lead deposition EQC	Pb(ClO4)2 0.01 M / HClO4 0.1 I		4-2-2009 11:04:15	13-6-2016 14:59:54	☆☆☆☆☆	Add
Demo 03 - Bipotentiostat measu	RRDE measurement	MAC80064#3	15-7-2013 13:45:21	13-6-2016 14:59:57	☆☆☆☆☆	Add
Demo 04 - Hydrodynamic linear	Fe2+/Fe3+, NaOH 0.2 M	AUT71848	31-8-2015 13:53:57	13-6-2016 14:59:59	☆☆☆☆☆	Add
Demo 05 - Fe(II) - Fe (III) on pcP	Fe2+/Fe3+ Reversibility Test - I	AUT71848	31-8-2015 14:40:14	13-6-2016 14:59:43	☆☆☆☆☆	Add
Demo 06 - Galvanostatic CV	Lead deposition on gold, galva	AUT71848	31-8-2015 11:27:11	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 07 - Chrono measurer	Example of fast options measu	AUT71848	1-9-2015 13:20:24	13-6-2016 14:59:32	☆☆☆☆☆	Add
Demo 08 - Supercapacitor cyclic	Supercapacitor, 3.3 F, CV, differ	AUT71848	1-9-2015 13:29:23	13-6-2016 14:59:37	☆☆☆☆☆	Add
Demo 09 - Supercapacitor impe	Supercapacitor, 3.3 F	AUT50229	1-9-2015 13:50:29	13-6-2016 14:59:38	☆☆☆☆☆	Add
Demo 10 - Differential pulse me	Differential pulse voltammetry:	AUT50477	18-8-2015 15:11:45	13-6-2016 14:59:38	☆☆☆☆☆	Add

Figure 241 Merging the selected items



NOTICE

The  button is only visible when two or items are selected in the **Library**.

A message will be displayed, showing the following information (see figure 242, page 199):

- **Name:** the name of the merged item. By default, NOVA will generate the [MERGED] 'Name of the first selected item' automatically as the name of the merged item.
- **Remarks:** the remarks for the merged item. By default, NOVA fills this input field with the remarks of all the selected items.




Figure 242 Default name and remarks are generated automatically

It is possible to specify a Name and Remarks for the merged file (see figure 243, page 199).




Figure 243 Specifying name and remarks

Click the **OK** button to merge the items. A new **Library** item will be added to the current location (see figure 244, page 199).




Figure 244 The merged item is added to the Library




Figure 246 Using the search function

After specifying the search string, press the **[Enter]** key to trigger the search. All Procedure, Data or Schedule items matching the search criteria will be displayed on a separate tab (see figure 247, page 201).




Figure 247 The procedures, data and schedule items matching the specified search string are shown in a dedicated tab

Figure 246 and Figure 247 illustrate how the search function can be used to find all items that contain the word **voltammetry** in the **Name** or **Remarks**.



NOTICE

The results are grouped by type in the results tab.



If needed, wildcards (*) can be used. A wildcard indicates that any word can replace it in the search string. In the example shown in *Figure 248*, two wildcards are used: *** in * with EtOH**. This search string format indicates any word can be replace the * when the search is executed.




Figure 248 Wildcards can be used in the search field

The results of the search string used in *Figure 248* are shown in *Figure 249*.




Figure 249 The results are shown in the dedicated tab

7 NOVA commands

NOVA is provided with an extensive set of commands which can be used to modify or create procedures. These commands can be arranged in sequence in order to match the experimental requirements. All the commands provided in NOVA are grouped into different sections:

- **Control:** this group contains commands for user interaction, flow control and external API interfacing. See *Chapter 7.1* for more details.
- **Measurement - general:** this group contains all the commands used to perform basic controls of the instrument. See *Chapter 7.2* for more details.
- **Measurement - cyclic and linear sweep voltammetry:** this group contains all the commands for cyclic and linear sweep voltammetry measurements. See *Chapter 7.3* for more details.
- **Measurement - voltammetric analysis:** this group contains all the commands for voltammetric analysis measurements. See *Chapter 7.4* for more details.
- **Measurement - chrono methods:** this group contains all the commands for time-resolved measurements. See *Chapter 7.5* for more details.
- **Measurement - impedance:** this group contains all the commands for impedance spectroscopy and electrochemical frequency modulation measurements. See *Chapter 7.6* for more details.
- **Data handling:** this group contains all commands designed to process the measured data. See *Chapter 7.7* for more details.
- **Analysis - general:** this group contains all general purpose data analysis commands. See *Chapter 7.8* for more details.
- **Analysis - impedance:** this group contains all the data analysis commands designed for impedance spectroscopy data. See *Chapter 7.9* for more details.
- **Metrohm devices:** this group contains commands that can be used to control supported Metrohm devices connected to the host computer. See *Chapter 7.10* for more details.
- **External devices:** this group contains commands that can be used to control supported external devices. See *Chapter 7.11* for more details.


The details of the command properties of the **Message** command are shown in *Figure 251*:

Figure 251 The properties of the **Message** command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Title:** the title of the message.
- **Message:** the contents of the message.
- **Use time limit:** a toggle provided to switch an automatic time-out of the message on or off.
- **Time limit (s):** the time limit, in s, after which the message is cleared if the *Use time limit* toggle is set to on (default 30 s).
- **Ask for input:** a toggle provided to specify if an input field should be shown in the message.
- **Value:** the default value to show in the input field if the *Ask for input* toggle is set to on.

7.1.2 Send email

	<p>This command can be used to send an email to the specified recipient during a procedure.</p>
---	---

The details of the command properties of the **Send email** command are shown in *Figure 252*:



Figure 253 Three modes are provided by the Repeat command

1. Repeat n times (default mode)
2. Repeat for multiple values
3. Timed repeat



NOTICE

The **Repeat** command description in the procedure editor is dynamically adjusted in function of the specified mode.

7.1.3.1 Repeat n times

The following properties are available when the command is used in the *Repeat n times* mode (see figure 254, page 207):




Figure 254 Repeat n times mode properties

- **Command name:** a user-defined name for the command.
- **Number of repetitions:** the number of repetitions in the repeat loop (default 10).

← Repeat for multiple values

Values ✕ ✎ +

Parameter 1

Add range

Begin value

End value

Number of values

Distribution ▼

Add range

Import Table

Column delimiter ▼

Decimal separator ▼

Rows to skip

Import Table

Figure 256 Overview of the Repeat for multiple values mode panel

The **Repeat for multiple values** panel can be used to build a table of values, in one or more columns. During the experiment, the repeat loop will cycle through each row of the table, and use the values of each column in the measurement.



NOTICE

All the columns specified in the table must have the same number of elements.

Repeat for multiple values

Values ✕ ✎ +

Scan rate (V/s)

Add range

Begin value

End value

Number of values

Distribution

Add range

Figure 258 Editing the column header



NOTICE

If the table contains more than one column, pressing the **[Tab]** key will validate the name of the selected column header and the column header of the next column will be selected for editing. Pressing the combination of **[Shift]** and **[Tab]** will do the same with the previous column header.



NOTICE

If a cell is selected when the  button is clicked, the column header of column containing the selected cell will be highlighted.

7.1.3.2.2 Manually add values to a table

To manually add values to a table, select the first available cell in a column and type the value (see figure 259, page 212).

The following properties are available:

- **Begin value:** the first value of the range.
- **End value:** the last value of the range.
- **Number of values/Number of values per decade:** the number of values in the range or the number of points per decade in the range.
- **Distribution:** the distribution used to calculate the range. Four distributions are available, selectable using the provided drop-down list:
 - **Linear:** the range is built using a linear distribution.
 - **Square root:** the range is built using a square root distribution.
 - **Logarithmic:** the range is built using a logarithmic distribution.
 - **Points per decade:** the range is built by calculating the number of decades in the range and by adding the specified number of points per calculated decade. This distribution is also logarithmic.

Table 6 provides an overview of the formulae used to calculate the distributions supported by the **Add range** option.

Table 6 The distributions used in the Add range option

Type	Increment, Δ	Distribution
Linear	$\Delta = \frac{\text{End} - \text{Start}}{N - 1}$	$V_i = (\text{Start} + (i - 1)\Delta)$
Square root	$\Delta = \frac{\sqrt{\text{End}} - \sqrt{\text{Start}}}{N - 1}$	$V_i = (\sqrt{\text{Start}} + (i - 1)\Delta)^2$
Logarithmic	$\frac{\text{LOG}(\text{End}) - \text{LOG}(\text{Start})}{N - 1}$	$\text{LOG}(V_i) = (\text{LOG}(\text{Start}) + (i - 1)\Delta)$
Points per decade	$\Delta = 10$	$V_i = (\text{Start} \cdot \Delta^{(i-1)})$

Click the **Add range** button to add a range of value to the table (see figure 261, page 214).

Figure 262 The **Import table** option can be used to import a table from a CSV file

The following properties are available:

- **Import table at runtime:** when this toggle is enabled, the table will be imported during the running procedure at the moment that the **Repeat for multiple values** command is executed. This toggle is disabled by default. Importing the table during a running measurement is recommended only for advanced NOVA users. Please read the section below entitled **Import table at runtime** for more information.
- **Filename:** name of the CSV file to be imported. This field will automatically be populated when the **Import Table** button is used to access the CSV file.
- **Column delimiter:** a drop-down list that provides the choice of column delimiter (Space, Tab, Comma (,), Semicolon (;) or Colon (:)).
- **Decimal separator:** a drop-down list that provides the choice of decimal separator (Dot (.) or Comma (,)).
- **Number of rows to skip:** defines the number of rows to skip when importing the data.

The **Import table** **Import Table** button is provided to access the CSV file through a Windows Explorer dialog.

3. Configure the file import properties according to the properties of the CSV file. Pay close attention to the **Column delimiter** and the **Decimal separator** properties; these must match the properties of the CSV file or the table may be imported with unexpected formatting. Testing the table import in a separate test procedure is highly recommended.
4. Set up the relevant links within the procedure. Verifying that the links are configured correctly is recommended before starting the procedure.
5. Start the procedure.



NOTICE

It is **not possible** to add columns to the table while the measurement is running. Columns must not be added to the CSV file after the table has been configured in the Repeat for multiple values command properties.



NOTICE

It **is possible** to add rows to the table while the measurement is running. When rows are added to the CSV file and the changes are saved before the Repeat for multiple values command tile is reached in the procedure sequence, the modified table will be repeated and the new rows will be handled as new Repeats. The number of Repeats will be increased according to the new number of rows.



NOTICE

When the table is imported during the running measurement **validation cannot be performed** on the contents of the table or on properties linked to the contents of the table. The user must be sure that the properties are within the capabilities of their instrument and of the NOVA software. Otherwise, an exception may occur.

7.1.3.2.5 Add additional columns

To add extra columns to the table, click the **+** button above the table. A drop-down list will be displayed offering a choice between two options (see figure 263, page 218):

- **A value column:** a column that contains only numbers.



NOTICE

A text column is identified by an uppercase T in the column header, as shown in *Figure 264*.

The new column can now be edited (*see figure 265, page 219*).

← Repeat for multiple values

Values ✕ ✎ | - +

	Scan rate (V/s)	Parameter 2 T
1	0,1	Scan rate 1
2	0,2	Scan rate 2
3	0,3	Scan rate 3
4	0,4	Scan rate 4
5	0,5	Scan rate 5

Add range

Begin value

End value

Number of values

Distribution

Figure 265 Editing the new column



NOTICE

The Add range option cannot be used when editing a *text* column. The cells of this type of column must be edited manually.

7.1.3.2.6 Moving columns in the table

When the table contains two or more columns, it is possible to rearrange the order of the column by clicking a column header and dragging it to another location, while holding the mouse button (*see figure 266, page 220*).




Figure 266 Moving the columns

When the mouse button is released, the column will be repositioned in the table (see figure 267, page 220).




Figure 267 The repositioned column

7.1.3.2.7 Delete values from the table

Values from the table can be removed in two ways:

- A single value from the table can be deleted. This only clears the content of the table cell.
- A complete row from the table can be deleted. This removes the complete row from the table.

To delete a single value from the table, select the cell containing the value and press the **[Delete]** key (see figure 268, page 221).




Figure 268 Deleting a value from the table

To delete a complete row of the table, click the index cell in front of the row to select the row and press the **[Delete]** key to delete the selected row of the table (see figure 269, page 221).




Figure 269 Deleting a complete row from the table

The selected row is completely removed from the table (see figure 270, page 221).




Figure 270 The selected row is removed from the table




Figure 273 The column contents are cleared

7.1.3.2.10 Remove columns from the table

It is possible to remove a column from the table by selecting the column and clicking the **+** button located above the table (see figure 274, page 223).




Figure 274 Removing a column from the table

The selected column is removed (see figure 275, page 223).





Figure 275 The selected column is removed



The **Increment** command can be used in two different modes, which can be selected using the provided drop-down list (see figure 277, page 225):




Figure 277 Two modes are provided by the Increment command

1. Increment with Value (default mode)
2. Increment with Signal



NOTICE

The **Increment** command description in the procedure editor is dynamically adjusted in function of the specified mode.

7.1.4.1 Increment with Value

The following properties are available when the command is used in the *Increment with Value* mode (see figure 278, page 225):




Figure 278 *Increment with Value* mode property

- **Command name:** a user-defined name for the command.




Figure 280 The signal used in the Increment with Signal mode

7.1.5 Play sound

 <p>Play sound</p>	<p>This command can be used to play a sound using a system or user defined source.</p>
---	--

The details of the command properties of the **Play sound** command are shown in *Figure 281*:

- **Command name:** a user-defined name for the command.
- **Format:** a string containing one or more *format items* (**{x}**, where **x** is an integer value, starting at 0). Each specified format item is linkable to another command property.



NOTICE

The numbers used to specify the format items in the **Build text** command must be unique and must be sequential, starting at 0.

7.1.7 .NET

 <p>.NET Call a static method</p>	<p>This command can be used to call functionality provided by a .NET API (Application Programming Interface).</p>
---	---



CAUTION

This command is intended for advanced users. The use of this command falls outside of the scope of this manual.

The **.NET** command can be used in six different modes, which can be selected using the provided drop-down list (*see figure 283, page 229*):




Figure 283 Six modes are provided by the .NET command

1. Call a static method (default mode)



NOTICE

The Class and Method fields are populated as soon as a valid Assembly file is specified.

7.1.7.2 Get or set a static field

The following properties are available when the command is used in the *Get or set a static field* mode (see figure 285, page 231):

The screenshot shows a 'Properties' dialog box with the title 'Get or set a static field'. It contains the following fields:

- Command name:** A text input field containing 'Get or set a static fiel'.
- Action:** A dropdown menu with 'Get or set a static fiel' selected.
- Assembly file:** A text input field with a 'Browse' button to its right.
- Class:** A text input field.
- Field:** A text input field.
- Direction:** A dropdown menu with 'Get' selected.

Figure 285 Get or set a static field mode properties

- **Command name:** a user-defined name for the command.
- **Assembly file:** specifies the path to the assembly file containing the functionality to call. A **Browse** button is provided to locate the file.
- **Class:** the class provided in the assembly file.
- **Field:** the field to get or set.
- **Direction (Get/Set):** a drop-down list that can be set to Get or Set.

7.1.7.3 Create object

The following properties are available when the command is used in the *Create an object* mode (see figure 286, page 232):

7.1.7.5 Call a method

The following properties are available when the command is used in the *Call a method* mode (see figure 288, page 233):

Figure 288 *Call a method* mode properties

- **Command name:** a user-defined name for the command.
- **Method:** the method to call.

7.1.7.6 Get or set a field

The following properties are available when the command is used in the *Get or set a field* mode (see figure 289, page 233):

Figure 289 *Get or set a field* mode properties

- **Command name:** a user-defined name for the command.
- **Field:** the field to get or set.
- **Direction (Get/Set):** a drop-down list that can be set to Get or Set.

7.1.8 Procedure Information

	<p>This command can be used to input information about the procedure, add remarks, and access the instrument serial number as linkable properties</p>
--	---

7.2 Measurement - general

Commands located in the **Measurement – general** group can be used to control the settings of the instrument during a procedure without performing an actual measurement.

The available commands are represented by a shortcut icon (see figure 292, page 235).




Figure 292 The Measurement - general commands

The following commands are available:

- **Autolab control:** a command which can be used to define the instrumental settings of the Autolab potentiostat/galvanostat and its optional modules *Autolab control* (see chapter 7.2.1, page 236)
- **Apply:** a command which can be used to define the applied potential or current *Apply* (see chapter 7.2.2, page 239)
- **Cell:** a command which can be used switch the cell on or off *Cell* (see chapter 7.2.3, page 240)
- **Wait:** a command which can be used to force the procedure to wait *Wait* (see chapter 7.2.4, page 240)
- **OCP measurement:** a command which can be used to measure the open circuit potential *OCP* (see chapter 7.2.5, page 246)
- **Set pH measurement temperature:** a command which can be used to set the pH measurement temperature *Set pH measurement temperature* (see chapter 7.2.6, page 247)
- **Reset EQCM delta frequency:** a command which can be used to reset the Δ Frequency signal measured by the EQCM module *Reset EQCM delta frequency* (see chapter 7.2.7, page 248)
- **Control Autolab R(R)DE:** a command used to control the rotation rate of the Autolab rotating disk electrode (RDE) or rotating ring disk electrode (RRDE)
- **MDE control:** a command used to control a mercury drop electrode stand connected to the Autolab using the IME663 or the IME303 module *MDE control* (see chapter 7.2.9, page 251)
- **Multi Autolab synchronization:** a command which can be used to create a synchronization point in a procedure for multi Autolab measurements *Synchronization* (see chapter 7.2.10, page 254)



7.2.1 Autolab control

	<p>This command can be used specify the hardware configuration of the instrument during a measurement. All the instrumental settings are configured using a dedicated dialog. The available settings depend on the hardware configuration.</p>
---	--

The details of the command properties of the **Autolab control** command are shown in *Figure 293*.




Figure 293 The properties of the Autolab control command

The following properties are available:

- **Command name:** a user-defined name for the command.

The More button can be used to edit the instrument settings for the **Autolab control** command. The Autolab control screen will be displayed (see *figure 294, page 236*).




Figure 294 The Autolab control editor



NOTICE

The settings provided in the **Autolab control** command dialog depend on the hardware setup. More details about the available settings are provided in the hardware description chapters located at the end of this document *Hardware description* (see chapter 16, page 873).

Whenever a setting, available in the Autolab control screen, is adjusted, this setting will be made available in the **Properties** panel (see figure 295, page 237).



Figure 295 Modified settings are automatically made available in the Properties panel

These settings can be directly edited in the **Properties** panel, without the need of opening the Autolab control screen. A button will be added to any modified setting in order to undo this modification and remove this setting from the **Properties** panel (see figure 296, page 238).




Figure 296 Modified settings can be directly edited or removed from the Properties panel

It is possible to add additional settings from the Autolab control screen to the **Properties** panel without modifying them, by clicking the  button (see figure 297, page 238).




Figure 297 Additional settings can be added to the Properties panel

These additional settings will become visible in the **Properties** panel (see figure 298, page 239).





Figure 298 Additional settings will be visible in the Properties panel

7.2.2 Apply



The details of the properties of the **Apply** command are shown in Figure 299:




Figure 299 The properties of the Apply command

The following properties are available:


- **Command name:** a user-defined name for the command.
- **Potential/Current:** the applied potential or current value, in V or A respectively.
- **With respect to drop-down list:** a drop-down list that provides the choice of the reference used to apply a potential value (only shown in potentiostatic mode). The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE

The **Apply** command description in the procedure editor is dynamically adjusted in function of the specified value.

7.2.3 Cell



This command can be used to switch the cell off or on.

The details of the command properties of the **Cell** command are shown in *Figure 300*:




Figure 300 The properties of the Cell command

The following properties are available:


- **Command name:** a user-defined name for the command.
- **Switch cell:** a  toggle control provided to switch the cell off or on.




NOTICE

The **Cell** command description in the procedure editor is dynamically adjusted in function of the toggle.

7.2.4 Wait

	<p>This command can be used to force the procedure to wait for a predefined amount of time (in seconds) or until a certain trigger signal is recorded a DIO port of the instrument or the input lines of a Metrohm 6.2148.010 Remote Box.</p>
---	---

The **Wait** command can be used in four different modes, which can be selected using the provided drop-down list (see figure 301, page 241):



The screenshot shows a 'Properties' window for a 'Wait 5 s' command. It includes a 'Command name' field with 'Wait 5 s' and a 'Duration' field with '5 s'. A 'Wait for' dropdown menu is open, displaying four options: 'Seconds', 'DIO', 'Remote inputs', and 'Metrohm device'. A mouse cursor is hovering over the 'DIO' option.

Figure 301 Four modes are provided by the Wait command

1. Wait for Seconds (default mode)
2. Wait for DIO
3. Wait for Remote inputs
4. Wait for Metrohm device



NOTICE

The **Wait** command description in the procedure editor is dynamically adjusted in function of the specified mode.

7.2.4.1 Wait for Seconds

The following properties are available when the **Wait** command is used in the *Wait for Seconds* mode (see figure 302, page 242):

- **Mask:** the trigger mask, specified as 8 or 4 bits. The expected bit sequence must be formatted using 1, 0 and X (1 indicates that the pin status must be 'high', 0 indicates that the pin status must be 'low', X indicates that the pin status may be both). In *Figure 303*, the command will force the procedure to wait until pins 8, 7, 2 and 1 are 'high', pins 6 and 5 are 'low'. The status of pins 4 and 3 is irrelevant (X status).
- **Use time limit:** a toggle is provided to enable or disable the time limit. When this option is enabled, the *Wait for DIO* command will stop waiting after the specified amount of time.
- **Time limit:** the time limit after which the command stops waiting (if the *Use time limit* toggle is on), in s.

7.2.4.3 Wait for Remote inputs



CAUTION

This mode requires a **Metrohm 6.2148.010 Remote Box** to be connected to the computer.

The following properties are available when the **Wait** command is used in the *Wait for Remote inputs* mode (see *figure 304*, page 243):

The screenshot shows a 'Properties' dialog box with the following fields:

- Command name: Wait for Remote input
- Wait for: Remote inputs (dropdown menu)
- Device name: DigitalIO_1
- Mask: 11001100
- Use time limit:

Figure 304 The properties of the *Wait for Remote inputs* mode

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Remote Box.
- **Mask:** the trigger mask, specified as byte. The expected byte must be formatted using 1 and 0 (1 indicates that the pin status must be 'high' and 0 indicates that the pin status must be 'low'). In *Figure 304*, the command will force the procedure to wait until IN7, IN6, IN3 and IN2 are set to 'high' state and IN5, IN4, IN1 and IN0 are set to 'low' state.
- **Use time limit:** a toggle is provided to enable or disable the time limit. When this option is enabled, the *Wait for Remote inputs* command will stop waiting after the specified amount of time.




Figure 306 Illustration of the parallel execution option for the Metrohm devices

- In *Figure 306, A*: Dosino 1 and Dosino 2 have parallel execution disabled. Both Dosino need to finish the Dose command before the CV staircase command can start.
- In *Figure 306, B*: parallel execution is enabled on Dosino 2 and disabled on Dosino 1. Dosino 2 starts dosing immediately after Dosino 1 is finished. The CV staircase command starts as soon as Dosino 2 starts dosing.
- In *Figure 306, C*: parallel execution is enabled on Dosino 1 and disabled on Dosino 2. Dosino 2 starts dosing at the same time as Dosino 1. Only when Dosino 2 is finished can the CV staircase command start.
- In *Figure 306, D*: parallel execution is enabled for both Dosino 1 and Dosino 2. All three commands start at the same time.

The *Wait for Metrohm device* mode can be used in a procedure to force the procedure to wait until the specified device finishes the command it is executing. This mode can thus be used to overrule the parallel execution of the device.



NOTICE

The *Wait for Metrohm device* mode has no effect on devices for which parallel execution is disabled.

- **Accept on dE/dt limit:** a  toggle control provided to specify if the OCP command should stop when the measured dE/dt value becomes equal or lower than the **dE/dt limit** value.
- **dE/dt limit:** the time derivative limit, in V/s. When this value is not 0, the recording of the OCP will stop when the time derivative of the potential is smaller or equal to the specified limit. This property is only available when the **Accept on dE/dt limit** toggle is enabled.
- **Estimated number of points:** a read-only field indicating the estimated number of points that will be collected. Calculated based on the **Duration** and the **Interval time**
- **Estimated duration:** a read-only field indicating the estimated duration of the measurement.

When the **OCP** command is executed, a two default plots will appear in the **Plots frame**:


- WEPotential vs. Time
- dWEPotential / dt vs. Time



NOTICE

The **OCP** command can be terminated manually using the skip  or stop  buttons. When the measurement is terminated manually, the OCP is determined by selecting the final OCP value, or applying the averaging properties defined by the user with the **Use average OCP toggle** and the **Time to average** property.

7.2.6 Set pH measurement temperature



This command can be used to specify the measurement temperature, for automatic pH correction (if the temperature is not measured through the T input of the **pX1000** module or if the pH is measured using a **pX** module).



CAUTION

This command requires a **pX1000** or **pX** module *pX1000 module* (see chapter 16.3.2.18, page 1162) installed in the Autolab.

This command performs the following mathematical adjustment:

$$\text{pH} = 7 - \frac{a}{b} + \frac{E}{b} \cdot \frac{T_{\text{cal}}}{T}$$



Where a and b are the intercept and the slope of the calibration curve, E is the measured potential, T is the specified temperature and T_{cal} is the calibration temperature.

The details of the properties of the **Set pH measurement temperature** command is shown in *Figure 308*:

Figure 308 The properties of the Set pH measurement temperature command

The following properties are available:


- **Command name:** a user-defined name for the command.
- **Temperature:** the measurement temperature, in °C.



NOTICE

The **Set pH temperature measurement** command description in the procedure editor is dynamically adjusted in function of the specified value.

7.2.7 Reset EQCM delta frequency



This command can be used to create a break-point in the procedure during which the signals from the **EQCM** module can be adjusted and zeroed, if necessary.



CAUTION

This command requires an **EQCM** module *EQCM module* (see chapter 16.3.2.10, page 1075) installed in the Autolab.

The details of the command properties of the **Reset EQCM delta frequency** command are shown in *Figure 309*.




Figure 309 The properties of the Autolab control command

The following properties are available:

- **Command name:** a user-defined name for the command.

When the **Reset EQCM delta frequency** command is executed, a dedicated window will be shown, providing additional control during the measurement (see figure 310, page 249).




Figure 310 The Determine EQCM zero frequency window displayed when the Reset EQCM delta frequency command is executed

The following information is shown in the window:



- **Top section:**
 - **EQCM(1). Δ Frequency:** displays the latest and average values of the measured EQCM(1).Frequency signal are displayed, in V.
 - **EQCM(1).Driving force:** displays the latest, average and minimum values of the measured EQCM(1).Driving force signal are displayed, in V.
 - **EQCM(1).Temperature:** displays the latest and average values of the measured EQCM(1).Temperature signal are displayed, in °C.
 - **Time to average:** a read-only field that indicates the duration, in seconds, used to determine the average values of the EQCM signals.
- **Middle section:** this section displays a real time plot of the EQCM(1).Driving force signal. Using the provided trimmer, it is possible to adjust the driving force as indicated in the EQCM User Manual.
- **Bottom section:**
 - **Zero Δf button:** this button can be used to reset the measured value of the EQCM(1). Δ Frequency signal. When this button is pressed, the value of the signal is recorded and then subtracted from the actual value, thus zeroing the signal. While the EQCM. (1). Δ Frequency signal is set to zero, the button is disabled.
 - **Clear plot button:** this button clears the plot of the EQCM(1).Driving force signal.
 - **Abort button:** this button can be used to force the procedure to stop.
 - **Accept button:** this button can be used to close the Determine EQCM zero frequency window. The procedure will continue using the last value of the EQCM(1). Δ Frequency signal.

7.2.8 Autolab R(R)DE control

	<p>This command can be used to control the Autolab rotating disk electrode (RDE) or rotating ring disk electrode (RRDE), connected to the Autolab and operated in remote control mode.</p>
---	--

The details of the command properties of the **Autolab R(R)DE control** command are shown in *Figure 311*:

Figure 311 The properties of the Autolab R(R)DE control command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Switch R(R)DE:** a  toggle control provided to switch the Autolab RDE or RRDE off or on.
- **Rotation rate:** specifies the rotation rate of the Autolab RDE or RRDE, in RPM.



NOTICE

The **Autolab R(R)DE control** command description in the procedure editor is dynamically adjusted in function of the toggle.

7.2.9 MDE control



This command can be used to control the Mercury Drop Electrode (MDE) using the **IME663** or the **IME303** interface and the **Metrohm 663 VA Stand**, the **Princeton Applied Research PAR303(A) Stand** or a compatible mercury drop electrode stand.



CAUTION

This command requires a **IME663** *IME663 module (see chapter 16.3.2.15, page 1130)* or **IME303** *IME303 module (see chapter 16.3.2.14, page 1124)* connected to the Autolab. When this command is used without a **IME663** or **IME303**, an **error** will be displayed for the command.

The **MDE control** command can be used in three different modes, which can be selected using the provided drop-down list (*see figure 312, page 252*):

Figure 313 Purge mode properties

- **Command name:** a user-defined name for the command.
- **Duration:** specifies the duration during which the N₂ purging functionality provided by the MDE is active, in s.

7.2.9.2 Set Stirrer

The following properties are available when the **MDE control** command is used in the *Set Stirrer* mode (see figure 314, page 253):

Figure 314 Set stirrer mode properties

- **Command name:** a user-defined name for the command.
- **Switch stirrer:** specifies the status of the built-in stirrer of the MDE through the provided  toggle. The specified status remains active until changed.

7.2.9.3 New drop

The following properties are available when the **MDE control** command is used in the *New drop* mode (see figure 315, page 254):



Properties ➔

5 new drops

Command name

Action

Number of new drops

Figure 315 New drop mode properties


- **Command name:** a user-defined name for the command.
- **Number of new drops:** specifies the number of drops to knock off the capillary of the MDE by activating the built-in tapper.



NOTICE

A 500 ms settling time is used each time the tapper is activated. This settling time can be adjusted in the hardware setup.

7.2.10 Synchronization

 <p>Software synchronization</p>	<p>This command can be used to create a synchronization point in the procedure.</p>
---	---

The **Synchronization** command can be used in two different modes, which can be selected using the provided drop-down list (see figure 316, page 254):

Properties ➔

Software synchronization

Command name

Synchronize

Number of instruments

Group name

Use time limit

Figure 316 Two modes are provided by the Synchronization

Two modes are provided by the Synchronization command

1. Software synchronization (default mode)
2. Hardware synchronization



NOTICE

The **Synchronization** command description in the procedure editor is dynamically adjusted in function of the specified mode.

7.2.10.1 Software synchronization

The following properties are available when the **Synchronization** command is used in the *Software synchronization* mode (see figure 317, page 255):

The screenshot shows a 'Properties' dialog box with a right-pointing arrow icon. The title is 'Software synchronization'. It contains the following fields:

- Command name: Software synchronizat
- Synchronize: Software (dropdown menu)
- Number of instruments: 1
- Group name: (empty text box)
- Use time limit: (toggle switch, currently off)

Figure 317 Software synchronization mode properties

- **Command name:** a user-defined name for the command.
- **Number of instruments:** the number Autolab instruments to synchronize. The synchronization is triggered as soon as this number is reached or when the optional time limit has been reached.
- **Group name:** defines a unique name for the synchronization command.
- **Use time limit:** a toggle that can be used to specify if a time limit should be used for the synchronization command.
- **Time limit:** specifies the time limit, in s.
- **Abort after time limit:** a toggle that can be used to specify if the measurement should be aborted if the time limit is reached.

7.3 Measurement - cyclic and linear sweep voltammetry commands

Commands located in the **Measurement – cyclic and linear sweep voltammetry** group can be used to perform programmed potential or current sweep measurements.

The available commands are represented by a shortcut icon (see figure 319, page 257).



Figure 319 The Measurement - cyclic and linear sweep voltammetry commands

The following commands are available:

- **CV staircase:** a command which can be used to perform staircase cyclic voltammetry measurements *CV staircase* (see chapter 7.3.1, page 257).
- **CV linear scan:** a command which can be used to perform linear scan cyclic voltammetry measurements *CV linear scan* (see chapter 7.3.2, page 260). This command requires the **SCAN250** or **SCANGEN** module *SCAN250 module* (see chapter 16.3.2.19, page 1169).
- **LSV staircase:** a command which can be used to perform staircase linear sweep voltammetry measurements *LSV staircase* (see chapter 7.3.3, page 262).

7.3.1 CV staircase

	<p>This command can be used to perform a staircase cyclic voltammetry measurement, in potentiostatic or galvanostatic conditions.</p>
---	---

The details of the properties of the **CV staircase** command are shown in Figure 320:

- **Scan rate:** the scan rate of the potential or current sweep, in V/s or A/s.

Four additional properties are shown as *read-only*:

- **Interval time:** the time interval between two consecutive points in the scan. This property is defined by the potential or current step and the scan rate.
- **Estimated number of points:** the estimated number of points in the scan. This property is defined by the start and stop potential or current and the potential or current step.
- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.
- **Number of stop crossings:** the number of times the potential or current scan will cross the stop potential or current value. This property is defined by the number of scans.



CAUTION

in order to properly identify the scans in the data, it is important to make sure that the following conditions are respected when defining the parameters of the **CV staircase** command:

- The stop value must be smaller than the upper vertex minus the step value.
- The stop value must be larger than the lower vertex plus the step value.



NOTICE

When this command is used in potentiostatic mode, a drop-down list provides the choice of the reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE

The **CV staircase** command provides access to additional options, through the More button *Additional measurement command properties* (see chapter 9, page 615).

Properties ➔

CV linear scan

Command name

Mode ▼

Start potential V_{REF} ▼

Upper vertex potential V_{REF} ▼

Lower vertex potential V_{REF} ▼

Number of scans

Scan rate V/s

Potential interval V

Interval time s

Estimated number of points

Estimated duration s

Number of vertex potential crossings

Figure 321 CV linear scan properties

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Mode:** specifies if the scan is performed in normal mode or high speed mode. This parameter is only shown when the optional **ADC10M** or **ADC750** module is present in the instrument. For high speed potential scans (more than 10 V/s), the high speed mode is recommended.
- **Start potential:** the start potential, in V.
- **Upper vertex potential:** the upper vertex potential, in V.
- **Lower vertex potential:** the lower vertex potential, in V.
- **Stop on:** specifies if the scan should stop on one of the vertices or on the start potential value, using the provided drop-down list. This parameter is only shown when the optional **SCAN250** module is present.
- **Number of scans:** the number of potential scans.
- **Potential interval:** the potential interval between two consecutive data points. The interval can be positive or negative. With a positive interval, the scan starts from the start potential towards the upper vertex potential. With a negative interval, the scan direction is reversed.
- **Scan rate:** the scan rate of the potential sweep, in V/s.



Four additional properties are shown as *read-only*:

- **Interval time:** the time interval between two consecutive points in the scan. This property is defined by the potential interval and the scan rate.
- **Estimated number of points:** the estimated number of points in the scan. This property is defined by the start and stop potential and the potential interval.
- **Estimated duration:** the estimated duration of the command, in s. This property is defined by the estimated number of points and the interval time as well as the duration of the underlying commands, if applicable.
- **Number of vertex/start potential crossings:** the number of times the potential scan will cross one of the potential vertices or the start potential. This property is defined by the number of scans.



NOTICE


For each potential value, a drop-down list provides the choice of the reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE

The **CV linear scan** command provides access to additional options, through the More button *Additional measurement command properties* (see chapter 9, page 615).

7.3.3 LSV staircase



This command can be used to perform a staircase linear sweep voltammetry measurement, in potentiostatic or galvanostatic conditions.

The details of the properties of the **LSV staircase** command are shown in *Figure 322*:

Properties

LSV staircase

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Scan rate V/s

Step V

Interval time s

Estimated number of points

Estimated duration s

[More](#)

Figure 322 LSV staircase properties

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Start potential/current:** the start potential or current value, in V or A respectively.
- **Stop potential/current:** the stop potential or current value, in V or A respectively.
- **Scan rate:** the scan rate of the potential or current sweep, in V/s or A/s.
- **Step:** the potential or current step, in V or A respectively.

Three additional properties are shown as *read-only*:

- **Interval time:** the time interval between two consecutive points in the scan. This property is defined by the potential or current step and the scan rate.
- **Estimated number of points:** the estimated number of points in the scan. This property is defined by the start and stop potential or current and the potential or current step.
- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.

- **Differential pulse voltammetry:** a command which can be used to perform a differential pulse voltammetry measurement. The differential current is calculated during the measurement *Differential pulse voltammetry* (see chapter 7.4.3, page 270).
- **Differential normal pulse voltammetry:** a command which can be used to perform a differential normal pulse voltammetry measurement. The differential current is calculated during the measurement *Differential normal pulse voltammetry* (see chapter 7.4.4, page 273).
- **Square wave voltammetry:** a command which can be used to perform a square wave voltammetry measurement. The differential current is calculated during the measurement *Square wave voltammetry* (see chapter 7.4.5, page 276).
- **Potentiometric stripping analysis:** a command which can be used to perform chemical and constant current potentiometric stripping analysis *PSA (Potentiometric stripping analysis)* (see chapter 7.4.6, page 279).
- **AC voltammetry:** a command which can be used to perform an AC voltammetry measurement *AC voltammetry* (see chapter 7.4.7, page 283).



CAUTION

The voltammetric analysis commands cannot be used in combination with the **ECD module**.

7.4.1 Sampled DC voltammetry



This command can be used to perform a sampled DC voltammetry measurement, in potentiostatic conditions.



CAUTION

This command can only be used in Potentiostatic mode.

The details of the properties of the **Sampled DC voltammetry** command are shown in *Figure 324*.




Figure 325 Overview of the measurement properties of the Sampled DC voltammetry command



NOTICE

When this command is used, a drop-down list provides the choice of the reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE


The **Sampled DC voltammetry** command provides access to additional options, through the More button *Additional measurement command properties* (see chapter 9, page 615).



NOTICE

The **Stop potential**, **Step** and **Interval time** properties can be modified in real time.

7.4.2 Normal pulse voltammetry



This command can be used to perform a normal pulse voltammetry measurement, in potentiostatic conditions.

- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.
- **Scan rate:** the calculated scan rate, in V/s, determined based on the step potential and the interval time.

Figure 327 represents the measurement properties of the **Normal pulse voltammetry** command, schematically.




Figure 327 Overview of the measurement properties of the Normal pulse voltammetry command



NOTICE

When this command is used, a drop-down list provides the choice of the reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE

The **Normal pulse voltammetry** command provides access to additional options, through the More button *Additional measurement command properties* (see chapter 9, page 615).



NOTICE

The **Stop potential**, **Base potential**, **Step**, **Normal pulse time** and **Interval time** properties can be modified in real time.

- **Modulation amplitude:** the amplitude of the potential modulation, in V.
- **Modulation time:** the duration of the potential modulation, in s.
- **Interval time:** the duration of the interval time, in s.

Three additional properties are shown as *read-only*:

- **Estimated number of points:** the estimated number of points. This property is defined by the start and stop potential and the potential step.
- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.
- **Scan rate:** the calculated scan rate, in V/s, determined based on the step potential and the interval time.



CAUTION

The definition of the **modulation amplitude** property is consistent with the definition used in the Metrohm Computrace and VIVA software packages. When this value is positive, the pulse will be applied in the same direction as the potential scan (positive pulse in the positive going direction and negative pulse in the negative going direction). When this value is negative, the pulse will be applied in the **reverse** direction as the potential scan (negative pulse in the positive going direction and positive pulse in the negative going direction). This definition differs from the definition used in NOVA 1.X. Procedures imported from NOVA 1.X are automatically converted to the new definition.

Figure 329 represents the measurement properties of the **Differential pulse voltammetry** command, schematically.




Figure 329 Overview of the measurement properties of the Differential pulse voltammetry command

In a **Differential pulse voltammetry** measurement, two consecutive current samples are collected for each step. The current value measured in the first part of the step corresponds to the WE(1).Base.Current signal while the current value measured at the end of the pulse corresponds to the WE(1).Pulse.Current signal. The differential value, corresponding to the WE(1). δ .Current signal is given by the difference of the pulse and the base current values



NOTICE

When this command is used, a drop-down list provides the choice of the reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE


The **Differential pulse voltammetry** command provides access to additional options, through the More button *Additional measurement command properties* (see chapter 9, page 615).



NOTICE

The **Stop potential**, **Step**, **Modulation amplitude**, **Modulation time** and **Interval time** properties can be modified in real time.

7.4.4 Differential normal pulse voltammetry



This command can be used to perform a differential normal pulse voltammetry measurement, in potentiostatic conditions.



CAUTION

This command can only be used in Potentiostatic mode.

The details of the properties of the **Differential normal pulse voltammetry** command are shown in *Figure 330*.

Properties ➔

Differential normal pulse

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Base potential V_{REF} ▼

Step V

Modulation amplitude V

Modulation time s

Normal pulse time s

Interval time s

Estimated number of points

Estimated duration s

Scan rate V/s

Figure 330 Differential normal pulse voltammetry properties

The following properties are available:

- **Command name:** a user-defined name for the command.



- **Start potential:** the start potential value, in V.
- **Stop potential:** the stop potential value, in V.
- **Base potential:** the base potential, in V.
- **Step:** the potential step, in V.
- **Modulation amplitude:** the amplitude of the potential modulation, in V.
- **Modulation time:** the duration of the potential modulation, in s.
- **Normal pulse time:** the duration of the normal pulse, in s.
- **Interval time:** the duration of the interval time, in s.



CAUTION

The definition of the **modulation amplitude** property is consistent with the definition used the Metrohm Computrace and VIVA software packages. When this value is positive, the pulse will be applied in the same direction as the potential scan (positive pulse in the positive going direction and negative pulse in the negative going direction). When this value is negative, the pulse will be applied in the **reverse** direction as the potential scan (negative pulse in the positive going direction and positive pulse in the negative going direction). This definition differs from the definition used in NOVA 1.X. Procedures imported from NOVA 1.X are automatically converted to the new definition.

Three additional properties are shown as *read-only*:

- **Estimated number of points:** the estimated number of points. This property is defined by the start and stop potential and the potential step.
- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.
- **Scan rate (V/s):** the calculated scan rate, determined based on the step potential and the interval time.

Figure 331 represents the measurement properties of the **Differential normal pulse voltammetry** command, schematically.




Figure 331 Overview of the measurement properties of the Differential normal pulse voltammetry command



CAUTION

The implementation of Differential normal pulse voltammetry is different from the description provided in *Electrochemistry* by C. M. A. Brett and A. M. Oliveira Brett, Oxford University Press, 1993.

In a **Differential normal pulse voltammetry** measurement, two consecutive current samples are collected for each step. The current value measured in pulse of the step corresponds to the WE(1).Pulse.Current signal while the current value measured at the end of the modulation corresponds to the WE(1).Modulation.Current signal. The difference between the modulation and the pulse current corresponds to the WE(1). δ .Current signal.



NOTICE

When this command is used, a drop-down list provides the choice of the reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE


The **Differential normal pulse voltammetry** command provides access to additional options, through the More button *Additional measurement command properties* (see chapter 9, page 615).



NOTICE

The **Stop potential**, **Base potential**, **Step**, **Modulation amplitude**, **Modulation time**, **Normal pulse time** and **Interval time** properties can be modified in real time.

7.4.5 Square wave voltammetry



This command can be used to perform a square wave voltammetry measurement, in potentiostatic conditions.



CAUTION

This command can only be used in Potentiostatic mode.

The details of the properties of the **Square wave voltammetry** command are shown in *Figure 332*.

Properties ➔

Square wave

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Step V

Modulation amplitude V

Frequency Hz

Estimated number of points

Interval time s

Estimated duration s

Scan rate V/s

More

Figure 332 Square wave voltammetry properties

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Start potential:** the start potential value, in V.
- **Stop potential:** the stop potential value, in V.
- **Step:** the potential step, in V.
- **Modulation amplitude:** the amplitude of the square wave, in V.
- **Frequency:** the frequency of the square wave, in Hz.

Four additional properties are shown as *read-only*:

- **Estimated number of points:** the estimated number of points. This property is defined by the start and stop potential and the potential step.
- **Interval time:** the calculated interval time, based on the value of the frequency.
- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.
- **Scan rate:** the calculated scan rate, in V/s, determined based on the step potential and the interval time.

Figure 333 represents the measurement properties of the **Square wave voltammetry** command, schematically.




Figure 333 Overview of the measurement properties of the Square wave voltammetry command

In a **Square wave voltammetry** measurement, two consecutive current samples are collected for each step. The current value measured in first half of the step corresponds to the WE(1).Forward.Current signal while the current value measured in the second half of the step corresponds to the WE(1).Backward.Current signal. The difference between the backward and forward currents corresponds to the WE(1). δ .Current signal



NOTICE

When this command is used, a drop-down list provides the choice of the reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE


The **Square wave voltammetry** command provides access to additional options, through the More button *Additional measurement command properties* (see chapter 9, page 615).



NOTICE

The **Stop potential**, **Step**, **Amplitude** and **Frequency** properties can be modified in real time.

7.4.6 PSA (Potentiometric stripping analysis)

	<p>This command can be used to perform potentiometric stripping analysis (PSA) measurements.</p>
---	--

The **PSA** command can be used in two different modes, which can be selected using the provided drop-down list (see figure 334, page 279):




Figure 334 Two modes are provided by the PSA command

1. Chemical (default mode)
2. Constant current



CAUTION

No additional signals can be measured by the **PSA** command. The sampling rate is set to the highest possible value during this type of measurement and the measured data cannot be displayed in real time. Options like cutoffs and counters cannot be used.

During a potentiometric stripping analysis measurement, the potential of the working electrode is recorded as a function of time while a chemical or electrochemical oxidation is taking place. As such, this method is the equivalent of a potentiometric titration, in which the titrant is added *in situ* at a constant rate. The measurement stops when the maximum time is reached or when the measured potential exceeds a user defined limit. A typical potentiometric stripping analysis potential profile is shown in Figure 335.




Figure 335 A typical potentiometric stripping analysis measurement

The voltage measurement E versus time is used to calculate the retention times dt/dE vs E . Figure 336 shows an example of the E vs t measurement and the resulting peak-shaped plot.




Figure 336 dt/dE versus potential curve

The peak voltage position is characteristic of the substance, the peak area is proportional to its concentration.

7.4.6.1 Chemical PSA

The following properties are available when the command is used in the *Chemical* mode (see figure 337, page 281):

Properties 

PSA

Command name

Mode

Potential limit V

Maximum time s

Estimated duration s

More

Figure 337 Chemical mode properties

- **Command name:** a user-defined name for the command.
- **Potential limit:** the maximum potential of the working electrode, in V. The measurement stops when the potential of the working electrode exceeds the specified value.
- **Maximum time:** the maximum duration of the stripping stage, in s. The measurement stops when this limit is reached.
- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the maximum time and the duration of underlying commands, if applicable.

Advanced properties are accessed by selecting the **More** button.

- **Filter:** specifies if a filter must be applied on the measured potential signal using the provided  toggle (Default: ON). The implemented filter is based on a moving average over the specified *Filter time* property.
- **Filter time:** the filter time, in s, used if the filter is On. This property is automatically set to 20 ms or 16.66 ms depending on the line frequency specified in the hardware (50 Hz or 60 Hz, respectively).



CAUTION


The filtering is based on a moving average over the specified *Filter time*. Some data at the onset of the measurement is required to establish the initial average. Stabilization of the filtering during the first approximately 20 ms of the measurement could lead to visible artifacts in the E vs. t curves at this timescale. If the data corresponding to the first 20 ms of the measurement is of interest, the measurement should be carried out with the **Filter** set to OFF.



CAUTION

The filtering is based on a moving average over the specified *Filter time*. Some data at the onset of the measurement is required to establish the initial average. Stabilization of the filtering during the first approximately 20 ms of the measurement could lead to visible artifacts in the E vs. t curves at this timescale. If the data corresponding to the first 20 ms of the measurement is of interest, the measurement should be carried out with the **Filter** set to OFF.

7.4.7 AC voltammetry



This command can be used to perform an AC voltammetry measurement, in potentiostatic conditions.



CAUTION

This command can only be used in Potentiostatic mode.

The details of the properties of the **AC voltammetry** command are shown in *Figure 339*.

- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.
- **Scan rate:** the calculated scan rate, in V/s, determined based on the step potential and the interval time.

Figure 340 represents the measurement properties of the **AC voltammetry** command, schematically.




Figure 340 Overview of the measurement properties of the AC voltammetry command

In a **AC voltammetry** measurement, the AC current is measured while the modulation is applied. This value is normally plotted against the applied potential. Alongside the AC current signal, the calculated impedance, admittance and phase shift are also available. Additional signals, provided by the Autolab instrument, can be sampled as well.



NOTICE

Frequency values in the range of 0.2 to 250 Hz are supported by the **AC voltammetry** command. Working with low frequency values (lower than 8 Hz) may require the adjustment of parameters such as the modulation time and the current range.



NOTICE

When this command is used, a drop-down list provides the choice of the reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE

The **Record signals** command does not apply any potential or current.

The details of the properties of the **Record signals** command are shown in *Figure 342*:

The screenshot shows a 'Properties' dialog box for the 'Record signals' command. The dialog has a title bar with a close button. Below the title bar, the text 'Record signals' is displayed. There are five input fields, each with a label and a value: 'Command name' with 'Record signals', 'Duration' with '5', 'Interval time' with '1', 'Estimated number of points' with '5', and 'Estimated duration' with '5'. Each field has a unit 's' to its right. A 'More' button is located at the bottom right of the dialog.

Figure 342 The properties of the *Record signals* command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Duration:** specifies the duration of the measurement, in s.
- **Interval time:** specifies the interval time used in the measurement, in s. The smallest interval time is 1.33 ms.
- **Estimated number of points:** this property shows the estimated number of points, determined from the specified Duration and Interval time.
- **Estimated duration:** this property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.

The **Record signals** command provides access to additional measurement options, through the More button *Additional measurement command properties* (see *chapter 9, page 615*).

The following additional properties are available (see *figure 343, page 288*):

7.5.1.2 Using the fast options

When the fast options are used, the same strategy as in the default mode used is for measuring the data points but the options are now verified after each user-defined fast interval time. This means that the options can be verified at a faster rate than the sampling rate.



NOTICE

The fast interval time must be smaller than the interval time and must be an integral fraction of the interval time.

When the fast options are used, the measurement options are decoupled from the sampling of the data. This is particularly useful for long measurements on a cell that requires the options to be tested with a short interval time.

Figure 344 shows an example of such a set of properties. The duration and interval time are set to 10 s and 0.2 seconds, respectively. This leads to an estimated number of points of 50. Setting the fast interval time to 20 ms, the options will be verified at a much faster rate than the sampling rate.




Figure 344 Example of fast options

7.5.1.3 Use the fast options and the time derivative threshold

When both the fast options and the time derivative threshold are used, the same strategy as in the previous mode is used. The time derivative value of one or more signals is determined using the fast interval time. For each time derivative signal, a threshold can be defined by the user. When the absolute value of a time derivative signal exceeds the specified threshold, the data points are measured using the fast interval time instead of the interval time. This means that the sampling rate can be modified depending on the derivative of one or more signals.



NOTICE

In order to use the time derivative threshold, at least one time derivative must be sampled.



NOTICE

The fast interval time must be smaller than the interval time and must be an integral fraction of the interval time.

Figure 345 shows an example of such a set of properties. The duration and interval time are set to 10 s and 0.2 seconds, respectively. This leads to an estimated number of points of 50. Setting the fast interval time to 20 ms, the options will be verified at a much faster rate than the sampling rate. Furthermore, if the dWE(1).Potential/dt signal is measured, a threshold value can be specified for this signal. In this example, the threshold is set to 1 V/s.





Figure 345 Example of time derivative threshold


Using these properties, the dWE(1).Potential/dt signal is calculated every 20 ms. If the value of this signal is larger (in absolute value) than the specified threshold, a data point is collected using the fast interval time instead of the interval time.

7.5.2 Chrono methods



This command can be used to record the signals during a user-defined sequence of potential or current steps.

The **Chrono methods** command can be used in two different modes, which can be selected using the provided drop-down list (see figure 346, page 291):



The screenshot shows a 'Properties' dialog box with a title bar and a close button. Under the heading 'Chrono methods', there are several input fields: 'Command name' (text box with 'Chrono methods'), 'Mode' (dropdown menu with 'Normal' selected and 'High speed' highlighted), 'Number of repeats' (text box with 'Normal'), 'Total duration' (text box), 'Estimated number of points' (text box with '0'), and 'Estimated duration' (text box with '0' and a 's' unit). A 'More' button is located at the bottom right.

Figure 346 Two modes are provided by the Chrono methods command

- Normal
- High speed



NOTICE

The **Mode** selection drop-down is only shown for instruments that are equipped with a fast sampling ADC module (**ADC10M** or **ADC750**). Please refer to *Chapter 16.3.2.1* and *Chapter 16.3.2.2* for more information. Instruments that are not fitted with a fast sampling ADC module can only use the Normal mode. For those instruments, the **Mode** drop-down control is not shown.

7.5.2.1 Chrono methods - Normal

The following properties are available when the **Chrono methods** command is used in **Normal** mode (see figure 347, page 292):




Figure 348 The Chrono methods Sequence editor in Normal mode



NOTICE

Additional settings can be adjusted in the Sequence editor. Please refer to *Chapter 9* for more information.

7.5.2.2 Chrono methods - High speed

The following properties are available when the **Chrono methods** command is used in **High speed** mode (see figure 349, page 293):

Figure 349 The properties of the Chrono methods command in High speed mode

The following properties are available:

- **Command name:** a user-defined name for the command.



- **Mode:** this property defines the measurement mode of the Chrono methods command. A drop-down control provides the choice between Normal (default) and High speed.
- **Interval time:** specifies the interval time, in s, used for the measurement.
- **Total duration:** indicates the expected duration of the chrono methods measurement, as read-only property, in s.
- **Estimated number of points:** this read-only property shows the estimated number of points.
- **Estimated duration:** this read-only property shows the estimated duration, in s. This property is defined by the estimated number of points and the interval time, as well as the duration of underlying commands, if applicable.



NOTICE

In **High speed** mode, the interval time used in the **Chrono methods** command is constant.

The sequence of steps used by the **Chrono methods** command can be edited by clicking the More button.



Figure 350 The Chrono methods Sequence editor in High speed mode



NOTICE

Additional settings can be adjusted in the Sequence editor. Please refer to *Chapter 9* for more information.

7.5.2.3 Using the Sequence editor

Using the buttons located in the top left corner of the Sequence editor panel, the sequence of steps used in the **Chrono methods** command can be constructed.

The following tasks can be carried out using the Sequence editor:

1. Add an item to the sequence (using the **+** button)
2. Remove an item from the sequence (using the **-** button)
3. Duplicate an item in the sequence (using the **⌘** button)
4. Move a sequence item up or down (using the **↑** button and the **↓** button)

Clicking the **+** button in the Chrono methods Sequence editor reveals a drop-down list that can be used to add one of the following items:

- **Step:** this item creates a step in the sequence. A step applies a potential or current value on the cell for the specified duration during which the data is recorded.
- **Level:** this item creates a level in the sequence. A level does not apply a potential or current value on the cell. The data is recorded during the specified duration.
- **Repeat:** this item creates a new sub-sequence in the main sequence, in which new items can be added. This sub-sequence can be repeated any number of times and the electrochemical response of the cell is sampled during the whole sub-sequence.
- **Repeat (unsampled):** this item creates a sub-sequence in the main sequence, in which new items can be added. This sub sequence can be repeated any number of times, but does not generate any data points, as the electrochemical response of the cell is unsampled for the whole sub-sequence. This can be useful for conditioning the electrode with a pulse sequence, or for reducing the number of data points recorded during long measurements.




Figure 352 The properties of the Step item

Depending on the mode in which the **Chrono methods** command is used, the Step item has the following properties:

- **Text:** a label that can be used to provide a name to the Step, for book-keeping purposes.
- **Duration:** the duration of the Step, in s.
- **Sample*:** a toggle that can be used to switch the sampling of data on or off. When the sampling is switched off, the specified potential or current value will be applied but no data points will be measured.
- **Interval time*:** the interval time, in s, used for sampling the data during the Step.
- **Estimated number of points*:** a read-only property that indicates the expected number of data points, based on the interval time and the duration.
- **Potential/Current:** the potential or current, in V or A, applied during the Step, respectively.



NOTICE

When the Chrono methods command is used in potentiostatic mode, a drop-down list provides the choice of reference used to apply a potential value. The potential can be specified with respect to the reference electrode potential (V_{REF}) or the open-circuit potential (V_{OCP}).



NOTICE

The properties indicated by a * are not available when the **Chrono methods** command is used in **High speed** mode.

Two additional advanced properties are available, through the Pre Autolab control and the Post Autolab control buttons. These buttons can be used to define a specific instrument setting before the Step is applied or after the step is applied, respectively. Both buttons provide access to the **Autolab control** properties.



NOTICE

For more information on the **Autolab control** command, please refer to *Chapter 7.2.1*.

7.5.2.3.2 Using the Level item

The Level item can be added to the Sequence editor by clicking the **+** button.

This item creates a Level in the sequence. The Level does not change the applied potential or current. The response of the cell is measured for the specified duration.

The properties of the added Level are shown in the dedicated sub-panel (*see figure 353, page 299*).




Figure 353 The properties of the Level item



NOTICE

The Level item is not available when the **Chrono methods** command is used in **High speed** mode.

The Level item has the following properties:

- **Text:** a label that can be used to provide a name to the Level, for bookkeeping purposes.
- **Duration:** the duration of the Level, in s.
- **Sample:** a toggle that can be used to switch the sampling of data on or off. When the sampling is switched off, the specified potential or current value will be applied but no data points will be measured.
- **Interval time:** the interval time, in s, used for sampling the data during the Level.
- **Estimated number of points:** a read-only property that indicates the expected number of data points, based on the interval time and the duration.

The Repeat item has the following properties:

- **Text:** a label that can be used to provide a name to the Repeat, for bookkeeping purposes.
- **Repeats:** the number of repetitions of the sub-sequence.

Two additional advanced properties are available, through the Pre Autolab control and the Post Autolab control buttons. These buttons can be used to define a specific instrument setting before the Repeat is started or after the Repeat is finished, respectively. Both buttons provide access to the **Autolab control** properties.



NOTICE

For more information on the **Autolab control** command, please refer to *Chapter 7.2.1*.

7.5.2.3.4 Using the Repeat (unsampled) item

The Repeat (unsampled) item can be added to the Sequence editor by clicking the **+** button.

This item creates a Repeat (unsampled) in the sequence. The Repeat (unsampled) item creates a sub-sequence to which new Steps or Levels can be added. The whole sub-sequence can be repeated, however no data will be sampled during this sub-sequence.

The properties of the added Repeat (unsampled) are shown in the dedicated sub-panel (see *figure 354, page 300*).




Figure 355 The properties of the Repeat (unsampled) item



NOTICE

The Repeat (unsampled) item is not available when the **Chrono methods** command is used in High speed mode.

The Repeat (unsampled) item has the following properties:

- **Text:** a label that can be used to provide a name to the Repeat (unsampled), for bookkeeping purposes.
- **Repeat (unsampled):** the number of repetitions of the sub-sequence.

Two additional advanced properties are available, through the Pre Autolab control and the Post Autolab control buttons. These buttons can be used to define a specific instrument setting before the Repeat (unsampled) is started or after the Repeat (unsampled) is finished, respectively. Both buttons provide access to the **Autolab control** properties.



NOTICE

For more information on the **Autolab control** command, please refer to *Chapter 7.2.1*.

7.6 Measurement - impedance commands

Commands located in the **Measurement – impedance** group can be used to perform impedance and impedance spectroscopy measurements or measurements involving sinewave modulations.



CAUTION

Impedance measurements require the optional **FRA32M** or **FRA2** module *FRA32M module (see chapter 16.3.2.13, page 1112)*.

The available commands are represented by a shortcut icon (*see figure 356, page 302*)





Figure 356 The Measurement - impedance commands

The following commands are available:

- **FRA measurement:** a command which can be used to perform a frequency scan *FRA measurement* (see chapter 7.6.1, page 303).
- **FRA single frequency:** a command which can be used to measure the impedance at a single frequency *FRA single frequency* (see chapter 7.6.2, page 305).
- **Electrochemical frequency modulation:** a command which can be used to perform Electrochemical Frequency Modulation measurements *Electrochemical Frequency Modulation* (see chapter 7.6.4, page 327).

7.6.1 FRA measurement



This command can be used to perform an impedance spectroscopy measurement through a frequency scan. This command requires the optional **FRA32M** or **FRA2** module.



CAUTION

The **FRA measurement** command requires the optional **FRA32M** or **FRA2** module *FRA32M module* (see chapter 16.3.2.13, page 1112).

The details of the properties of the **FRA measurement** command are shown in *Figure 357*:

- **Amplitude:** specifies the amplitude to apply during the frequency scan. The units depend on the specified mode and on the *Input connection* property (V, A, or external units).
- **Use RMS amplitude:** a  toggle control provided to specify if the amplitude value is the root mean squared (RMS) or top value.
- **Wave type:** specifies the type of signal used during the frequency scan. The choice is provided between the default single sine or the multi sine wave types (Single sine, 5 sines or 15 sines).
- **Input connection:** specifies if the measurement should be carried out internally (through the PGSTAT) or externally, using the external inputs provided on the front panel of the **FRA32M** or **FRA2** module.
- **Estimated duration:** the estimated duration of the command, in s. This property is defined by the frequencies and the acquisition settings and the duration of underlying commands, if applicable.

Table 7 provides an overview of the formulae used to calculate the distributions supported by the FRA measurement command.

Table 7 The distributions used in the FRA measurement command

Type	Increment, Δ	Distribution
Linear	$\Delta = \frac{\text{End} - \text{Start}}{N - 1}$	$V_i = (\text{Start} + (i - 1)\Delta)$
Square root	$\Delta = \frac{\sqrt{\text{End}} - \sqrt{\text{Start}}}{N - 1}$	$V_i = (\sqrt{\text{Start}} + (i - 1)\Delta)^2$
Logarithmic	$\frac{\text{LOG}(\text{End}) - \text{LOG}(\text{Start})}{N - 1}$	$\text{LOG}(V_i) = (\text{LOG}(\text{Start}) + (i - 1)\Delta)$
Points per decade	$\Delta = 10^{\frac{\text{FLOOR}\left(\left(\frac{\text{LOG}(\text{Start})}{\text{LOG}(\text{End})}\right)\right) N_{\text{Decade}}}{\text{LOG}(\text{Start}) - \text{LOG}(\text{End})}}$	$V_i = (\text{Start} \cdot \Delta^{(i-1)})$



NOTICE

The **FRA measurement** command provides access to additional options, through the  button *Additional properties* (see chapter 7.6.3, page 307).

7.6.2 FRA single frequency

- **Input connection:** specifies if the measurement should be carried out internally (through the PGSTAT) or externally, using the external inputs provided on the front panel of the **FRA32M** or **FRA2** module.
- **Estimated duration:** the estimated duration of the command, in s. This property is defined by the frequency and the acquisition settings and the duration of underlying commands, if applicable.



NOTICE

The **FRA single frequency** command provides access to additional options, through the More button *Additional properties* (see chapter 7.6.3, page 307).

7.6.3 Additional properties

Unlike the other measurement command which common additional properties, described in detail in *Chapter 9*, the **FRA measurement** and the **FRA single frequency** commands have specific additional properties, which can be accessed by clicking the More button or by double-clicking the command tile in the procedure editor (see figure 359, page 307).

The screenshot shows a 'Properties' dialog box for the 'FRA measurement' command. The dialog contains the following fields and controls:

- Command name: FRA measurement
- First applied frequency: 1E+05 Hz
- Last applied frequency: 0,1 Hz
- Number of frequencies: 10 per decade
- Frequency step type: Points per decade (dropdown menu)
- Amplitude: 0,01 V_{RMS}
- Use RMS amplitude:
- Wave type: Sine (dropdown menu)
- Input connection: Internal (dropdown menu)
- Estimated duration: 200,12 s
- A 'More' button is located at the bottom right of the dialog, with a mouse cursor pointing to it.

Figure 359 Accessing the additional properties of the FRA measurement command

7.6.3.1 Advanced properties - sampler

The **Sampler** section provides sampling settings used by the FRA measurement (see figure 361, page 309).

The screenshot shows the 'FRA measurement' settings window. On the left, there is a navigation menu with 'Sampler' selected. The main area is divided into two sections: 'Basic' and 'Advanced'.

Basic section:

- Maximum integration time: 0,125 s
- Minimum number of integration cycles: 1
- Sample time domain: [checkbox]
- Sample frequency domain: [checkbox]
- Sample DC: [checkbox]
- Calculate admittance: [checkbox]

Advanced section:

- Transfer function: Re - j Im
- Lowest bandwidth: High stability
- Number of cycles to reach steady state: 10
- Maximum time to reach steady state: 1 s
- With a minimum fraction of a cycle: 0
- Automatic amplitude correction: [checkbox]
- Iterative: [checkbox]
- Amplitude threshold percentage: 5 %
- Automatic resolution correction: [checkbox]
- Iterative: [checkbox]
- Minimum resolution: 32 %
- Maximum amount of re-measurements: 25

Figure 361 The Sampler settings

The available properties are divided in two sub-sections:

- **Basic:** these are basic acquisition properties that define how FRA measurements are carried out.
- **Advanced:** these are advanced settings, predefined to an optimal value for most measurements. Please read the following section of the manual carefully before adjusting these properties.

The following **Basic** properties are available:

- **Maximum integration time:** the longest time, in seconds, during which the AC response of the cell is recorded for data analysis. Long integration time values increase the duration of the measurement but improve the signal to noise ratio. The default value is 0.125 s.
- **Minimum integration cycles:** defines the minimum cycles of the AC response to record for data analysis (this value overrides the previous property at low frequencies). The minimum value is one cycle and the maximum number is 16. Integrating over a large number of cycles increases the duration of the measurement but improves the signal to noise ratio. This value must be an integer. The default value is 1.



- **Sample time domain:** defines if the time domain information should be sampled during a FRA measurement, using the provided toggle. The time domain information consists of the raw potential and current sine waves. This information can be used to build a Lissajous plot, or to evaluate the signal to noise ratio and to verify the linearity of the cell response. The Time domain, Potential (AC) and Current (AC) and Potential resolution and Current resolution signals are added to the data for each individual frequency when this option is on.
- **Sample frequency domain:** defines if the frequency domain information should be sampled during a FRA measurement, using the provided toggle. The frequency domain information consists of the calculated FFT results obtained from the measured time domain. The frequency domain information can be used to evaluate the measured frequency contributions.
- **Sample DC:** defines if the DC component of the two input signals (Potential and Current, or external signals) information should be sampled during a FRA measurement, using the provided toggle. This option is active by default.
- **Calculate admittance:** specifies if the admittance values must be calculated during the measurement (Y' , $-Y''$), using the provided toggle.

The **Integration time** and **Minimum number of cycles to integrate** define the duration of the data acquisition segment for each frequency. These two properties are competing against one another and depending on the frequency of the applied signal: the measurement duration will be defined by one of these two properties:

- If the *Frequency* is larger than $(1/\text{Integration time})$, the acquisition time will be defined by the **Integration time** property.
- If the *Frequency* is smaller than $(1/\text{Integration time})$, the acquisition time will be defined by the **Minimum number of cycles to integrate** property.

The following **Advanced** properties are available:

- **Transfer function:** this property defines how the impedance is expressed in terms of its real (Re, or Z') and imaginary (Im or Z'') components, using the provided drop-down list. By default, the Re-jIm convention is used. Using this toggle it is possible to acquire impedance data using the alternative convention.

- **Lowest bandwidth:** defines the lowest bandwidth setting used by the Autolab during the measurement (High stability, High speed and Ultra high speed). In High stability, the bandwidth will automatically be set to High stability for frequencies below 10 kHz and to High speed for frequencies below 100 kHz. The Ultra high speed mode is used for frequencies above 100 kHz. When this property is set to High speed, the High stability setting is not used. When this property is set to Ultra high speed, only this bandwidth setting is used. High stability is set by default and is recommended for most measurements.
- **Number of cycles to reach steady state:** defines the number of cycles to apply in between two consecutive frequencies before resuming data acquisition (0 – 30000 cycles). The default value is 10.
- **Maximum time to reach steady state:** defines the maximum amount of time to wait between two consecutive frequencies before resuming data acquisition (this value overrides the previous setting at low frequencies). The maximum value is 30000 s. The default value is 1.
- **With a minimum fraction of a cycle:** defines the minimal fraction of a cycle to wait before the response can be recorded (overrides the previous setting at very low frequencies). This value can be set between 0 and 1. The default value is 0.
- **Automatic amplitude correction:** defines if the automatic amplitude correction algorithm is used during FRA measurements, using the provided  toggle. Two automatic amplitude correction modes are available. The selection of the correction algorithm is defined by the position of the **Iterative** toggle:
 - **Normal correction:** this correction mode is enabled when the toggle of the **Iterative** property is set to . In this mode, the amplitude is measured at each frequency during the measurement and the applied amplitude is adjusted for the next frequency point. This means that each frequency value is only adjusted once.
 - **Iterative correction:** this correction mode is enabled when the toggle of the **Iterative** property is set to . In this mode, the amplitude is measured at each frequency and adjusted until the applied amplitude is equal to the expected amplitude within the tolerances specified by the **amplitude threshold percentage property**. This means that each frequency value can be adjusted multiple times before the correct amplitude is measured.
- **Amplitude threshold percentage:** defines the threshold value used to control the applied amplitude. When the **Automatic amplitude correction** property is on and set to *iterative*, the applied amplitude will be considered to be equal to the required amplitude if it fits within the specified **Amplitude threshold percentage**. This value is defined as a percentile value. The default value is 5 %.



- **Automatic resolution correction:** defines if the automatic resolution correction algorithm is used during FRA measurements, using the provided  toggle. Two automatic resolution correction modes are available. The selection of the correction algorithm is defined by the position of the **Iterative** toggle:
 - **Normal correction:** this correction mode is enabled when the toggle of the **Iterative** property is set to . In this mode, the resolution is measured at each frequency during the measurement and the gain factors and current ranges are adjusted, if applicable, for the next frequency point. This means that each frequency value is only adjusted once.
 - **Iterative correction:** this correction mode is enabled when the toggle of the **Iterative** property is set to . In this mode, the resolution is measured at each frequency and adjusted until the measured resolution is higher or equal to the specified **minimum resolution**. This means that each frequency value can be adjusted multiple times before the correct resolution is measured.
- **Minimum resolution:** defines the minimum resolution value to reach on both input channels of the impedance analyzer module. When the **Automatic resolution correction** property is on and set to *iterative*, the data will be remeasured until the resolution of both input signals is higher or equal to the specified **Minimum resolution**. This value is defined as a percentile value. The default value is 32 %.
- **Maximum amount of re-measurements:** defines the maximum number of re-measurements allowed if the **Automatic amplitude correction** property is on. The default value is 25.

The **Number of cycles to reach steady state**, **Maximum time to reach steady state**, **With a minimum fraction of a cycle** properties define the duration of the stabilization segment between two consecutive frequencies. These properties are competing against one another, and depending on the frequency of the applied signal, the timing will be defined by one of these three properties:

- If the (Number of cycles to reach steady state/Frequency) is smaller than Maximum time to reach steady state, the settling time will be defined by the **Number of cycles to reach steady state** property.
- If the (Number of cycles to reach steady state/Frequency) is larger than Maximum time to reach steady state, the settling time will be defined by the **Maximum time to reach steady state** property.
- For very low frequencies, if (1/Frequency) is larger than the **Maximum time to reach steady state**, the settling time will be defined by the **With the minimum fraction of a cycle** property.

Figure 362 provides a schematic overview of the five timing properties. At low frequency, the integration time is overruled by the minimum number

of cycles to integrate. At low frequency, the number of cycles to reach steady state property is overruled by the maximum time to reach steady state property or the by the with a minimum fraction of cycle property, in a similar way.




Figure 362 Overview of the timing properties (Integration time: 2 s, Minimum integration cycles: 4, Number of cycles to wait for steady state: 4, Maximum time to reach steady state: 4 s, Minimal cycle fraction of to wait for steady state: 0.5)

7.6.3.2 Advanced properties - options

The **Options** section provides automatic current ranging used by the FRA measurement (see figure 377, page 325).




Figure 363 The Options settings

The following properties are available:

- **Automatic current ranging:** sets the automatic current ranging option on or off, using the provided toggle.
- **Highest current range:** defines the highest allowed current range, using the provided drop-down list.
- **Lowest current range:** defines the lowest allowed current range, using the provided drop-down list.



In FRA measurements, using the **FRA measurement** or the **FRA single frequency** commands, it is possible to use the **Booster10A** or **Booster20A** current range in the automatic current ranging option (see figure 364, page 314).




Figure 364 The current range of the Booster10A and Booster20A is available in the Automatic current ranging option

7.6.3.3 Advanced properties - plots

The **Plots** section provides plot settings used by the FRA measurement (see figure 365, page 314).




Figure 365 The Plots settings



NOTICE

The available plots depend on the properties specified in the **Sampler** section.

The following plots are available:

- **Nyquist impedance:** plots the measured $-Z''$ values versus the measured Z' values, using isometric axes.
- **Nyquist admittance:** plots the calculated $-Y''$ values versus the calculated Y' values, using isometric axes.
- **Bode:** plots the phase (in opposed values) and the logarithm of the measured impedance (Z), versus the logarithm of the frequency.
- **AC vs t:** plots the raw sinewave amplitudes for the potential and the current signals versus the time.
- **Resolution vs t:** plots the instrumental resolution for the measured potential and current signals versus the time.
- **Lissajous:** plots the raw AC current versus the raw AC potential.

7.6.3.4 Advanced properties - summary

The **Summary** section provides a complete summary of the frequency scan parameters defined in the **FRA measurement** command (see figure 366, page 315).

FRA measurement

Sampler
Options
Plots
Summary

	Frequency (Hz)	Amplitude (V_{RMS})	Wave type	Integration time (s)	Integration cycles
1	100000	0,01	Sine	0,125	1
2	10000	0,01	Sine	0,125	1
3	1000	0,01	Sine	0,125	1
4	100	0,01	Sine	0,125	1
5	10	0,01	Sine	0,125	1
6	1	0,01	Sine	0,125	1
7	0,1	0,01	Sine	0,125	1

Add range

First applied frequency Hz
 Last applied frequency Hz
 Number of frequencies
 Frequency step ▾
 Amplitude V
 Wave type ▾

Figure 366 The Summary section of the FRA measurement command



NOTICE

This section is only available for the **FRA measurement** command.

Each frequency in the scan is displayed in the table, along with the amplitude values, wave type, minimum integration time and maximum number of cycles to integrate (see figure 366, page 315).

A tooltip displays information about each individual frequency in the scan (see figure 367, page 316).

FRA measurement

Sampler
Options
Plots
Summary

	Frequency (Hz)	Amplitude (V _{RMS})	Wave type	Integration time (s)	Integration cycles
1	100000	0,01	Sine	0,125	1
2	10000	0,01	Sine	0,125	1
3	1000	0,01	Sine	0,125	1
4	100	0,01	Sine	0,125	1
5	10	0,01	Sine	0,125	1
6	1	0,01	Sine	0,125	1
7	0,1	0,01	Sine	0,125	1

Tooltip: TOP: 0.014142 V

Add range

First applied frequency: 100000 Hz
 Last applied frequency: 0,1 Hz
 Number of frequencies: 50
 Frequency step: Logarithmic
 Amplitude: 0,01 V
 Wave type: Sine

Add range

Figure 367 Information on each individual frequency is displayed in a tooltip

The summary section can be used to fine tune the frequency scan in three different ways:

1. The properties of one or more frequencies in the table can be adjusted manually *Manual modification of the properties* (see chapter 7.6.3.4.1, page 317).
2. Additional frequencies can be added manually to the table *Manual modification of the properties* (see chapter 7.6.3.4.1, page 317).
3. The frequencies can be sorted ascending or descending *Sorting the table* (see chapter 7.6.3.4.3, page 321).

7.6.3.4.1 Manual modification of the properties

If needed, one of more properties can be overruled in the summary table shown in the **Summary** section. Double click the field to be edited in the table and modify the value (see figure 368, page 317).

The screenshot shows the 'FRA measurement' interface. On the left, there is a sidebar with 'Summary' selected. The main area contains a table with the following data:

	Frequency (Hz)	Amplitude (V _{RMS})	Wave type	Integration time (s)	Integration cycles
1	100000	0,02	Sine	0,125	1
2	10000	0,01	Sine	0,125	1
3	1000	0,01	Sine	0,125	1
4	100	0,01	Sine	0,125	1
5	10	0,01	Sine	0,125	1
6	1	0,01	Sine	0,125	1
7	0,1	0,01	Sine	0,125	1

Below the table is the 'Add range' section with the following fields:

- First applied frequency: 100000 Hz
- Last applied frequency: 0,1 Hz
- Number of frequencies: 50
- Frequency step: Logarithmic
- Amplitude: 0,01 V
- Wave type: Sine

An 'Add range' button is located at the bottom right of the 'Add range' section.

Figure 368 Modifying individual values in the Summary section

Press the **[Enter]** key to validate the new value. The new value will be used instead of the initial value specified in the **FRA measurement** command properties (see figure 369, page 318).



FRA measurement

Sampler
Options
Plots
Summary

	Frequency (Hz)	Amplitude (V _{RMS})	Wave type	Integration time (s)	Integration cycles
1	100000	0,02	Sine	0,125	1
2	10000	0,01	Sine	0,125	1
3	1000	0,01	Sine	0,125	1
4	100	0,01	Sine	0,125	1
5	10	0,01	Sine	0,125	1
6	1	0,01	Sine	0,125	1
7	0,1	0,01	Sine	0,125	1

Add range

First applied frequency: 100000 Hz
 Last applied frequency: 0,1 Hz
 Number of frequencies: 50
 Frequency step: Logarithmic
 Amplitude: 0,01 V
 Wave type: Sine

Add range

Figure 369 The new value is updated in the FRA editor Summary section

7.6.3.4.2 Adding frequencies to the table

If needed, additional frequencies can be added to the table presented in the **Summary** section, by editing the properties located in the **Add range** sub-panel and clicking the **Add range** button (see figure 370, page 318).

FRA measurement

Sampler
Options
Plots
Summary

	Frequency (Hz)	Amplitude (V _{RMS})	Wave type	Integration time (s)	Integration cycles
1	100000	0,02	Sine	0,125	1
2	10000	0,01	Sine	0,125	1
3	1000	0,01	Sine	0,125	1
4	100	0,01	Sine	0,125	1
5	10	0,01	Sine	0,125	1
6	1	0,01	Sine	0,125	1
7	0,1	0,01	Sine	0,125	1

Add range

First applied frequency: 0,1 Hz
 Last applied frequency: 0,001 Hz
 Number of frequencies: 1 per decade
 Frequency step: Points per decade
 Amplitude: 0,01 V
 Wave type: Sine

Add range

Figure 370 Adding additional frequencies to the range

The **Summary** section will be updated and the new frequencies will be added to the table (see figure 371, page 319).

The screenshot shows the 'FRA measurement' interface. On the left, there is a navigation menu with 'Summary' highlighted. The main area contains a table with the following data:

	Frequency (Hz)	Amplitude (V_{RMS})	Wave type	Integration time (s)	Integration cycles
1	100000	0,02	Sine	0,125	1
2	10000	0,01	Sine	0,125	1
3	1000	0,01	Sine	0,125	1
4	100	0,01	Sine	0,125	1
5	10	0,01	Sine	0,125	1
6	1	0,01	Sine	0,125	1
7	0,1	0,01	Sine	0,125	1
8	0,1	0,01	Sine	0,125	1
9	0,01	0,01	Sine	0,125	1
10	0,001	0,01	Sine	0,125	1

Below the table is the 'Add range' configuration panel with the following settings:

- First applied frequency: 0,1 Hz
- Last applied frequency: 0,001 Hz
- Number of frequencies: 1 per decade
- Frequency step: Points per decade
- Amplitude: 0,01 V
- Wave type: Sine

An 'Add range' button is located at the bottom right of the configuration panel.

Figure 371 The updated summary

It is also possible to delete frequencies from the table by selecting one or more rows of the table and pressing the **[Delete]** button (see figure 372, page 320).



FRA measurement

Sampler
Options
Plots
Summary

	Frequency (Hz)	Amplitude (V_{RMS})	Wave type	Integration time (s)	Integration cycles
1	100000	0,02	Sine	0,125	1
2	10000	0,01	Sine	0,125	1
3	1000	0,01	Sine	0,125	1
4	100	0,01	Sine	0,125	1
5	10	0,01	Sine	0,125	1
6	1	0,01	Sine	0,125	1
7	0,1	0,01	Sine	0,125	1
8	0,1	0,01	Sine	0,125	1
9	0,01	0,01	Sine	0,125	1
10	0,001	0,01	Sine	0,125	1

Add range

First applied frequency Hz
 Last applied frequency Hz
 Number of frequencies
 Frequency step
 Amplitude V
 Wave type

Figure 372 It is possible to delete frequencies from the table
 The selected frequencies will be deleted from the table (see figure 373, page 321).

← FRA measurement

Sampler
Options
Plots
Summary

	Frequency (Hz)	Amplitude (V _{RMS})	Wave type	Integration time (s)	Integration cycles
1	10000	0,01	Sine	0,125	1
2	1000	0,01	Sine	0,125	1
3	100	0,01	Sine	0,125	1
4	10	0,01	Sine	0,125	1
5	1	0,01	Sine	0,125	1
6	0,1	0,01	Sine	0,125	1
7	0,1	0,01	Sine	0,125	1
8	0,01	0,01	Sine	0,125	1
9	0,001	0,01	Sine	0,125	1

Add range

First applied frequency Hz
 Last applied frequency Hz
 Number of frequencies
 Frequency step ▾
 Amplitude V
 Wave type ▾

Figure 373 The selected frequencies are removed from the table

7.6.3.4.3 Sorting the table

It is possible to sort the frequencies ascending or descending in the **Summary** section by clicking the Frequency column header. A small arrow will be displayed in the column header indicating the sorting direction. Clicking the header again will cycle between ascending and descending (see figure 374, page 322).



← FRA measurement

Sampler

Options

Plots

Summary

	Frequency (Hz)	Amplitude (V _{RMS})	Wave type	Integration time (s)	Integration cycles
1	0,001	0,01	Sine	0,125	1
2	0,01	0,01	Sine	0,125	1
3	0,1	0,01	Sine	0,125	1
4	0,1	0,01	Sine	0,125	1
5	1	0,01	Sine	0,125	1
6	10	0,01	Sine	0,125	1
7	100	0,01	Sine	0,125	1
8	1000	0,01	Sine	0,125	1
9	10000	0,01	Sine	0,125	1
10	100000	0,02	Sine	0,125	1

Add range

First applied frequency Hz

Last applied frequency Hz

Number of frequencies per decade

Frequency step

Amplitude V

Wave type

Figure 374 Sorting the frequencies in ascending or descending direction



NOTICE

It is only possible to the frequency values in the table.

7.6.3.5 Multi sine measurements

For very low frequency measurements, the multi sine option can be used to increase the acquisition speed. This option allows you to create a frequency scan in which low frequency single sine signals are replaced by a linear combination of five or fifteen frequencies. Each multi sine signal generates a number of data points equal to the number of components in the linear combination.

The multi sine option is **only** available for low frequencies. The maximum frequencies for multi sine measurements are listed in *Table 8*.

Table 8 The frequency limit for multi sine measurement for the FRA32M and FRA2 modules


Module	5 sine frequency limit	15 sine frequency limit
FRA32M	320 Hz	32 Hz
FRA2	3472 Hz	315.2 Hz

To create a multi sine frequency scan, the **Wave type** property can be set in the **Properties** panel (see figure 375, page 323).

The screenshot shows the 'Properties' panel for an 'FRA measurement'. The 'Wave type' dropdown menu is expanded, showing three options: 'Sine', '5 sines', and '15 sines'. A mouse cursor is hovering over the '5 sines' option. Other visible settings include: Command name: FRA measurement; First applied frequency: 1E+05 Hz; Last applied frequency: 0,1 Hz; Number of frequencies: 10 per decade; Frequency step type: Points per decade; Amplitude: 0,01 V_{RMS}; Use RMS amplitude: checked; Input connection: Sine; Estimated duration: 5 sines, 15 sines. A 'More' button is located at the bottom right of the panel.

Figure 375 The Wave type property can be used to create a multi sine measurement

Depending on the available hardware (**FRA32M** or **FRA2**), the frequency scan will be generated and the multi sine signals will be used for the low frequencies in the range. The **Summary** section displays the wave type in the table (see figure 376, page 324).



← FRA measurement

External

Sampler
Plots
Summary

←V

Use conversion factor

Conversion factor 1

Unit V

←X

Channel name X Channel

Use conversion factor

Conversion factor 1

Unit V

←Y

Channel name Y Channel

Use conversion factor

Conversion factor 1

Unit V

Transfer function

Definition X/Y

Name H

Unit

Phase name -Phase


Real part name H'

Imaginary part name -H''

Figure 377 The External settings

The settings shown in the **External** section can be used to specify how external impedance measurements are carried out. Four sub-sections are available in the External section:

- ← **V**: these settings define the conversion from generated AC voltage to the AC signal used for external impedance measurements. The following properties are available:
 - **Use conversion factor**: defines if a conversion factor must be used, using the provided toggle.
 - **Conversion factor**: defines the conversion factor, if applicable.
 - **Unit**: defines the units of the external AC signal, if applicable.



←V

Use conversion factor

Conversion factor 0,001

Unit RPM



- → **X**: these settings define the conversion used for the X input of the impedance analyzer module. The following properties are available:
 - **Channel name**: defines the name of the input signal.
 - **Use conversion factor**: defines if a conversion factor must be used, using the toggle.
 - **Conversion factor**: defines the conversion factor used for the input signal.
 - **Unit**: defines the units of the input signal.

–X

Channel name	Rotation rate
Use conversion factor	<input checked="" type="checkbox"/>
Conversion factor	1000
Unit	RPM

- → **Y**: these settings define the conversion used for the X input of the impedance analyzer module. The following properties are available:
 - **Channel name**: defines the name of the input signal.
 - **Use conversion factor**: defines if a conversion factor must be used, using the toggle.
 - **Conversion factor**: defines the conversion factor used for the input signal.
 - **Unit**: defines the units of the input signal.


–Y

Channel name	Current
Use conversion factor	<input checked="" type="checkbox"/>
Conversion factor	0,001
Unit	mA

- **Transfer function**: these settings define how the external transfer function is calculated. The following properties are available:
 - **Definition**: defines how the transfer function is calculated, using the provided drop-down list (X/Y or Y/X).
 - **Name**: defines the name of the transfer function.
 - **Unit**: defines the unit of the transfer function.
 - **Phase name**: defines the name of the phase angle.
 - **Real part name**: defines the name of the real part of the transfer function.
 - **Imaginary part name**: defines the name of the imaginary part to the transfer function.

Transfer function	
Definition	X/Y
Name	H(EHD)
Unit	
Phase name	-Phase
Real part name	H'(EHD)
Imaginary part name	-H''(EHD)

7.6.4 Electrochemical Frequency Modulation



This command can be used to perform an electrochemical frequency modulation (EFM) measurement. This command requires the optional **FRA32M** module.



CAUTION

The **Electrochemical Frequency Modulation** command requires the optional **FRA32M** module *FRA32M module (see chapter 16.3.2.13, page 1112)*.



CAUTION

The **Electrochemical Frequency Modulation** command can only be used in **potentiostatic**. Automatic current ranging and other options like cutoffs and counters are not supported by this command.

The details of the properties of the **Electrochemical Frequency Modulation** command are shown in *Figure 378*:

- **Model:** specifies the model used to analyze the data, using the provided drop-down list. Three different models are provided:
 - **Activation control:** a model that provides a general purpose model of a corroding system.
 - **Diffusion control:** a model that can be used for systems for which no cathodic current is observed.
 - **Passivating:** a model that can be used for systems for which no anodic current is observed.
- **Density:** specifies the density of the sample in g/cm³.
- **Equivalent weight:** defines the equivalent weight of the sample in g/mol of exchanged electrons.
- **Surface area:** defines the area of the sample, in cm².
- **Estimated duration:** the estimated duration of the command, in s. This property is defined by the base frequency and the acquisition settings and the duration of underlying commands, if applicable.



NOTICE

The **Electrochemical Frequency Modulation** command provides access to additional options, through the More button *Additional properties* (see chapter 7.6.3, page 307).

During a measurement, the **Electrochemical Frequency Modulation** command will apply a multi sine perturbation in potentiostatic conditions using a signal containing two decoupled frequencies. These two frequencies are obtained using the specified **Base frequency**, **Multiplier 1** and **Multiplier 2**. The two frequencies used in the measurement are reported in the **Properties** panel (**Frequency 1** and **Frequency 2**).



CAUTION

Multiplier 1 and **Multiplier 2** must not be integral multiples one another.

The current resulting from the application of the two sinewaves is recorded and converted from the time domain to the frequency domain. The second and third harmonic of the base frequencies are determined as well as the intermodulated frequencies. *Table 7.6.4* provides a summary of the frequencies analyzed in the measured signal.

Table 9 The frequencies used by the Electrochemical Frequency Modulation command



Current spectrum component	Determined from
$i_{2\omega_1}$	Second harmonic of ω_1 (or ω_2)
$i_{3\omega_1}$	Third harmonic of ω_1 (or ω_2)
i_{ω_1, ω_2}	ω_1 (or ω_2)
$i_{2\omega_2 \pm \omega_1}$	Second order of ω_1 or ω_2 intermodulated with ω_2 or ω_1
$i_{\omega_2 \pm \omega_1}^2$	ω_1 or ω_2 intermodulated with ω_2 or ω_1
i_{ω_1, ω_2}^2	ω_1 (or ω_2)

Using the measured values, the **Electrochemical Frequency Modulation** command can calculate the following corrosion indicators:

- **b_a** : the anodic Tafel slope, in V/dec.
- **b_c** : the cathodic Tafel slope, in V/dec.
- **i_{corr}** : the corrosion current, in A.

These values are obtained using mathematical expression that depend on the selected model.

1. **Activation control model**

$$i_{corr} = \frac{i_{\omega_1, \omega_2}^2}{2\sqrt{8i_{\omega_1, \omega_2} \cdot i_{2\omega_2 \pm \omega_1} - 3i_{\omega_2 \pm \omega_1}^2}}$$

$$b_a = 2,303 \frac{i_{\omega_1, \omega_2} U_0}{i_{\omega_2 \pm \omega_1} + \sqrt{8i_{\omega_1, \omega_2} \cdot i_{2\omega_2 \pm \omega_1} - 3i_{\omega_2 \pm \omega_1}^2}}$$

$$b_c = 2,303 \frac{i_{\omega_1, \omega_2} U_0}{-i_{\omega_2 \pm \omega_1} + \sqrt{8i_{\omega_1, \omega_2} \cdot i_{2\omega_2 \pm \omega_1} - 3i_{\omega_2 \pm \omega_1}^2}}$$

2. **Diffusion control model**

$$i_{corr} = \frac{i_{\omega_1, \omega_2}^2}{2i_{\omega_2 \pm \omega_1}}$$

$$b_a = 2,303 \frac{i_{\omega_1, \omega_2} U_0}{2i_{\omega_2 \pm \omega_1}}$$

$$b_c = \infty$$

3. **Passivating model**

$$i_{corr} = \frac{i_{\omega_1, \omega_2}^2}{2i_{\omega_2 \pm \omega_1}}$$

$$b_a = \infty$$

$$b_c = 2,303 \frac{i_{\omega_1, \omega_2} U_0}{2i_{\omega_2 \pm \omega_1}}$$

U_0 is the applied potential amplitude, in V, in all three models.

Using the calculated values, additional corrosion indicators can be calculated:

- **j_{corr}** : the corrosion current density, in A/cm².
- **Corrosion rate**: the corrosion rate, in mm of material/year.
- **Polarization resistance**: the converted polarization resistance, in Ohm.

Finally, two causality factors can be calculated, regardless of the model:

- **Causality factor (2)**: a second order consistency factor. This value should be as close as possible to 2.
- **Causality factor (3)**: a third order consistency factor. This value should be as close as possible to 3.

These factors are calculated according to:

$$CF_2 = \frac{i_{\omega_2 \pm \omega_1}}{i_{2\omega_1}}$$

$$CF_3 = \frac{i_{2\omega_2 \pm \omega_1}}{i_{3\omega_1}}$$



NOTICE

The values of b_a and b_c are reported in absolute value.



NOTICE

The **Electrochemical Frequency Modulation** technique is based on *Electrochemical Frequency Modulation: A New Electrochemical Technique for Online Corrosion Monitoring*, R.W. Bosch, J. Hubrecht, W. F. Bogaerts, and B.C. Syrett, Corrosion, 57 (2001).

 <p>Windower</p>	<p>This command can be used to extract a subset of signal values from measured or calculated signals, to be used in the data analysis.</p>
---	--

The details of the properties of the **Windower** command are shown in *Figure 380*:




Figure 380 The properties of the Windower command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Source:** this is the source signal used to window the data. The source signal is one of the available signals available in the data. Only one source signal can be selected.
- **Begin** the start value used by the windower for the specified source signal.
- **End:** the end value used by the windower for the specified source signal.



NOTICE

The **Begin** and **End** values can be fine-tuned in the additional properties panel, available using the More button.

The boundaries can be fine-tuned in the additional properties panel of the **Windower** command. In this panel more than one row can be used in order to add more boundaries to the table shown in the panel *Figure 382*.




Figure 382 Adding additional boundaries to the Windower command

When the **Windower** command is used on data, the *Use distinct values* toggle is available. When this toggle is activated, the table is replaced by a list of checkboxes of the available value for the selected source signal (see *figure 383, page 335*).




Figure 383 Using the distinct values selection




Figure 385 The boundaries can be converted from table view to checkboxes view

7.7.2 Build signal

	<p>This command can be used to create signals based on measured or calculated signals and procedure properties, to be used in the data analysis.</p>
--	--

The details of the command properties of the **Build signal** command are shown in *Figure 386*:




Figure 386 The properties of the Build signal command

The following properties are available:

- **Command name:** a user-defined name for the command.



NOTICE

By default, the **Build signal** command does not have pre-defined search criteria and it will list all available properties and signals in the **Select** panel.

7.7.2.1 Using the Filter

The **Filter** sub-panel can be used to specify search criteria in order to filter the available signals and properties. To filter the available signals and properties, a filter can be created using the table located in the **Filter** sub-panel. First a *Filter type* can be selected from the first available drop-down list (see figure 388, page 339).




Figure 388 Selecting the Filter type

The drop-down list provides the choice between two possible filters:

- **Command type:** this filter provides the possibility to filter the properties and signals based on the type (or name) of a command.
- **Options:** this filter provides the possibility to filter the properties and signals based on a command option.

After specifying the *Filter type* property in the table, it is possible to select an available argument in first available cell of the second column of the table. The possible arguments are provided in a drop-down list (see figure 389, page 339).




Figure 389 Selecting the Command type



NOTICE

The available arguments provided in the drop-down list depend on the commands used in the procedure.

Once the argument is specified, the list of available signals and properties in the **Select** sub-panel will be updated. Only the signals and properties that match the filter criteria specified in the **Filter** sub-panel will be displayed (see figure 390, page 340).




Figure 390 The filtered list of signals and properties

It is possible to create multiple conditions for filtering the properties and signals. For example, it is possible to add a first *Filter type* and *Command type*, shown in Figure 391.




Figure 391 Combining multiple filter conditions

It is then possible to add a second *Filter type* and *Command type*, as shown in Figure 392.






Figure 392 Adding a second condition to the filter

This filter condition will show only the signals and properties provided by all the **Record signals (> 1 ms)** commands located inside a **Repeat n times** command. These signals and properties will be shown in the **Select** sub-panel (see figure 393, page 341).




Figure 393 The updated list of signals and properties

7.7.2.2 Selecting the signals

The signals and properties shown in the **Select** panel can be added to the **Build signal** command by using the provided toggle (see figure 394, page 342).



Figure 394 Selecting signals or properties in the Build signal editor

If a filter is used, only the signals and properties provided by the commands or options that fit the selection criteria specified by the filter will be shown (see figure 395, page 342).




Figure 395 Selecting filtered signals or properties in the Build signal editor

For the signals listed in the **Select** sub-panel, it is also possible to define a range of values, by specifying an index range in the input fields provided next to each of the selected signals (see figure 396, page 343).

Build signal

Filter

	Filter type	Command type
First filter on	Command type	Record signals (>1 ms)

Select

Parameter	Select	From	To	Signal name
Corrected time	<input type="checkbox"/>			Corrected time
Index	<input checked="" type="checkbox"/>	1	100	Index
Interval time (s)	<input type="checkbox"/>			Interval time (s)
Time	<input checked="" type="checkbox"/>	1	100	Time
WE(1).Current	<input checked="" type="checkbox"/>	1	100	WE(1).Current
WE(1).Potential	<input checked="" type="checkbox"/>	1	100	WE(1).Potential

Figure 396 Specifying a range of values for the selected signals



CAUTION

When a range is specified it is necessary to specify two values (From and To).



NOTICE

It is only possible to specify a range for signals.

If needed, a specific name can be specified in the Signal name field for each selected signal. By default, the name of the signal is shown in the Signal name column. However, it is possible to specify a custom name by typing it into the provided input field, as shown in *Figure 397*.




Figure 397 A custom name can be specified if needed



CAUTION

Each user-defined Signal name must be unique.

7.7.2.3 Remove a filter

To remove one of the filter conditions specified in the **Filter** sub-panel, click on the cell to select the complete row (see figure 398, page 344).




Figure 398 Click the row to select it

With the row selected, press the **[Delete]** key to remove the row from the table. The list of available signals and properties will be updated (see figure 399, page 345).





Figure 399 Deleting the row will trigger the filter to be updated

7.7.3 Calculate signal



This command can be used to calculate signals based on measured or calculated signals and procedure properties, to be used in the data analysis.

The details of the command properties of the **Calculate signal** command are shown in *Figure 400*:




Figure 400 The properties of the Calculate signal command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Signal name:** the name of the calculated signal.
- **Unit:** the unit of the calculated signal. The unit can either be typed in the input field or it can be picked from the drop-down list.
- **Expression:** the mathematical expression used to calculate the signal.
- **Calculate single value:** a toggle provided to force the **Calculate signal** command to return a single value.




Figure 402 Additional functions can be accessed through using the provided button

The **Calculate signal** command automatically parses the mathematical expression in order to discriminate between mathematical operators and the arguments of the operators. In Figure 403, the *Current* string is identified as argument and the **10LOG** and **ABS** strings are identified as mathematical operators.




Figure 403 The Calculate signal automatically differentiates between operators and argument



If an error is made in the mathematical expression, the **Expression** field will be highlighted in red and a tooltip will indicate the nature of the error (see figure 404, page 348).




Figure 404 The expression is automatically tested for errors

After the expression is parsed, the identified arguments are listed in the **Link parameters** sub-panel (see figure 405, page 348).




Figure 405 The identified arguments are listed in the Link parameters sub-panel

Depending on how the **Calculate signal** command is used, the arguments of the mathematical expression can be specified in two different ways:

1. Using the provided dropdown list *Linking arguments using the drop-down list* (see chapter 7.7.3.1, page 349).
2. Using a link *Linking arguments with links* (see chapter 7.7.3.2, page 350).

7.7.3.1 Linking arguments using the drop-down list

If the **Calculate signal** is stacked onto a command that provides linkable properties, as shown in *Figure 406*, the arguments using the mathematical expression of the **Calculate signal** command can be directly linked using the drop-down list provided in the **Link parameter** sub-panel.




Figure 406 Using the Calculate signal when stacked onto another command

The available linkable properties provided by the command on which the **Calculate signal** is stacked are automatically populated in the drop-down list (see *figure 407*, page 350).




Figure 408 Using the Calculate signal in an arbitrary location in the procedure

In the **Edit links** screen, it is possible to link the argument of the **Calculate signal** command to another property in the procedure (see figure 409, page 351).




Figure 409 Linking the argument of the Calculate signal command



NOTICE

More information on linking commands, please refer to *Chapter 10.13*.

7.7.3.3 Mathematical operators

Table 10 provides an overview of the mathematical or logical operators available using a dedicated **button** in the **Calculate signal** editor.

Table 10 *Mathematical and logical operators provided in the Calculate signal editor*



Mathematical operator	Button	Explanation
$\text{abs}(x)$		Determines the absolute value of the argument x
$\text{asinh}(x)$		Determines the inverse hyperbolic sine of the argument x
$\text{asin}(x)$		Determines the inverse sine of the argument x
$\text{sinh}(x)$		Determines the hyperbolic sine of the argument x
$\text{sin}(x)$		Determines the sine of the argument x
$\text{sqrt}(x)$		Determines the square root of the argument x
$\text{fac}(x)$		Determines the factorial of the argument x
$\text{min}(x)$		Determines the minimum value of the argument x
$\text{acosh}(x)$		Determines the inverse hyperbolic cosine of the argument x
$\text{acos}(x)$		Determines the inverse cosine of the argument x
$\text{cosh}(x)$		Determines the hyperbolic cosine of the argument x
$\text{cos}(x)$		Determines the cosine of the argument x
x^y		Raises the argument x to the power of y

Mathematical operator	Button	Explanation
$\max(x)$		Determines the maximum value of the argument x
$\operatorname{atanh}(x)$		Determines the inverse hyperbolic tangent of the argument x
$\operatorname{atan}(x)$		Determines the inverse tangent of the argument x
$\operatorname{tanh}(x)$		Determines the hyperbolic tangent of the argument x
$\tan(x)$		Determines the tangent of the argument x
$\exp(x)$		Determines the exponential function of the argument x
$\operatorname{mean}(x)$		Determines the average value of the argument x
pi		The constant number π
$(x)==(y)$		Determines if the argument x is equal to the argument y
$(x)<(y)$		Determines if the argument x is smaller than the argument y
$(x)<=(y)$		Determines if the argument x is smaller or equal to the argument y
$\ln(x)$		Determines the natural logarithm of the argument x



Mathematical operator	Button	Explanation
stddev(x)		Determines the standard deviation of the argument x
e(x)		Determines the exponential function of the argument x
$(x) \neq (y)$		Determines if the argument x not equal to the argument y
$(x) > (y)$		Determines if the argument x is larger than the argument y
$(x) \geq (y)$		Determines if the argument x is larger or equal to the argument y
log10(x)		Determines the 10 base logarithm of the argument x

7.7.3.4 Additional functions

Table 11 provides an overview of the mathematical or logical operators available using the  button in the **Calculate signal** editor.



NOTICE

When the operators use more than one argument, the arguments need to be separated by a semi-colon (;).

Table 11 Additional functions provided in the Calculate signal editor

Function	Explanation
ACOSEC(x)	Returns the inverse cosecant of the argument x
ACOSECH(x)	Returns the hyperbolic inverse cosecant of the argument x
ACOT(x)	Returns the inverse cotangent of the argument x

Function	Explanation
ACOTH(x)	Returns the hyperbolic inverse cotangent of the argument x
ASEC(x)	Returns the inverse secant of the argument x
ASECH(x)	Returns the hyperbolic inverse secant of the argument x
CEIL(x)	Rounds the argument x to the next available integer
COMPLEXDIV ARG(x1;i1;x2;i2)	Determines the argument of the complex division of $(x1-j\mathbf{i}1)$ by $(x2-j\mathbf{i}2)$
COMPLEXDIV IMAG(x1;i1;x2;i2)	Determines the imaginary part of the complex division of $(x1-j\mathbf{i}1)$ by $(x2-j\mathbf{i}2)$
COMPLEXDIV MOD(x1;i1;x2;i2)	Determines the modulus of the complex division of $(x1-j\mathbf{i}1)$ by $(x2-j\mathbf{i}2)$
COMPLEXDIV REAL(x1;i1;x2;i2)	Determines the real part of the complex division of $(x1-j\mathbf{i}1)$ by $(x2-j\mathbf{i}2)$
COMPLEXMULT ARG(x1;i1;x2;i2)	Determines the argument of the complex multiplication of $(x1-j\mathbf{i}1)$ by $(x2-j\mathbf{i}2)$
COMPLEXMULT IMAG(x1;i1;x2;i2)	Determines the imaginary part of the complex multiplication of $(x1-j\mathbf{i}1)$ by $(x2-j\mathbf{i}2)$
COMPLEXMULT MOD(x1;i1;x2;i2)	Determines the modulus of the complex multiplication of $(x1-j\mathbf{i}1)$ by $(x2-j\mathbf{i}2)$
COMPLEXMULT REAL(x1;i1;x2;i2)	Determines the real part of the complex multiplication of $(x1-j\mathbf{i}1)$ by $(x2-j\mathbf{i}2)$
COSEC(x)	Returns the cosecant of the argument x
COSECH(x)	Returns the hyperbolic cosecant of the argument x



Function	Explanation
COT(x)	Returns the cotangent of the argument x
COTH(x)	Returns the hyperbolic cotangent of the argument x
DEGTORAG(x)	Converts the angle x from degrees to radians
DERIVATIVE(x;y;z)	Returns the z^{th} derivative of the argument x against the argument y
EXPAND(x;y)	Expands argument x by a factor of y
FFT FREQUENCY(x)	Returns the frequency of the Fast Fourier Transform of the argument x
FFT IMAG(x;bool)	Returns the real component of the Fast Fourier Transform of the argument x determined using a normal FFT (bool = 0) or normalized FFT (bool = 1)
FFT REAL(x;bool)	Returns the imaginary component of the Fast Fourier Transform of the argument x determined using a normal FFT (bool = 0) or normalized FFT (bool = 1)
FLOOR(x)	Rounds the argument x to the previous available integer
FPART(x)	Returns the fractional part of the argument x
INDEXER(x)	Indexes the argument x starting at a value of 1
INTEGRATE(x;y)	Returns the integral of the argument x against the argument y
ITEM(x;y)	Returns the y^{th} item of the argument x
LENGTH(x)	Return the length of the argument x

Function	Explanation
NLOG(x)	Returns the natural logarithm of the argument x
RADTODEG(x)	Converts the angle x from radians to degrees
ROUND(x)	Returns the rounded value of the argument x
SAVITZKY GOLAY(x;left;right;order)	Applies the Savitzky Golay smoothing on the argument x , using the specified left and right points and the specified polynomial order
SEC(x)	Returns the secant of the argument x
SECH(x)	Returns the hyperbolic secant of the argument x
SHRINK STANDARD(x;y)	Shrinks the size of the argument x to a new size of y keeping one value every n , where n is x/y
SHRINK DIFFERENTIAL (x;y;z)	Shrinks the size of the argument x to a new size of z using the derivative of the source arguments dy/dx , keeping the z highest derivative values
SHRINK DIFFERENTIAL ORDINATE(x;y;z)	Shrinks the size of the argument y to a new size of z using the derivative of the source arguments dy/dx , keeping the z highest derivative values
SHRINK MEAN(x;y)	Shrinks the size of the argument x to a new size of y using the average value of n values, where n is x/y
SIGNIFICANTS(x;y)	Formats the argument x to y significant digits
SPIKEREJECT(x)	Applies a spike rejection algorithm to the argument x



NOTICE

This command can be used without an Autolab connected to the computer.

The details of the properties of the **Import data** command are shown in *Figure 411*:

The screenshot shows a 'Properties' dialog box for the 'Import data' command. It contains the following fields and controls:

- Command name:** A text input field containing 'Import data'.
- File name:** A text input field that is currently empty, followed by a 'Browse' button.
- Column delimiter:** A drop-down menu with 'Space' selected.
- Decimal separator:** A drop-down menu with '.' selected.
- Number of rows to skip:** A text input field containing '2'.
- More:** A button located at the bottom right of the dialog.

Figure 411 The properties of the *Import data* command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **File name:** the path to the file containing the data to import. A Browse button is provided in order to specify the location of the file through a Windows Explorer dialog.
- **Column delimiter:** a drop-down list that provides the choice of column delimiter (Space, Tab, Comma (,), Semicolon (;) or Colon (:)). This property only applies to ASCII files.
- **Decimal separator:** a drop-down list that provides the choice of decimal separator (Dot (.) or Comma (,)). This property only applies to ASCII files.
- **Number of rows to skip:** defines the number of rows to skip when importing the data. This only applies to ASCII files.

The file name and location can be specified directly or using Browse button. The Windows Explorer dialog provides the means to import GPES, FRA or any type of file (see *figure 412*, page 360).




Figure 412 The file type can be adjusted in the Windows Explorer dialog



NOTICE

The **Import data** command automatically adjusts the properties displayed based on the extension of the specified file.

Additional settings are available by clicking the **More** in the **Properties** panel. A new screen will be displayed as shown as *Figure 411*.




Figure 413 Additional settings are available for the Import data command

Depending on the type of file imported, the following additional settings are available:

- **Columns:** the Columns editor can be used when importing ASCII file to specify the number of columns in the source file, assign a name to each column and specify units for the data in each column, if applicable.

- **Plots:** the Plot editor can be used to specify how the data imported by the file should be displayed. This editor is available for all the file types and the use of the plot editor is explained in *Chapter 9.5*.

To specify the number of columns in the ASCII file import using the **Import Data** command, the Columns editor can be used. By default, two rows are specified in the table and additional rows can be added for additional columns in the table. A signal name can be provided in each text cell in the table and units can be added as well. For example, if the imported file has three columns, with the first one being Time, the second Current and the third one Potential, the Columns editor can be adjusted as shown in *Figure 414*.




Figure 414 Using the Columns editor

Clicking a cell in the table shows a drop-down list with a number of pre-defined signal names or units. If needed, a custom name and unit can be specified by typing directly in the selected cell of the table (see *figure 415*, page 362).





Figure 415 Specifying a custom name



NOTICE

GPES and FRA data files can also be imported directly using the  button located in the **Actions** panel of the **Dashboard Actions** (see chapter 4.1, page 85).

7.7.6 Export data



This command can be used to export data to certain file formats: ASCII, Excel, ZView and RelaxIS file formats are supported.



NOTICE

This command can be used without an Autolab connected to the computer.

The details of the properties of the **Export data** command are shown in Figure 416:

Properties ✕

Export data

Command name

File name

Working folder

File format ▼

Number of columns

Column delimiter ▼

Decimal separator ▼

File mode ▼

Remarks


Write column headers

Figure 416 The properties of the Export data command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **File name:** the name of the exported data file. A button is provided in order to specify the name of the file through a Windows Explorer dialog.
- **Working folder:** the folder where the exported data file will be saved. A button is provided in order to specify the location through a Windows Explorer dialog.
- **File format:** specifies the type of output file (ASCII, Excel, RelaxIS or ZView) using the provided drop-down list.
- **Number of columns:** specifies the number of columns in the output file. This property only applies to ASCII files.
- **Column delimiter:** a drop-down list that provides the choice of column delimiter (Space, Tab, Comma (,), Semicolon (;) or Colon (:)). This property only applies to ASCII files.
- **Decimal separator:** a drop-down list that provides the choice of decimal separator (Dot (.) or Comma (,)). This property only applies to ASCII files.

The details of the command properties of the **Generate index** command are shown in *Figure 417*.



Properties →

Generate index

Command name

Figure 417 The properties of the *Generate index* command

The following properties are available:

- **Command name:** a user-defined name for the command.




NOTICE

The **Generate index** command can not be used stand alone. This command is designed to work in conjunction with another command providing the data to index. The **Generate index** command can be *stacked* onto the command providing the source data *Stacking commands* (see chapter 10.12, page 674).

7.7.8 Shrink data

 <p>Shrink data</p>	<p>This data can be used to shrink the source data to a smaller set.</p>
--	--

The details of the properties of the **Shrink data** command are shown in *Figure 418*:



Properties →

Shrink data

Command name

New number of data points

Shrink method ▼

Figure 418 The properties of the *Shrink data* command



The following properties are available:

- **Command name:** a user-defined name for the command.
- **New number of data points:** specifies the number of data points generated by the shrink algorithm. This value must be smaller or equal than the number of points in the source data. The number of data points in the source data divided by the specified *New number of data points* provides the reduction factor, n .
- **Shrink method** specifies the shrink method, using the provided drop-down list. Three shrink methods are available:
 - **Standard:** using this method, points are removed from the source data without a specific selection argument. This method keep a data point out of every n data points.
 - **Mean:** using this method, the average value of n data points in the source data is determined and stored in the shrunked data.
 - **Differential:** using this method, the selection is based on the differential of the source data, dY/dX . The points with the highest differential are kept while the other points are discarded.



NOTICE

The **Shrink data** command can not be used stand alone. This command is designed to work in conjunction with another command providing source data used by the **Shrink data** command. The **Shrink data** command can *linked* to the command providing the source data *Links* (see chapter 10.13, page 678).



NOTICE

The **Shrink data** needs a *X source* signal and a *Y source* signal.

Additional properties are available by clicking the More in the **Properties** panel. A new screen will be displayed (see figure 419, page 367).




Figure 419 Additional properties are available for Shrink data command

The following additional settings are available:

- **Custom plots:** the Plot editor can be used to specify how the shrunk data should be displayed. The use of the custom plot editor is explained in *Chapter 9.5*.

7.8 Analysis - general commands

Analysis - general commands can be used to perform data analysis on measured data or to integrate data analysis steps in a procedure.

The available commands are represented by a shortcut icon (see figure 420, page 367).




Figure 420 The Analysis - general commands

The following commands are available:

- **Smooth:** a command which can be used to smooth measured data to remove noise or spikes *Smooth* (see chapter 7.8.1, page 368).



- **Peak search:** a command which can be used to find peaks in measured data *Peak search* (see chapter 7.8.2, page 372).
- **Regression:** a command which can be used to perform a regression on measured data *Regression* (see chapter 7.8.3, page 373).
- **Derivative:** a command which can be used to calculate the first derivative of the provided data *Derivative* (see chapter 7.8.4, page 376).
- **Integrate:** a command which can be used to calculate the integral of the provided data *Integrate* (see chapter 7.8.5, page 378).
- **Interpolate:** a command which can be used to determine data points by linear interpolation of measured data *Interpolate* (see chapter 7.8.6, page 379).
- **FFT analysis:** a command which can be used to transform time domain data into frequency domain data by applying a Fast Fourier Transform on the provided data *FFT analysis* (see chapter 7.8.7, page 379).
- **Convolution:** a command which can be used to perform a convolution analysis on the provided data *Convolution* (see chapter 7.8.8, page 381).
- **Calculate charge:** a command which can be used to determine the charge from the measured current *Calculate charge* (see chapter 7.8.9, page 385).
- **Hydrodynamic analysis:** a command which can be used to perform a Levich and Koutecký-Levich analysis on measured data recorded using forced convection, using the Autolab rotating disk electrode (RDE) or the Autolab rotating ring disk electrode (RRDE) *Hydrodynamic analysis* (see chapter 7.8.10, page 386).
- **ECN spectral noise analysis:** a command which can be used to analyze electrochemical noise (ECN) data *ECN spectral noise analysis* (see chapter 7.8.11, page 387).
- **iR drop correction:** a command which can be used to correct measured data for ohmic losses *iR drop correction* (see chapter 7.8.12, page 391).
- **Baseline correction:** a command which can be used to subtract a baseline from the measured data *Baseline correction* (see chapter 7.8.13, page 392).
- **Corrosion rate analysis:** a command which can be used to analyze linear polarization data and determine the corrosion rate *Corrosion rate analysis* (see chapter 7.8.14, page 398).

7.8.1 Smooth

 <p style="text-align: center;">SG Smooth</p>	<p>This command can be used to smooth data and remove spikes.</p>
--	---



The **Smooth** command can be used in two different modes, which can be selected using the provided drop-down list (see figure 421, page 369):

The screenshot shows a 'Properties' dialog box for the 'Smooth' command. The 'Mode' dropdown menu is open, showing two options: 'Savitzky Golay' (the current selection) and 'FFT'. Other settings include 'Command name' (Smooth), 'Spike rejection' (Savitzky Golay), 'Polynomial order' (FFT), 'Smooth level' (Level 2), and 'Number of points left/right' (4). A 'More' button is visible at the bottom right.

Figure 421 Two modes are provided by the Smooth command

1. Savitzky-Golay (SG) smooth (default mode)
2. FFT smooth



NOTICE

The Savitzky-Golay (SG) smoothing method is described in Anal. Chem.,36, 1627 (1964). It involves a polynomial fit through the experimental data. This method is also called weighted moving averaging.



NOTICE

The **Smooth** command description in the procedure editor is dynamically adjusted in function of the specified mode.

7.8.1.1 SG Smooth

The following properties are available when the command is used in the *SG Smooth* mode (see figure 422, page 370):



NOTICE

The **Number of points left/right** defines the size of the weighted moving average.

7.8.1.2 FFT Smooth

The following properties are available when the command is used in the *FFT Smooth* mode (see figure 423, page 371):

The screenshot shows a 'Properties' dialog box for the 'Smooth' command. The title bar says 'Properties' with a close button. Below the title, the word 'Smooth' is displayed. The dialog contains the following fields:

- Command name: A text input field containing 'Smooth'.
- Mode: A dropdown menu currently set to 'FFT'.
- Filter type: A dropdown menu currently set to 'Band pass'.
- Frequency 1: A text input field followed by 'Hz'.
- Frequency 2: A text input field followed by 'Hz'.

A 'More' button is located at the bottom right of the dialog.

Figure 423 *FFT Smooth mode properties*

- **Command name:** a user-defined name for the command.
- **Filter type:** defines the type of FFT filter used by the command, using the provided drop-down list. Four different filter types are available for the *FFT Smooth* mode:
 - **Low pass:** all the contributions from frequencies higher than the user-selected cutoff frequency are rejected. This method can be used to remove high frequency noise from a measurement.
 - **High pass:** all the contributions from frequencies lower than the user-selected cutoff frequency are rejected. This method can be used to remove low frequency noise from a measurement.
 - **Band pass:** only the contributions from frequencies within a user-defined frequency range are kept. All frequencies that fall outside of the user defined range are rejected.
 - **Band stop:** all the contributions from frequencies within a user-defined frequency range are rejected. Only the frequencies that fall outside of the user defined range are kept.
- **Frequency 1:** defines the first frequency limit used by the command, in Hz.
- **Frequency 2:** defines the second frequency limit used by the command, in Hz.


- **Number of points in search window:** this property defines the number of points that must be located above and below a zero crossing of the first derivative of the signal (dY/dX), in order to qualify as a peak. This setting is useful to discriminate between noise and real peaks. The default value is 6.
- **Peak type:** defines the type of peaks to search, forward or reverse, using the provided dropdown list. Using the forward setting, NOVA will search for regular peaks (anodic peak during the positive going scan or cathodic peak in the opposite direction). The reverse setting allows NOVA to search for peaks in the opposite direction.
- **Start X:** defines the initial abscissa used for the peak search.
- **End X:** defines the final abscissa used for the peak search.



NOTICE

When the **Start X** and **End X** properties are not defined, the peaks will be searched in the whole range of X values provided in the source data.

7.8.3 Regression

 <p>Regression</p>	<p>This command can be used to perform a regression on the source data. The source data contains X and Y values.</p>
---	--

The **Regression** command provides three different modes, which can be selected using the provided drop-down list (see figure 425, page 373):

Properties ➔

Regression

Command name

Mode ▼

Use offset

Direction

Start X

End X

Figure 425 The Regression command provides three regression modes


7.8.3.2 Polynomial regression

The following properties are available when the **Regression** command is used in *Polynomial* mode (see figure 427, page 375):

The image shows a 'Properties' dialog box for the 'Regression' command. The title bar says 'Properties' with a close button. Below the title, the word 'Regression' is displayed. The dialog contains several controls: a text input field for 'Command name' containing 'Regression'; a dropdown menu for 'Mode' set to 'Polynomial'; a text input field for 'Polynomial order' containing '1'; a 'Best fit' toggle switch which is currently turned off; a dropdown menu for 'Direction' set to 'All'; and two empty text input fields for 'Start X' and 'End X'.

Figure 427 The properties of the Polynomial mode of the Regression command

- **Command name:** a user-defined name for the command.
- **Polynomial order:** specifies the polynomial order used by the regression.
- **Best fit:** specifies if the best fit should be used in the regression, using the provided toggle. Depending on this toggle, the following equations are used:
 - **Best fit off:** performs a regression using the specified polynomial order.
 - **Best fit on:** performs a regression using all the polynomial functions up to the maximum polynomial order. The regression providing the smallest χ^2 (Chi-squared) is automatically selected by the software.
- **Direction:** specifies the direction to use in the calculation, using the provided dropdown list. Three directions are available:
 - **All:** all the data provided in the source data is used for the regression. This is the default direction.
 - **Forward:** only the data values in the positive going direction is used for the regression.
 - **Reverse:** only the data values in the negative going direction is used for the regression.
- **Start X:** defines the initial abscissa used for the regression.
- **End X:** defines the final abscissa used for the regression.



The details of the command properties of the **Derivative** command are shown in *Figure 429*.




Figure 429 The properties of the Derivative command

The following properties are available:

- **Command name:** a user-defined name for the command.

This command needs to be *linked* to sourced data *Links* (see *chapter 10.13, page 678*). The **Derivative** command provides two input anchoring points and four output anchoring points (see *figure 430, page 377*).





Figure 430 The anchoring points for linking the Derivative command

The command uses the two input signals to calculate the derivative of the second signal versus the first signal.



7.8.5 Integrate

	<p>This command can be used to integrate a curve and determine the area.</p>
---	--

The details of the command properties of the **Integrate** command are shown in *Figure 431*.




Figure 431 The properties of the Integrate command

The following properties are available:

- **Command name:** a user-defined name for the command.

This command needs to be *linked* to sourced data *Links* (see chapter 10.13, page 678). The **Integrate** command provides two input anchoring points and four output anchoring points (see figure 432, page 378).





Figure 432 The anchoring points for linking the Integrate command

The command uses the two input signals to calculate the integral of the second signal versus the first signal.

7.8.6 Interpolate



The details of the command properties of the **Interpolate** command are shown in *Figure 433*.





Figure 433 The properties of the Interpolate command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Search:** a drop-down list allowing the selection of the value to search for (Y value or X value). By default, the command searches for a Y value.
- **Search at X value/Search at Y value:** the location used by the **Interpolate** command. This property is automatically adjusted depending on the **Search** drop-down list.

7.8.7 FFT analysis



The details of the command properties of the **FFT analysis** command are shown in *Figure 434*.




Figure 434 The properties of the FFT analysis command

The following properties are available:

- **Command name:** a user-defined name for the command.

This command needs to be *linked* to sourced data *Links* (see chapter 10.13, page 678). The **FFT analysis** command provides two input anchoring points and six output anchoring points (see figure 435, page 380).




Figure 435 The anchoring points for linking the FFT analysis command

The command uses the two input signals to transform the time domain data to frequency domain data. The frequency, amplitude as well as the real and imaginary parts of the amplitude are returned.



CAUTION

The **FFT analysis** is intended to be used on source data formatted with the *Time* signal the X data. When another signal is used, the **FFT analysis** command will be executed but the *Frequency* signal calculated by the command will no longer be an actual frequency.

7.8.8 Convolution

 <p>Convolution</p>	<p>This command can be used to apply the convolution analysis on the measured data or to integrate this analysis technique in the procedure.</p>
--	--



NOTICE

This command can only be used on measurements containing the *Time* and *WE(1).Current* signals.

The **Convolution** command can be used in six different modes, which can be selected using the provided drop-down list (see figure 436, page 381):




Figure 436 Six modes are provided by the Convolution command

1. Time semi-derivative (default mode)
2. Time semi-integral
3. GO differintegration
4. FRLT differintegration
5. Spherical convolution
6. Kinetic convolution

The time semi-integral algorithm uses a semi-integral transformation of a time dependent function, $f(t)$, according to:

$$\frac{d^{-1/2}}{dt^{-1/2}} f(t)$$

7.8.8.3 G0 differintegration (Grünwald-0)

The following properties are available when the command is used in the *G0 differintegration* mode (see figure 439, page 383):

The image shows a software interface titled 'Properties' with a close button in the top right. Under the heading 'Convolution', there are three input fields: 'Command name' containing 'Convolution', 'Type' set to 'G0 differintegration' with a dropdown arrow, and 'Order' containing '-0,5'. A 'More' button is located at the bottom right of the dialog.

Figure 439 *G0 differintegration mode properties*

- **Command name:** a user-defined name for the command.
- **Order:** the order used in the G0 differintegration algorithm (default: -0.5).

The G0 differintegration algorithm can be used to carry out differintegration to any user-defined order. Specific *Order* values provide a mathematical equivalence with other transformations:

- For an *Order* value of 1, the operation is the equivalent of a derivative.
- For an *Order* value of -1, the operation is the equivalent of an integration.
- For an *Order* value of 0.5, the operation is the equivalent of a time semi-derivative method.
- For an *Order* value of -0.5, the operation is the equivalent of a time semi-integral method.

Error in results increases with the length of the interval and accumulates, i.e. error in latter points is larger than in earlier ones. Important advantage is that this algorithm does not require the value of the function for $t = 0$, which makes it very well suited for transformation of chronoamperometric data, where $i_{t \rightarrow 0} = 0$. The disadvantage of the algorithm is that the total number of operations is proportional to the square of the number of data points, so calculation time grows fast with the length of the data set. The fundamentals of this algorithm are described in Oldham KB, J. *Electroanal. Chem.* 121 (1981) 341-342.

- **Electrode radius:** the radius of the electrode, in cm.
- **Diffusion coefficient:** the diffusion coefficient, in cm^2/s .

is used to carry out convolution of the data measured using a spherical electrode and staircase potential waveform. Values of the diffusion coefficient and the electrode radius are necessary. Details of the algorithm can be found in S.O. Engblom, K.B. Oldham, *Anal. Chem.* 62 (1990) 625-630.

7.8.8.6 Kinetic

The following properties are available when the command is used in the *Kinetic convolution mode* (see figure 441, page 384):

The screenshot shows a 'Properties' dialog box with a 'Convolution' section. It contains three input fields: 'Command name' with the value 'Convolution', 'Type' with a dropdown menu set to 'Kinetic', and 'Rate constant' with the value '0,001' and a unit 's⁻¹' to its right. A 'More' button is located at the bottom right of the dialog.

Figure 442 Kinetic convolution mode properties

- **Command name:** a user-defined name for the command.
- **Rate constant:** the rate constant of the chemical reaction, in s^{-1} .

The kinetic convolution algorithm carries out kinetic convolution according to F.E. Woodard, R.D. Goodin, P.J. Kinlen, *Anal. Chem.* 56 (1984) 1920-1923. This convolution requires the value of the rate constant of irreversible homogeneous follow-up reaction (EC_i mechanism).

7.8.9 Calculate charge

	<p>This command can be used to calculate the charge by integrating the measured current against time. The total charge is reported in Coulomb (C).</p>
--	--

The details of the command properties of the **Calculate charge** command are shown in Figure 443.





Figure 443 The properties of the Calculate charge command



NOTICE

This command can only be used on measurements containing the *Time* and *WE(1).Current* signals.

7.8.10 Hydrodynamic analysis



This command can be used to perform a Levich and Koutecký-Levich analysis on hydrodynamic data.



NOTICE

This command is intended to be used in combination with a rotating disk or rotating ring disk electrode, controlled by NOVA.

The details of the command properties of the **Hydrodynamic analysis** command are shown in *Figure 444*:




Figure 444 The property of the Hydrodynamic analysis command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Current index:** the index of the current value used in the by the command. The index of a current value located in the mass transport-limited region should be specified for the Levich analysis and the index of a current value located in the mixed kinetic-mass transported region should be specified for the Koutecký-Levich analysis. For all rotation rates used in the procedure, the current value at the specified index will be used.

The rotation of the electrode creates a convective drag from the bulk of the solution towards the surface of the electrode, resulting in a mixed control of mass transport, involving a convective part which depends on the square root of the angular frequency of the electrode and diffusion layer which also depends on this property. Under these experimental conditions, the limiting current values, i_l and kinetic current i_k , are related to the rotation rate of the working electrode according to the Levich equation and Koutecký-Levich equation:

$$i_l = 0.62 \cdot nFD^{2/3} \nu^{-1/6} C^\infty \sqrt{\omega}$$

$$\frac{1}{i} = \frac{1}{i_k} + \frac{1}{0.62 \cdot nFD^{2/3} \nu^{-1/6} C^\infty \sqrt{\omega}}$$


Where A is the geometric area of the electrode, in cm^2 , n is the number of electrons involved in the electrochemical reaction, F is the Faraday constant, D is the diffusion coefficient of the electroactive species, in cm^2/s , ν is the kinematic viscosity in cm^2/s and $\sqrt{\omega}$ is the square root of the angular frequency of the rotating electrode, in $(\text{rad/s})^{1/2}$.



NOTICE

The **Hydrodynamic analysis** command automatically carries out two linear regressions using the **Regression** command.

7.8.11 ECN spectral noise analysis



This command can be used to analyze electrochemical noise measurements (ECN).

The **ECN spectral noise analysis** command can be used in two different modes, which can be selected using the provided drop-down list (see figure 445, page 388):




Figure 445 Two modes are provided by the ECN spectral noise analysis command

- **FFT:** a spectral noise analysis that uses the *Fast Fourier Transform* method.
- **MEM:** a spectral noise analysis that uses the *Maximum Entropy* method.



NOTICE

This command can only be used on measurements containing the *Time*, *ECN(1).Potential* and *WE(1).Current* signals.

Electrochemical noise data is generally analyzed by computing the spectral density of the measured data. This can be achieved by transforming the time domain information to a frequency domain spectrum, using the **Fast Fourier Transformation (FFT)** or the **Maximum Entropy Method (MEM)**.

Traditional time domain to frequency domain transformation assumes that the data outside of the measured time segment is either zero or that the data in this segment repeats periodically. This hypothesis is not valid for electrochemical noise data. In order to satisfy these requirements and to avoid edge effects in the data, it is common practice to apply a **window function** on the time domain data. This calculation involves the multiplication of the time domain data by a function which is zero at the extremes of the time domain data and rises smoothly to unity value in its center.

Alongside the power spectra determined by the transformation of the data into the frequency domain, the ECN spectral noise analysis also calculates the following statistical indicators:

- Noise resistance, R_n .
- Pitting index (or localization index), PI

- Current and potential skewness
- Current and potential kurtosis

The noise resistance, R_n , is given by:

$$R_n = \frac{\sigma_v}{\sigma_i}$$

Where σ_v and σ_i are the standard deviations of the measured potential and current, respectively. The value of the noise resistance is reported in Ohm.

The pitting index, or localization index, PI, is given by:

$$PI = \frac{\sigma_i}{i_{RMS}} = \sqrt{\frac{\sum_{j=1}^N (i_j - \bar{i})^2}{\sum_{j=1}^N i_j^2}}$$

Where i_{RMS} is the root mean squared value of the measured current. The pitting index can be between 0 and 1. A value close to 0 is observed for systems in which the measured current values show only small deviation with respect to the average current value. On the other hand, the pitting index will be close to 1 when the individual current values are significantly deviating from the average current value. This value is therefore an indication of the distribution of the current values recorded during an electrochemical noise experiment.

Skewness and kurtosis are additional indicators calculated according to the following equations, respectively:

$$\frac{1}{N} \sum_{j=1}^N \left(\frac{X_j - \bar{X}}{\sigma} \right)^3$$

$$\frac{1}{N} \sum_{j=1}^N \left(\frac{X_j - \bar{X}}{\sigma} \right)^4$$

7.8.11.1 Fast Fourier Transform (FFT)

The following properties are available when the command is used in the FFT mode (see figure 446, page 390):

Properties 

ECN spectral noise analysis

Command name

Mode

MEM coefficients

Subtract baseline

Window function

Figure 447 MEM mode properties

- **Command name:** a user-defined name for the command.
- **MEM coefficients:** specifies the number of coefficients to be used in the MEM algorithm. The default value is 20.
- **Subtract baseline:** a toggle which can be used to enable or disable baseline subtraction. When this property is enabled, a linear regression is used to subtract the baseline from the measured potential or current values.
- **Window function:** defines the type of windowing function used for the MEM algorithm, using the provided drop-down list. The default function is the Square function.



NOTICE

For a detailed description of the *Window functions* used in NOVA, the reader is invited to refer to W. H. Press, S. A. Teukolsky, W. T. Vetterling, B. P. Flannery, Numerical Recipes – The Art of Scientific Computing, 3rd edition, Cambridge University Press, 2007.

7.8.12 iR drop correction

 <p>iR drop correction</p>	<p>This command can be used to correct measured data for ohmic drop losses.</p>
---	---

The command properties of the **iR drop correction** command are shown in *Figure 448*:




Figure 449 Four modes are provided by the Baseline correction command

1. Linear (default mode)
2. Polynomial
3. Exponential
4. Moving average



NOTICE

The **Baseline correction** command description in the procedure editor is dynamically adjusted in function of the specified mode.

For the *Linear*, *Polynomial* and *Exponential* mode, the points defining the location of the baseline used in the correction can be specified using the More button (see figure 450, page 393).




Figure 450 The points defining the baseline are specified in a dedicated editor



A table will be displayed in a new screen (see figure 451, page 394).




Figure 451 The selected point table

Using the provided editor, it is possible to define the location of the two or more points. The location of each point is defined by specifying a X and Y coordinate (see figure 452, page 394).




Figure 452 Specifying the baseline points

Depending on mode of the **Baseline correction**, the following number of points need to be defined:

- **Linear mode:** two or more points are required to define a linear baseline.
- **Polynomial:** $n+1$ or more points are required to define a polynomial based of order n .
- **Exponential:** two or more points are required to define an exponential baseline.

The *Snap to data* toggle can be used to force the baseline to be snapped to the nearest data point in the source data. If this toggle is on, then the Y coordinate will be ignored and the data point nearest to the specified X abscissa will be used.

Clicking the button clears the whole table.

7.8.13.1 Linear

The following properties are available when the command is used in the *Linear* mode (see figure 453, page 395):

The screenshot shows a 'Properties' dialog box with a close button (X) in the top right corner. Below the title, the text 'Baseline correction' is displayed. There are two input fields: 'Command name' with the text 'Baseline correction' and 'Mode' with a dropdown menu showing 'Linear'. A 'More' button is located at the bottom right of the dialog.

Figure 453 The Linear properties

- **Command name:** a user-defined name for the command.

The points defining the location of the baseline used in the correction can be specified using the button (see figure 450, page 393).

7.8.13.2 Polynomial

The following properties are available when the command is used in the *Polynomial* mode (see figure 454, page 396):

- **Use offset:** defines if an offset should be used in the exponential baseline correction, using the provided toggle. Depending on this toggle, the following equations are used:
 - **Use offset off:** performs a baseline correction using the equation $y = be^{cx}$.
 - **Use offset on:** performs a baseline correction using the equation $y = a + be^{cx}$.

The points defining the location of the baseline used in the correction can be specified using the button (see figure 450, page 393).

7.8.13.4 Moving average

The following properties are available when the command is used in the *Moving average* mode Figure 456:

Figure 456 The Moving average properties

- **Command name:** a user-defined name for the command.
- **Window size:** defines the number of points in the moving average window (default: 2).

The moving average baseline correction performs the following steps:

1. The source data is grouped into segments of n points; where n corresponds to the **Window size** property.
2. The average value of each segment is calculated.
3. The source data is reduced from m data points to m/n averages.
4. Each i^{th} average value is compared to the average value of its immediate neighboring values, at $i-1$ and $i+1$.
 - a. For positive going sweeps, if the i^{th} average value is higher than the average value of the averages at $i-1$ and $i+1$, then the i^{th} average value is replaced by the average value of the averages at $i-1$ and $i+1$.
 - b. For negative going sweeps, if the i^{th} average value is lower than the average value of the averages at $i-1$ and $i+1$, then the i^{th} average value is replaced by the average value of the averages at $i-1$ and $i+1$.

7.8.14.1 Tafel Analysis

The following properties are available when the command is used in the *Tafel Analysis* mode (see figure 458, page 399):

The screenshot shows a 'Properties' dialog box with the following fields and values:

- Command name: Corrosion rate analysis
- Mode: Tafel Analysis
- Density: 7,86 g/cm³
- Equivalent weight: 27,925 g/mol
- Surface area: 1 cm²
- Perform fit: (off)

A 'More' button is located at the bottom right of the dialog.

Figure 458 Tafel Analysis mode properties


- **Command name:** a user-defined name for the command.
- **Density:** specifies the density of the sample in g/cm³.
- **Equivalent weight:** defines the equivalent weight of the sample in g/mol of exchanged electrons.
- **Surface area:** defines the area of the sample, in cm².
- **Perform fit:** specifies if the data should be fitted, using the provided toggle. Depending on this toggle, the corrosion rate analysis command carries out the following analysis:
 - **Perform fit off:** only the Tafel slopes are determined and the corrosion rate analysis is carried out using the results of the Tafel slope analysis.
 - **Perform fit on:** the Tafel slope analysis is carried out and the data is fitted using the **Butler-Volmer** equation. The corrosion rate analysis is carried out based on the results of the fit.

The Butler-Volmer equation is given by:

$$i = i_{\text{corr}} \left(e^{\frac{2,303(E-E_{\text{corr}})}{b_a}} - e^{-\frac{2,303(E-E_{\text{corr}})}{b_c}} \right)$$

Where i is the measured current, i_{corr} is the corrosion exchange current, E is the applied potential, E_{corr} is the corrosion potential and b_a and b_c are the Tafel slopes, in V/decade, respectively.

The **Corrosion rate analysis** command requires the definition of four points, defining the location of the linear parts of the anodic and cathodic



← Corrosion rate analysis

Selected points

X	Y
-0,15	0
-0,25	0
-0,4	0
-0,5	0

Reset

Figure 461 Specifying the four points



NOTICE

If needed, the location of these points can be finetuned after the measurement is finished *Corrosion rate analysis* (see chapter 12.8, page 796).



NOTICE

NOVA will automatically select the data points closest to the specified points in the table when the command is executed.



NOTICE

The Y coordinates can be set to 0.

Clicking the  button clears the whole table.

7.8.14.2 Polarization Resistance

The following properties are available when the command is used in the *Polarization Resistance* mode (see figure 462, page 402):

Properties ➔

Corrosion rate analysis

Command name

Mode ▼

Density g/cm³

Equivalent weight g/mol

Surface area cm²

|ba| V/dec

|bc| V/dec

Range mV

Figure 462 Polarization resistance mode properties

- **Command name:** a user-defined name for the command.
- **Density:** specifies the density of the sample in g/cm³.
- **Equivalent weight:** defines the equivalent weight of the sample in g/mol of exchanged electrons.
- **Surface area:** defines the area of the sample, in cm².
- **|ba|:** defines the absolute value of the anodic Tafel slope value, in V/decade of current.
- **|bc|:** defines the absolute value of the cathodic Tafel slope value, in V/decade of current.
- **Range:** defines the potential range, in mV, around the observed corrosion potential, in which the analysis is carried out. The specified value will be used on both sides of the corrosion potential.

The *Polarization Resistance* mode uses the Stern-Geary equation to determine the corrosion current, i_{corr} , according to:

$$i_{corr} = \frac{|b_a| \cdot |b_c|}{2.303(|b_a| + |b_c|) R_p}$$

Where b_a and b_c are the specified Tafel slopes, in absolute value, and R_p is the inverted slope of the linear regression carried out in the specified Range around the observed corrosion potential.



NOTICE

The Polarization Resistance analysis method is based on M. Stern, A. L. Geary, *J ECS* Vol. 104, No. 1, 56-63, 1957.

7.9 Analysis - impedance

Analysis - impedance commands can be used to perform data analysis on measured impedance data or to integrate data analysis steps in a procedure.

The available commands are represented by a shortcut icon (see figure 463, page 403).




Figure 463 The Analysis - impedance commands

The following commands are available:

- **Electrochemical circle fit:** a command which can be used to quickly fit a semi-circle in a Nyquist plot using a $R(RQ)$ equivalent circuit *Electrochemical circle fit* (see chapter 7.9.1, page 403).
- **Fit and simulation:** a command which can be used to fit measured impedance data with a user-defined equivalent circuit *Fit and simulation* (see chapter 7.9.2, page 405).
- **Kramer-Kronig test:** a command which can be used to perform the Kramer-Kronig test on measured impedance data *Kronig-Kramers test* (see chapter 7.9.3, page 449).
- **Include all FRA data:** a command that can be used to calculate additional values from the measured impedance *Include all FRA data* (see chapter 7.9.4, page 452).
- **Potential scan FRA data:** a command that can be used to calculate values from measured potential scan FRA data to perform a Mott-Schottky analysis *Potential scan FRA data* (see chapter 7.9.5, page 453).

7.9.1 Electrochemical circle fit



 <p>Electrochemical circle fit</p>	<p>This tool can be used to fit a semi-circle in a Nyquist plot with a R(RQ) equivalent circuit using three user defined points.</p>
---	--

The details of the command properties of the **Electrochemical circle fit** command are shown in *Figure 464*.

Properties ↔

Electrochemical circle fit

Command name

Figure 464 The properties of the Electrochemical circle fit command

The following properties are available:

- **Command name:** a user-defined name for the command.

In order to use this command, three or more point defining the location of the semi-circle in the Nyquist plot need to be defined. These points can be specified using the button (see *figure 464, page 404*).

A table will be displayed in a new screen (see *figure 465, page 404*).

← **Electrochemical circle fit**

Selected points

Snap to data

X	Y

Figure 465 The selected point table

Using the provided editor, it is possible to define the location of three or more points. The location of each point is defined by specifying a X and Y coordinate (see *figure 466, page 405*).





Figure 466 Specifying the points to define the semi-circle

The *Snap to data* toggle can be used to force the points defining the semi-circle to be snapped to the nearest data point in the source data. If this toggle is on, then the Y coordinate will be ignored and the data point nearest to the specified X coordinate.

Clicking the button clears the whole table.

7.9.2 Fit and simulation

	<p>This command allows to use the fit and simulation analysis tool. Measured data can be fitted (or simulated) using a pre-defined equivalent circuit. The equivalent circuit is defined using the dedicated Circuit editor.</p>
---	--

The details of the properties of the **Fit and simulation** command are shown in *Figure 467*:

- **Use weight factor:** a toggle which can be used to define whether a weight factor should be used during the calculation. If weight factors are used, each point is multiplied by a weight factor equal to the inverse of the square of the impedance modulus. If this option is not used, the weight factor is the same for each point, i.e. the inverse of the square root of the average of the impedance modulus.
- **Fit or simulation:** defines the calculation method using the provided drop-down list. Using the Fit method, the software will try to find the most suitable values for the parameters of each element defined in the equivalent circuit, starting with initial, user-defined values. The simulation method simply calculates the impedance values for the equivalent circuit, as it is defined by the user.
- **Measurement data format:** defines the type of data of the data, using the provided drop-down list.

The equivalent circuit can be specified either by typing a string in the **Circuit description** field, as shown in *Figure 468*.

The screenshot shows a 'Properties' dialog box with a 'Fit and simulation' section. The 'Circuit description' field is active and contains the text 'R(RQ)'. To the right of this field is an 'Edit' button. Below the 'Circuit description' field are several other input fields: 'Command name' (Fit and simulation), 'Maximum number of iterations' (300), 'Maximum change in χ^2 (scaled)' (0,001), 'Max iterations without improvement' (50), 'Fitting style' (Impedance), 'Use weight factor' (checked), 'Fit or Simulation' (Fit), and 'Measurement data format' (Impedance). A 'More' button is located at the bottom right of the dialog box.

Figure 468 Editing the circuit description

Alternatively, it is possible to define the equivalent circuit in a dedicated editor, by clicking the button (see *figure 469*, page 408).




Figure 470 The Equivalent circuit editor

Detailed analysis of the data obtained during an electrochemical impedance measurement is usually performed by fitting the experimental data with an equivalent circuit, based on the Boukamp model. Many circuit elements can be used to fit the experimental data with a model. However, the equivalent circuit must be constructed carefully, since a given experimental data set can be fitted with more than one unique equivalent circuit.

The following tasks can be carried out in the Equivalent circuit editor:

1. Drawing the equivalent circuit using individual circuit elements
2. Generate an equivalent circuit from a CDC string
3. Loading a pre-defined equivalent circuit
4. Importing and exporting equivalent circuits
5. Advanced editing
6. Edit element properties
7. Creating linkable properties
8. Save circuit to Library



7.9.2.1 Circuit elements

The **Fit and simulation** command allows the definition of an equivalent circuit using the elements shown in *Table 12*.

Table 12 Overview of the available equivalent circuit elements


Element	Symbol
R, resistance	
C, capacitance	
L, inductance	
Q, constant phase element	
W, Warburg impedance	
O, Warburg - short circuit terminus	
T, Warburg - open circuit terminus	
G, Gerischer impedance	
B2, Bisquert #2	

All of the circuit elements are fitted with one input connection and one output connection.

These circuit elements can be arranged in series or in parallel, using the connectors shown in *Table 13*.

Table 13 Overview of the available connectors

Connector	Symbol
Serial	
Parallel split	

Connector	Symbol
Parallel join	

All of the connectors are fitted with one of or more input connections and one or more output connections.

7.9.2.1.1 Resistance, R

The resistance circuit element is represented by the letter **R** and identified by the following symbol:



This element is used to typically represent solution resistance or charge transfer resistance.

The impedance of the resistance is provided by:

$$Z_R = R$$

The properties of the R element are shown in *Figure 471*.




Figure 471 The properties of the R element

The following properties are available:

- **Start:** the start value of the resistance, in Ohm. The default value is 500 Ohm.
- **Fitted:** the fitted value of the resistance, in Ohm. This value is only available after the **Fit and simulation** has been executed.
- **Fixed:** specifies if the value can be modified by the **Fit and simulation** command, using the provided checkbox. Fixed properties are shown in red in the Equivalent circuit editor.
- **Min:** specifies the minimum value for the resistance, in Ohm. The default value is -1 TOhm.



- **Max:** specifies the maximum value for the resistance, in Ohm. The default value is 1 TOhm.
- **Apply limits:** specifies if the Min. and Max. limits should be used by the **Fit and simulation** command, using the provided checkbox. When this property is on, the value of the resistance will be kept between the specified Min. and Max., otherwise, the value will be allowed to take any possible value.
- **Input:** creates an input anchoring point for linking purposes, using the specified checkbox.
- **Output:** creates an output anchoring point for linking purposes, using the specified checkbox.



NOTICE

In order to create input and output anchoring points for linking, a **unique** name must be specified for the element. Please refer to *Chapter 7.9.2.7* for more information.

7.9.2.1.2 Capacitance, C

The capacitance circuit element is represented by the letter **C** and identified by the following symbol:



This element is used to typically represent double layer capacitance of the electrochemical interface.

The impedance of the capacitance is provided by:

$$Z_c = \frac{-j}{\omega C}$$

The properties of the C element are shown in *Figure 472*.

C (F)

Start

Fitted

Fixed

Min

Max

Apply limits

Input

Output

Figure 472 The properties of the C element

The following properties are available:

- **Start:** the start value of the capacitance, in F. The default value is 1 μ F.
- **Fitted:** the fitted value of the capacitance, in F. This value is only available after the **Fit and simulation** has been executed.
- **Fixed:** specifies if the value can be modified by the **Fit and simulation** command, using the provided checkbox. Fixed properties are shown in red in the Equivalent circuit editor.
- **Min:** specifies the minimum value for the capacitance, in F. The default value is 1 pF.
- **Max:** specifies the maximum value for the capacitance, in F. The default value is 100 kF.
- **Apply limits:** specifies if the Min. and Max. limits should be used by the **Fit and simulation** command, using the provided checkbox. When this property is on, the value of the capacitance will be kept between the specified Min. and Max., otherwise, the value will be allowed to take any possible value.
- **Input:** creates an input anchoring point for linking purposes, using the specified checkbox.
- **Output:** creates an output anchoring point for linking purposes, using the specified checkbox.



NOTICE

In order to create input and output anchoring points for linking, a **unique** name must be specified for the element. Please refer to *Chapter 7.9.2.7* for more information.

7.9.2.1.3 Inductance, L

The inductance circuit element is represented by the letter **L** and identified by the following symbol:



This element is used to typically represent adsorption process on the electrochemical interface.


The impedance of the inductance is provided by:

$$Z_L = j\omega L$$

The properties of the L element are shown in *Figure 473*.

7.9.2.1.4 Constant phase element, Q

The constant phase element circuit element is represented by the letter **Q** and identified by the following symbol:



This element is used to typically represent the non-ideal behavior of the electrochemical double layer.

The impedance of the constant phase element is provided by:

$$Z_Q = \frac{1}{Y_0(j\omega)^n}$$

The properties of the Q element are shown in *Figure 474*.

Y0 (Mho)

Start: 1E-06

Fitted: [] [↕]

Fixed

Min: 1E-15

Max: 100000

Apply limits

Input

Output

N

Start: 1

Fitted: [] [↕]

Fixed

Min: 0

Max: 1

Apply limits

Input

Output

Figure 474 The properties of the Q element

The Q element is defined by two values:

- **Y0**: the admittance value, in Mho.
- **n**: the exponent used in the expression of the constant phase element.

The following **specific** properties are available for the **Y0** value:



- **Start:** the start value of the admittance, **Y0**, of the constant phase element, in Mho. The default value is 1 μ Mho.
- **Fitted:** the fitted value of the admittance, **Y0**, constant phase element, in Mho. This value is only available after the **Fit and simulation** has been executed.
- **Min:** specifies the minimum value for the admittance, **Y0**, of the constant phase element, in Mho. The default value is 1 fMho.
- **Max:** specifies the maximum value for the admittance, **Y0**, of the constant phase element, in Mho. The default value is 100 kMho.

The following **specific** properties are available for the **n** value:

- **Start:** the start value of the exponent, **n**, of the constant phase element. The default value is 1.
- **Fitted:** the fitted value of the exponent, **n**, of the constant phase element. This value is only available after the **Fit and simulation** has been executed.
- **Min:** specifies the minimum value for the exponent, **n**, of the constant phase element. The default value is 0.
- **Max:** specifies the maximum value for the exponent, **n**, of the constant phase element. The default value is 1.

The following **common** properties are available:

- **Fixed:** specifies if the value can be modified by the **Fit and simulation** command, using the provided checkbox. Fixed properties are shown in red in the Equivalent circuit editor.
- **Apply limits:** specifies if the Min. and Max. limits should be used by the **Fit and simulation** command, using the provided checkbox. When this property is on, the value will be kept between the specified Min. and Max., otherwise, the value will be allowed to take any possible value.
- **Input:** creates an input anchoring point for linking purposes, using the specified checkbox.
- **Output:** creates an output anchoring point for linking purposes, using the specified checkbox.



NOTICE

In order to create input and output anchoring points for linking, a **unique** name must be specified for the element. Please refer to *Chapter 7.9.2.7* for more information.

7.9.2.1.5 Warburg, W

The Warburg circuit element is represented by the letter **W** and identified by the following symbol:



This element is used to typically represent the semi-infinite diffusion of electroactive species.

The impedance of the Warburg is provided by:

$$Z_w = \frac{1}{Y_0 \sqrt{j\omega}}$$

The properties of the W element are shown in *Figure 475*.

Figure 475 The properties of the W element

The following properties are available:

- **Start:** the start value of the admittance, in Mho. The default value is 100 mMho.
- **Fitted:** the fitted value of the admittance, in Mho. This value is only available after the **Fit and simulation** has been executed.
- **Fixed:** specifies if the value can be modified by the **Fit and simulation** command, using the provided checkbox. Fixed properties are shown in red in the Equivalent circuit editor.
- **Min:** specifies the minimum value for the admittance, in Mho. The default value is 1 pMho.
- **Max:** specifies the maximum value for the admittance, in Mho. The default value is 1 TMho.
- **Apply limits:** specifies if the Min. and Max. limits should be used by the **Fit and simulation** command, using the provided checkbox. When this property is on, the value of the admittance will be kept between the specified Min. and Max., otherwise, the value will be allowed to take any possible value.

The figure shows two screenshots of software property dialog boxes. The top dialog is titled "Y0 (Mho)" and contains the following fields and options: "Start" (0,001), "Fitted" (empty), "Fixed" (checkbox), "Min" (1E-15), "Max" (1000), "Apply limits" (checked), "Input" (checkbox), and "Output" (checkbox). The bottom dialog is titled "B" and contains the following fields and options: "Start" (0,1), "Fitted" (empty), "Fixed" (checkbox), "Min" (1E-06), "Max" (1000), "Apply limits" (checked), "Input" (checkbox), and "Output" (checkbox).

Figure 476 The properties of the O element

The O element is defined by two values:

- **Y0**: the admittance value, in Mho.
- **B**: the factor associated with the thickness of the diffusion layer.

The following **specific** properties are available for the **Y0** value:

- **Start**: the start value of the admittance, **Y0**, of the O element, in Mho. The default value is 1 mMho.
- **Fitted**: the fitted value of the admittance, **Y0**, of the O element, in Mho. This value is only available after the **Fit and simulation** has been executed.
- **Min**: specifies the minimum value for the admittance, **Y0**, of the O element, in Mho. The default value is 1 fMho.
- **Max**: specifies the maximum value for the admittance, **Y0**, of the O element, in Mho. The default value is 1 kMho.

The following **specific** properties are available for the **B** value:

- **Start**: the start value of the thickness factor, **B**, of the O element. The default value is 0.1.

Figure 477 The properties of the T element

The T element is defined by two values:

- **Y0**: the admittance value, in Mho.
- **B**: the factor associated with the thickness of the diffusion layer.

The following **specific** properties are available for the **Y0** value:

- **Start**: the start value of the admittance, **Y0**, of the T element, in Mho. The default value is 1 mMho.
- **Fitted**: the fitted value of the admittance, **Y0**, T element, in Mho. This value is only available after the **Fit and simulation** has been executed.
- **Min**: specifies the minimum value for the admittance, **Y0**, of the T element, in Mho. The default value is 1 fMho.
- **Max**: specifies the maximum value for the admittance, **Y0**, of the T element, in Mho. The default value is 1 kMho.

The following **specific** properties are available for the **B** value:

- **Start**: the start value of the thickness factor, **B**, of the T element. The default value is 0.1.

The image shows two panels from the NOVA software interface. The top panel is titled 'Ka' and contains the following controls: a 'Start' field with the value '0,5', a 'Fitted' field with an empty box and a right-pointing arrow button, a 'Fixed' checkbox (unchecked), a 'Min' field with '1E-06', a 'Max' field with '1000', an 'Apply limits' checkbox (checked), an 'Input' checkbox (unchecked), and an 'Output' checkbox (unchecked). The bottom panel is titled 'Y0 (Mho)' and contains: a 'Start' field with '0,001', a 'Fitted' field with an empty box and a right-pointing arrow button, a 'Fixed' checkbox (unchecked), a 'Min' field with '1E-15', a 'Max' field with '1000', an 'Apply limits' checkbox (checked), an 'Input' checkbox (unchecked), and an 'Output' checkbox (unchecked).

Figure 478 The properties of the G element

The G element is defined by two values:

- **Ka**: the kinetic constant of the chemical reaction.
- **Y0**: the admittance value, in Mho.

The following **specific** properties are available for the **Ka** value:

- **Start**: the start value of the kinetic constant, **Ka**, of the Gerischer. The default value is 0.5.
- **Fitted**: the fitted value of the kinetic constant, **Ka**, of the Gerischer. This value is only available after the **Fit and simulation** has been executed.
- **Min**: specifies the minimum value for the kinetic constant, **Ka**, of the Gerischer. The default value is 1 μ .
- **Max**: specifies the maximum value for the kinetic constant, **Ka**, of the Gerischer. The default value is 1000.

The following **specific** properties are available for the **Y0** value:

- **Start**: the start value of the admittance, **Y0**, of the Gerischer, in Mho. The default value is 1 mMho.



- **Fitted:** the fitted value of the admittance, **Y0**, Gerischer, in Mho. This value is only available after the **Fit and simulation** has been executed.
- **Min:** specifies the minimum value for the admittance, **Y0**, of the Gerischer, in Mho. The default value is 1 fMho.
- **Max:** specifies the maximum value for the admittance, **Y0**, of the Gerischer, in Mho. The default value is 1 kMho.

The following **common** properties are available:

- **Fixed:** specifies if the value can be modified by the **Fit and simulation** command, using the provided checkbox. Fixed properties are shown in red in the Equivalent circuit editor.
- **Apply limits:** specifies if the Min. and Max. limits should be used by the **Fit and simulation** command, using the provided checkbox. When this property is on, the value will be kept between the specified Min. and Max., otherwise, the value will be allowed to take any possible value.
- **Input:** creates an input anchoring point for linking purposes, using the specified checkbox.
- **Output:** creates an output anchoring point for linking purposes, using the specified checkbox.




NOTICE

In order to create input and output anchoring points for linking, a **unique** name must be specified for the element. Please refer to *Chapter 7.9.2.7* for more information.

7.9.2.1.9 Bisquert #2, B2

The Bisquert #2 circuit element is represented by the letter **B2** and identified by the following symbol:



This element is a transmission line element derived from the classical model for a porous or mixed-phase electrode of thickness L. The model is represented in *Figure 479*.




Figure 479 Overview of the general transmission line model used in the B2 element

In the **B2** element, the **X** element used in the transmission line is represented by a parallel combination of a resistor (R) and a constant phase element (Q).

This transmission line is often used in the world of dye-sensitized solar cells (DSC) and in general systems that analyze the combination of charge transport, accumulation and recombination.

The impedance of this equivalent circuit element may be written as:

$$Z_{B2} = \frac{X_1 X_2}{X_1 + X_2} \left(L + \frac{2\lambda}{\sinh\left(\frac{L}{\lambda}\right)} \right) + \lambda \frac{X_1^2 + X_2^2}{X_1 + X_2} \cot \operatorname{anh}\left(\frac{L}{\lambda}\right)$$

Where λ is given by:

$$\lambda = \sqrt{\frac{X_3}{X_1 + X_2}}$$

This element is a composite element, consisting of three types of parallel (RQ) element combinations. For the properties of the R element and the Q element, please refer to *Chapter 7.9.2.1.1* and *Chapter 7.9.2.1.4*, respectively.

The **B2** element provides one additional property, L, representing the length of the transmission line, shown in *Figure 480*.



NOTICE

For more information on the Bisquert #2 transmission line model, please refer to J. Bisquert, G. Garcia-Belmonte, F. Fabregat-Santiago, A. Compte, *Electrochemistry Communications* 1999, 1:9:429-435 and J. Bisquert; *Phys. Chem. Chem. Phys.*, Vol. 2 (2000), pp. 4185-4192.

7.9.2.1.10 Serial connection

The **Serial connection** can be used to place two circuit elements in **series**. The **Serial connection** is represented by the following symbol:



The **Serial connection** has one input connection and one output connection.

7.9.2.1.11 Parallel split connection

The **Parallel split connection** can be used to place two or more circuit elements in **parallel**. The **Parallel split** connection creates a parallel arrangement. The **Parallel split** connection is represented by the following symbol:



The **Parallel split** connection has one input connection and two output connections.

If needed additional output connections can be created, by right-clicking the connection and selecting the *Add output* option from the context menu (see figure 481, page 427).




Figure 481 Adding additional output connections

The additional output will be added to the element (see figure 482, page 428).




Figure 482 The additional output is added to the Parallel split connection




NOTICE

The same menu can be used to remove outputs from the element.

7.9.2.1.12 Parallel join connection

The **Parallel join connection** can be used to place two or more circuit elements in **parallel**. The **Parallel join** connection closes a parallel arrangement. The **Parallel join** connection is represented by the following symbol:



The **Parallel join** connection has two input connections and one output connection.

If needed additional input connections can be created, by right-clicking the connection and selecting the *Add input* option from the context menu (see figure 483, page 428).




Figure 483 Adding additional input connections

The additional input be added to the element (see figure 484, page 429).




Figure 484 The additional input is added to the Parallel join connection



NOTICE

The same menu can be used to remove inputs from the element.

7.9.2.2 Build a custom equivalent circuit

The Equivalent Circuit Editor window can also be used to draw the equivalent circuit by connecting individual element to one another, graphically. It is possible to add a circuit element to the editor from the *Insert* option available in the **Edit** menu (see figure 485, page 429).




Figure 485 Adding a circuit element from the Edit menu.

It is also possible to add a circuit element by right-clicking the editor window and using the context menu (see figure 486, page 430).




Figure 486 Adding a circuit element from the right-click menu
The selected circuit element will be added to the Equivalent circuit editor window (see figure 487, page 431).




Figure 487 The circuit element is added to the editor

Once two or more circuit elements or connectors are added to the Equivalent circuit editor, they can be linked to one another.

Each circuit element is fitted with one input connection (•) and one output connection (⊖). The connectors can have one or more input connections and one or more output connections *Circuit elements* (see chapter 7.9.2.1, page 410).

The following rules are used when creating custom equivalent circuits:

- It is only possible to connect an output connection of one element or connector to the input connection of an adjacent element or connector.
- A valid equivalent circuit can only have one free input connection and one free output connection.
- A valid link between an input connection and an output connection is represented by a closed loop symbol ⊙.

To create a link between two items, click one connection, and while holding the mouse button, drag this connection close to the connection of the next item (see figure 488, page 432).




Figure 488 Linking two items in the Equivalent circuit editor

When the two ends are close enough in the editor, the software will automatically create a link.

Using this method, any equivalent circuit respecting the rules detailed above can be created (see figure 489, page 432).




Figure 489 The custom made equivalent circuit

When the circuit is ready, it is possible to verify if there are errors by selecting the *Generate CDC from circuit* option from the Tools menu *Build a custom equivalent circuit* (see chapter 7.9.2.2, page 429).






Figure 490 Generating a CDC string from the equivalent circuit

If no errors are detected in the equivalent circuit, a valid CDC string (Circuit Description Code) will be displayed. If errors are detected, an error message will be shown (see figure 491, page 433).




Figure 491 An error message is shown when the circuit is invalid

7.9.2.3 Generate an equivalent circuit from a CDC string

To define the equivalent circuit from a CDC string (Circuit Description Code), select the *Generate Circuit from CDC* option from Tools menu (see figure 492, page 434).




Figure 492 Select the Generate Circuit from CDC option to manually enter a CDC string

A new window that can be used to input the CDC string will be displayed (see figure 493, page 434).




Figure 493 The CDC string can be entered using the proper formatting

To define the CDC string, the following syntax rules must be followed:

- Any of the nine element symbols defined in Table 12 can be used.
- Element placed in parallel must be written between ().
- Element placed in series must be written between [].

Once the CDC string is defined, click the OK button to create the circuit. The equivalent circuit will be drawn in the Equivalent Circuit Editor win-

dow, displaying the default initial values of the circuit elements (see figure 494, page 435).




Figure 494 The equivalent circuit is generated from the CDC string

If the CDC string is invalid, an error message will be displayed (see figure 495, page 435).




Figure 495 An error message is displayed if the CDC string is invalid

7.9.2.4 Load pre-defined circuit from a list

It is possible to choose an equivalent circuit from a pre-defined list of typical or user-defined circuits. To do this, select the *Open Circuit* option from the Circuit menu (see figure 496, page 436).




Figure 496 Opening the Circuit library

A new window will be displayed, showing two tabs (see figure 497, page 436).




Figure 497 The library provides two lists of pre-defined equivalent circuits



- **Pre-defined circuits:** this list contains a number of typical equivalent circuits.
- **User-defined circuits:** this list contains user-defined circuits.

Select the required equivalent circuit from either list and click the OK button. The selected equivalent circuit will be drawn in the Equivalent Circuit Editor window. The default or user-defined initial values will be displayed in blue (see figure 498, page 437).




Figure 498 The equivalent circuit is loaded from the circuit library



NOTICE

If the *Insert Circuit* option is selected instead of the *Open circuit* option, the selected circuit will be added to the Equivalent Circuit Editor without clearing the editor first.

7.9.2.5 Importing and Exporting equivalent circuits

It is possible to export equivalent circuits or to import equivalent circuits using the **File** menu (see figure 499, page 438).

- Copy/Cut and Paste element(s): select one or more elements in the equivalent editor by dragging a box around the circuit element and selecting the *Copy* or *Cut* option from the **Edit** menu (or the right-click menu or the **[CTRL] + [C]** and **[CTRL] + [X]** keyboard shortcuts) to copy them to the clipboard (see figure 500, page 439). The copied elements can then be pasted into the equivalent circuit editor, using the *Paste* option from the **Edit** menu (or the right-click menu or the **[CTRL] + [V]** keyboard shortcut), as shown in Figure 7.9.2.6.




Figure 500 Selected circuit elements can be copied/pasted directly in the editor




Figure 501 Elements pasted into the equivalent circuit editor have the same parameter values as the source elements

- **Convert Q element to pseudo capacitance:** this option can be used to convert a constant phase element **Q** element placed in parallel with a resistance **R** element to be converted to a pseudocapacitance, **C**. The conversion is performed according to:

$$C_{\text{pseudo}} = Y_0^{\frac{1}{n}} \cdot R^{\left(\frac{1}{n}-1\right)}$$

Where C_{pseudo} is the resulting pseudo capacitance, in F, Y_0 is the admittance value of the constant phase element, R is the resistance value and n is the exponent of the constant phase element. This formula is applicable to systems that can be accurately described with an RC circuit, as *blocking electrodes*. For more information and to reference this equation please refer to the following book: M. E. Orazem and B. Tribollet *Electrochemical Impedance Spectroscopy* Wiley, 2008, page 236.

To use this conversion tool, right click a **Q** element in parallel with a **R** element and select the *Convert to pseudo capacitance option* from the context menu as shown in *Figure 503*. The Q element will be converted to an equivalent capacitance value (see *figure 504, page 442*).




Figure 503 Converting a Q element to a pseudo capacitance




Figure 504 The Q element is converted to a C

7.9.2.7 Editing equivalent circuit properties

When the equivalent circuit is ready, it is possible to edit the properties of each of the circuit elements. To edit the properties of one of the element, click the element to select it. The selected element will be highlighted, as shown in Figure 505.




Figure 505 Selecting the equivalent circuit element

With the element selected, move the mouse pointer over the **Properties** tab on the right-hand side. The properties panel will be expanded, revealing the properties of the selected element (see figure 506, page 443).




Figure 506 Displaying the properties panel

The properties panel shows one or more containers for each element, which can be expanded or collapsed to reveal or to hide advanced variables (see figure 507, page 443).




Figure 507 The basic and advanced properties




Figure 508 Specifying the name of the element



NOTICE

The element name **must** be unique!

Once the name is specified, it is possible to check the **Input** and/or **Output** checkboxes in the **Properties** panel (see figure 509, page 446).




Figure 509 Specifying the linking behavior for the circuit element

This will create an input and output anchoring point for the element property, allowing it to be linked to another command properties (see figure 510, page 446).




Figure 510 The element property can now be linked

7.9.2.9 Saving equivalent circuits to the Library

Any equivalent circuit can be saved to the **Library**. Saved equivalent circuits will become available to the user from the *Open/Insert Circuit* option as described in *Chapter 7.9.2.4*.

To save the circuit to the Library, select the *Save Circuit* option from the **Circuit** menu (see figure 511, page 447).




Figure 511 Using the Save circuit option

A new window will be displayed, prompting for the name of the equivalent circuit (see figure 512, page 447).




Figure 512 Saving the equivalent circuit

Specify a name for the circuit and press the OK button to save it in the database.




Figure 514 Right-clicking the saved circuit allows renaming, deleting or quick access to the file location

7.9.3 Kronig-Kramers test

 <p>Kronig-Kramers test</p>	<p>This command allows testing a set of measured data using the Kronig-Kramers equations. This test provides an estimation of the 'goodness to fit' of the data set.</p>
---	--

The details of the properties of the **Kronig-Kramers test** command are shown in Figure 515:




Figure 515 The properties of the Kronig-Kramers test command

The following properties are available:

special circuit used in the test is a series of RC circuits (for impedance representation). This circuit are shown *Figure 516*.




Figure 516 Circuit used in for Kramers-Kronig test on impedance presentation

By default, the number of (RC) subcircuits is equal to the number of data points. If there is a chance that the measured signal was very noisy, the number of subcircuits may be reduced to avoid over-fitting and, consequently, including the noise in the model.

The result of the test is the value of pseudo, χ_{ps}^2 , the sum of squares of the relative residuals. In each case the χ^2 for the real and the imaginary part is reported (overall χ^2 is a sum of real and imaginary χ^2). Large χ^2 values indicate that the data quality is low. A small value, on the other hand, usually indicates a good fit.

The equations used in the Kramers-Kronig test are provided below:

$$\chi_{ps}^2 = \sum_{i=1}^n \frac{[Z_{re,i} - Z_{re}(\omega_i)]^2 + [Z_{im,i} - Z_{im}(\omega_i)]^2}{|Z(\omega_i)|^2}$$

$$\chi_{re}^2 = \sum_{i=1}^n \frac{[Z_{re,i} - Z_{re}(\omega_i)]^2}{|Z(\omega_i)|^2}$$

$$\chi_{im}^2 = \sum_{i=1}^n \frac{[Z_{im,i} - Z_{im}(\omega_i)]^2}{|Z(\omega_i)|^2}$$

What is actually large and small depends on the number and the value of data points. As a rule of thumb, values lower than 10^{-6} usually means an excellent fit, reasonable between 10^{-5} and 10^{-6} , marginal between 10^{-4} and 10^{-5} and bad for even higher values. Moreover, the residuals should be small and randomly distributed around zero.

The test can be carried out on real part, imaginary part or both part of admittance/impedance (complex fit). In the case of fit on one part only, the second part of the measured data set is generated using Kramers-Kronig transformation (using the assumption that the system obeys Kramers-Kronig criteria) and then χ^2 for the second part is computed.

In addition to χ^2 , the serial or parallel (depending on representation) R, L and C values are computed. These values do not have any special meaning and they simply belong to the set of results of Kramers-Kronig test. In



particular, they should not be associated with any serial or parallel elements present in the system or its equivalent circuit representation.



NOTICE

The detailed discussion of the Kramers-Kronig test, the theory underlying the choice of properties, and a refined interpretation of the outcomes can be found in B.A. Boukamp, J. Electrochem. Soc. 142, 1885 (1995). It is advised to read this article before this command is used.

7.9.4 Include all FRA data



This command automatically calculates the admittance data based on measured impedance data.

The details of the properties of the **Include all FRA data** command are shown in *Figure 517*:

Figure 517 The properties of the *Include all FRA data* command

The following properties are available:

- **Command name:** a user-defined name for the command.
- **Geometric capacitance:** specifies the value of the geometric capacitance, ϵ , in F (default: 1). This value is used in the calculation of the permittivity.

The **Include all FRA data** command can be used to automatically calculate and display additional information that can be derived mathematically from impedance data. This command calculates the following additional values:

- **Real admittance, Y'** : the real part of the admittance. This value is calculated according to:

$$Y' = \frac{Z'}{|Z|^2}$$

- **Imaginary admittance, $-Y''$** : the imaginary part of the admittance. This value is calculated according to:

$$-Y'' = \frac{-Z''}{|Z|^2}$$

- **Angular frequency, ω** : the angular frequency, in rad/s. This value is calculated according to:

$$\omega = 2\pi f$$

- **Real permittivity, $Y\varepsilon'$** : the real part of the permittivity. This value is calculated according to:

$$Y\varepsilon' = \Re\left(\frac{Y' - jY''}{j\omega\varepsilon}\right)$$


- **Imaginary permittivity, $-Y\varepsilon''$** : the imaginary part of the permittivity. This value is calculated according to:

$$-Y\varepsilon'' = \Im\left(\frac{Y' - jY''}{j\omega\varepsilon}\right)$$

- **Series capacitance, C_s** : the series capacitance. This value is calculated according to:

$$C_s = -\frac{1}{\omega Z''}$$

7.9.5 Potential scan FRA data

 <p>Potential scan FRA data</p>	<p>This command automatically calculates the values required to create a Mott-Schottky plot.</p>
--	--

The details of the properties of the **Potential scan FRA data** command are shown in *Figure 518*:

Properties ➔

Potential scan FRA data

Command name

Rs Ω

Figure 518 The properties of the Potential scan FRA data command

- **Dosino:** a command that can be used to control a Metrohm 800 Dosino connected to the host computer *Dosino* (see chapter 7.10.1, page 455).
- **Sample Processor:** a command that can be used to control a Metrohm 814, 815 or 858 Sample Processor connected to the host computer *Sample Processor* (see chapter 7.10.2, page 460)
- **Stirrer:** a command that can be used to control a Metrohm 801 Magnetic Stirrer or a Metrohm 802 Rod Stirrer or Metrohm 741 Magnetic Stirrer connected to a 804 Titration Stand connected to the host computer *Stirrer* (see chapter 7.10.3, page 468).
- **Remote I/O:** a command that can be used to control a Metrohm 6.2148.010 Remote Box connected to the host computer *Remote* (see chapter 7.10.4, page 469).

7.10.1 Dosino

	<p>This command can be used to control the Metrohm 800 Dosino connected to the computer.</p>
---	---

The **Dosino** command can be used in six different modes, which can be selected using the provided drop-down list (see figure 520, page 455):




Figure 520 Six modes are provided by the *Dosino* command

1. Prepare
2. Dose (default mode)
3. Empty
4. Fill
5. To end

7.10.1.2 Prepare

The *Prepare* mode of the **Dosino** command can be used to prepare the Dosino by rinsing and filling the connected tubes and the dosing cylinder. The tubes of the Dosino should be freed from air bubbles at least once a day by carrying out a full prepare cycle. This process will take time depending on the length of the tubes.

During the preparation process, the dosing cylinder as well as the connected tubings are completely filled. The volume required to fill the tubings is determined based on the parameters specified in the hardware setup of the Dosino.



NOTICE

The Fill port, defined in the Dosino hardware setup, is used to refill the Dosino during the preparation process. See *Chapter 5.5.1.2* for more information.



NOTICE

Ports that are set to *inactive* in the Dosino hardware setup are not used in the preparation process. See *Chapter 5.5.1.2* for more information.

The following properties are available when the **Dosino** command is used in the *Prepare* mode (see *figure 522*, page 457):

Figure 522 Prepare mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Dosino.
- **Cycles:** the number of cycles used to prepare the Dosino.

7.10.1.4 Fill

The *Fill* mode of the **Dosino** command can be used to completely refill the dosing cylinder of the specified Dosino. The liquid is aspirated through the Fill port defined in the Dosino hardware setup *Dosino hardware setup* (see chapter 5.5.1.2, page 151).

The following properties are available when the **Dosino** command is used in the *Fill* mode (see figure 524, page 459):

Figure 524 Fill mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Dosino.

7.10.1.5 To end

The *To end* mode of the **Dosino** command can be used to eject the contents of the dosing cylinder through the specified port. The piston stops at the specified end volume. This is useful for pipetting functions or for removing air bubbles from the dosing cylinder.

The following properties are available when the **Dosino** command is used in the *To end* mode (see figure 525, page 459):

Figure 525 Exchange mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Dosino.
- **Port:** the port used to perform the to end action.



7.10.1.6 Exchange


The *Exchange* mode of the **Dosino** command can be used to prepare a dosing cylinder for exchange. The dosing cylinder is filled and the stop-cock is moved to the exchange position. The cylinder is filled by aspirating the necessary volume via the Fill port specified in the Dosino hardware setup *Dosino hardware setup* (see chapter 5.5.1.2, page 151).

The following properties are available when the **Dosino** command is used in the *Exchange* mode (see figure 526, page 460):

Figure 526 Exchange mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Dosino.

7.10.2 Sample Processor

 <p>Move</p>	<p>This command can be used to control the Metrohm 814, 815 or 858 Sample Processor connected to the computer.</p>
---	---

The **Sample Processor** command can be used in eight different modes, which can be selected using the provided drop-down list (see figure 527, page 461):




Figure 527 Eight modes are provided by the Sample Processor command

1. Move (default mode)
2. Lift
3. Valve
4. Pump
5. Swing
6. Stirrer
7. Inject
8. Peristaltic pump



NOTICE

The last two modes are only available when using the **Metrohm 858 Professional Sample Processor**.




NOTICE

The **Sample Processor** command description in the procedure editor is dynamically adjusted in function of the specified mode.

7.10.2.2 Lift

The *Lift* mode of the **Sample Processor** command can be used to set the position of the lift on the specified Sample Processor tower. The position of the lift can be specified between 0 mm (top of the tower) and the maximum position defined in the Sample Processor hardware setup *Sample Processor hardware setup* (see chapter 5.5.2.2, page 158).

The following properties are available when the **Sample Processor** command is used in the *Lift* mode (see figure 529, page 463):



The screenshot shows a 'Properties' dialog box with a close button in the top right corner. The title is 'Lift 0 mm'. Below the title are five input fields: 'Command name' with the value 'Lift 0 mm', 'Action' with a dropdown menu showing 'Lift', 'Device name' with the value '4589', 'Tower' with the value '1', and 'Position' with the value '0' and a unit label 'mm' to its right.

Figure 529 Lift mode properties

- **Device name:** the identifying name of the Sample Processor.
- **Tower:** specifies which tower is used by the command.
- **Position:** specifies the position of the lift, in mm, with respect to the top of the tower.

7.10.2.3 Valve

The *Valve* mode of the **Sample Processor** command can be used to activate or deactivate valves mounted on the back plane of a tower.



NOTICE

Valves remain on or off until modified by the procedure or through the Sample processor manual control panel *Metrohm Sample Processor control panel* (see chapter 5.5.2, page 155).

The following properties are available when the **Sample Processor** command is used in the *Valve* mode (see figure 530, page 464):




Figure 531 Pump mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Sample Processor.
- **Tower:** specifies which tower is used by the command.
- **Pump 1:** specifies the state of pump 1 through a dedicated toggle.
- **Pump 2:** specifies the state of pump 2 through a dedicated toggle.

7.10.2.5 Swing

The *Swing* mode of the **Sample Processor** command can be used to change the position of the swing head installed on the specified Sample Processor tower.

The following properties are available when the **Sample Processor** command is used in the *Swing* mode (see figure 532, page 465):




Figure 532 Swing mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Sample Processor.



- **Tower:** specifies which tower is used by the command.
- **Angle:** specifies the angle of the swing arm with respect to the tower, in °. The range of value depends on the type of swing head mounted on the swing arm.

7.10.2.6 Stir

The *Stirrer* mode of the **Sample Processor** command can be used to control the rotation rate of a **Metrohm 802 Rod Stirrer** or **Metrohm 741 Magnetic Stirrer** connected to the Sample Processor tower.

The following properties are available when the **Sample Processor** command is used in the *Stirrer* mode (see figure 533, page 466):




Figure 533 *Stirrer mode properties*

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Sample Processor.
- **Tower:** specifies which tower is used by the command.
- **Rotation rate:** the rotation rate, specified between -15 and 15. A value of 0 will stop the stirrer.

7.10.2.7 Inject

The *Inject* mode of the **Sample Processor** command can be used to set the position of the injection valve. The connections to the injection valve can be toggled between the **Fill** position and the **Inject** position (see figure 534, page 466).




Figure 534 *The positions of the injection valve*

1 Fill position

2 Inject position



NOTICE

This mode can only be used in combination with the **Metrohm 858 Professional Sample Processor** fitted with the injection valve.

The following properties are available when the **Sample Processor** command is used in the *Inject valve* mode (see figure 534, page 466):

Figure 535 *Inject valve mode properties*

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Sample Processor.
- **Position:** specifies the position of the inject valve, using the provided drop-down list.

7.10.2.8 Peristaltic pump

The *Peristaltic pump* mode of the **Sample Processor** command can be used to control the peristaltic pump installed on the Sample Processor.



NOTICE

This mode can only be used in combination with the **Metrohm 858 Professional Sample Processor** fitted with the peristaltic pump.


The following properties are available when the **Sample Processor** command is used in the *Peristaltic pump* mode (see figure 536, page 468):



Figure 536 Peristaltic pump mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Sample Processor.
- **Rotation rate:** the rotation rate, specified between -15 and 15. A value of 0 will stop the pump.

7.10.3 Stirrer

	<p>This command can be used to control the Metrohm 801 Magnetic Stirrer or Metrohm 804 Titration Stand in combination with the Metrohm 802 Rod Stirrer connected to the computer.</p>
--	--

The details of the command properties of the **Stirrer** command are shown in *Figure 537*:

Figure 537 The properties of the Stirrer command

The following properties are available:


- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Stirrer.
- **Rotation rate:** the rotation rate, specified between -15 and 15. A value of 0 will stop the stirrer.



NOTICE

The **Stirrer** command description in the procedure editor is dynamically adjusted in function of the specified value.

7.10.4 Remote



This command can be used to control the **Met-rohm 6.2148.010 Remote Box** connected to the computer.

The **Remote** command can be used in two different modes, which can be selected using the provided drop-down list (see figure 538, page 469):




Figure 538 Two modes are provided by the Remote command

1. Remote inputs (default mode)
2. Remote outputs



NOTICE

The **Remote** command description in the procedure editor is dynamically adjusted in function of the specified mode.



CAUTION

The **Metrohm 6.2148.010 Remote Box** can also be used in combination with the **Wait** command *Wait* (see chapter 7.2.4, page 240). When the Remote Box is connected to the computer, the **Wait** command provides one additional mode, *Wait for Remote Inputs*, which uses the eight input lines provided by the Remote Box.

7.10.4.1 Remote inputs

The *Remote inputs* mode of the **Remote** command can be used to read the state of the 8 input lines (numbered IN7 to IN0). The state of each input line can be either 'low' or 'high' state, represented by a 0 or a 1, respectively.

The following properties are available when the **Remote** command is used in the *Remote inputs* mode (see figure 539, page 470):

Figure 539 Remote inputs mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Remote Box.
- **Inputs:** specifies the state of the 8 input lines, during the execution of the procedure. The state is returned as a string of 8 characters, consisting of '0' and '1', representing the state of the input lines, from IN7 to IN0.



NOTICE

The state of the 8 input lines of the Remote Box is determined when the command is executed.

7.10.4.2 Remote outputs

The *Remote outputs* mode of the **Remote** command can be used to set the state of the 14 output lines (numbered OUT13 to OUT0). The state of each output line can be set to either 'low' or 'high' state, represented by a 0 or a 1, respectively.

The following properties are available when the **Remote** command is used in the *Remote outputs* mode (see figure 540, page 471):

Figure 540 Remote outputs mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying name of the Remote Box.
- **Outputs:** specifies the state of the 14 output lines is specified as a 14 character string, consisting of '0' and '1', representing the state of the output lines, from OUT13 to OUT0.



NOTICE

The state of the 14 output lines of the Remote Box is persistent until changed or until the Remote Box is powered down.

The screenshot shows a 'Properties' dialog box for a 'Software trigger' command. The 'Action' dropdown menu is open, showing two options: 'Software trigger' (which is highlighted in green) and 'DIO trigger'. Other fields in the dialog include 'Command name' (Software trigger), 'Device name' (Software trigger), 'Start wavelength' (1333 nm), 'Stop wavelength' (1333 nm), 'Integration time' (500 ms), and 'Number of averages' (1). There is also a checkbox for 'Enable light source shutter control' which is currently unchecked. A 'More' button is located at the bottom right of the dialog.

Figure 542 Two modes are provided by the Spectroscopy command

1. Software trigger (default mode)
2. DIO trigger



NOTICE

The **Spectroscopy** command description in the procedure editor is dynamically adjusted in function of the specified mode.

7.11.1.1 Software trigger

The following properties are available when the command is used in the *Software trigger* mode (see figure 543, page 474):

- **Shutter open:** specifies the state of the light source shutter, using the provided  toggle. This property is only visible if the *Enable light source shutter control* property is set to on.



NOTICE

The light source shutter will remain in the specified state until changed.




CAUTION

This mode requires a physical connection between the light source and the Autolab DIO connector if the *Enable light source shutter control* property is set to on. Please refer to the Spectrophotometer User Manual for more information.

7.11.1.2 DIO trigger

The following properties are available when the command is used in the *DIO trigger* mode (see figure 544, page 475):



Properties ✕

DIO trigger

Command name

Action

Device name

Start wavelength nm

Stop wavelength nm

Integration time ms

Number of averages

Get spectrum counter

Calculate Absorbance and Transmittance

More

Figure 544 DIO trigger mode properties

- **Command name:** a user-defined name for the command.



- **Device name:** specifies the name of the spectrophotometer used in the measurement.
- **Start wavelength:** specifies the start wavelength used by the spectrophotometer, in nm.
- **Stop wavelength:** specifies the start wavelength used by the spectrophotometer, in nm.
- **Integration time:** specifies the integration time used by the spectrophotometer, in ms.
- **Number of averages:** specify the number of averages used by the spectrophotometer.
- **Get spectrum counter:** the counter value used by the triggering command.
- **Calculate Absorbance and Transmittance:** specifies if the measured values should be converted to absorbance and transmittance using values of a dark spectrum and reference spectrum using the provided toggle.

If the **Calculate Absorbance and Transmittance** property is on, it is necessary to link two single spectra to the **Spectroscopy** command. Two input anchoring points will be added to the command (see figure 545, page 476).




Figure 545 Dark and reference spectra can be linked to the Spectroscopy command

Using these two anchoring points, a dark spectrum and a reference spectrum can be linked to the **Spectroscopy** command in order to convert the measured values to absorbance and transmittance.

These values are calculated using the measured values (S_{Measured}), the linked *Dark* spectrum values (S_{Dark}) and the linked *Reference* spectrum values ($S_{\text{Reference}}$) according to:

- **Absorbance:**

$$A = -\log\left(\frac{S_{\text{Measured}} - S_{\text{Dark}}}{S_{\text{Reference}} - S_{\text{Dark}}}\right)$$

- **Transmittance:**

$$T = 100 \cdot \left(\frac{S_{\text{Measured}} - S_{\text{Dark}}}{S_{\text{Reference}} - S_{\text{Dark}}}\right)$$



CAUTION

The linked dark and reference spectra must be measured in the same conditions as those of the **Spectroscopy** command they are linked to.



CAUTION

The linked dark and reference spectra may be imported from saved data files using the **Import data** command. The spectroscopic measured values (S_{measured}) column from the imported data file can be linked to the appropriate anchoring point, either the **Input dark spectrum** or **Input reference spectrum**. The imported data must have been measured with an Autolab spectrophotometer using the same conditions as those of the **Spectroscopy** command they are linked to.



CAUTION

This mode requires a physical connection between the spectrophotometer and the Autolab DIO connector. Please refer to the Spectrophotometer User Manual for more information.

To use the **Spectroscopy** command in *DIO trigger* mode in a NOVA procedure, the command needs to be stacked onto the electrochemical



measurement command that it is used with. *Figure 546* provides an example.




Figure 546 Stacking the spectroscopy command on a measurement command

Using this configuration, the **Spectroscopy** command used in *DIO trigger* mode will be executed whenever the parent measurement command (**LSV staircase** in *Figure 546*) will send a DIO trigger to the spectrophotometer.



NOTICE

More information on the stacking of commands can be found in *Chapter 10.12*.

At the end of a measurement, the electrochemical data will be provided by the parent measurement command and the spectroelectrochemical data will be provided by the **Spectroscopy** command (see *figure 547*, page 479).






Figure 547 The spectroscopy and electrochemistry data is available in the Spectroscopy command

7.11.2 External device control



Send

This command can be used to interface to external instruments connected through the RS-232 standard.


The RS-232 standard describes a communication method where information is sent bit by bit on a physical channel. The information must be broken up in data words. The length of a data word is variable (usually between 5 and 8 bits). For proper transfer additional bits are added for synchronization and error checking purposes.



NOTICE

Interfacing to external devices through the RS-232 standard requires a properly configured COM port on the computer.

The **External device control** command can be used in four different modes, which can be selected using the provided drop-down list (see figure 548, page 480):



Properties [Close]

Initialize

Command name:

Action:

Device name:

Device type:

Port name:

Baud rate:

Data bits:

New line:

Parity:

Stop bits:

Hand shake:

Parallel operation:

Figure 549 Initialize mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** specifies the name of the device to initialize at the beginning of the measurement. The device name must be unique and will be used to identify the connected device in NOVA.
- **Port name:** specifies the COM port used to control the external device (replace *x* with the COM port number).
- **Baud rate:** specifies the baud rate used to communicate with the external device.
- **Data bits:** specifies the number of data bits (8 by default).
- **New line:** specify the character used to create a new line (`\n` by default).
- **Parity:** specifies the parity, using the provided drop-down list (None, Odd, Even, Mark, Space).
- **Stop bits:** specifies the number of stop bits, using the provided drop down list (0, 1, 2, 1.5).
- **Handshake:** specifies the handshaking mode for the communication, using the provided drop down list (None, X on X off, Request to send, Request to send X on X off)

- **Command:** a string defining the format of the expected string from the external device with placeholders ({0}, {1}, ...) for (variable) parameters in the string.



NOTICE

To receive a data string from the external device, the *Send* mode is first used to send a specific data string to the external device. The *Receive* is then added to the procedure to read the reply string from the external device.

7.11.2.4 Close

The following properties are available when the command is used in the *Close* mode (see figure 552, page 483):

Figure 552 Close mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the name of the device to close at the end of the measurement.



CAUTION

Always use the *Close* mode to close the communication to an external device at the end of each measurement.


7.11.3 RHD control



This command can be used to control the **Autolab RHD Microcell HC** connected to the computer.

7.11.3.2 Set temperature

The following properties are available when the **RHD control** command is used in the *Set temperature* mode (see figure 555, page 485):



The screenshot shows a 'Properties' dialog box with a close button in the top right corner. Below the title bar, the text 'Set 25 °C' is displayed. There are four input fields: 'Command name' containing 'Set 25 °C', 'Mode' with a dropdown menu showing 'Set temperature', 'Device name' which is empty, and 'Temperature' containing '25' with a '°C' unit label to its right.

Figure 555 Set temperature mode properties

- **Command name:** a user-defined name for the command.
- **Device name:** the identifying device name of the Autolab RHD Microcell HC controller.
- **Temperature:** the target temperature to set on the Autolab RHD Microcell HC controller.



NOTICE

When the **RHD control** command is used in the *Set temperature* mode, the command will be executed and will hold until the temperature stabilization conditions, defined in the hardware setup of the Autolab RHD Microcell HC controller are reached *Autolab RHD Microcell HC control panel* (see chapter 5.3, page 134). This command cannot be skipped or stopped.

- Impedance spectroscopy
 - FRA impedance potentiostatic
 - FRA impedance galvanostatic
 - FRA potential scan
 - FRA current scan
 - FRA time scan potentiostatic
 - FRA time scan galvanostatic

Electrochemical Frequency Modulation

The rest of this chapter provides a detailed description of each procedure provided as a default in NOVA.

8.1 Cyclic voltammetry

NOVA provides five default procedures for cyclic voltammetry. These procedures can be used to perform a cyclic potential or current scan and record the response of the cell. Some of these procedures require optional hardware extensions.

The following procedures are available:

- Cyclic voltammetry potentiostatic
- Cyclic voltammetry galvanostatic
- Cyclic voltammetry current integration (requires the **FI20** or **on-board integrator**, please refer to *Chapter 16.3.2.11* for more information)
- Cyclic voltammetry linear scan (requires the **SCAN250** or **SCANGEN** module, please refer to *Chapter 16.3.2.19* for more information)
- Cyclic voltammetry linear scan high speed (requires the **SCAN250** or **SCANGEN** module and **ADC10M** or **ADC750** module, please refer to *Chapter 16.3.2.19* and *Chapter 16.3.2.1* for more information).

8.1.1 Cyclic voltammetry potentiostatic

The default **Cyclic voltammetry potentiostatic** procedure provides an example of a typical *staircase* cyclic voltammetry procedure in potentiostatic mode (see *figure 556*, page 487).




Figure 556 The default Cyclic voltammetry potentiostatic procedure

The procedure samples the following signals (see figure 558, page 489):




Figure 558 The sampler of the CV staircase command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 559, page 489):




Figure 559 The options of the CV staircase command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 560, page 490):




Figure 560 The plots of the CV staircase command

- i vs E: WE(1).Current versus Potential applied

The procedure also has the value of alpha property available in the Advanced section. This value is set to the default value of 1 (see figure 561, page 490).




Figure 561 The advanced settings of the CV staircase command

8.1.2 Cyclic voltammetry galvanostatic

The default **Cyclic voltammetry galvanostatic** procedure provides an example of a typical *staircase* cyclic voltammetry procedure in galvanostatic mode (see figure 562, page 490).




Figure 562 The default Cyclic voltammetry galvanostatic procedure



NOTICE

The galvanostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure has the following measurement properties, specified for the **CV staircase** command (see figure 563, page 491):

Properties ➔

CV staircase

Command name	<input type="text" value="CV staircase"/>
Start current	<input type="text" value="0"/> A
Upper vertex current	<input type="text" value="0,001"/> A
Lower vertex current	<input type="text" value="-0,001"/> A
Stop current	<input type="text" value="0"/> A
Number of scans	<input type="text" value="1"/>
Scan rate	<input type="text" value="0,0001"/> A/s
Step	<input type="text" value="2,44E-06"/> A
Interval time	<input type="text" value="0,024414"/> s
Estimated number of points	<input type="text" value="1641"/>
Estimated duration	<input type="text" value="40,063"/> s
Number of stop crossings	<input type="text" value="2"/>

More

Figure 563 The measurement properties of the CV staircase command

- **CV staircase**

- Start current: 0 A
- Upper vertex current: 0.001 A
- Lower vertex current: -0.001 A
- Stop current: 0 A
- Number of scans: 1
- Step: 2.44 μ A
- Scan rate: 0.0001 A/s



The procedure samples the following signals (see figure 564, page 492):




Figure 564 The sampler of the CV staircase command

- WE(1).Current
- WE(1).Potential (averaged)
- Time

The procedure plots the following data (see figure 565, page 492):




Figure 565 The plots of the CV staircase command

- E vs i: WE(1).Potential versus Current applied

The procedure also has the value of alpha property available in the Advanced section. This value is set to the default value of 1 (see figure 566, page 493).




Figure 566 The advanced settings of the CV staircase command

8.1.3 Cyclic voltammetry potentiostatic current integration



CAUTION

This procedure requires the optional **FI20** module or the **on-board integrator FI20 module** (see chapter 16.3.2.11, page 1083).

The default **Cyclic voltammetry potentiostatic** current integration procedure provides an example of a typical *staircase* cyclic voltammetry procedure in potentiostatic mode, using the optional **FI20** module or **on-board integrator** (see figure 567, page 493).




Figure 567 The default Cyclic voltammetry potentiostatic current integration procedure

The charge determined during each step is used to recalculate the total current.



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure samples the following signals (see figure 569, page 495):




Figure 569 The sampler of the CV staircase command

- WE(1).Current
- WE(1).Potential
- Integrator(1).Integrated Current
- Time

The procedure plots the following data (see figure 570, page 495):




Figure 570 The plots of the CV staircase command

Properties	
CV linear scan	
Command name	CV linear scan
Mode	Normal
Start potential	0 V _{REF}
Upper vertex potential	1 V _{REF}
Lower vertex potential	-1 V _{REF}
Stop on	Vertex
Number of scans	1,25
Scan rate	0,1 V/s
Potential interval	0,00244 V
Interval time	0,0244 s
Estimated number of points	2050
Estimated duration	50,02 s
Number of vertex potential crossings	3
More	

Figure 573 The properties of the CV linear scan command

▪ **CV linear scan**

- Start potential: 0 V, versus reference electrode
- Upper vertex potential: 1 V, versus reference electrode
- Lower vertex potential: -1 V, versus reference electrode
- Number of scans: 1,25
- Potential interval: 0.00244 V
- Scan rate: 100 mV/s

The procedure samples the following signals (see figure 574, page 498):




Figure 574 The sampler of the CV linear scan command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 575, page 498):




Figure 575 The options of the CV linear scan command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 576, page 499):




Figure 576 The plots of the CV linear scan command

- i vs E: WE(1).Current versus Potential applied

8.1.5 Cyclic voltammetry potentiostatic linear scan high speed



CAUTION

This procedure requires the optional **SCAN250** or **SCANGEN** module in combination with the optional **ADC10M** or **ADC750** module. *SCAN250 module (see chapter 16.3.2.19, page 1169) and ADC10M module (see chapter 16.3.2.1, page 998).*

The default **Cyclic voltammetry potentiostatic linear scan high speed** procedure provides an example of a typical **linear scan** cyclic voltammetry procedure at very high scan rate in potentiostatic mode, using the optional **SCAN250** or **SCANGEN** module in combination with the optional **ADC10M** or **ADC750** module (*see figure 577, page 500*).

Properties ➔

CV linear scan

Command name

Mode ▼

Start potential V_{REF} ▼

Upper vertex potential V_{REF} ▼

Lower vertex potential V_{REF} ▼

Stop on ▼

Number of scans

Scan rate V/s

Potential interval V

Interval time s

Estimated number of points

Estimated duration s

Number of vertex potential crossings

Figure 578 The properties of the CV linear scan command

▪ **CV linear scan**

- Start potential: 0 V, versus reference electrode
- Upper vertex potential: 1 V, versus reference electrode
- Lower vertex potential: -1 V, versus reference electrode
- Number of scans: 1,25
- Potential interval: 0.00056 V
- Scan rate: 100 V/s

The procedure samples the following **ADC10M** or **ADC750** settings (see figure 579, page 502):




Figure 579 The ADC10M or ADC750 settings of the CV linear scan command

- Channel 1: WE(1).Potential, Gain 1, unfiltered
- Channel 2: WE(1).Current, Gain 1, unfiltered

The procedure plots the following data (see figure 580, page 502):




Figure 580 The plots of the CV linear scan command

- i vs E : WE(1).Current versus Potential applied



NOTICE

The measured data cannot be displayed in real-time. The data is only available at the end of the measurement.

8.2 Linear sweep voltammetry

NOVA provides four default procedures for linear sweep voltammetry. These procedures can be used to perform a potential or current sweep and record the response of the cell. Some of these procedure require optional hardware extensions.

The following procedures are available:

- Linear sweep voltammetry potentiostatic
- Linear sweep voltammetry galvanostatic
- Linear polarization
- Hydrodynamic linear sweep (requires the Autolab rotating disk electrode (RDE) or Autolab rotating ring disk electrode (RRDE), please refer to the Autolab RDE/RRDE User Manual for more information)
- Hydrodynamic linear sweep with RRDE (requires the Autolab rotating ring disk electrode (RRDE), please refer to the Autolab RDE/RRDE User Manual for more information)
- Spectroelectrochemical linear sweep voltammetry (requires an Autolab or Avantes spectrophotometer)

8.2.1 Linear sweep voltammetry potentiostatic

The default **Linear sweep voltammetry potentiostatic** procedure provides an example of a typical *staircase* Linear sweep voltammetry procedure in potentiostatic mode (see figure 581, page 503).




Figure 581 The default Linear sweep voltammetry potentiostatic procedure




Figure 583 The sampler of the LSV staircase command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 584, page 505):




Figure 584 The options of the LSV staircase command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 585, page 506):




Figure 585 The plots of the LSV staircase command

- i vs E: WE(1).Current versus Potential applied

The procedure also has the value of alpha property available in the Advanced section. This value is set to the default value of 1 (see figure 586, page 506).




Figure 586 The advanced settings of the LSV staircase command

8.2.2 Linear sweep voltammetry galvanostatic

The default **Linear sweep voltammetry galvanostatic** procedure provides an example of a typical *staircase* Linear sweep voltammetry procedure in galvanostatic mode (see figure 587, page 506).




Figure 587 The default Linear sweep voltammetry galvanostatic procedure



NOTICE

The galvanostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure has the following measurement properties, specified for the **LSV staircase** command (see figure 588, page 507):

Properties	
LSV staircase	
Command name	LSV staircase
Start current	0 A
Stop current	0,001 A
Scan rate	0,0001 A/s
Step	2,44E-06 A
Interval time	0,024414 s
Estimated number of points	410
Estimated duration	10,01 s
More	

Figure 588 The measurement properties of the LSV staircase command

- **LSV staircase**

- Start current: 0 A
- Stop current: 0.001 A
- Step: 2.44 μ A
- Scan rate: 0.0001 A/s

The procedure samples the following signals (see figure 589, page 508):




Figure 589 The measurement properties of the LSV staircase command

- WE(1).Current
- WE(1).Potential (averaged)
- Time

The procedure plots the following data (see figure 590, page 508):




Figure 590 The measurement properties of the LSV staircase command

- E vs i: WE(1).Potential versus Current applied

The procedure also has the value of alpha property available in the Advanced section. This value is set to the default value of 1 (see figure 591, page 509).




Figure 591 The advanced settings of the LSV staircase command

8.2.3 Linear polarization

The default **Linear polarization** procedure provides an example of a typical *staircase* corrosion measurement according to *ASTM G5-14* in potentiostatic mode (see figure 592, page 509).




Figure 592 The default Linear polarization procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure has the following measurement properties, specified for the **LSV staircase** command (see figure 593, page 510):




Figure 593 The measurement properties of the LSV staircase command

▪ **LSV staircase**

- Start potential: -0.1 V, versus open circuit potential
- Stop potential: 0.1 V, versus open circuit potential
- Step: 0.001 V
- Scan rate: 1 mV/s

The procedure samples the following signals (see figure 594, page 510):




Figure 594 The sampler of the LSV staircase command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 595, page 511):




Figure 595 The options of the LSV staircase command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 596, page 511):




Figure 596 The plots of the LSV staircase command

- Log(i) vs E: Log(WE(1).Current) versus Potential applied

The procedure also has the value of alpha property available in the Advanced section. This value is set to the default value of 1 (see figure 597, page 511).




Figure 597 The advanced settings of the LSV staircase command



NOTICE

The open circuit potential is measured by the **OCP determination** command located before the **LSV staircase** command. Please refer to *Chapter 7.2.5* for more information.



NOTICE

The procedure includes a **Corrosion rate analysis** command to automatically analyze the measured data. Please refer to *Chapter 7.8.14* for more information.

8.2.4 Hydrodynamic linear sweep



CAUTION

This procedure requires the optional **Autolab rotating disk electrode (RDE)** or **Autolab rotating ring disk electrode (RRDE)** connected to the Autolab using the motor controller. The procedure is designed to remotely control the rotation rate. For more information, please refer to the Autolab RDE/RRDE User Manual.

The default **Hydrodynamic linear sweep** procedure provides an example of a typical *staircase* linear sweep voltammetry procedure in potentiostatic mode in combination with the **Autolab rotating disk electrode (RDE)** or **Autolab rotating ring disk electrode (RRDE)** (see *figure 598, page 512*).




Figure 598 The default Hydrodynamic linear sweep procedure





NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The Hydrodynamic linear sweep voltammetry procedure performs a linear sweep voltammetry using the Autolab RDE or Autolab RRDE, with six different rotation rates. The rotation rate of the Autolab RDE or Autolab RRDE is set using the **R(R)DE** command linked to the values of a **Repeat** command.

The **Repeat** command is used in the *Repeat for multiple values* mode and is preconfigured to cycle through six rotation rates, starting at 500 RPM until 3000 RPM, using a square root distribution (see figure 599, page 513).

Repeat for multiple values

Values ✕ ✎ +

	Rotation rate (RPM)
1	500
2	831,92
3	1247,9
4	1747,9
5	2331,9
6	3000

Add range

Begin value

End value

Number of values

Distribution

Figure 599 The repeat loop used in the default Hydrodynamic linear sweep procedure

The **Rotation rate** parameter, created by **Repeat** command, is linked to the **R(R)DE** command included in the repeat loop (see figure 600, page 514).




Figure 600 The link used to control the rotation rate of the R(R)DE

This procedure is intended to be used with the **Remote** switch of the Autolab motor controller engaged (on the back plane of the controller) and with a BNC cable connected between the DAC164 ←1 connector (Vout for the μ Autolab type II, μ Autolab type III, PGSTAT101, M101, PGSTAT204 and M204) and the Remote input plug on the back plane of the Autolab RDE motor controller (see figure 601, page 514).




Figure 601 The back plane of the Autolab motor controller

The procedure has the following measurement properties, specified for the **LSV staircase** command (see figure 602, page 515):

Properties 

LSV staircase

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Scan rate V/s

Step V

Interval time s

Estimated number of points

Estimated duration s

Figure 602 The measurement properties of the LSV staircase command

▪ **LSV staircase**

- Start potential: 1 V, versus reference electrode
- Stop potential: 0 V, versus reference electrode
- Step: -0.00244 V
- Scan rate: 100 mV/s



NOTICE

The *Step potential* value is negative because the potential scan is performed in the negative going direction.

The procedure samples the following signals (see figure 603, page 516):




Figure 605 The plots of the LSV staircase command

- i vs E: WE(1).Current versus Potential applied

The procedure also has the value of alpha property available in the Advanced section. This value is set to the default value of 1 (see figure 606, page 517).




Figure 606 The advanced settings of the LSV staircase command



NOTICE

The procedure includes a **Hydrodynamic analysis** command to automatically analyze the measured data. Please refer to *Chapter 7.8.10* for more information.



8.2.5 Hydrodynamic linear sweep with RRDE



CAUTION

This procedure requires the **BA** module *BA module* (see chapter 16.3.2.3, page 1011).



CAUTION

This procedure requires the optional **Autolab rotating ring disk electrode (RRDE)** connected to the Autolab using the motor controller. The procedure is designed to remotely control the rotation rate. For more information, please refer to the Autolab RDE/RRDE User Manual.

The default **Hydrodynamic linear sweep with RRDE** procedure provides an example of a typical *staircase* linear sweep voltammetry procedure in potentiostatic mode in combination with the **Autolab rotating ring disk electrode (RRDE)** (see figure 607, page 518).




Figure 607 The default Hydrodynamic linear sweep with RRDE procedure




NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The Hydrodynamic linear sweep with RRDE procedure performs a linear sweep voltammetry using the Autolab RRDE, with six different rotation rates. The rotation rate of the Autolab RRDE is set using the **R(R)DE** command linked to the values of a **Repeat** command.

The **Repeat** command is used in the *Repeat for multiple values* mode and is preconfigured to cycle through six rotation rates, starting at 500 RPM until 3000 RPM, using a square root distribution (see figure 599, page 513).



Repeat for multiple values

Values ✕ ✎ +

	Rotation rate (RPM)
1	500
2	831,92
3	1247,9
4	1747,9
5	2331,9
6	3000

Add range

Begin value

End value

Number of values

Distribution

Figure 608 The repeat loop used in the default Hydrodynamic linear sweep procedure

The **Rotation rate** parameter, created by **Repeat** command, is linked to the **R(R)DE** command included in the repeat loop (see figure 600, page 514).




Figure 609 The link used to control the rotation rate of the R(R)DE

This procedure is intended to be used with the **Remote** switch of the Autolab motor controller engaged (on the back plane of the controller) and with a BNC cable connected between the DAC164 ←1 connector (Vout for the PGSTAT204 and M204) and the Remote input plug on the back plane of the Autolab RDE motor controller (see figure 601, page 514).



Figure 610 The back plane of the Autolab motor controller

The procedure has the following measurement properties, specified for the **LSV staircase** command (see figure 611, page 521):

Properties

LSV staircase

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Scan rate V/s

Step V

Interval time s

Estimated number of points

Estimated duration s

More

Figure 611 The measurement properties of the LSV staircase command

▪ **LSV staircase**

- Start potential: 0.6 V, versus reference electrode
- Stop potential: -0.4 V, versus reference electrode
- Step: -0.00244 V
- Scan rate: 100 mV/s



NOTICE

The *Step potential* value is negative because the potential scan is performed in the negative going direction.

The settings of the **BA** module, used to control the ring, are defined using the **Autolab control** command located at the beginning of the procedure (see figure 612, page 522).




Figure 613 The sampler of the LSV staircase command

- WE(1).Current (averaged)
- WE(1).Potential
- WE(2).Current (averaged)
- Time

The procedure uses the following options (see figure 614, page 523):




Figure 614 The options of the LSV staircase command

- Automatic current ranging
 - Highest current range: 10 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 605, page 517):




Figure 615 The plots of the LSV staircase command

- i vs E: WE(1).Current versus Potential applied
- i(WE2) vs E: WE(2).Current versus Potential applied

The procedure also has the value of alpha property available in the Advanced section. This value is set to the default value of 1 (see figure 606, page 517).




Figure 616 The advanced settings of the LSV staircase command



NOTICE

The procedure includes a **Hydrodynamic analysis** command to automatically analyze the measured data. Please refer to *Chapter 7.8.10* for more information.

8.2.6 Spectroelectrochemical linear sweep



CAUTION

This procedure requires an optional **Autolab spectrophotometer** or supported **Avantes spectrophotometer** connected to the Autolab using the required trigger cable.

The default **Spectroelectrochemical linear sweep** procedure provides an example of a typical *staircase* linear sweep voltammetry procedure in potentiostatic mode in combination with the **Autolab spectrophotometer** or supported **Avantes spectrophotometer** (see figure 617, page 525).




Figure 617 The default Spectroelectrochemical linear sweep procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The **Spectroelectrochemical linear sweep voltammetry** procedure performs a linear sweep voltammetry using the spectrophotometer connected to the computer. The **Spectroscopy** command, included three times in this procedure, is used to measure the dark and reference spectra of the sample, before the linear sweep voltammetry measurement starts and the sample spectra during the execution of the LSV staircase commands, synchronized using a dedicated counter.

The procedure has the following measurement properties, specified for the **LSV staircase** command (see figure 618, page 526):



Properties ↗

LSV staircase

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Scan rate V/s

Step V

Interval time s

Estimated number of points

Estimated duration s

Figure 618 The measurement properties of the LSV staircase command

■ **LSV staircase**

- Start potential: 0 V, versus reference electrode
- Stop potential: 1 V, versus reference electrode
- Step: 0.00244 V
- Scan rate: 100 mV/s

The procedure samples the following signals (see figure 619, page 526):

← LSV staircase

Sampler

Options

Plots

Advanced

Signal	Sample	Average	d/dt
WE(1).Current	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
WE(1).Potential	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Power	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Resistance	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Charge	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Integrator(1).Charge	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Integrator(1).Integrated Current	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Time	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Sample alternating

Figure 619 The sampler of the LSV staircase command

- WE(1).Current (averaged)
- WE(1).Potential

- Time

The procedure uses the following options (see figure 620, page 527):




Figure 620 The options of the LSV staircase command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA
- Counters
 - Get spectrum when counter = 50, reset option on



NOTICE

The counter option specified in the options of the **LSV staircase** command is used to trigger the acquisition of a spectrum on the connected Autolab or Avantes spectrophotometer. This counter is repeated every 50 points.

The procedure plots the following data (see figure 621, page 528):




Figure 621 The plots of the LSV staircase command

- i vs E: WE(1).Current versus Potential applied

The procedure also has the value of alpha property available in the Advanced section. This value is set to the default value of 1 (see figure 622, page 528).




Figure 622 The advanced settings of the LSV staircase command

The **Spectroscopy** command stacked on the **LSV staircase** command is used to acquire the spectroscopy data during the measurement and collect all the of the measured data at the end of the measurement. This command has a number of additional pre-defined plots (see figure 623, page 529):

← DIO trigger

Default plots

	Enabled	Plot number	Options
Sample	<input checked="" type="checkbox"/>	4	Edit

Custom plots - +

Text	X	Y	Z	Enabled	Plot number	Options
Absorbance vs λ	Wavelength	Absorbance	Index	<input checked="" type="checkbox"/>	5	Edit
Transmittance vs λ	Wavelength	Transmittance	Index	<input checked="" type="checkbox"/>	6	Edit

Figure 623 Additional plots defined in the Spectroscopy command

- Sample: measured spectroscopy data versus wavelength
- Absorbance vs λ : calculated absorbance versus wavelength
- Transmittance vs λ : calculated transmittance versus wavelength



NOTICE

The absorbance and transmittance values are calculated using the dark and reference data collected by the two **Spectroscopy** commands located before the **LSV staircase** command in the procedure.

8.3 Voltammetric analysis

NOVA provides six default procedures for voltammetric analysis. These procedures can be used to perform a potential sweep with optional pulses or sinewaves and record the response of the cell.



CAUTION

All the procedures included in this group require the optional **IME663** or the optional **IME303**. Please refer to *Chapter 16.3.2.15* and *Chapter 16.3.2.14* for more information.

The following procedures are available:

- Sampled DC polarography
- Normal pulse voltammetry
- Differential pulse voltammetry



- Differential normal pulse voltammetry
- Square wave voltammetry
- AC voltammetry

8.3.1 Sampled DC polarography



CAUTION

This procedure requires a **IME663** *IME663 module* (see chapter 16.3.2.15, page 1130) or **IME303** *IME303 module* (see chapter 16.3.2.14, page 1124) connected to the Autolab. When this procedure is used without a **IME663** or **IME303**, an **error** will be displayed for the command.



NOTICE

To use this procedure without the optional **IME663** or the **IME303**, please delete the **Electrode preconditioning** command group and the **Equilibration** command group.

The default **Sampled DC polarography** procedure provides an example of a typical measurement using the *Sampled DC* method (see figure 624, page 530).




Figure 624 The default Sampled DC polarography procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

This procedure include two command groups, used to the control the mercury drop electrode.

- **Electrode preconditioning:** this command group is used to create new mercury drops at the beginning of the procedure. The commands in this group are used to purge the solution for the specified duration, create the specified number of new drops and switch the stirrer on (see figure 625, page 531).




Figure 625 The Electrode preconditioning group

- **Equilibration:** this command group is used to create an equilibration step in the procedure. The commands in this group are used to switch the stirrer off and wait for the specified amount of time (see figure 626, page 531).




Figure 626 The Equilibration group

The procedure has the following measurement properties, specified for the **Sampled DC** command (see figure 627, page 532):




Figure 627 The measurement properties of the Sampled DC command

- **Sampled DC**
 - Start potential: -1.2 V, versus reference electrode
 - Stop potential: 0.05 V, versus reference electrode
 - Step: 0.005 V
 - Interval time: 0.5 s

The procedure samples the following signals (see figure 628, page 532):



Figure 628 The sampler of the Sampled DC command

- WE(1).Current (averaged)
- WE(1).Potential

- Time

The procedure uses the following options (see figure 629, page 533):




Figure 629 The options of the Sampled DC command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA
- Counters
 - When counter = 1, Autolab control, Reset

The Counters option, using the procedure, is used to create a new drop with every potential step. The details of the Autolab control action are shown in Figure 630.



Figure 630 The Autolab control option triggered with the Counter option

The procedure plots the following data (see figure 631, page 534):




Figure 631 The plots of the Sampled DC command

- i vs E: WE(1).Current versus Potential applied

8.3.2 Normal pulse voltammetry



CAUTION

This procedure requires a **IME663** *IME663 module* (see chapter 16.3.2.15, page 1130) or **IME303** *IME303 module* (see chapter 16.3.2.14, page 1124) connected to the Autolab. When this procedure is used without a **IME663** or **IME303**, an **error** will be displayed for the command.



NOTICE

To use this procedure without the optional **IME663** or the **IME303**, please delete the **Electrode preconditioning** command group and the **Equilibration** command group.

The default **Normal pulse voltammetry** procedure provides an example of a typical measurement using the *Normal pulse* method (see figure 632, page 535).




Figure 632 The default Normal pulse voltammetry procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

This procedure include two command groups, used to the control the mercury drop electrode.

- **Electrode preconditioning:** this command group is used to create new mercury drops at the beginning of the procedure. The commands in this group are used to purge the solution for the specified duration, create the specified number of new drops and switch the stirrer on (see figure 633, page 536).




Figure 633 The Electrode preconditioning group

- Equilibration:** this command group is used to create an equilibration step in the procedure. The commands in this group are used to switch the stirrer off and wait for the specified amount of time (see figure 634, page 536).




Figure 634 The Equilibration group

The procedure has the following measurement properties, specified for the **Normal pulse** command (see figure 635, page 537):

Properties ➔

Normal pulse

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Base potential V_{REF} ▼

Step V

Normal pulse time s

Interval time s

Estimated number of points

Estimated duration s

Scan rate V/s

More

Figure 635 The measurement properties of the Normal pulse command

▪ **Normal pulse**

- Start potential: -1.2 V, versus reference electrode
- Stop potential: 0.07 V, versus reference electrode
- Step: 0.005 V
- Base potential: 0 V, versus reference electrode
- Normal pulse time: 0.07 s
- Interval time: 0.5 V

The procedure samples the following signals (see figure 636, page 538):




Figure 636 The sampler of the Normal pulse command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 637, page 538):




Figure 637 The options of the Normal pulse command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 638, page 538):




Figure 638 The plots of the Normal pulse command

- i vs E : WE(1).Current versus Potential applied

8.3.3 Differential pulse voltammetry



CAUTION

This procedure requires a **IME663** *IME663 module* (see chapter 16.3.2.15, page 1130) or **IME303** *IME303 module* (see chapter 16.3.2.14, page 1124) connected to the Autolab. When this procedure is used without a **IME663** or **IME303**, an **error** will be displayed for the command.



NOTICE

To use this procedure without the optional **IME663** or the **IME303**, please delete the **Electrode preconditioning** command group and the **Equilibration** command group.

The default **Differential pulse voltammetry** procedure provides an example of a typical measurement using the *Differential pulse* method (see figure 639, page 539).




Figure 639 The default Differential pulse voltammetry procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

This procedure include two command groups, used to the control the mercury drop electrode.

- **Electrode preconditioning**: this command group is used to create new mercury drops at the beginning of the procedure. The commands in this group are used to purge the solution for the specified duration, create the specified number of new drops and switch the stirrer on (see figure 640, page 540).




Figure 640 The Electrode preconditioning group

- **Equilibration:** this command group is used to create an equilibration step in the procedure. The commands in this group are used to switch the stirrer off and wait for the specified amount of time (see figure 641, page 540).




Figure 641 The Equilibration group

The procedure has the following measurement properties, specified for the **Differential pulse** command (see figure 642, page 541):



Properties ➔

Differential pulse

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Step V

Modulation amplitude V

Modulation time s

Interval time s

Estimated number of points

Estimated duration s

Scan rate V/s

More

Figure 642 The measurement properties of the Differential pulse command

▪ **Differential pulse**

- Start potential: -1.2 V, versus reference electrode
- Stop potential: 0.05 V, versus reference electrode
- Step: 0.005 V
- Modulation amplitude: 0.025 V
- Modulation time: 0.05 s
- Interval time: 0.5 s

The procedure samples the following signals (see figure 643, page 542):




Figure 643 The sampler of the Differential pulse command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 644, page 542):




Figure 644 The options of the Differential pulse command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 645, page 542):




Figure 645 The plots of the Differential pulse command

- δi vs E: $\delta[WE(1).Current]$ versus Potential applied

8.3.4 Differential normal pulse voltammetry



CAUTION

This procedure requires a **IME663** *IME663 module* (see chapter 16.3.2.15, page 1130) or **IME303** *IME303 module* (see chapter 16.3.2.14, page 1124) connected to the Autolab. When this procedure is used without a **IME663** or **IME303**, an **error** will be displayed for the command.



NOTICE

To use this procedure without the optional **IME663** or the **IME303**, please delete the **Electrode preconditioning** command group and the **Equilibration** command group.

The default **Differential normal pulse voltammetry** procedure provides an example of a typical measurement using the *Differential normal pulse* method (see figure 646, page 543).




Figure 646 The default Differential normal pulse voltammetry procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

This procedure include two command groups, used to the control the mercury drop electrode.

- **Electrode preconditioning:** this command group is used to create new mercury drops at the beginning of the procedure. The commands in this group are used to purge the solution for the specified duration, create the specified number of new drops and switch the stirrer on (see figure 647, page 544).




Figure 647 The Electrode preconditioning group

- Equilibration:** this command group is used to create an equilibration step in the procedure. The commands in this group are used to switch the stirrer off and wait for the specified amount of time (see figure 648, page 544).




Figure 648 The Equilibration group

The procedure has the following measurement properties, specified for the **Differential normal pulse** command (see figure 649, page 545):

Properties ➔

Differential normal pulse

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Base potential V_{REF} ▼

Step V

Modulation amplitude V

Modulation time s

Normal pulse time s

Interval time s

Estimated number of points

Estimated duration s

Scan rate V/s

More

Figure 649 The measurement properties of the Differential normal pulse command

- **Differential normal pulse**
 - Start potential: -1.2 V, versus reference electrode
 - Stop potential: 0.07 V, versus reference electrode
 - Step: 0.005 V
 - Base potential: 0 V, versus reference electrode
 - Modulation amplitude: 0.025 V
 - Normal pulse time: 0.025 s
 - Interval time: 0.5 V
 - Modulation time: 0.025 s

The procedure samples the following signals (see figure 650, page 546):



Figure 650 The sampler of the Differential normal pulse command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 651, page 546):




Figure 651 The options of the Differential normal pulse command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 652, page 546):




Figure 652 The plots of the Differential normal pulse command

- δi vs E: $\delta[WE(1).Current]$ versus Potential applied

8.3.5 Square wave voltammetry



CAUTION

This procedure requires a **IME663** *IME663 module* (see chapter 16.3.2.15, page 1130) or **IME303** *IME303 module* (see chapter 16.3.2.14, page 1124) connected to the Autolab. When this procedure is used without a **IME663** or **IME303**, an **error** will be displayed for the command.



NOTICE

To use this procedure without the optional **IME663** or the **IME303**, please delete the **Electrode preconditioning** command group and the **Equilibration** command group.

The default **Square wave voltammetry** procedure provides an example of a typical measurement using the *Square wave* method (see figure 653, page 547).




Figure 653 The default Square wave voltammetry procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

This procedure include two command groups, used to the control the mercury drop electrode.

- **Electrode preconditioning:** this command group is used to create new mercury drops at the beginning of the procedure. The commands in this group are used to purge the solution for the specified duration, create the specified number of new drops and switch the stirrer on (see figure 654, page 548).




Figure 654 The Electrode preconditioning group

- **Equilibration:** this command group is used to create an equilibration step in the procedure. The commands in this group are used to switch the stirrer off and wait for the specified amount of time (see figure 655, page 548).





Figure 655 The Equilibration group

The procedure has the following measurement properties, specified for the **Square wave** command (see figure 656, page 549):



Properties →

Square wave

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Step V

Modulation amplitude V

Frequency Hz

Estimated number of points

Interval time s

Estimated duration s

Scan rate V/s

More

Figure 656 The measurement properties of the Square wave command

■ **Square wave**

- Start potential: -1.2 V, versus reference electrode
- Stop potential: 0.07 V, versus reference electrode
- Step: 0.005 V
- Amplitude: 0.02 V
- Frequency: 25 Hz

The procedure samples the following signals (see figure 657, page 550):




Figure 657 The sampler of the Square wave command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 658, page 550):




Figure 658 The options of the Square wave command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 659, page 550):




Figure 659 The plots of the Square wave command

- δi vs E: $\delta[WE(1).Current]$ versus Potential applied

8.3.6 AC voltammetry



CAUTION

This procedure requires a **IME663** *IME663 module* (see chapter 16.3.2.15, page 1130) or **IME303** *IME303 module* (see chapter 16.3.2.14, page 1124) connected to the Autolab. When this procedure is used without a **IME663** or **IME303**, an **error** will be displayed for the command.



NOTICE

To use this procedure without the optional **IME663** or the **IME303**, please delete the **Electrode preconditioning** command group and the **Equilibration** command group.

The default **AC voltammetry** procedure provides an example of a typical measurement using the *AC voltammetry* method (see figure 660, page 551).



Figure 660 The default AC voltammetry procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

This procedure include two command groups, used to the control the mercury drop electrode.

- **Electrode preconditioning:** this command group is used to create new mercury drops at the beginning of the procedure. The commands in this group are used to purge the solution for the specified duration, create the specified number of new drops and switch the stirrer on (see figure 661, page 552).




Figure 661 The Electrode preconditioning group

- **Equilibration:** this command group is used to create an equilibration step in the procedure. The commands in this group are used to switch the stirrer off and wait for the specified amount of time (*see figure 662, page 552*).




Figure 662 The Equilibration group

The procedure has the following measurement properties, specified for the **AC voltammetry** command (*see figure 663, page 553*):

Properties ➔

AC voltammetry

Command name

Start potential V_{REF} ▼

Stop potential V_{REF} ▼

Step V

Modulation amplitude V_{RMS}

Modulation time s

Frequency Hz

Interval time s

Harmonic

Estimated number of points

Estimated duration s

Scan rate V/s

More

Figure 663 The measurement properties of the AC voltammetry command

▪ **AC voltammetry**

- Start potential: -1.2 V, versus reference electrode
- Stop potential: 0.05 V, versus reference electrode
- Step: 0.005 V
- Modulation amplitude: 0.025 V RMS
- Modulation time: 0.2 s
- Frequency: 37 Hz
- Interval time: 0.6 s
- Harmonic: 1

The procedure samples the following signals (see figure 664, page 554):



Figure 664 The sampler of the AC voltammetry command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 665, page 554):




Figure 665 The options of the AC voltammetry command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 666, page 555):




Figure 666 The plots of the AC voltammetry command

- i(AC) vs E: AC current versus Potential applied

8.4 Chrono methods

NOVA provides nine default procedures for chrono methods. These procedures can be used to perform time resolved measurements. Some of these procedures require optional hardware extensions.

The following procedures are available:

- Chrono amperometry ($\Delta t > 1$ ms)
- Chrono coulometry ($\Delta t > 1$ ms) (requires the **FI20** module or the **on-board integrator**, please refer to *Chapter 16.3.2.11* for more information)
- Chrono potentiometry ($\Delta t > 1$ ms)
- Chrono amperometry fast
- Chrono coulometry fast (requires the **FI20** module or the **on-board integrator**, please refer to *Chapter 16.3.2.11* for more information)
- Chrono potentiometry fast
- Chrono amperometry high speed (requires the **ADC10M** or **ADC750** module, please refer to *Chapter 16.3.2.1* for more information)
- Chrono potentiometry high speed (requires the **ADC10M** or **ADC750** module, please refer to *Chapter 16.3.2.1* for more information)
- Chrono charge discharge

Properties

Record signals

Command name

Duration s

Interval time s

Estimated number of points

Estimated duration s

More

Figure 668 The measurement properties of the Record signals command

- **Record signals**

- Duration: 5 s
- Interval time: 0.01 s

The procedure samples the following signals (see figure 669, page 557):

Record signals

Fast options

Sampler

Options

Plots

Signal	Sample	Average	d/dt
WE(1).Current	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
WE(1).Potential	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Power	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Resistance	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Charge	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Time	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Sample alternating

Figure 669 The sampler of the Record signals command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 670, page 558):




Figure 670 The options of the Record signals command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 1 μA

The procedure plots the following data (see figure 671, page 558):




Figure 671 The plots of the Record signals command

- i vs t: WE(1).Current versus time

8.4.2 Chrono coulometry ($\Delta t > 1$ ms)



CAUTION

This procedure requires the optional **FI20** module or the **on-board integrator FI20 module** (see chapter 16.3.2.11, page 1083).

The default **Chrono coulometry ($\Delta t > 1$ ms)** procedure provides an example of a typical chrono coulometric measurement using a sequence of potential steps (see figure 672, page 559).




Figure 672 The default Chrono coulometry ($\Delta t > 1$ ms) procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure uses a sequence of three potential values (specified through the **Apply** command) followed by three **Record signals** commands.



NOTICE

The smallest possible interval time for the **Record signals** command is 1.3 ms.

The potential values applied are 0 V, 0.5 V and -0.5 V. The **Record signals** commands have the following measurement properties (see figure 673, page 560):




Figure 675 The plots of the Record signals command

- Q vs t: Integrator(1).Charge versus time

8.4.3 Chrono potentiometry ($\Delta t > 1$ ms)

The default **Chrono potentiometry ($\Delta t > 1$ ms)** procedure provides an example of a typical chrono potentiometric measurement using a sequence of current steps (see figure 676, page 561).




Figure 676 The default Chrono potentiometry ($\Delta t > 1$ ms) procedure



NOTICE

The galvanostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure uses a sequence of three current values (specified through the **Apply** command) followed by three **Record signals** commands.




Figure 678 The sampler of the Record signals command

- WE(1).Potential (averaged)
- WE(1).Current
- Time

The procedure plots the following data (see figure 679, page 563):




Figure 679 The plots of the Record signals command

- E vs t: WE(1)Potential versus time

Properties

Chrono methods

Command name

Mode

Number of repeats

Total duration s

Estimated number of points

Estimated duration s

Figure 681 The measurement properties of the Chrono methods command

The procedure uses a sequence of four potential values applied and measured through the **Chrono methods** command (see figure 682, page 565).

Chrono methods

Sequence

Sampler

Plots

Step

Step

Step

Basic

Text

Duration s

Sample

Interval time s

Estimated number of points

Potential V_{REF}

Advanced

Figure 682 The details of the Chrono methods command

The potential values applied are 0 V, 0.3 V, -0.3 V and 0 V, applied versus the reference potential. The **Chrono methods** command has the following measurement properties (see figure 682, page 565):



- **Chrono methods**
 - Duration: 0.01 s
 - Interval time: 0.0001 s

The procedure samples the following signals (see figure 683, page 566):




Figure 683 The sampler of the Chrono methods command

- WE(1).Current (averaged)
- Time

The procedure plots the following data (see figure 684, page 566):




Figure 684 The plots of the Chrono methods command

- i vs t: WE(1).Current versus time

8.4.5 Chrono coulometry fast



CAUTION

This procedure requires the optional **FI20** module or the **on-board integrator FI20 module** (see chapter 16.3.2.11, page 1083).

The default **Chrono coulometry fast** procedure provides an example of a typical chrono coulometric measurement using a sequence of potential steps, with a short interval time (see figure 685, page 567).




Figure 685 The default Chrono coulometry fast procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure uses a sequence of two potential values (specified in the **Chrono methods** command).



NOTICE

The smallest possible interval time for the **Chrono methods** command is 100 μs .

The **Chrono methods** command has the following measurement properties (see figure 686, page 568):

- **Chrono methods**
 - Duration: 0.5 s
 - Interval time: 0.001 s

The procedure samples the following signals (see figure 688, page 569):




Figure 688 The sampler of the Chrono methods command

- WE(1).Current
- Integrator(1).Charge
- Time

The procedure plots the following data (see figure 689, page 569):




Figure 689 The plots of the Chrono methods command

- Q vs t: Integrator(1).Charge versus time

Properties

Chrono methods

Command name

Mode

Number of repeats

Total duration s

Estimated number of points

Estimated duration s

Figure 691 The measurement properties of the Chrono methods command

The procedure uses a sequence of four current values applied and measured through the **Chrono methods** command (see figure 692, page 571).

Chrono methods

Sequence

Sampler

Plots

Step

Basic

Text

Duration s

Sample

Interval time s

Estimated number of points

Current A

Advanced

Figure 692 The details of the Chrono methods command

The current values applied are 0 A, 0.003 A, -0.003 A and 0 A. The **Chrono methods** command has the following measurement properties (see figure 692, page 571):



- **Chrono methods**
 - Duration: 0.01 s
 - Interval time: 0.0001 s

The procedure samples the following signals (see figure 693, page 572):




Figure 693 The sampler of the Chrono methods command

- WE(1).Potential (averaged)
- Time

The procedure plots the following data (see figure 694, page 572):




Figure 694 The plots of the Chrono methods command

- E vs t: WE(1).Potential versus time

8.4.7 Chrono amperometry high speed



CAUTION

This procedure requires the optional **ADC10M** or **ADC750** module *ADC10M module (see chapter 16.3.2.1, page 998)*.

The default **Chrono amperometry high speed** procedure provides an example of a typical chrono amperometric measurement using a sequence of potential steps, with a short interval time (*see figure 695, page 573*).




Figure 695 The default Chrono amperometry high speed procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control (see chapter 7.2.1, page 236)*.

The procedure uses a sequence of four potential values (specified in the **Chrono methods** command).



NOTICE

The smallest possible interval time for the **Chrono methods** command in *high speed* mode is 100 ns with the **ADC10M** and 1.33 μ s with the **ADC750**.

The **Chrono methods** command has the following measurement properties (*see figure 696, page 574*):

The procedure samples the following **ADC10M** or **ADC750** settings (see figure 698, page 575):




Figure 698 The ADC10M or ADC750 settings of the Chrono methods command

- Channel 1: WE(1).Potential, Gain 1, unfiltered
- Channel 2: WE(1).Current, Gain 1, unfiltered

The procedure plots the following data (see figure 699, page 576):




Figure 700 The default Chrono potentiometry high speed procedure



NOTICE

The galvanostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure uses a sequence of four current values (specified in the **Chrono methods** command).



NOTICE

The smallest possible interval time for the **Chrono methods** command in *high speed* mode is 100 ns with the **ADC10M** and 1.33 μ s with the **ADC750**.

The **Chrono methods** command has the following measurement properties (see figure 701, page 578):

The procedure samples the following **ADC10M** or **ADC750** settings (see figure 703, page 579):




Figure 703 The ADC10M or ADC750 settings of the Chrono methods command

- Channel 1: WE(1).Potential, Gain 1, unfiltered
- Channel 2: WE(1).Current, Gain 1, unfiltered

The procedure plots the following data (see figure 704, page 580):



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure uses a **Repeat** command, used in the *Repeat n times* mode, containing a sequence of two potential values (specified through the **Apply** command) followed by two **Record signals** commands.



NOTICE

The smallest possible interval time for the **Record signals** command is 1.3 ms.

The potential values applied are 1.2 V and -0.5 V. The **Record signals** commands have the following measurement properties (see figure 706, page 581):

Properties	
Record signals	
Command name	Record signals
Duration	2,5 s
Interval time	0,01 s
Estimated number of points	250
Estimated duration	2,5 s
More	

Figure 706 The measurement properties of the Record signals command

- **Record signals**
 - Duration: 2.5 s
 - Interval time: 0.01 s

The procedure samples the following signals (see figure 707, page 582):




Figure 707 The sampler of the Record signals command

- WE(1).Current (averaged)
- WE(1).Potential
- Time

The procedure uses the following options (see figure 708, page 582):




Figure 708 The options of the Record signals command

- Automatic current ranging
 - Highest current range: 1 mA
 - Lowest current range: 1 μ A

The procedure plots the following data (see figure 709, page 583):




Figure 709 The plots of the Record signals command

- i vs t: WE(1).Current versus time

8.4.10 Chrono charge discharge galvanostatic

The default **Chrono charge discharge galvanostatic** procedure provides an example of a typical galvanostatic charge and discharge measurement using a sequence of applied current steps (see figure 710, page 583).




Figure 710 The default Chrono charge discharge galvanostatic procedure



NOTICE

The duration of the **Record signals** commands are deliberately long because there is a **Cutoff** on the Potential signal set in the **Options** tab of the **Record signals** command tiles. In a typical charge discharge galvanostatic measurement, the charge step continues until a desired potential maximum is reached, then the procedure moves on to the discharge step. The discharge step continues until a desired potential minimum is reached, then the procedure moves on to the next charge step, or is completed. For more information about the potential **cutoff** please continue to the **Options** of the **Record signals** command. *Chrono charge discharge galvanostatic (see chapter 8.4.10, page 583)*

The **Record signals** commands sample the following signals:




Figure 712 The sampler of the Record signals command

- WE(1).Potential
- Time

The first **Record signals** command tile (charge step) uses the following options:




Figure 715 The plots of the Record signals commands

- E vs t: WE(1).Potential versus time

8.5 Potentiometric stripping analysis

NOVA provides two default procedures for potentiometric stripping analysis (PSA).

The following procedures are available:

- Potentiometric stripping analysis
- Potentiometric stripping analysis (Constant current)

8.5.1 Potentiometric stripping analysis

The default **Potentiometric stripping analysis** procedure provides an example of a typical measurement using the **PSA** command (see figure 716, page 587).




Figure 716 The default Potentiometric stripping analysis procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).



The procedure has the following measurement properties, specified for the **PSA** command (see figure 717, page 588):

Figure 717 The measurement properties of the PSA command

- **PSA**
 - Potential limit: -0.001 V, versus reference electrode
 - Maximum time: 10 s
 - Filter: on
 - Filter time: 0.020 s or 0.0166 s

The procedure plots the following data (see figure 718, page 588):

Text	X	Y	Z	Plot number	Options
$\delta t/\delta E$ vs E	WE(1).Potential	$\delta t/\delta E$	$\delta t/\delta E$	1	Edit
E vs t	Time	WE(1).Potential	WE(1).Potential	2	Edit

Figure 718 The plots of the PSA command

- $\delta t/\delta E$ vs E: $\delta t/\delta WE(1).Potential$ versus WE(1).Potential
- E vs t: WE(1).Potential versus time

8.5.2 Potentiometric stripping analysis constant current

The default **Potentiometric stripping analysis** procedure provides an example of a typical measurement using the **PSA constant current** command (see figure 719, page 589).




Figure 719 The default Potentiometric stripping analysis (Constant current) procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure has the following measurement properties, specified for the **PSA constant current** command (see figure 720, page 589):

The 'Properties' dialog box for 'PSA Constant current' contains the following fields:

- Command name: PSA Constant current
- Mode: Constant current (dropdown menu)
- Constant current: 1E-06 A
- Potential limit: 0,8 V
- Maximum time: 10 s
- Filter: [Progress bar]
- Filter time: 0,02 s
- Estimated duration: 10 s

A 'More' button is located at the bottom right of the dialog.

Figure 720 The measurement properties of the PSA constant current command



- **PSA**
 - Constant current: 1 μA
 - Potential limit: 0.8 V, versus reference electrode
 - Maximum time: 10 s
 - Filter: on
 - Filter time: 0.020 s or 0.0166 s

The procedure plots the following data (see figure 721, page 590):

PSA Constant current					
Custom plots					+
Text	X	Y	Z	Plot number	Options
$\delta t/\delta E$ vs E	WE(1).Potential	$\delta t/\delta E$	$\delta t/\delta E$	1	Edit
E vs t	Time	WE(1).Potential	WE(1).Potential	2	Edit

Figure 721 The plots of the PSA constant current command

- $\delta t/\delta E$ vs E: $\delta t/\delta \text{WE}(1).\text{Potential}$ versus WE(1).Potential
- E vs t: WE(1).Potential versus time

8.6 Impedance spectroscopy

NOVA provides six default procedures for impedance spectroscopy. These procedures can be used to perform a cyclic potential or current scan and record the response of the cell.



CAUTION

These procedures require the optional **FRA32M** or **FRA2** module *FRA32M module (see chapter 16.3.2.13, page 1112)*.

The following procedures are available:

- FRA impedance potentiostatic
- FRA impedance galvanostatic
- FRA potential scan
- FRA current scan
- FRA time scan potentiostatic
- FRA time scan galvanostatic

Additionally, a default procedure for Electrochemical Frequency Modulation measurement is also included in this group.

8.6.1 FRA impedance potentiostatic

The default FRA impedance potentiostatic procedure provides an example of an electrochemical impedance spectroscopy measurement in potentiostatic conditions (see figure 722, page 591).




Figure 722 The default FRA impedance potentiostatic procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure has the following measurement properties, specified for the **FRA measurement** command (see figure 723, page 592):

FRA measurement

Sampler

Options
Plots
Summary

Basic

Maximum integration time 0,125 s

Minimum number of integration cycles 1

Sample time domain

Sample frequency domain

Sample DC

Calculate admittance

Accuracy

Advanced

Transfer function Re - j Im

Lowest bandwidth High stability

Number of cycles to reach steady state 10

Maximum time to reach steady state 1 s

With a minimum fraction of a cycle 0

Automatic amplitude correction

Iterative

Amplitude threshold percentage 5 %

Automatic resolution correction

Iterative

Minimum resolution 32 %

Maximum amount of re-measurements 25

Figure 724 The sampler of the FRA measurement command

- DC signals

The procedure uses the following options (see figure 725, page 593):

FRA measurement

Sampler
Options
Plots
Summary

Automatic current ranging

Highest current range 10 mA

Lowest current range 100 nA

Figure 725 The options of the FRA measurement command

- Automatic current ranging
 - Highest current range: 10 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 726, page 594):

Properties ➔

FRA measurement

Command name

First applied frequency Hz

Last applied frequency Hz

Number of frequencies per decade

Frequency step type ▼

Amplitude A_{RMS}

Use RMS amplitude

Wave type ▼

Input connection ▼

Estimated duration s

More

Figure 728 The measurement properties of the FRA measurement command

- **FRA measurement**
 - Start frequency: 100 kHz
 - Stop frequency: 0.1 Hz
 - Number of frequencies per decade: 10
 - Amplitude: 10 μ A
 - Use RMS amplitude: yes
 - Wave type: sine
 - Input connection: internal

The procedure samples the following signals (see figure 729, page 596):




Figure 729 The sampler of the FRA measurement command

- DC signals

The procedure plots the following data (see figure 730, page 596):




Figure 730 The plots of the FRA measurement command

- Nyquist impedance
- Bode

8.6.3 FRA potential scan

The default FRA potential scan procedure provides an example of an electrochemical impedance spectroscopy measured repeated for a pre-defined series of DC potentials (see figure 731, page 597).




Figure 731 The default FRA potential scan procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The FRA potential scan procedure performs an impedance measurement at twelve different potential values. The potential values are set using a **Repeat** command.

The **Repeat** command is used in the *Repeat for multiple values* mode and is preconfigured to cycle through twelve potentials values, starting at 1.2 V until 0.1 V, using a linear distribution (see figure 732, page 598).



← Repeat for multiple values

Values ✕ ✎ +

	Potential (V)
1	1,2
2	1,1
3	1
4	0,9
5	0,8
6	0,7
7	0,6
8	0,5
9	0,4
10	0,3
11	0,2
12	0,1

Add range

Begin value

End value

Number of values

Distribution

Add range

Figure 732 The repeat loop using the default FRA potential scan procedure

The *Potential* parameter, created by the **Repeat** command, is linked to the **Apply** command included in the repeat loop (see figure 733, page 598).

← Edit links

Repeat for multiple values

Repetition number

Potential (V)

Apply 1,2 V

Potential

Figure 733 The link used to set the potential values

The procedure has the following measurement properties, specified for the **FRA measurement** command (see figure 734, page 599):

The screenshot shows a 'Properties' dialog box for the 'FRA measurement' command. The settings are as follows:

Property	Value	Unit
Command name	FRA measurement	
First applied frequency	1E+05	Hz
Last applied frequency	1	Hz
Number of frequencies	1	per decade
Frequency step type	Points per decade	
Amplitude	0,01	V _{RMS}
Use RMS amplitude	<input checked="" type="checkbox"/>	
Wave type	Sine	
Input connection	Internal	
Estimated duration	15,736	s

A 'More' button is located at the bottom right of the dialog.

Figure 734 The measurement properties of the FRA measurement command

- **FRA measurement**
 - Start frequency: 100 kHz
 - Stop frequency: 1 Hz
 - Number of frequencies per decade: 1
 - Amplitude: 0.010 V
 - Use RMS amplitude: yes
 - Wave type: sine
 - Input connection: internal

The procedure samples the following signals (see figure 735, page 600):




Figure 735 The sampler of the FRA measurement command

- DC signals

The procedure uses the following options (see figure 736, page 600):




Figure 736 The options of the FRA measurement command

- Automatic current ranging
 - Highest current range: 100 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 737, page 601):




Figure 737 The plots of the FRA measurement command

- Nyquist impedance
- Bode

The **Potential scan FRA data** command located at the end of the procedure is used to automatically generate Mott-Schottky plots. For more information on this command, please refer to *Chapter 7.9.5*.

8.6.4 FRA current scan

The default FRA current scan procedure provides an example of an electrochemical impedance spectroscopy measured repeated for a pre-defined series of DC currents (see figure 738, page 601).




Figure 738 The default FRA current scan procedure



NOTICE

The galvanostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The FRA current scan procedure performs an impedance measurement at seven different current values. The current values are set using a **Repeat** command.

The **Repeat** command is used in the *Repeat for multiple values* mode and is preconfigured to cycle through seven current values, starting at - 3 mA until 3 mA, using a linear distribution (see figure 739, page 602).




Figure 739 The repeat loop using the default FRA current scan procedure

The *Current* parameter, created by the **Repeat** command, is linked to the **Apply** command included in the repeat loop (see figure 740, page 603).






Figure 740 The link used to set the current values

The procedure has the following measurement properties, specified for the **FRA measurement** command (see figure 741, page 603):

Properties +

FRA measurement

Command name

First applied frequency Hz

Last applied frequency Hz

Number of frequencies per decade

Frequency step type

Amplitude A_{RMS}

Use RMS amplitude

Wave type

Input connection

Estimated duration s

Figure 741 The measurement properties of the FRA measurement command



- **FRA measurement**
 - Start frequency: 100 kHz
 - Stop frequency: 1 Hz
 - Number of frequencies per decade: 1
 - Amplitude: 10 μ A
 - Use RMS amplitude: yes
 - Wave type: sine
 - Input connection: internal

The procedure samples the following signals (see figure 742, page 604):




Figure 742 The sampler of the FRA measurement command

- DC signals

The procedure plots the following data (see figure 743, page 605):




Figure 743 The plots of the FRA measurement command

- Nyquist impedance
- Bode

8.6.5 FRA time scan potentiostatic

The default FRA time scan procedure provides an example of an electrochemical impedance spectroscopy measured at fixed time intervals, in potentiostatic mode (see figure 744, page 605).




Figure 744 The default FRA time scan potentiostatic procedure



NOTICE

The potentiostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure uses a **Repeat** command, used in *Timed repeat* mode, for a pre-defined duration of 200 s and interval time of 20 s (see figure 745, page 606).



Properties ➔

Repeat timed

Command name

Repeat

Number of repetitions

Duration s

Interval time s

Figure 745 The properties of the Repeat command

The procedure has the following measurement properties, specified for the **FRA measurement** command (see figure 746, page 606):

Properties ➔

FRA measurement

Command name

First applied frequency Hz

Last applied frequency Hz

Number of frequencies per decade

Frequency step type

Amplitude V_{RMS}

Use RMS amplitude

Wave type

Input connection

Estimated duration s

Figure 746 The measurement properties of the FRA measurement command

- **FRA measurement**

- Start frequency: 100 kHz
- Stop frequency: 1 kHz
- Number of frequencies per decade: 1
- Amplitude: 0.010 V
- Use RMS amplitude: yes
- Wave type: sine
- Input connection: internal

The procedure samples the following signals (see figure 747, page 607):

FRA measurement

Sampler

Basic

Maximum integration time 0,125 s

Minimum number of integration cycles 1

Sample time domain

Sample frequency domain

Sample DC

Calculate admittance

Accuracy

Advanced

Transfer function Re - j Im

Lowest bandwidth High stability

Number of cycles to reach steady state 10

Maximum time to reach steady state 1 s

With a minimum fraction of a cycle 0

Automatic amplitude correction

Iterative

Amplitude threshold percentage 5 %

Automatic resolution correction

Iterative

Minimum resolution 32 %

Maximum amount of re-measurements 25

Figure 747 The sampler of the FRA measurement command

- DC signals

The procedure uses the following options (see figure 748, page 607):

FRA measurement

Options

Automatic current ranging

Highest current range 10 mA

Lowest current range 100 nA

Figure 748 The options of the FRA measurement command



- Automatic current ranging
 - Highest current range: 10 mA
 - Lowest current range: 100 nA

The procedure plots the following data (see figure 749, page 608):




Figure 749 The plots of the FRA measurement command

- Nyquist impedance
- Bode

Additionally, the procedure gathers all the measured data points and plots the following time resolved data:

- Z vs t
- -phase vs t

8.6.6 FRA time scan galvanostatic

The default FRA time scan procedure provides an example of an electro-chemical impedance spectroscopy measured at fixed time intervals, in galvanostatic mode (see figure 750, page 608).




Figure 750 The default FRA time scan galvanostatic procedure



NOTICE

The galvanostatic mode is selected at the beginning of the procedure using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

The procedure uses a **Repeat** command, used in *Timed repeat* mode, for a pre-defined duration of 200 s and interval time of 20 s (see figure 751, page 609).

Properties ☰

Repeat timed

Command name

Repeat

Number of repetitions

Duration s

Interval time s

Figure 751 The properties of the Repeat command

The procedure has the following measurement properties, specified for the **FRA measurement** command (see figure 752, page 610):




Figure 753 The sampler of the FRA measurement command

- DC signals

The procedure plots the following data (see figure 754, page 611):




Figure 754 The plots of the FRA measurement command

- Nyquist impedance
- Bode

Properties ☰

Electrochemical Frequency Modulation

Command name

Base frequency Hz

Multiplier 1

Multiplier 2

Frequency 1 Hz

Frequency 2 Hz

Amplitude V_{TOP}

Number of cycles

Model ▼

Density g/cm³

Equivalent weight g/mol

Surface area cm²

Estimated duration s

Figure 756 The measurement properties of the Electrochemical Frequency Modulation command

▪ **Electrochemical Frequency Modulation**

- Base frequency: 0,1 kHz
- Multiplier 1: 2
- Multiplier 2: 5
- Amplitude: 0.010 V
- Number of cycles: 4
- Model: Activation Control
- Density: 7,87 g/cm³
- Equivalent weight: 27,92 g/mol
- Surface area: 1 cm²

The procedure plots the following data (see figure 757, page 614):



← Electrochemical Frequency Modulation

Default plots

	Enabled	Plot number	Options
Nyquist $-Z''$ vs Z'	<input type="checkbox"/>	<input type="text"/>	
Bode modulus	<input type="checkbox"/>	<input type="text"/>	
Bode phase	<input type="checkbox"/>	<input type="text"/>	
E(AC) vs t	<input checked="" type="checkbox"/>	<input type="text"/>	Edit
i(AC) vs t	<input checked="" type="checkbox"/>	<input type="text"/>	Edit
E(resolution) vs t	<input type="checkbox"/>	<input type="text"/>	
i(resolution) vs t	<input type="checkbox"/>	<input type="text"/>	
Lissajous	<input type="checkbox"/>	<input type="text"/>	

Custom plots +

Text	X	Y	Z	Enabled	Plot number	Options
E vs f	Frequency domain	Potential frequency domain		<input checked="" type="checkbox"/>	<input type="text"/>	Edit
i vs f	Frequency domain	Current density		<input checked="" type="checkbox"/>	<input type="text"/>	Edit

Figure 757 The plots of the Electrochemical Frequency Modulation command

- E(AC) vs t
- i(AC) vs t
- E vs f
- i vs f

9 Additional measurement command properties

Most of the measurement commands in NOVA have additional properties which can be accessed through the [More](#) button, as shown *Figure 758*.

The screenshot shows a 'Properties' dialog box for a 'CV staircase' measurement command. The parameters are as follows:

Property	Value	Unit
Command name	CV staircase	
Start potential	0	V _{REF}
Upper vertex potential	1	V _{REF}
Lower vertex potential	-1	V _{REF}
Stop potential	0	V _{REF}
Number of scans	1	
Scan rate	0,1	V/s
Step	0,00244	V
Interval time	0,0244	s
Estimated number of points	1640	
Estimated duration	40,016	s
Number of stop crossings	2	

A 'More' button is located at the bottom right of the dialog box.

Figure 758 Additional properties are provided by most measurement commands



NOTICE

Not all measurement command provide additional options and the provided option may change, depending on the measurement command.

This section provides information on the following additional properties:

- **Sampler:** the sampler defines which signals to sample during the measurement.



- **Options:** the options are additional measurement settings that affect how the data is measured.
- **Plots:** the plots define how the measured data should be plotted.
- **Advanced:** advanced acquisition properties used during the measurement.



NOTICE

The **Advanced** properties are only available for the **CV staircase** and **LSV staircase** commands.

9.1 Sampler

The sampler defines which signals are measured or calculated by the command and how these signals should be measured. For each available signal,  toggles are provided to control the sampler settings (see figure 759, page 616).




Figure 759 The sampler defines which signals are measured or calculated by the command



NOTICE

The list of signals displayed in the sampler depends on the hardware setup.

Depending on the type of signal, the following settings can be defined in the sampler:

- **Sample:** this setting defines that the signal is sampled by the command. A single analog-to-digital conversion is performed for a sampled signal.
- **Average:** this setting defines that a sampled signal must be averaged. When a sampled signal is averaged, as many analog-to-digital conversions are performed and an averaged value is stored. Averaging a signal significantly improves the signal-to-noise ratio.
- **d/dt:** this setting defines that the time derivative of a sampled signal must be calculated.



NOTICE

The **average** setting is only available for signals that can be sampled. Some of the signals provided in the sampler are calculated (Power, Resistance and Charge) while other signals are digitized by a dedicated optional module (EQCM signals).



NOTICE

Up to six signals can be averaged during a measurement.

Additionally, a toggle is provided for the **Sample alternating** setting, below the sampler table (*see figure 759, page 616*). This setting defines how averaged signals are sampled by the command:

- **Sample alternating off:** when this setting is off, all averaged signals are sampled, sequentially. The WE(1).Current signal is always sampled last.
- **Sample alternating on:** when this setting is on, all averaged signals are sampled at the same time, alongside the WE(1).Current signal.

Clicking the button opens a new screen that provides additional information on the exact timing of the sampler, in μs . The signals are provided in a table, as shown in *Figure 760*.



← Sampler

Signal	Start time (μs)	Duration (μs)
Time	18000	0
WE(1).Potential	18000	200
WE(1).Current	18200	6100
WE(1).Power WE(1).Resistance	24300	100
Total duration		6400

Figure 760 Detailed view of timing used by the sampler

The signals are listed in the table in chronological order to sampling. The Start time column provides the time, in μs , after which the sampling of the signal starts, with respect to the beginning of the interval time. The Duration column provides the duration, in μs , during which each signal is sampled. Depending on the type of signal and on the sampling method, the following durations are used:

- **Time:** the duration of the sampling of the Time signal is always 0 μs .
- **Sampled signals:** the duration of the sampling of signals that are not averaged is at most 200 μs .
- **Calculated signals:** the duration of the calculations carried out for the determination of calculated signals is at most 100 μs .

9.2 Automatic current ranging

The **Automatic current ranging** option specifies which of the available current ranges can be used by the measurement command. When this option is used, the instrument will automatically select the most suitable current range available. The instrument will also change the current range in the following cases:

- **Current overload:** the measured current exceeds the current overload threshold. The active current range is adjusted to the next available higher range.
- **Current underload:** the measured current exceeds the current underload threshold. The active current range is adjusted to the next available lower range.



NOTICE

Five consecutive overload or underload detections are required to trigger a change in current range.

Automatic current ranging settings are defined in the dedicated table, in the **Automatic current ranging** sub-panel (see figure 761, page 619).




Figure 761 The Automatic current ranging option is defined in a dedicated sub-panel

Three properties can be specified for each working electrode:

- **Enabled:** a toggle used to set the Automatic current ranging on or off.
- **Highest current range:** defines the highest possible current range, using the provided drop-down list.
- **Lowest current range:** defines the lowest possible current range, using the provided drop-down list.

Additionally, the **Optimize current range** toggle is located below the table. When this setting is on, the instrument will automatically adjust the current range of each electrode for which the **Automatic current ranging** option is enabled to the most suitable current range **before** the command starts measuring.



NOTICE

It is highly recommended to use the **Optimize current range** option whenever using the **Automatic current ranging** option in order to ensure that each measurement starts in the most suitable current range.

- **Only once:** specifies if the cutoff action should be executed only once or each time the cutoff condition is met, using the provided toggle .
- **Detections:** defines the number of consecutive detections required to trigger the cutoff.
- **Link as:** defines a unique name for the cutoff **Value** that can be used to link to other command parameters in the procedure editor.



NOTICE


The **Reverse scan** direction action is only available for the **LSV staircase** and the **CV staircase** commands.

9.3.1 Cutoff configuration

The following steps describe how to add and configure a cutoff.

1 Add a cutoff to the list

Click on the **+** button to add a cutoff to the table.




NOTICE

A cutoff on the WE(1).Current is automatically generated.

2 Specify the signal

Click on the cell of the **Signal** column and select the signal to use in the cutoff using the provided drop-down list.





3 Specify the inequality

Click on the cell of the **When** column and select the inequality to use in the cutoff using the provided drop-down list (< or >).

Cutoffs							+
Signal	When	Value	Action	Only once	Detections	Link as	
WE(1).Current	>	1	Stop command	<input type="checkbox"/>	4		

4 Specify the value

Specify the threshold value for the signal used in the cutoff in the corresponding cell of the **Value** column.

Cutoffs							+
Signal	When	Value	Action	Only once	Detections	Link as	
WE(1).Current	>	0,001	Stop command	<input type="checkbox"/>	4		

5 Specify the action

Click on the cell of the **Action** column and select the action to use in the cutoff using the provided drop-down list.

Cutoffs							+
Signal	When	Value	Action	Only once	Detections	Link as	
WE(1).Current	>	1	Stop command	<input type="checkbox"/>	4		

6 Set the only once property

Use the provided toggle to define if the cutoff should be triggered only once or continuously.

Cutoffs							+
Signal	When	Value	Action	Only once	Detections	Link as	
WE(1).Current	>	1	Stop command	<input checked="" type="checkbox"/>	4		

7 Set the number of detections

Specify the number of detections value for the signal used in the cutoff in the corresponding cell of the **Detections** column.



Cutoffs							+
Signal	When	Value	Action	Only once	Detections	Link as	
WE(1).Current	>	1	Stop command	<input type="checkbox"/>	10		

8 Specify a unique linkable name

If required, a *unique* linkable name can be specified in the **Link as** column. If a name is specified, the threshold specified in the Value column can be linked to another command parameter in the procedure. Using this link, the actual threshold value can be modified during the execution of the procedure.

Cutoffs							+
Signal	When	Value	Action	Only once	Detections	Link as	
WE(1).Current	>	1	Stop command	<input type="checkbox"/>	10	Link	



NOTICE

The **Link as** property is optional and can be left empty if no link is required.



NOTICE

To remove a cutoff from the table, select the row of the cutoff and click the  button above the **Cutoffs** table.

9.3.2 Combining cutoffs

It is possible to define more than one cutoff condition in the **Cutoffs** table. Depending on how the cutoffs conditions are defined, it is possible to arrange two or more cutoffs in two different ways:

- **OR arrangement:** each cutoff condition is defined as a standalone cutoff. The action defined for each of them is triggered whenever the corresponding threshold value is reached. This corresponds to a *OR* logical operator. The measurement command will be affected by each individual cutoff separately.
- **AND arrangement:** the two or more cutoff conditions can be joined with a *AND* action in order to trigger a single action when each of the involved cutoffs is triggered.

Figure 763 shows an example of three cutoff conditions. The first cutoff monitors the value of the WE(1).Current signal and forces the command to stop if this signal exceeds 1 mA. The second cutoff monitors the WE(1).Potential signal. When the value of this signal exceeds 1.2 V, the



third cutoff will be monitored. When the third cutoff, specified on the WE(1).Charge signal is triggered, the complete procedure will be stopped.

Cutoffs						
Signal	When	Value	Action	Only once	Detections	Link as
WE(1).Current	>	0,001	Stop command	<input type="checkbox"/>	10	
WE(1).Potential	>	1,2	And	<input type="checkbox"/>	4	
WE(1).Charge	>	3	Stop complete procedure	<input type="checkbox"/>	4	

Figure 763 Multiple cutoffs



NOTICE

The second and third cutoff shown in *Figure 763* are connected by a grey line on the left-hand side of the table, indicating that both cutoffs have a **AND** relationship.

9.4 Counters

Counters can be used during a measurement to perform dedicated actions whenever a condition associated with the counter is triggered. Each counter accumulates during a measurement, and it is possible to assign a specific instrumental action when a counter reaches a user defined value.

Since the counters are intrinsically linked to the measured data, the events triggered by the counters are directly correlated to the data points.

Counters are defined in the dedicated table, in the **Counters** sub-panel (see *figure 764, page 624*).

Counters				
When	Value	Action	Reset	Properties

Figure 764 The Counters are defined in a dedicated sub-panel

Five properties are defined per counter:

- **When:** defines the equality or inequality for the counter.
- **Value:** defines the counter threshold value.

- **Action:** defines the action taken when the counter is triggered. Three actions are available:
 - **And:** no action is taken when the counter is triggered. Instead, this counter is joined to one or more counters conditions. When all the counters joined with the **And** action are triggered, the collective action is executed.
 - **Pulse:** a user-defined TTL pulse is generated at the DIO connector.
 - **Autolab control:** an instance of the Autolab control command is executed.
 - **Shutter control:** defines the state of the shutter of a connected Autolab or Avantes light source with TTL control.
 - **Get spectrum:** triggers the acquisition of a spectrum on a connected Autolab or Avantes spectrophotometer.
- **Reset:** specifies if the counter should be reset when it is triggered, using the provided toggle .
- **Properties:** defines the properties of the Action defined in the Action column.

9.4.1 Counter configuration

The following steps describe how to configure a counter.

1 Add a counter to the list

Click on the **+** button to add a counter to the table.



NOTICE

A counter is automatically generated.

2 Specify the counter (in)equality

Click on the cell of the **When** column and select the equality or inequality to use in the counter using the provided drop-down list (< , = or >).




3 Specify the value

Specify the threshold value for the counter in the corresponding cell of the **Value** column.




4 Specify the action

Click on the cell of the **Action** column and select the action to use in the counter using the provided drop-down list.




5 Set the reset property

Use the provided toggle to define if the counter should be reset after it is triggered.



6 Define the properties of the specified action

Use the provided properties frame to define the properties of the Action defined for the counter.



NOTICE

To remove a counter from the table, select the row of the counter and click the  button above the **Counters** table.

9.4.2 Counter action - Pulse

The **Pulse** action can be used to send a TTL (Transistor-Transistor Logic) pulse to an external device when the condition defined for the counter is met. This pulse can be used to trigger the external device to perform a specific action.

The properties of the Pulse are defined in the dedicated frame, on the right-hand side of the **Counters** table (see figure 765, page 628).



Properties		
DIO connector	P1	▼
DIO port	Port A	▼
Pulse value	0	⁰¹ ₁₀
End value	00000000	¹ ₂₃
Duration	4000	µs

Figure 765 The Pulse properties are defined in the frame on the right-hand side of the Counters table

The following properties are available:

- **DIO connector (P1 or P2):** defines the DIO connector used to send the pulse.
- **Port (A, B or C):** defines the DIO port used to send the pulse.
- **Pulse value:** the decimal or binary expression of the 8 bit pulse state of the specified DIO port.
- **End value:** the decimal or binary expression of the 8 bit end state of the specified DIO port.
- **Duration (µs):** the duration of the pulse, in µs.



NOTICE

It is possible to switch from binary expression to decimal expression and from decimal expression to binary expression by clicking the ¹₂₃ and ⁰¹₁₀ buttons located next to the Pulse value and End value fields, respectively.



NOTICE

More information on the DIO ports and connectors can be found in *Chapter 16.3.1.3*.

For example, using the settings specified in *Figure 766*, the following pulse will be generated from DIO connector P1, port B:

1. **From initial state to Pulse value:** the pulse will start from the initial state of the DIO port. It will then go to the **Pulse** value defined in the **Properties** frame. In this example, the **Pulse** value is *10000000* (or 128 in decimal).



2. **From Pulse value to End value:** after 10 ms, the DIO port will transition from the **Pulse** value to the **End** value. In this example, the **End** value is *00000000*.

Properties	
DIO connector	P1 ▼
DIO port	Port B ▼
Pulse value	10000000 ¹ ₂₃
End value	00000000 ¹ ₂₃
Duration	10000 μs

Figure 766 Example of a Pulse

9.4.3 Counter action - Autolab control

The **Autolab control** action can be used to set the properties of the instrument using an instance of the **Autolab control** command.

The properties of the Pulse are defined in the dedicated frame, on the right-hand side of the **Counters** table (see figure 767, page 629).

Properties
<div style="border: 1px solid gray; padding: 5px; display: inline-block;"> Edit Autolab control </div> 

Figure 767 The Autolab control properties are defined in the frame on the right-hand side of the Counters table

Clicking the Edit Autolab control button opens the Autolab control editor (see figure 768, page 630).




Figure 768 The Autolab control editor



NOTICE

For more information on the **Autolab control** command, please refer to *Chapter 7.2.1*.

9.4.4 Counter action - Shutter control

The **Shutter control** action can be used to open or close the shutter of a connected Autolab or Avantes light source by setting the required DIO value on the specified connector.

The properties of the **Shutter control** action are defined in the dedicated frame, on the right-hand side of the **Counters** table (see *figure 769, page 630*).




Figure 769 The properties of the Shutter control action

The following properties are available:

- **DIO connector (P1 or P2):** defines the DIO connector used to control the light source shutter.
- **Shutter open:** defined the state of the shutter, using the provided  toggle. When the shutter is off, no light comes out of the light source. When the shutter is on, light can come out of the light source.



NOTICE

The light source shutter will remain in the specified state until changed.




CAUTION

The **Shutter control** action only works with light sources that support TTL control that are used with this mode enabled.

9.4.5 Counter action - Get spectrum

The **Get spectrum** action can be used to synchronize the acquisition of a spectrum on a connected Autolab or Avantes spectrophotometer by sending a TTL pulse of required length.

The properties of the **Get spectrum** action are defined in the dedicated frame, on the right-hand side of the **Counters** table (see *figure 770, page 631*).



The screenshot shows a window titled "Properties". Inside, there is a label "DIO connector" followed by a dropdown menu. The dropdown menu is currently displaying "P1" and has a small downward-pointing arrow on its right side.

Figure 770 The property of the Get spectrum action

The following property are available:

- **DIO connector (P1 or P2):** defines the DIO connector used to send the trigger.



NOTICE

The third and fourth counter shown in *Figure 771* are connected by a grey line on the left-hand side of the table, indicating that both counters have a **AND** relationship.

9.5 Plots

The plots define which how measured or calculated signals are plotted during a measurement. Two plots groups are shown in this section (see *figure 772, page 633*):

- **Default plots:** a table containing a list of preconfigured plots.
- **Custom plots:** a table that can be used to define custom plots.

CV staircase

Sampler
Options
Plots
Advanced

Default plots

	Enabled	Plot number	Options
i vs E	<input checked="" type="checkbox"/>	<input type="text"/>	Edit
i vs t	<input type="checkbox"/>	<input type="text"/>	
Log(i) vs E	<input type="checkbox"/>	<input type="text"/>	
Log(i) vs Log(t)	<input type="checkbox"/>	<input type="text"/>	
E vs i	<input type="checkbox"/>	<input type="text"/>	
E vs t	<input type="checkbox"/>	<input type="text"/>	
E vs Log(i)	<input type="checkbox"/>	<input type="text"/>	
E vs Log(t)	<input type="checkbox"/>	<input type="text"/>	

Custom plots +

Text	X	Y	Z	Enabled	Plot number	Options

Figure 772 The Plots define how the measured data is displayed during a measurement



NOTICE

The plots listed in the **Default** plots table depend on the measurement command and on the signal defined in the **Sampler**.

Custom plots								+
Text	X	Y	Z	Enabled	Plot number	Options		

Figure 774 The Custom plots table

1 Add a plot to the list

Click on the **+** button to add a plot to the table.

Custom plots								+
Text	X	Y	Z	Enabled	Plot number	Options		
								+

Adds a new plot.

2 Specify the name of the plot

Specify the name of the custom plot by typing the name in the first cell of the **Text** column.

Custom plots								-	+
Text	X	Y	Z	Enabled	Plot number	Options			
Example								Edit	

3 Specify the signal for the X axis


Click the first available cell in the **X** column and select the signal to plot on the X axis using the provided drop-down list.

Custom plots								-	+
Text	X	Y	Z	Enabled	Plot number	Options			
Example	Time							Edit	

- Potential applied
- Time
- WE(1).Current
- WE(1).Potential
- Scan
- Index


4 Specify the signal for the Y axis

Click the first available cell in the **Y** column and select the signal to plot on the Y axis using the provided drop-down list.



5 Specify the signal for the Z axis

Click the first available cell in the **Z** column and select the signal to plot on the Z axis using the provided drop-down list.



NOTICE

To remove a plot from the table, select the row of the plot and click the  button above the **Custom plots** table.

9.5.3 Plot options

To edit the plot options, the button  located next to each enabled plot or each custom plot is provided (see figure 775, page 637).




Figure 775 Editing the plot options

Clicking this button displays the plot options editor screen. The controls on this screen can be used to define the plot settings of the corresponding plot (see figure 776, page 638).

9.5.3.1 Data option

The **Data** sub-panel can be used to defined general properties of the data (see figure 777, page 639).




Figure 777 The Data sub-panel

The following properties can be edited in the **Data** sub-panel:

- **Point style:** defines the point style, color and size used by the plot, using dedicated drop-down lists.
- **Line style:** defines the line style, color and size used by the plot, using dedicated drop-down lists.
- **Y-axis placement:** specify the location of the Y axis using the provided drop-down list. The choice is provided between left and right.

9.5.3.2 Axes option

The **Axes** sub-panel can be used to defined general properties of the plot axes (see figure 778, page 639).




Figure 778 The Axes sub-panel

The following properties can be edited in the **Axes** sub-panel for each individual axis:

- **Label:** defines the label of the axis. When this field is left empty, the name of the signal plotted on this axis will be used instead.
- **Scale type:** defines the scale type of the axis, using the provided drop-down list. The choice is provide between linear and logarithmic.

- **Title font:** the font used for the title. The font type and size can be specified using dedicated drop-down lists. The format of the title can be edited by toggling the bold formatting or italic formatting on or off using the dedicated buttons.
- **Show grid:** enables or disables the chart grid, using the provided  toggle.

9.6 Automatic integration time

The **Automatic integration time** option specifies which of the available integration time constants can be used by the measurement command.



NOTICE

This option is only available for instruments that are fitted with the optional **FI20** module or the **on-board integrator**. For more information, please refer to *Chapter 16.3.2.11*.

When this option is used, the instrument will automatically select the most suitable integration time constant.

Automatic integration time settings are defined in the dedicated table, in the **Automatic integration time** sub-panel (see figure 780, page 641).





Figure 780 The Automatic integration time option is defined in a dedicated sub-panel

Three properties can be specified for each working electrode:

- **Enabled:** a  toggle used to set the Automatic integration time on or off.
- **Highest time constant:** defines the highest possible time constant, using the provided drop-down list.
- **Lowest time constant:** defines the lowest possible time constant, using the provided drop-down list.

of the interval time. The value of the applied delay, in s, is updated in the field below the input field for the *Alpha value* (see figure 784, page 643).



The screenshot shows a software interface for a 'CV staircase' measurement. On the left, there is a vertical menu with options: 'Sampler', 'Options', 'Plots', and 'Advanced' (which is highlighted with a green arrow). The main area contains two input fields: 'Alpha value' with a value of '0.5' and 'Delay' with a value of '0.010818' followed by a unit 's'. The 'Delay' value is shown to be automatically updated based on the 'Alpha value'.

Figure 784 The Delay value is automatically updated when the Alpha value is modified



NOTICE

The actual delay depends on the interval time.



NOTICE

For more information on the *Alpha value* property, please refer to M. Saralathan, R.A. Osteryoung, J. Electroanal. Chem. 222, 69 (1987).

10.1 Creating a new procedure

To create a new procedure, click the **New procedure** button in the **Actions** panel on the dashboard (see figure 785, page 645).




Figure 785 Click the New procedure button to create a new procedure

A new tab will be created, providing an empty procedure editor that can be used to create a customized procedure (see figure 786, page 645).




Figure 786 A new tab is created

The new procedure editor displays the three main panels:

The screenshot shows a 'Properties' dialog box for a procedure named 'Chrono amperometry (Δt > 1 ms)'. The dialog has a title bar with a maximize button. Below the title, the procedure name is displayed. There are three input fields: 'Procedure name' containing 'Chrono amperometry (...)', 'Remarks' containing 'Chrono amperometry (Δt > 1 ms)', and 'Estimated duration' containing '15' with a unit 's' to the right. Below these fields are two buttons: 'Edit End status' and 'More'. A 'Tags' section follows, featuring a 'Rating' of five stars (all filled) and a 'Tags' button with an 'Add' label.

Figure 788 Procedure properties of the default Chrono amperometry procedure

Additional buttons are available to edit the of the procedure options or the End status of the instrument. Please refer to *Chapter 10.2* and *Chapter 10.3* for more information.

The controls provided in the **Tags** sub-panel can be used to assign a rating and tags to the procedure. This provides further options for bookkeeping purposes. More information on the use of the rating and tags can be found in *Chapter 6.8*.

10.2 Global options and global sampler

For all procedures, it is possible to define global options and global sampler settings. If defined, these settings will be used for all the commands unless overruled for a specific command in the procedure. To define these settings, click the **More** button in the **Properties** panel (see figure 789, page 648).




Figure 789 Click the More button to open define the global options and global sampler

A new screen will be displayed, as shown in Figure 790, showing two different sections:

- **Sampler:** the settings in this section define the global sampler settings (see figure 790, page 648).
- **Options:** the settings in this section define the global options settings (see figure 791, page 649).




Figure 790 The global sampler settings




Figure 791 The global options settings

The settings in **Sampler** section can be used to specify which signals have to be sampled during the procedure. More information about the **Sampler** can be found in *Chapter 9.1*.

If commands are already located in the procedure, these commands will be displayed in the Apply global sampler to subsection. Using the provided checkboxes, it is possible to define on which the global sampler settings need to be applied (see *figure 792, page 650*).




Figure 792 It is possible to define on which commands the global sampler needs to be applied

The same applies to the global options. Using the provided checkboxes, it is possible to define on which the global options settings need to be applied (see figure 793, page 650).




Figure 793 It is possible to define on which commands the global options needs to be applied

10.3 End status Autolab

The End status Autolab is an option that can be defined for all procedures to define the how the Autolab system should behave when the procedure stops. The settings defined in the End status Autolab will be activated when the following events occur:

- The procedure stops normally
- The procedure is stopped by the user
- The procedure is stopped by a cutoff

To edit the End status Autolab, click the  button (see figure 794, page 651).




Figure 794 Click the Edit End status button to open the End status Autolab editor

The End status Autolab screen will be displayed. In this screen all of the instrumental settings and module settings can be specified (see figure 795, page 652).




Figure 795 The End status Autolab editor



NOTICE

The End status Autolab editor is the same as the one for the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

10.4 Procedure tracks

Procedures in NOVA consist of at least one main **Track** of commands, which are executed in sequence. Each command can be used to create a sub-track in which additional commands can be located. Commands located in each sub-track are executed sequentially when the parent command located in the main track is executed.

A simple example is provided by the default *Hydrodynamic linear sweep* procedure, available from the **Default** procedures, in the **Library** (please refer to *Chapter 8.2.4* for a complete description of this procedure). This procedure contains one main track in which a **Repeat** command is located. This command is configured in the *Repeat for multiple values* mode. The rotation rates required for this measurement are pre-defined in the **Repeat** command (see figure 796, page 653).




Figure 796 The Hydrodynamic linear sweep voltammetry procedure

Additional commands are located in the sub-track of the **Repeat** command. These commands are visible when the **Repeat** command is selected (see figure 796, page 653).

When the procedure is executed, the four first commands located in the main track are executed in sequence. When the **Repeat** command is executed, the four commands located in the sub-track are executed in sequence and repeated six times as defined by the **Repeat** command. When all six repetitions are completed, the procedure resumes the main track of the procedure.

10.5 Procedure wrapping

The procedure editor frame has a limited width. When a procedure track has more commands than can be displayed in a single line, the software will wrap the track and display the commands on multiple lines. In the example below, the procedure has a single track, wrapped on two lines. The last **Cell** command is located on the second line (see figure 797, page 653).




Figure 797 Long tracks are wrapped on several lines if needed

If the NOVA window is resized, or if the **Properties** panel on the right-hand side or the **Command** panel on the left hand side are resized or collapsed, the procedure editor will be readjusted, and if possible, the commands will be displayed on a single line (see figure 798, page 654).




Figure 798 Commands are relocated to a single line when enough room is available

10.6 Procedure zooming

The procedure editor frame has a limited width. If needed, the size of the items in the procedure editor frame can be adjusted with the controls located in the top right corner of the frame (see figure 799, page 654).




Figure 799 Zoom controls are provided in the procedure editor

Using this function will either scale the size of the items and the text up or down (between 200 % and 50 % of the original size), as shown in Figure 800.




Figure 800 Zooming the procedure editor out

The following zooming controls are available:

- **Zoom out:** decreases the scaling of the items and text shown on screen. The  button or **[CTRL] + [-]** keyboard shortcut can be used to do this.
- **Zoom to 100%:** resets the scaling of the items and text shown on screen to the default size. The  button or **[F4]** keyboard shortcut can be used to do this.
- **Zoom in:** increases the scaling of the items and text shown on screen. The  button or **[CTRL] + [=]** keyboard shortcut can be used to do this.

10.7 Command groups

Command groups can be created in a procedure in order to structure a complex procedure. Grouping commands allows to group commands in a container command in the main track (or a sub-track). This is useful for hiding commands that are not directly important to the measurement, or it can simply be used to group commands that have similar functionality. A group also provide the possibility to relocate several commands in the procedure at once.

10.7.1 Grouping commands

It is possible to group commands located in a procedure. Creating groups provides the benefit of locating command into groups which can then be moved in the procedure editor. Grouping commands can also be used to create a clear procedure structure, which is especially useful for complex procedures. Grouped commands are replaced by a **Group** item in the procedure editor in which the grouped commands are relocated in the original order.

To group commands, select the commands in the procedure editor and click the  button located in the top right corner of the procedure editor. It is also possible to use the keyboard shortcut consisting of the **[CTRL] + [G]** combination (see figure 801, page 655).




Figure 801 Grouping commands

The grouped commands will be replaced in the procedure track by the Group item. Clicking this item in the procedure editor reveals the grouped command in the track below (see figure 802, page 656).




Figure 802 The grouped commands are visible when the group is selected

10.7.2 Ungrouping commands

It is possible to ungroup commands that are located in a group. Ungrouping commands removes the group from the procedure editor and promotes the involved commands to the next available procedure track above the group.

To ungroup grouped commands, select the group in the procedure editor and click the  button located in the top right corner of the procedure editor. It is also possible to use the keyboard shortcut consisting of the **[SHIFT] + [CTRL] + [G]** combination (see figure 803, page 656).




Figure 803 Ungrouping grouped commands

The group will be removed from the procedure editor and the grouped commands will be restored in the track above the former group in the order in which these commands were located in the group (see figure 804, page 657).




Figure 804 The ungrouped commands are restored in the procedure editor

10.7.3 Renaming groups

For bookkeeping purposes, it is possible to change the name of a group of commands. This is useful when creating complex procedures as each group can be given a relevant name. To change the name of a group in the procedure, select the group and change the name in the **Properties** panel located on the right-hand side of the screen (see figure 805, page 657).




Figure 805 Changing the name of a group

After validation of the new name, the procedure editor will be updated, displaying the name of the group (see figure 806, page 657).




Figure 806 The group is renamed



NOTICE

It is not possible to change the status of a disabled command while the procedure is running.

If the disabled command is a **Group** command or if the disabled command has additional commands stacked below it, all the commands located in this group or stack are disabled (see figure 809, page 659).




Figure 809 Disabling a Group command in the procedure editor

10.8.2 Enabling commands

Disabled commands can be enabled again in the procedure editor. Enabling a disabled command restores this command to its previous state in the procedure.

To enable one or more *disabled* commands in the procedure editor, select the command or commands and click the ✓ button located in the top-right corner of the procedure editor (see figure 810, page 659).




Figure 810 Enabling disabled commands

When a disabled **Group** command is enabled, all the commands contained in the group are enabled. This also applies to command located in a disabled stack.





Figure 811 The original procedure

The **Wait** command will be added between the **Cell** command and the **CV staircase** command.


1 Select the command to add

Click the command to add to the procedure.




2 Drag the command in the procedure editor

While holding the mouse button, drag the command to the procedure editor frame.




3 Place the command in the procedure

Place the new command in the procedure.



4 Finalize the insertion of the command

Release the mouse button to validate the position of the new command in the procedure.



10.9.1.1.2 Using the drag and drop method to add commands to a command group or a sub-track

The following steps illustrate how to use the drag and drop method for adding commands to a command group or to a sub-track. This method is used to add a **Message** command to the following procedure (see figure 812, page 663).





Figure 812 The original procedure

The **Message** command will be added at the beginning of the Group track, before the **Apply** command.


1 Select the Group command

Click the **Group** command in the procedure editor to display the commands located in the group below the main procedure track.




2 Select the command to add

Click the command to add to the procedure.




3 Drag the command in the procedure editor

While holding the mouse button, drag the command to the procedure editor frame.




4 Place the command in the procedure

Place the new command in the procedure.



5 Finalize the insertion of the command

Release the mouse button to validate the position of the new command in the procedure.



10.9.1.2 Adding commands using the double click method

The *double click* method provides the means to quickly add commands to a procedure. Double clicking a command in the commands browser adds the command to the end of the active track in the procedure editor.

For example, double clicking the **Play sound** command adds this command after the **Cell** command in the *Cyclic voltammetry potentiostatic* procedure (see figure 813, page 666).




Figure 813 Adding a command to a procedure using the double click method

The command is added to the procedure editor (see figure 814, page 667).




Figure 814 The command is added to the procedure

If the procedure contains more than one track, the double-clicked command is added at the very end of the active track. In the example shown below, when the contents of the **Group** are visible, double-clicking the **Play sound** command adds the command at the end of the sequence in the **Group** (see figure 815, page 667).




Figure 815 The command is added to the active track

If the contents of the **Group** command are not displayed in the procedure editor, double-clicking the **Play sound** command leads to the same result as in Figure 814.

10.9.2 Removing commands

It is possible to remove commands from a NOVA procedure, by selecting the command or commands and clicking the  button or pressing the **[Delete]** key (see figure 816, page 668).



Figure 816 Select the command or commands to delete

It is also possible to use the *Delete* option available from the **Edit** menu. Deleted commands are removed from the procedure and all commands located after the deleted commands are shifted leftwards to replace the deleted commands (see figure 817, page 668).




Figure 817 The commands located on the right of the deleted command are shifted leftwards



NOTICE

Removing a **Group** command or a command stack removes all the commands in the group or in the stack from the procedure.

10.10 Moving commands

Commands located in a procedure can be moved and relocated anywhere in the procedure using the drag and drop method. Whenever a command is moved, the other commands located in the procedure are shifted leftwards or rightwards in order to create room for the moved command.



NOTICE

It is only possible to move one command at a time.

10.10.1 Moving commands using the drag and drop method

The following steps illustrate how to use the drag and drop method for moving commands in the procedure editor. This method is used to move a **Message** command, located in a **Group** command, to the beginning of the following procedure (see figure 818, page 669).





Figure 818 The original procedure


1 Select the command to move

Click the command to move in the procedure editor.




2 Drag the command to move

While holding the mouse button, move the selected command to a new location in the procedure editor. The other commands located in the procedure editor will be shifted in order to make room for the moved command.



3 Release the mouse button

Release the mouse button to validate the new position of the command in the procedure.



10.10.2 Using the drag and drop method to move commands to a command group or a sub-track

The following steps illustrate how to use the drag and drop method for moving commands to a command group or to a sub-track . This method is used to move a **Message** command, located at the beginning of the procedure, to the **Group** command (see figure 819, page 670).





Figure 819 The original procedure


1 Select the destination Group command

Click the **Group** command in the procedure editor to display the commands located in the group below the main procedure track.




2 Select the command to move

Click the command to move in the procedure editor.



3 Drag the command to move

While holding the mouse button, move the selected command to a new location in the procedure editor. The other commands located in the procedure editor will be shifted in order to make room for the moved command.






Figure 820 Selecting multiple commands in the procedure editor

Once the commands are selected, it is possible to use the editing tools like cut or copy and paste or drag and drop to edit the procedure. When the one of the commands in the selection is dragged through the procedure editor, an indicator will be shown in the top left corner of the command, indicating the number of additional commands moved at the same time (see figure 821, page 673).




Figure 821 The number of commands dragged at the same time is indicated in the procedure editor

In the case of Figure 821, one additional command is included in the selection. A **+1** indicator is therefore shown in the top left corner of the **CV staircase** command.

The selected commands can be repositioned in the procedure editor using the drag and drop method (see figure 822, page 673).



Figure 822 Relocating commands using the drag and drop method

The releasing the mouse confirms the new position of the commands in the procedure editor (see figure 823, page 674).

10.12.1 Creating command stacks

The following steps illustrate how to create a command stack using the drag and drop method. This method is used to add a **Calculate signal** command to the following procedure (see figure 825, page 675).





Figure 825 The original procedure


1 Select the command to add

Click the command to add to the procedure.




2 Drag the command in the procedure editor

While holding the mouse button, drag the command to the procedure editor frame.



3 Drag the command onto the parent command

Drag the command onto the parent command.




NOTICE

The command shrinks when located onto the parent command.

4 Finalize the command stack

Release the mouse button, validating the location of the command and the command stack.



10.12.2 Remove commands from stacks

To remove a command from a stack, simply select the command in the stack and click the  button or press the **[Delete]** key. The command will be removed from the stack (see figure 826, page 677).




Figure 826 Select the command to remove from the stack




Figure 829 The Repeat command has a list of predefined rotation rates

The rotation rate of a rotating disk electrode (RDE) or rotating ring-disk electrode (RRDE) can therefore be automatically adjusted during a measurement. In this procedure, the **Repeat** command is linked to a R(R)DE command in the repeat loop.

This chapter explains how the linking mechanism works in NOVA and how to perform the following link-related actions:

- View links
- Create links
- Edit links

10.13.1 Viewing links

It is possible to show the links of any linked command in NOVA by selecting such a command anywhere in the procedure and clicking the  button located in the top right corner of the procedure editor (see figure 830, page 680) or by using the keyboard shortcut **[CTRL] + [L]**.




Figure 830 Viewing procedure links

The links are displayed in the dedicated **Edit link** screen. All the links involving the properties of the selected command are shown (see figure 831, page 680).




Figure 831 The Edit links screen

The following information is represented in the **Edit links** screen:

- **Commands:** all commands linked to the selected command are represented in the **Edit links** screen. Commands located before the selected command are represented on the left of the screen and command located after the selected command are represented on the right of the screen.
- **Linkable properties:** all linkable properties of the commands represented in the **Edit links** screen are represented. These properties are identified with a name and one or more anchoring points.
- **Anchoring points:** one or more anchoring points, identified by a symbol, are represented for each linkable property. Anchoring points locating on the left of a linkable property are *output* points. Anchoring points located on the right of a linkable property are *input* points.
- **Link line:** one or more grey lines connecting two or more anchoring points.

The grey link line connects the anchoring points located on the right of the *Rotation rate (RPM)* property of the **Repeat** command with the anchoring point located on the left of the *Rotation rate* property of the **R(R)DE** command.

This link indicates that all the values defined in the **Repeat** command will be used to change the *Rotation rate* property of the **R(R)DE** command during the measurement. The *Rotation rate (RPM)* property of the **Repeat** command is used as an output and the *Rotation rate* property of the **R(R)DE** command is used as an input, indicating the direction of the link.

10.13.2 Creating links

It is possible to create links in any NOVA procedure containing linkable commands. To illustrate this option, the following procedure template will be used (see figure 832, page 681).




Figure 832 The procedure used to illustrate the creation of links in NOVA



NOTICE

The procedure used in this example is created by adding a **Message** command, used as an Input, at the beginning of the default Cyclic voltammetry potentiostatic procedure.

Links can be created between two or more commands in the **Edit links** screen. To create links between two or more command, it is necessary to select these command in the procedure editor and click the button, located in the top right corner for the procedure editor or use the **[CTRL] + [L]** keyboard shortcut.



NOTICE


The button is only visible when two or more commands are selected in the procedure editor.

Two procedures are available for creating links:

1. Creating links between two commands
2. Creating links between more than two commands

3 Set the input anchoring point

While still holding the mouse button, move the line on top of the *Potential* input anchoring point of the **Apply** command and release the mouse button. The link line will be drawn between the two properties.



The link is now created.

4 Close the Edit link screen

Click the  button located in the top left corner to close the **Edit links** screen.

The properties are now linked. The created link will force both properties to be the same at any point during the measurement or whenever either one is modified by the user.



NOTICE


The **Apply 0 V** command is dynamically changed to **Apply 0,1 V** after the link is created.

10.13.2.2 Creating a link between more than two commands



NOTICE


The steps detailed in this section apply to an example procedure. These steps can be repeated for any procedure containing two or more linkable commands.



The first link is now created.

4 Set the second output anchoring point

Click the output anchoring point of the *Potential* property of the **Apply** command and, while holding the mouse button, drag a line towards the input anchoring point of the *Start potential* property of the **CV staircase** command.



5 Set the second input anchoring point

While still holding the mouse button, move the line on top of the *Start potential* input anchoring point of the **CV staircase** command and release the mouse button. The link line will be drawn between the two properties.




Figure 833 Selecting the Message command first

Both the **Apply** and **CV staircase** commands are shown at the right-hand side of the **Message** command because both commands are located after the **Message** command in the procedure. The **Apply** and **CV staircase** commands are shown above one another, with only one command in focus and the other out of focus (greyed out). The links can be edited between the **Message** command and the command in focus (**CV staircase** in Figure 833).

If the **CV staircase** command is selected first, the **Message** and **Apply** commands will be displayed in a different way in the **Edit links** screen (see figure 834, page 688).




Figure 834 Selecting the CV staircase command first

Both the **Message** and **Apply** commands are shown at the left-hand side of the **CV staircase** command because both commands are located before the **CV staircase** command in the procedure. The **Message** and **Apply** commands are shown above one another, with only one command in focus and the other out of focus (greyed out). The links can be edited between the **CV staircase** command and the command in focus (**Apply** in Figure 833).

In both cases, it is possible to click the greyed out command to switch the focus in the **Edit links** screen and view or edit the links of the other command (see figure 835, page 689).




Figure 835 Switching the command focus in the Edit links screen



NOTICE

It is also possible to use the mouse wheel to quickly scroll through the out of focus commands in the **Edit links** screen.

10.13.3 Editing links

It is possible to edit or remove links in a procedure at any time. To edit or remove links, it is necessary to open the Edit links screen by selecting one or more linked command and clicking the  button or using the **[CTRL] + [L]** keyboard shortcut (see figure 836, page 689).




Figure 836 Opening the Edit links screen

This will open the Edit links screen (see figure 837, page 690).




Figure 837 The Edit links screen

In the Edit links screen, it is possible to reroute an existing link by clicking one of the ends of the link and moving this to another anchoring point. In Figure 838 the link is moved from the Start potential anchoring point to the Stop potential anchoring point.




Figure 838 Changing the link to the Start potential property to the Stop potential property

It is also possible to delete a link, by clicking one end of a link and pulling it away from the anchoring point, as shown in Figure 839.




Figure 839 Removing a link

If the mouse button is released, the link will be removed (see figure 840, page 691).




Figure 840 The link is removed

10.14.1 Saving a My command

To illustrate how to create a My command, the following starting procedure will be used (see figure 842, page 693). The modified **CV staircase** command will be saved in this example.





Figure 842 The initial procedure used to create a My command

1 Modify the source command

Modify the source command that will be saved as a My command.


2 Save the command

Select the command to save and click the  button in the top-right corner of the **Procedure editor** panel.



3 Specify name and remarks

An input dialog will be displayed.



Name CV staircase

Remarks

4 Validate the name and remarks


Click the button to validate the name and remarks and save the command in My commands.

Command name

Remarks

5 The command is added to the group

The saved command is added to the My commands group.



10.14.2 Editing My commands

Commands that have been saved as My commands can be edited, exported or removed. All editing actions can be accessed by right-clicking the My command in the **Commands** panel and by selecting the required action from the context menu (*see figure 843, page 696*).




Figure 843 My commands can be edited, exported or deleted

Selecting the **Edit** option displays the name and remarks editor (see figure 844, page 696).




Figure 844 It is possible to adjust command name and remarks

Selecting the **Export** option displays a Windows Explorer dialog which can be used to specify a location and a file name for the My command. This command will be exported to the specified location with the specified name (see figure 845, page 697).




Figure 845 Specifying the name of the command file



NOTICE

The extension used by My commands is `.noi`.

Selecting the **Delete** option removes the My command from the computer. A validation message will be displayed (see figure 846, page 697).

Remove command

Are you sure you want to remove "My CV staircase" from my commands?



Figure 846 A confirmation message is shown

Click the `Yes` button to delete the command.




Figure 848 Starting a procedure from the Measurement menu

These three options will start the procedure using the **Default** instrument.



NOTICE

The serial number of the **Default** instrument is shown in the top left corner of the procedure editor.



NOTICE

It is possible to change the **Default** instrument *Change the default instrument (see chapter 5.1, page 95)*.

It is also possible to start a measurement on any available instrument by specifying on which instrument to run the procedure, using the Measurement menu (see figure 849, page 700).




Figure 849 Specifying the instrument on which to run the measurement

When a procedure is started, the following tasks are carried out:

1. The procedure is tested for warnings or errors *Procedure validation* (see chapter 11.2, page 701).
2. A new tab opens, with the same name of the source procedure. A clone of the procedure is created in the new tab. The new tab will be used to record and display the measured data while the source procedure remains unchanged *Procedure cloning* (see chapter 11.3, page 702).
3. A **Plots** frame will appear at the bottom of the screen. This frame will display all the measured data according to the properties defined in the procedure *Plots frame* (see chapter 11.4, page 704).

During a measurement, it is also possible to carry out a number of actions:

1. It is possible to modify some of the measurement properties *Real-time properties modification* (see chapter 11.5.1, page 708).
2. It is possible to hold or stop the procedure and it is possible to skip the command being executed *Procedure control* (see chapter 11.5.2, page 711).
3. It is possible to reserve the scan direction, if applicable *Reverse scan direction* (see chapter 11.5.3, page 712).
4. It is possible to display the instrument **Manual control** panel *Display the Manual control panel* (see chapter 11.5.4, page 713).
5. It is possible to enable or disable plots *Enable and disable plots* (see chapter 11.5.5, page 715).

At the end of the measurement, the following tasks are carried out:

1. At the end of the measurement, the information displayed in the new tab will be time stamped for bookkeeping purposes.

- Post validation is carried out at the very end and information or warning messages are shown, if applicable, indicating possible improvements of the procedure.

11.2 Procedure validation

Whenever a procedure is started, NOVA will verify the properties defined for each command in the procedure and test if these are compatible with the instrument the procedure is started on.

If a **Warning** is detected, the procedure will not start immediately and instead a message will be displayed to the user, providing information about the encountered **Warning** (see figure 850, page 701).




Figure 850 A Warning is detected

It is possible to click the **OK** button and ignore the **Warning** or click the **Cancel** button and return to the procedure editor to adjust the procedure.



NOTICE

Ignoring a **Warning** is **not** recommended!

If an **Error** is detected, then the procedure will not be allowed to continue and a message will be displayed providing information on the **Error** (see figure 851, page 702).




Figure 851 An Error is detected

It is then only possible to click the to close the message and return to the procedure editor.



NOTICE

If **Errors** and **Warnings** are detected, the **Errors** are listed before the **Warnings** in the validation message, as shown in Figure 851.

If no **Warnings** or **Errors** are detected, the procedure is started and the measurement begins.

11.3 Procedure cloning

After validation, the procedure starts. A **clone** of the source procedure is created in a new **tab**. The new tab will have the same name as the tab contained the source procedure. Cloning the source procedure is convenient because it creates a new version of the original procedure that can be modified during the experiment. The source procedure remains unchanged in the original tab.

The procedure then starts in the new tab (see figure 852, page 703).




Figure 852 The procedure starts in a new tab



NOTICE

The serial number of the instrument on which the procedure is started is reported below the name of the procedure.

The running state is indicated by the spinning wheel symbol, , shown in the tab as well as for the running command (see figure 852, page 703).

The buttons located in the top right corner of the procedure editor of the running procedure can be used to either skip to the next command in the procedure () , pause the running command () or stop the whole procedure () . More information can be found in *Chapter 11.5.2*.



11.4 Plots frame

When a procedure starts, an additional **Plots** frame is opened at the bottom of the screen (see figure 853, page 704).




Figure 853 The Plots frame is created at the bottom of the screen

This frame is used for displaying plots during a measurement. All the plots defined in the procedure are created in the **Plots** frame. During the measurement, whenever data becomes available, the plots are populated with measured data points.



NOTICE

Plots for which no data is available are shown in the **Plots** frame slightly greyed out. Whenever data becomes available for a plot, the plot will be shown normally.

The **Plots** frame can be resized to increase or decrease the size of the plots shown in the frame. It is also possible to undock the **Plots** frame, by clicking the button in the top right corner of the frame (see figure 854, page 705).




Figure 854 Undocking the Plots frame

A new window will be created, displaying the contents of the **Plots** frame. Zooming in and out buttons are provided in the top right corner to increase or decrease the size of the plots in the window (see figure 855, page 706).




Figure 855 The undocked Plots frame



NOTICE

The undocked **Plots** frame can be closed at any time.

11.4.1 Displaying multiple plots

When a procedure generates multiple plots, these plots can be arranged in two different ways:

- **Sequence arrangement:** all the plots defined in the procedure are shown in sequence in the **Plots** frame scaled to the largest available space. If more plots are defined than can be arranged in the Plots frame, a scrollbar will be added to the frame. Using this scrollbar, it is possible to change the plots shown in the frame (see figure 856, page 707).




Figure 856 Plots shown in sequence arrangement

- Tiled arrangement:** all the plots defined in the procedure are shown in the **Plots** frame and are shrunk to size required to show each plot in the frame. No scrollbar is added to the Plots frame in this case (see figure 857, page 707).




Figure 857 Plots displayed in tiled arrangement




Figure 858 Selecting a command during a measurement shows the available properties



NOTICE

The greyed out properties cannot be modified in real-time.



NOTICE

Changing the properties of a command that has already been executed is possible but will not have any effect on the running procedure. This modification may however become active when the measured data set is converted to a new procedure, as explained in *Convert data to procedure* (see chapter 11.10, page 745).

It is possible to specify a new value for one or more of the available properties (see figure 859, page 710).




Figure 859 Modifying the number of scans

A new value will be validated by pressing the **[Enter]** key or by clicking away from property value being edited. The new value will be validated before becoming active. If the new value is not acceptable for the edited property, it will be displayed with a red frame around it, indicating that it is invalid (see figure 860, page 710).




Figure 860 New properties are validated before becoming active



If the new value is valid, it will be updated in the running procedure and used in the applicable command instead of the original value (*see figure 861, page 711*).




Figure 861 The new property is used during the measurement



NOTICE

Modifying procedure properties in real-time does not affect the source procedure from which the procedure was started.



NOTICE

All real time modifications of measurement properties are logged into the data grid and stored in the data file.

11.5.2 Procedure control

The buttons located in the top right corner of the procedure editor of the running procedure can be used to either skip to the next command in the procedure (▶), pause the running command (⏸) or stop the whole procedure (■). The procedure editor will update the status of a command affected by these controls, if applicable (*see figure 862, page 712*).




Figure 862 Holding the CV staircase command



NOTICE

All interactions with the procedure controls buttons are logged into the data grid and stored in the data file.

11.5.3 Reverse scan direction



NOTICE

This option is only available for the **CV staircase** command and the **LSV staircase** command.

During a measurement, it may be possible to modify the scan direction by clicking the ↺ button in the top right corner of the procedure editor. This button is only shown while the command that supports this option is running (see figure 863, page 713).




Figure 863 Reversing the scan direction

If the scan direction is reversed, the command will continue running until the requirements for that command specified by the user are fulfilled.



NOTICE

All interactions with the reverse scan button are logged into the data grid and stored in the data file.

11.5.4 Display the Manual control panel

At any time during a measurement, it is possible to display the **Manual control** panel of the instrument involved in the measurement. This can be done by selecting the *Manual control* option from the **View** menu or by pressing the **[F10]** shortcut key (see figure 864, page 714).




Figure 864 Displaying the instrument manual control

The **Manual control** panel can be used to modify some of the hardware controls during a measurement (see figure 865, page 714).




Figure 865 The Manual control panel can be used to modify instrument settings during a measurement



11.5.5 Enable and disable plots

While a procedure is running, it is possible to click the **More** button in the command **Properties** panel or to double click a measurement command in the procedure to adjust the plot settings (see figure 866, page 715).




Figure 866 Double click a measurement command to adjust the plots shown in the Plots frame

A new screen will be shown, presenting controls that can be used to adjust the plots visibility (see figure 867, page 716).




Figure 867 Plots can be enabled or disabled at any time during a measurement



NOTICE

The screen shown in *Figure 867* is the same as the one shown in *Figure 772*, without the **Sampler** and the **Options**, which cannot be modified in real-time.

In this screen, it is possible to disable pre-defined plots or to enable new plots, if needed, using the provided toggles. It is possible to disable a pre-defined plot in the **Plots** frame directly by right-clicking a plot to disable and selecting the corresponding option from the context menu (*see figure 868, page 717*).




Figure 868 Quickly disabling a plot in the Plots frame

The plot will be removed from the **Plots** frame (see figure 869, page 717).




Figure 869 The selected plot is disabled



NOTICE

Disabled plots can be enabled again using the method described at the beginning of this Section (see figure 867, page 716).

11.5.6 Q+ and Q- determination



NOTICE

This option is only available for the **CV staircase** command.

During the execution of the **CV staircase** command, after each scan is completed, the anodic and cathodic charge (Q_+ and Q_-) is automatically determined from each cyclic voltammogram and reported in the **Properties** panel (see figure 870, page 718).




Figure 870 The values of Q_+ and Q_- are automatically added to the Properties panel

The Q_+ and Q_- values are determined at the end of each scan. These values are reported in C.





NOTICE

The values of Q+ and Q- are also saved alongside the other electrochemical signals sampled during the measurement.

11.6 End of measurement

When a measurement finishes, the measured data becomes available for evaluation and analysis. Depending on the settings defined in the **NOVA Options**, the data may or may not be saved automatically *Options* (see chapter 1.9, page 13).

NOVA will also carry out the following activities at the end of each measurement:

1. **Time stamping:** the measured data is time stamped using the time and date of the beginning of the measurement *Procedure time stamp* (see chapter 11.6.1, page 719).
2. **Post validation:** the measured data is evaluated and information or warnings are provided, if applicable *Post validation* (see chapter 11.6.2, page 720).

Data files are differentiated from procedure files in the **Procedure editor** panel by green shading of the command tiles.

11.6.1 Procedure time stamp

At the end of measurement, the procedure is issued a **time stamp**. The time stamp corresponds to the starting time of the measurement (see figure 871, page 719).




Figure 871 The procedure is time stamped at the end of each measurement




Figure 873 A post validation warning message

Post validation messages generally provide indications which can be used for finetune the measurement conditions.

11.7 Specify plot preview

Whenever a data set is saved in the **Library**, a plot preview is created. This plot preview can be displayed in a tooltip in the **Library** to provide a preview of the data as shown in *Chapter 6.9*.

By default, the first plot in the **Plot** frame is used as a plot preview, however it is possible to specify another plot as the preview plot at any time. To change the plot preview, right-click the plot to use and select the *Set as preview plot* from the context menu (see figure 874, page 722).




Figure 874 Specifying the plot preview



NOTICE

The new plot preview will be updated when the modifications to the data file are saved.

11.8 Detailed plot view

It is possible to double click a plot shown in the **Plots** frame to obtain a larger view of the plot, change some of the plot properties or toggle to a 3D view of the plot, if available. The detailed view of the plot replaces the procedure editor view (see figure 875, page 723).




Figure 875 Detailed view of a plot

The detailed plot view provides the following controls:

- **Plot panel:** a large panel showing the selected plot. A number of buttons are located in the top right corner of this frame to add a data analysis command, view the data marker or toggle the 3D view on or off.
- **Properties panel:** a panel that can be used to change the plot properties during the measurement. This panel can be collapsed if necessary, by clicking the  button.

Clicking the  button closes the detailed plot view and returns to the procedure editor.

11.8.1 Plot properties

The **Properties** panel, shown in the right hand side of the screen, can be used to modify the plot properties of the active plot at any time (see figure 876, page 724).

11.8.2 Toggle the 3D view

Clicking the  button in the top right corner of the **Plot** panel toggles the 3D view on or off (see figure 877, page 725).




Figure 877 Toggling the 3D view on or off

The 3D view shows the same data using one additional Z axis. The plot can be rotated using by clicking and dragging the mouse.



NOTICE

It is only possible to display the data in 3D when a signal that can be plotted in real-time has been assigned to the Z axis of the plot.

11.8.3 Toggle the step through data mode

Clicking the  button in the top right corner of the **Plot** panel toggles the *Step through data* mode on or off (see figure 878, page 726).




Figure 878 Toggling the step through data mode on or off

When the *Step through data* mode is on, an additional indicator is added to the plot, showing the X and Y coordinates of the point indicated by the arrow, in the case of a 2D plot, and the X, Y and Z coordinates of the point indicated by the arrow, in the case of a 3D plot.



NOTICE

The indicator is always shown for the first data point of the plot.

Using the mouse, it is possible to perform the following action (2D plot):

- Click anywhere in the plot area: the indicator is relocated to the closest data point of the plot.

Using the keyboard it is possible to perform the following actions (2D and 3D plot):

- [**←**]/[**→**]: the indicator can be moved by 1 point at a time.
- [**←**]/[**→**] and [**CTRL**]: the indicator can be moved by 10 points at a time.
- [**←**]/[**→**] and [**CTRL**] and [**SHIFT**]: the indicator can be moved by 100 points at a time.

11.8.4 Add an analysis command

Clicking the  button in the top right corner of the **Plot** panel displays a popout menu from which an analysis command can be selected (see figure 879, page 727).




Figure 879 Adding an analysis command

The selected analysis command will be added to the procedure and will be applied on the active plot.



NOTICE

The analysis commands displayed in the popout menu depend on the type of data shown in the active plot.



NOTICE

More information on data analysis is provided in *Chapter 12*.



11.8.5 Zooming options

The controls located above the plot frame provide the means to zoom in and out on the plot and provide the means to rescale the plot for optimal display (see figure 880, page 728).




Figure 880 Zooming options are located above the plot

The following zooming options are available:

- **Zoom out:** increases the scaling of the X and Y axis on 2D plots and X, Y and Z axis on 3D plots. The  button or **[CTRL] + [-]** keyboard shortcut can be used to do this.
- **Fit view:** adjusts the scaling of the X and Y axis on 2D plots and X, Y and Z axis on 3D plots. The  button or **[F4]** keyboard shortcut can be used to do this.
- **Zoom in:** decreases the scaling of the X and Y axis on 2D plots and X, Y and Z axis on 3D plots. The  button or **[CTRL] + [=]** keyboard shortcut can be used to do this.



NOTICE

It is also possible possible to manipulate the scaling of the plot by using the **View** menu and by using the mouse directly on the plot.

11.8.6 Print plot

NOVA support the printing of plots to a printer connected to the computer. It is possible to print the visible plot, by clicking the  button, located above the plot (see figure 881, page 729).




Figure 881 Printing the visible plot

A Print Settings/Preview window will be displayed (see figure 882, page 730).




Figure 882 The Print Settings/Preview window

The following settings can be edited:

- **Printer:** specifies the printer used to print the plot. The printer can be selected using the provided drop-down list and the settings of the printer can be adjusted using the dedicated **Setup...** button. The **Print** button can be used to print the plot on the selected printer using the specified settings.
- **Header:** specifies an optional header. The font can be specified using the dedicated **Font...** button.
- **Footer:** specifies an optional footer. The font can be specified using the dedicated **Font...** button.
- **Orientation:** specifies the orientation of the plot. Radio buttons provide the choice between Portrait and Landscape.
- **Margins:** specifies the margin settings (top, bottom, left and right).

- **Content:** specifies additional options for the printing output. The following additional controls are available:
 - **Bitmap/Vector:** specifies the rendering of the plot in the preview. Radio buttons provide the choice between Bitmap (pixel) output or Vector output.
 - **Keep aspect ratio:** a checkbox that can be used to specify if the aspect ration of the plot should be maintained or not.
 - **Use background fill:** a checkbox that can be used to specify if the background of the plot should be visible or not.
 - **Use graph fill:** a checkbox that can be used to specify if the plot background should be visible or not.
 - **Refresh:** a button that can be used to refresh the preview.



NOTICE

The Use graph fill checkbox has no effect in the current version of NOVA.

11.8.7 Export plot to image file

NOVA support the exporting of plots to an image file, which can be used in third party applications. Two types of image types can be used when exporting plots:

- **Pixel based output:** the data is exported to a pixel based file format, with or without compression (*.bmp, *.png, *.jpg, *.tiff, *.gif).
- **Vector based output:** the data is exported to a vector based file format (*.emf, *.svg, *.wmf).

It is possible to export the visible plot to an image file, by clicking the  button, located above the plot (*see figure 883, page 732*).




Figure 883 Exporting the plot to an image file

A popout menu will be displayed, as shown in Figure 884, providing the means of specifying the size of the image to export in pixels (in the case of a pixel based output file) or in arbitrary units (in the case of a vector based output file).




Figure 884 Specifying the size of the exported image

Clicking the **Export** button displays a Windows explorer dialog which can be used to specify the path, name and file type used to create the output image file (see figure 885, page 733).




Figure 885 Specifying the name, location and type of output file

11.8.8 Relocate plots

It is possible, when the measurement is finished, to change the location of the plots using the drag and drop method directly in the **Plots** frame (see figure 886, page 734).




Figure 886 Plot positions can be adjusted after a measurement is finished

Click and drag a plot in the **Plots** frame to adjust its position. A grey line will be shown, indicating the new position of the dragged plot (see figure 887, page 734).




Figure 887 A grey line shows the new position of the plot

Releasing the mouse button confirms the new position of the plot in the **Plots** frame (see figure 888, page 735).




Figure 888 The plots are rearranged when the mouse button is released

If the selected plot is dragged over another plot in the frame and the mouse button is released, the selected plot will be added to the existing plot as an overlay (see figure 889, page 736).




Figure 889 Dragging a plot onto another plot

The two plots will now be displayed in the same location (see figure 890, page 736).




Figure 890 The plots are now assigned the same location

11.9 Viewing the data grid

When a measurement is finished, it is possible to inspect the details of all the data and events recorded by each measurement command in the procedure in the **data grid**. The data grid can be accessed by selecting a command and clicking the  button, located in the top right corner of the procedure editor (see figure 891, page 737).




Figure 891 Opening the data grid from the procedure editor

It is also possible to display the data grid directly from the detailed view of a plot (see figure 892, page 738).




Figure 892 Opening the data grid from the detailed plot view

The data grid will be displayed. The data grid contains all the data and events recorded by the selected measurement command (see figure 893, page 738).

CV staircase

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range
0,00244141	7,32305	4,28467E-9	0,0037323	1	1	1,51173E-6	0	10 nA
0,00488281	7,34745	5,66406E-9	0,00618896	1	2	1,51173E-6	0	10 nA
0,00732422	7,37185	9,45435E-9	0,00868225	1	3	1,51173E-6	0	10 nA
0,00976563	7,39625	1,27686E-8	0,0110535	1	4	1,51173E-6	0	10 nA
0,012207	7,42065	1,42822E-8	0,0134888	1	5	1,51173E-6	0	10 nA
0,0146484	7,44505	1,53625E-8	0,0159668	1	6	1,51173E-6	0	10 nA
0,0170898	7,46945	1,72028E-8	0,018335	1	7	1,51173E-6	0	10 nA
0,0195313	7,49385	2,23083E-8	0,0207825	1	8	1,51173E-6	0	10 nA
0,0219727	7,51825	2,52716E-8	0,0233002	1	9	1,51173E-6	0	10 nA
0,0244141	7,54265	2,56836E-8	0,0256683	1	10	1,51173E-6	0	10 nA
0,0268555	7,56705	2,79755E-8	0,0281464	1	11	1,51173E-6	0	10 nA
0,0292969	7,59145	3,23334E-8	0,0306854	1	12	1,51173E-6	0	10 nA
0,0317383	7,61585	3,59009E-8	0,0330444	1	13	1,51173E-6	0	10 nA

Figure 893 The data grid shows all the measured data for the measurement command

11.9.1 Current range logged in the data grid

The current range used for recording the WE(1).Current in any electrochemical measurement is always reported in the data grid (see figure 894, page 739).

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range
0,0512695	7,81105	5,492249E-8	0,0527039	1	21	1,51173E-6	0	10 nA
0,0537109	7,83545	5,770264E-8	0,0549011	1	22	1,51173E-6	0	10 nA
0,0561523	7,85985	5,867615E-8	0,0575256	1	23	1,51173E-6	0	10 nA
0,0585938	7,88425	6,03302E-8	0,0599976	1	24	1,51173E-6	0	10 nA
0,0610352	7,90865	6,378174E-8	0,0622559	1	25	1,51173E-6	0	100 nA
0,0634766	7,93305	6,396484E-8	0,0646667	1	26	1,51173E-6	0	100 nA
0,065918	7,95745	6,863403E-8	0,0671692	1	27	1,51173E-6	0	100 nA
0,0683594	7,98185	6,988525E-8	0,0696716	1	28	1,51173E-6	0	100 nA
0,0708008	8,00625	7,150269E-8	0,071991	1	29	1,51173E-6	0	100 nA
0,0732422	8,03065	7,65686E-8	0,0745544	1	30	1,51173E-6	0	100 nA
0,0756836	8,05505	7,888794E-8	0,0769348	1	31	1,51173E-6	0	100 nA
0,078125	8,07945	8,178711E-8	0,0794067	1	32	1,51173E-6	0	100 nA
0,0805664	8,10385	8,21228E-8	0,0818176	1	33	1,51173E-6	0	100 nA

Figure 894 The current range by the instrument is logged in the data grid

If the current range was modified by the procedure during a measurement, this will be visible in the Current range column of the data grid.



NOTICE

Only the current range of the Autolab PGSTAT instrument is logged during a measurement.

11.9.2 Events logged in the data grid

Events taking place during a measurement are logged in the data grid. The following columns may become visible in the data grid if an applicable event was detected during a measurement:

- **Overloads:** these events correspond to situations where a current, voltage or temperature overload was detected during a measurement.
- **Cutoffs:** these events correspond to situations where a cutoff condition is met.
- **Counters:** these events correspond to situations where a counter is activated.
- **User events:** these events correspond to situations where the user changed a measurement property during a measurement or used a flow control option (stop, pause, reverse scan direction) provided by NOVA.

Figure 895 shows an example of events logged in the data grid.



LSV staircase

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Index	Current range	Overload	Cutoffs	Counters	User events
1,04401	12,075	1,04309E-6	1,04492	18	1 µA				
1,04645	12,3191	1,04553E-6	1,04706	19	1 µA				
1,04889	12,5633	1,04828E-6	1,0495	20	1 µA			Pulse	
1,05133	12,8074	1,05042E-6	1,05194	21	1 µA				
1,05377	13,0516	1,05253E-6	1,05408	22	1 µA				
1,05621	13,2957	1,05499E-6	1,05652	23	1 µA				
1,05865	13,5399	1,05743E-6	1,05927	24	1 µA				
1,0611	13,784	1,06018E-6	1,0611	25	1 µA				
1,06354	14,0281	1,06262E-6	1,06445	26	1 µA				
1,06598	14,2723	1,06537E-6	1,06689	27	1 µA				
1,06842	14,5164	1,06781E-6	1,06903	28	1 µA				Scan rate from 0.01 to 0.1 V/s
1,07086	14,7606	1,06964E-6	1,07117	29	1 µA				
1,0733	15,0047	1,06903E-6	1,0733	30	1 µA			Pulse	
1,07574	15,2488	1,07422E-6	1,07605	31	1 µA				
1,07819	15,493	1,0791E-6	1,0791	32	1 µA				
1,08063	15,7371	1,08185E-6	1,08093	33	1 µA				
1,08307	15,9813	1,08093E-6	1,08337	34	1 µA				
1,08551	16,2254	1,08093E-6	1,08582	35	1 µA				

Figure 895 Events are logged in the data grid

11.9.3 Formatting the data grid

The formatting of the columns can be modified by right-clicking one of the column headers and selecting the required number formatting from the context menu (see figure 896, page 740).

CV staircase

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range
0,00244141	7,32305	4,28467E		1	1	1,51173E-6	0	10 nA
0,00488281	7,34745	5,66406E		1	2	1,51173E-6	0	10 nA
0,00732422	7,37185	9,45435E		1	3	1,51173E-6	0	10 nA
0,00976563	7,39625	1,27686E		1	4	1,51173E-6	0	10 nA
0,012207	7,42065	1,42822E		1	5	1,51173E-6	0	10 nA
0,0146484	7,44505	1,53625E-8	0,0159668	1	6	1,51173E-6	0	10 nA
0,0170898	7,46945	1,72028E-8	0,018335	1	7	1,51173E-6	0	10 nA
0,0195313	7,49385	2,23083E-8	0,0207825	1	8	1,51173E-6	0	10 nA
0,0219727	7,51825	2,52716E-8	0,0233002	1	9	1,51173E-6	0	10 nA
0,0244141	7,54265	2,56836E-8	0,0256683	1	10	1,51173E-6	0	10 nA
0,0268555	7,56705	2,79755E-8	0,0281464	1	11	1,51173E-6	0	10 nA
0,0292969	7,59145	3,23334E-8	0,0306854	1	12	1,51173E-6	0	10 nA
0,0317383	7,61585	3,59009E-8	0,0330444	1	13	1,51173E-6	0	10 nA

Figure 896 The formatting used in the data grid can be specified

The number of significant digits or decimals can also be specified for each signal, by extending the context menu and specifying the required precision (see figure 897, page 741).

CV staircase

Potential applied (V)	Time (s)	WE(1)	Potential (V)	Scan	Index	Q+	Q-	Current range
0,00244141	7,32305	4,2846	3	1	1	1,51173E-6	0	10 nA
0,00488281	7,34745	5,6640	96	1	2	1,51173E-6	0	10 nA
0,00732422	7,37185	9,4543	25	1	3	1,51173E-6	0	10 nA
0,00976563	7,39625	1,2768	5	1	4	1,51173E-6	0	10 nA
0,012207	7,42065	1,4282	1	1	5	1,51173E-6	0	10 nA
0,0146484	7,44505	1,53625E-8	0,015966	1	6	1,51173E-6	0	10 nA
0,0170898	7,46945	1,72028E-8	0,018335	1	7	1,51173E-6	0	10 nA
0,0195313	7,49385	2,23083E-8	0,020782	1	8	1,51173E-6	0	10 nA
0,0219727	7,51825	2,52716E-8	0,023300	1	9	1,51173E-6	0	10 nA
0,0244141	7,54265	2,56836E-8	0,025668	1	10	1,51173E-6	0	10 nA
0,0268555	7,56705	2,79755E-8	0,028146	1	11	1,51173E-6	0	10 nA
0,0292969	7,59145	3,23334E-8	0,030685	1	12	1,51173E-6	0	10 nA
0,0317383	7,61585	3,59009E-8	0,033044	1	13	1,51173E-6	0	10 nA
					14			
					15			
					16			
					17			
					18			
					19			
					20			

Figure 897 The number of significant digits or decimals can be specified



NOTICE

The formatting of the columns in the data grid is saved when the file is saved.

11.9.4 Sorting the data grid

It is possible to sort the contents of the data grid by clicking one of the column header. This will sort the content of the column ascending or descending and the other columns of the data grid will be sorted based on the new order of the sorted column. Clicking the column header cycles from ascending sorting to descending sorting (see figure 898, page 742).



CV staircase

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index ▼	Q+	Q-	Current range
0,546875	12,7643	5,488586E-7	0,547699	1	224	1,51173E-6	0	100 nA
0,544434	12,7399	5,465698E-7	0,54541	1	223	1,51173E-6	0	100 nA
0,541992	12,7155	5,436401E-7	0,542877	1	222	1,51173E-6	0	100 nA
0,539551	12,6911	5,424194E-7	0,540375	1	221	1,51173E-6	0	100 nA
0,537109	12,6667	5,388184E-7	0,537903	1	220	1,51173E-6	0	100 nA
0,534668	12,6423	5,362244E-7	0,535675	1	219	1,51173E-6	0	100 nA
0,532227	12,6179	5,370178E-7	0,533112	1	218	1,51173E-6	0	100 nA
0,529785	12,5935	5,339966E-7	0,53064	1	217	1,51173E-6	0	100 nA
0,527344	12,5691	5,283508E-7	0,52832	1	216	1,51173E-6	0	100 nA
0,524902	12,5447	5,267029E-7	0,525909	1	215	1,51173E-6	0	100 nA
0,522461	12,5203	5,250549E-7	0,523285	1	214	1,51173E-6	0	100 nA
0,52002	12,4959	5,237122E-7	0,520996	1	213	1,51173E-6	0	100 nA
0,517578	12,4715	5,19165E-7	0,518555	1	212	1,51173E-6	0	100 nA

Figure 898 Sorting the contents of the data grid



NOTICE

A column sorted in ascending mode is indicated by the ▼ symbol. A column sorted in descending mode is indicated by the ▲ symbol.

11.9.5 Changing the order of the columns in the data grid

It is possible to change the order of the columns. To move a column in the data grid, click one of the column headers and drag the mouse left or right in the grid, while holding the mouse button (see figure 899, page 742).

CV staircase

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range
0,00244141	7,32305	4,284668E-9	0,0037323	1	1	1,51173E-6	0	10 nA
0,00488281	7,34745	5,664063E-9	0,00618896	1	2	1,51173E-6	0	10 nA
0,00732422	7,37185	9,454346E-9	0,00868225	1	3	1,51173E-6	0	10 nA
0,00976563	7,39625	1,276855E-8	0,0110535	1	4	1,51173E-6	0	10 nA
0,012207	7,42065	1,428223E-8	0,0134888	1	5	1,51173E-6	0	10 nA
0,0146484	7,44505	1,536255E-8	0,0159668	1	6	1,51173E-6	0	10 nA
0,0170898	7,46945	1,720276E-8	0,018335	1	7	1,51173E-6	0	10 nA
0,0195313	7,49385	2,230835E-8	0,0207825	1	8	1,51173E-6	0	10 nA
0,0219727	7,51825	2,527161E-8	0,0233002	1	9	1,51173E-6	0	10 nA
0,0244141	7,54265	2,568359E-8	0,0256683	1	10	1,51173E-6	0	10 nA
0,0268555	7,56705	2,797546E-8	0,0281464	1	11	1,51173E-6	0	10 nA
0,0292969	7,59145	3,233337E-8	0,0306854	1	12	1,51173E-6	0	10 nA
0,0317383	7,61585	3,590088E-8	0,0330444	1	13	1,51173E-6	0	10 nA

Figure 899 The order of the columns can be modified

Release the mouse button validate the new location of the column (see figure 900, page 743).



Index	Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Q+	Q-	Current range
1	0,00244141	7,32305	4,284668E-9	0,0037323	1	1,51173E-6	0	10 nA
2	0,00488281	7,34745	5,664063E-9	0,00618896	1	1,51173E-6	0	10 nA
3	0,00732422	7,37185	9,454346E-9	0,00868225	1	1,51173E-6	0	10 nA
4	0,00976563	7,39625	1,276855E-8	0,0110535	1	1,51173E-6	0	10 nA
5	0,012207	7,42065	1,428223E-8	0,0134888	1	1,51173E-6	0	10 nA
6	0,0146484	7,44505	1,536255E-8	0,0159668	1	1,51173E-6	0	10 nA
7	0,0170898	7,46945	1,720276E-8	0,018335	1	1,51173E-6	0	10 nA
8	0,0195313	7,49385	2,230835E-8	0,0207825	1	1,51173E-6	0	10 nA
9	0,0219727	7,51825	2,527161E-8	0,0233002	1	1,51173E-6	0	10 nA
10	0,0244141	7,54265	2,568359E-8	0,0256683	1	1,51173E-6	0	10 nA
11	0,0268555	7,56705	2,797546E-8	0,0281464	1	1,51173E-6	0	10 nA
12	0,0292969	7,59145	3,233337E-8	0,0306854	1	1,51173E-6	0	10 nA
13	0,0317383	7,61585	3,590088E-8	0,0330444	1	1,51173E-6	0	10 nA

Figure 900 The new order of columns in the data grid




NOTICE

The order of the columns in the data grid is saved when the file is saved.


11.9.6 Exporting the data from the data grid

Finally, the data grid can also be used to export the data in the grid to an ASCII file or an Excel file. To export the data, click the  button in the top right corner (see figure 901, page 743).



Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range
0,00244141	7,32305	4,284668E-9	0,0037323	1	1	1,51173E-6	0	10 nA
0,00488281	7,34745	5,664063E-9	0,00618896	1	2	1,51173E-6	0	10 nA
0,00732422	7,37185	9,454346E-9	0,00868225	1	3	1,51173E-6	0	10 nA
0,00976563	7,39625	1,276855E-8	0,0110535	1	4	1,51173E-6	0	10 nA
0,012207	7,42065	1,428223E-8	0,0134888	1	5	1,51173E-6	0	10 nA
0,0146484	7,44505	1,536255E-8	0,0159668	1	6	1,51173E-6	0	10 nA
0,0170898	7,46945	1,720276E-8	0,018335	1	7	1,51173E-6	0	10 nA
0,0195313	7,49385	2,230835E-8	0,0207825	1	8	1,51173E-6	0	10 nA
0,0219727	7,51825	2,527161E-8	0,0233002	1	9	1,51173E-6	0	10 nA
0,0244141	7,54265	2,568359E-8	0,0256683	1	10	1,51173E-6	0	10 nA
0,0268555	7,56705	2,797546E-8	0,0281464	1	11	1,51173E-6	0	10 nA
0,0292969	7,59145	3,233337E-8	0,0306854	1	12	1,51173E-6	0	10 nA
0,0317383	7,61585	3,590088E-8	0,0330444	1	13	1,51173E-6	0	10 nA

Figure 901 The data points can be exported using the provided button
A menu will pop-out, as shown in Figure 902.




The screenshot shows a data grid titled "CV staircase" with columns: Potential applied (V), Time (s), WE(1).Current (A), WE(1).Potential (V), Scan, Index, Q+, Q-, and Current range. The data is presented in a table with 13 rows. A pop-up menu is open on the right, showing export settings: File format (ASCII), Write column headers (checked), Column delimiter (Semicolon (;)), and Decimal separator (Comma (,)). An "Export" button is visible at the bottom right of the menu.

Potential applied (V)	Time (s)	WE(1).Current (A)	WE(1).Potential (V)	Scan	Index	Q+	Q-	Current range
0,00244141	7,32305	4,284668E-9	0,0037323	1	1	1,51173E-6	0	10 nA
0,00488281	7,34745	5,664063E-9	0,00618896	1	2	1,51173E-6	0	10 nA
0,00732422	7,37185	9,454346E-9	0,00868225	1	3	1,51173E-6	0	10 nA
0,00976563	7,39625	1,276855E-8	0,0110535	1	4	1,51173E-6	0	10 nA
0,012207	7,42065	1,428223E-8	0,0134888	1	5	1,51173E-6	0	10 nA
0,0146484	7,44505	1,536255E-8	0,0159668	1	6	1,51173E-6	0	10 nA
0,0170898	7,46945	1,720270E-8	0,018335	1	7	1,51173E-6	0	10 nA
0,0195313	7,49385	2,230835E-8	0,0207825	1	8	1,51173E-6	0	10 nA
0,0219727	7,51825	2,527161E-8	0,0233002	1	9	1,51173E-6	0	10 nA
0,0244141	7,54265	2,568359E-8	0,0256683	1	10	1,51173E-6	0	10 nA
0,0268555	7,56705	2,797546E-8	0,0281464	1	11	1,51173E-6	0	10 nA
0,0292969	7,59145	3,233337E-8	0,0306854	1	12	1,51173E-6	0	10 nA
0,0317383	7,61585	3,590088E-8	0,0330444	1	13	1,51173E-6	0	10 nA

Figure 902 The export settings are specified in the a dedicated pop-out menu


The menu can be used to specify the file format, using the provided drop-down list (see figure 903, page 744).



This close-up shows the export settings menu. The "File format" dropdown is set to "ASCII". The "Write column headers" checkbox is checked. The "Column delimiter" dropdown is set to "Excel". The "Decimal separator" dropdown is set to "Comma (,)". An "Export" button is located at the bottom right.

Figure 903 The data can be exported as ASCII or Excel

When the data is exported as ASCII (Comma Separated Values), additional settings can be specified (see figure 904, page 744). These settings depend on the required output format of the data.



This close-up shows the export settings menu for ASCII format. The "File format" dropdown is set to "ASCII". The "Write column headers" checkbox is checked. The "Column delimiter" dropdown is set to "Semicolon (;)". The "Decimal separator" dropdown is set to "Comma (,)". An "Export" button is located at the bottom right.

Figure 904 The settings used to export data to a ASCII file

When the data is exported as Excel, the file is automatically formatted (see figure 905, page 745).




Figure 905 The settings used to export data to an Excel file

The following settings can be specified:

- **File format:** specifies the format of the output file (ASCII or Excel), using the provided drop-down list.
- **Write column headers:** a toggle that can be used to indicate if the names of the signals need to be added to the output file.
- **Column delimiter:** specifies the symbol used as a column separator, using the provided drop-down list. This property is only available for ASCII output.
- **Decimal separator:** specifies the decimal separator symbol used in the output file, using the provided drop-down list. This property is only available for ASCII output.

Clicking the button displays a save dialog window which can be used to specify the filename and location (see figure 906, page 745).




Figure 906 Specifying the filename and location

11.10 Convert data to procedure

At the end of measurement, the measured data can be converted to a new procedure. In order to do this, click the  button in the top right corner of the procedure editor (see figure 907, page 746).




Figure 907 Converting data to procedure

If the procedure was modified during the measurement or if data analysis tools were added to the data, a message will be displayed when the button is clicked, providing the means to define how the data should be converted to a new procedure (see figure 908, page 746).




Figure 908 The changes can be kept or discarded

- Clicking the **Yes** button will convert the modified data to a new procedure. All modifications and changes that were carried out during and after the measurement will be added to the source procedure used to generate this data.
- Clicking the **No** button will convert the data to a new procedure but will discard all the changes that were carried out during and after the measurement.
- Clicking the **Cancel** button will cancel this action.

12 Data analysis

When data has been measured, it is possible to use the data analysis commands provided in NOVA to analyze the data. To analyze acquired data in NOVA, it is necessary to add the required command to the measured procedure and apply the function of these commands on the measure data.



NOTICE

Data analysis commands can be added to the initial procedure or the procedure after the measurement is finished.

To add a data analysis command to a measured procedure, two methods can be used:

- Drag and drop the analysis command in the procedure
- Use the contextual shortcut  button, located in the top right corner of the procedure editor

The functionality of the data analysis commands is explained in the previous chapters and will not be detailed again in this chapter. This chapter focuses on the use of these commands on **measured** data. Only the commands that provide controls that are used in a specific way on analyzed data are detailed in this chapter.

The following commands are detailed:

- Smooth
- Peak search
- Regression
- Integrate
- Interpolate
- Baseline correction
- Corrosion rate
- Hydrodynamic analysis
- Electrochemical circle fit
- Fit and simulation




Figure 910 The Smooth command is added to the procedure



NOTICE

For more information on the properties of the **Smooth** command, please refer to *Chapter 7.8.1*.

Clicking the **More** button opens a new screen in which the additional controls of the **Smooth** command are shown for the scope of data analysis. The plot on the left hand side shows the source data and the curve drawn by the **Smooth** command. The properties of the **Smooth** command are all set to their default values (see figure 911, page 750).




Figure 911 The additional controls of the Smooth command

Since the **Smooth** command has two different modes, each mode provides dedicated controls. The *Mode* drop-down list can be used to change the mode of the command.

12.1.1 SG mode

In *Savitzky-Golay* mode (SG), the **Smooth** command applies a smoothing algorithm based on the Savitzky-Golay algorithm. The properties can be adjusted in the panel on the right hand side (see figure 912, page 751).




Figure 912 The Smooth command with the default properties of the Savitzky-Golay mode

A preview of the smoothed curve, obtained by using the properties defined on the right hand side, is shown in green, overlaid on the source data. Using this preview, it is possible to fine tune the properties and see the effect on the expected result of the smoothing (see figure 913, page 751).




Figure 913 The preview curve is automatically adjusted when the properties are changed in the panel on the right hand side

12.2 Peak search

The **Peak search** command provides additional controls that can be used when the command is used to analyze data. To use the **Peak search** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands and possible plots on which these command can be applied (see figure 917, page 754).




Figure 917 Adding a Peak search command to the *i vs E* plot

The **Peak search** command is added to the procedure editor. Clicking the command shows the properties in the dedicated panel on the right hand side (see figure 918, page 755).




Figure 918 The Peak search command is added to the procedure



NOTICE

For more information on the properties of the **Peak search** command, please refer to *Chapter 7.8.2*.

Clicking the **More** button opens a new screen in which the additional controls of the **Peak search** command are shown for the scope of data analysis. The plot on the left hand side shows the source data and the peaks identified by the **Peak search** command. The properties of the **Peak search** command are all set to their default values (see figure 919, page 756).




Figure 919 The additional controls of the Peak search command

Since the **Peak search** command has two different search modes (see figure 920, page 756):

- **Automatic peak:** the peaks are automatically found based on the properties defined in the panel on the right hand side.
- **Manual:** peaks are identified by specifying a baseline manually.




Figure 920 The search mode can be set to Automatic or Manual



CAUTION

Switching search mode will clear all previously found peaks.

12.2.1 Automatic search mode

When the *Search mode* is set to Automatic, the peaks are automatically identified by the **Peak search** command. If the search properties are modified, the command is automatically refreshed.



NOTICE

The Automatic Search mode uses a Linear tangent type of baseline. This baseline is only available for the Automatic Search mode.

12.2.2 Manual peak search

When the *Search mode* is set to Manual, the type of baseline used to find the peaks can be specified using the provided *Baseline mode* dropdown list (see figure 921, page 757).




Figure 921 The base line mode can be specified using the provided dropdown list

The following baseline modes are available:

- Exponential
- Zero base
- Polynomial
- Linear curve cursor
- Linear free cursor
- Linear front
- Linear rear
- Linear front tangent

- Linear rear tangent

To carry out a manual peak search, it is necessary to draw a baseline using the mouse pointer. To do this, click the plot near the beginning of a peak and drag the mouse pointer across the plot area to draw the region in which the baseline should be located (see figure 922, page 758).




Figure 922 Defining a manual baseline

A line will be drawn as the mouse is dragged across the plot area. Click the mouse button is again to define the end of the baseline. When the baseline is defined, the peaks will be identified based on the properties defined in the **Properties** panel (see figure 923, page 758).




Figure 923 The baseline is drawn on the plot



NOTICE

The points used to define the baseline are identified by small vertical lines on the plot.

At any point it is possible to click the **Reset** button to remove all the peaks found (see figure 924, page 759).




Figure 924 Resetting the peak search results

12.2.2.1 Exponential

This option uses an exponential baseline in the determination of the peaks.

To define the baseline, click on the plot area to define the start point of the baseline and drag the mouse across the plot area to define the baseline. Click the mouse button again to define the end point of the baseline. When the start point and end point of the baseline have been defined, the X coordinates of these points will be used to find the closest points on the curve and the exponential baseline will be drawn between these points (see figure 925, page 760).




Figure 925 Manual peak search using the exponential baseline

12.2.2.2 Zero base

Using the **zero base** no baseline is used in the determination of the peak.

To define the baseline, click on the plot area to define the start point of the baseline and drag the mouse across the plot area to define the baseline. Click the mouse button again to define the end point of the baseline. When the start point and end point of the baseline have been defined the data point on Y axis, with the highest absolute value, located within the range defined by the start point and end point of the baseline is identified as a peak (see figure 926, page 761).




Figure 926 Manual peak search using the zero base baseline



NOTICE

The zero base search method locates the absolute maximum value of the curve in the curve segment closest to the first point defining the **search window**.

12.2.2.3 Polynomial

This baseline uses a polynomial function in the determination of the peaks.

To define the baseline, click on the plot area to define the start point of the baseline. Drag the mouse across the plot and click again to add way-points for the polynomial function. This can be repeated as many times as required. To define the end point of the baseline, press the **[Enter]** key on the keyboard. The end point will be set to the last location of the mouse pointer.

When the end point has been defined, the X coordinates of the start point and end point will be used to find the closest points on the curve and the polynomial baseline will be drawn between these points (see figure 925, page 760).




Figure 927 Manual peak search using the polynomial baseline

12.2.2.4 Linear curve cursor

This option uses a linear baseline in the determination of the peaks.

To define the baseline, click on the plot area to define the start point of the baseline and drag the mouse across the plot area to define the baseline. Click the mouse button again to define the end point of the baseline. When the start point and end point of the baseline have been defined, the X coordinates of these points will be used to find the closest points on the curve and the linear baseline will be drawn between these points (see figure 928, page 763).




Figure 928 Manual peak search using the linear curve cursor baseline

12.2.2.5 Linear free cursor

This option uses a linear baseline in the determination of the peaks.

To define the baseline, click on the plot area to define the start point of the baseline and drag the mouse across the plot area to define the baseline. Click the mouse button again to define the end point of the baseline. This baseline is not connected to the nearest data points on the curve (see figure 929, page 764).




Figure 929 Manual peak search using the linear free cursor baseline

12.2.2.6 Linear front

This option finds peaks by extending a tangent baseline located in front of the peak.

To define the baseline, click on the plot area to define the start point of the baseline and drag the mouse across the plot area to define the baseline. Click the mouse button again to define the end point of the baseline. When the start point and end point of the baseline have been defined, the tangent is extended frontwards and the peak is located (*see figure 930, page 765*).




Figure 930 Manual peak search using the linear front baseline

12.2.2.7 Linear rear

This option finds peaks by extending the baseline located after the peak.

To define the baseline, click on the plot area to define the start point of the baseline and drag the mouse across the plot area to define the baseline. Click the mouse button again to define the end point of the baseline. When the start point and end point of the baseline have been defined, the tangent is extended backwards and the peak is located (*see figure 931, page 766*).




Figure 931 Manual peak search using the linear rear baseline

12.2.2.8 Linear front tangent

This option finds peaks by extending the baseline located in front of the peak.

To define the baseline, click on the plot area to define the start point of the baseline and drag the mouse across the plot area to define the baseline. Click the mouse button again to define the end point of the baseline. When the start point and end point of the baseline have been defined, the software automatically connects the baseline to the curve at the data point for which the first derivative is the closest to the slope of drawn baseline (see figure 932, page 767).




Figure 932 Manual peak search using the linear front tangent baseline

12.2.2.9 Linear rear tangent

This option finds peaks by extending the baseline located after the peak.

To define the baseline, click on the plot area to define the start point of the baseline and drag the mouse across the plot area to define the baseline. Click the mouse button again to define the end point of the baseline. When the start point and end point of the baseline have been defined, the software automatically connects the baseline to the curve at the data point for which the first derivative is the closest to the slope of drawn baseline (see figure 933, page 768).




Figure 933 Manual peak search using the linear rear tangent baseline

12.2.3 Manual adjustments

Peaks found using the manual search mode can be finetuned at any time by moving one of the small vertical lines defining the location of the baseline on the plot (see figure 934, page 768).




Figure 934 It is possible to modify the baseline points

To modify one of the baseline markers, click the marker and while holding the mouse button, move the marker left or right along the curve. While the mouse button is held, the coordinates of the mouse button on the curve are shown (see figure 935, page 769).




Figure 935 Moving a baseline marker shows the coordinates

When the mouse button is released, a new baseline is determined based on the new location of the displaced marker and the peak will be redetermined automatically. The results will also be automatically updated (see figure 936, page 769).




Figure 936 The new baseline is calculated and the peak is updated

12.3 Regression analysis

The **Regression** command provides additional controls that can be used when the command is used to analyze data. To use the **Regression** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands and possible plots on which these command can be applied (see figure 938, page 771).




Figure 938 Adding a Regression command to the Nyquist plot

The **Regression** command is added to the procedure editor. Clicking the command shows the properties in the dedicated panel on the right hand side (see figure 939, page 772).




Figure 939 The Regression command is added to the procedure



NOTICE

For more information on the properties of the **Regression** command, please refer to *Chapter 7.8.3*.

Clicking the **More** button opens a new screen in which the additional controls of the **Regression** command are shown for the scope of data analysis. The plot on the left hand side shows the source data and the curve drawn by the **Regression** command. The properties of the **Regression** command are all set to their default values (see figure 940, page 773).




Figure 940 The additional controls of the Regression command

Using the mouse, it is possible to manually draw the area of the plot on which the **Regression** command should be executed. By clicking and holding the mouse button, a specific area can be drawn. This area will be delimited by two vertical lines and will be shown with a light green background (see figure 941, page 773).




Figure 941 Manually defining boundaries for the Regression command



NOTICE

The results of the **Regression** command are automatically recalculated each time one of the properties is modified or each time the **Regression** command is used on a specific area of the plot.



Once the boundaries of the **Regression** command have been specified, it is possible to manually adjust these boundaries by clicking either one of the boundary lines and dragging the line left or right (see figure 942, page 774).




Figure 942 Adjusting the boundaries of the Regression command

It is also possible to fine tune the properties and the boundaries manually in the **Properties** panel on the right hand side (see figure 943, page 774).




Figure 943 Finetuning the properties of the Regression command

Finally, clicking the **Reset** button resets all the properties of the **Regression** command back to the default values (see figure 944, page 774).




Figure 944 Resetting the properties of the Regression command

12.4 Integrate

The **Integrate** command provides additional controls that can be used when the command is used to analyze data. To use the **Integrate** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands and possible plots on which these command can be applied (see figure 945, page 775).




Figure 945 Adding a *Integrate* command to the *i vs t* plot

The **Integrate** command is added to the procedure editor. Clicking the command shows the properties in the dedicated panel on the right hand side (see figure 946, page 776).




Figure 946 The Integrate command is added to the procedure



NOTICE

For more information on the properties of the **Integrate** command, please refer to *Chapter 7.8.5*.

Clicking the **More** button opens a new screen in which the additional controls of the **Integrate** command are shown for the scope of data analysis. The plot on the left hand side shows the source data and the area calculated by the **Integrate** command is shown on the right-hand side of the plot (see figure 947, page 777).




Figure 947 The additional controls of the Integrate command

By default, the whole plot is integrated. The boundaries used for the integration of the data are represented by vertical lines on either side of the plot (see figure 947, page 777). Using the mouse, it is possible to manually adjust these boundaries by clicking either one of the boundary lines and dragging the line left or right (see figure 948, page 778).




Figure 948 Adjusting the boundaries of the Integrate command



NOTICE

The results of the **Integrate** command are automatically recalculated each time one of the properties is modified or each time the **Integrate** command is used on a specific area of the plot.

Finally, clicking the **Reset** button resets all the properties of the **Integrate** command back to the default values (see figure 949, page 779).




Figure 949 Resetting the properties of the Integrate command

12.5 Interpolate

The **Interpolate** command provides additional controls that can be used when the command is used to analyze data. To use the **Interpolate** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands and possible plots on which these command can be applied (see figure 950, page 780).




Figure 950 Adding a Interpolate command to the δi vs E plot

The **Interpolate** command is added to the procedure editor. Clicking the command shows the properties in the dedicated panel on the right hand side (see figure 951, page 780).




Figure 951 The Interpolate command is added to the procedure



NOTICE

For more information on the properties of the **Interpolate** command, please refer to *Chapter 7.8.6*.

Clicking the **More** button opens a new screen in which the additional controls of the **Interpolate** command are shown for the scope of data analysis. The plot on the left hand side shows the source data and the properties of the **Interpolate** command are shown on the right-hand side of the plot (see figure 952, page 781).




Figure 952 The additional controls of the Interpolate command

By default, the **Interpolate** command is executed, searching for Y value at a X position of 0. The results that match this search criteria are listed in the **Results** panel and indicated by the lines on the plot (see figure 953, page 782).




Figure 953 It is possible to adjust the properties of the Interpolate command

The dark green line indicates the location of the position at which the **Interpolate** command is carried out. The light green line indicates the position of the value(s) found by the **Interpolate** command. It is possible to change the position at which the Interpolate command is carried out by changing the value in the provided field in the **Properties** panel (see figure 954, page 782).




Figure 954 The Interpolate command is updated when the properties are changed

The command will be updated and the new results will be displayed graphically and in the **Results** panel. It is also possible to move the dark green line indicating the position at which the **Interpolate** command is carried out using the mouse and dragging the line across the plot (see figure 955, page 783).






Figure 955 It is possible to move the line indicating the position at which the Interpolate command is carried out

The command will be updated when the mouse button is released and the results will be updated as indicated above (see figure 956, page 783).




Figure 956 The Interpolate command is updated when the properties are changed

When the **Interpolate** command is able to find more than one value, as shown in Figure 957, each of the values found will be listed in the **Results** panel and will be indicated graphically by light green lines on the plot.




Figure 957 More than one value can be found by the Interpolate command

12.6 Hydrodynamic analysis

The **Hydrodynamic analysis** command provides additional controls that can be used when the command is used to analyze data. To use the **Hydrodynamic analysis** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands and possible plots on which these command can be applied (see figure 958, page 785).




Figure 958 Adding a Hydrodynamic analysis command to the measurement

The **Hydrodynamic analysis** command is added to the procedure editor. Clicking the command shows the properties in the dedicated panel on the right hand side (see figure 959, page 786).




Figure 959 The Hydrodynamic analysis command is added to the procedure



NOTICE

For more information on the properties of the **Hydrodynamic analysis** command, please refer to *Chapter 7.8.10*.

Clicking the **More** button opens a new screen in which the additional controls of the **Hydrodynamic analysis** command are shown for the scope of data analysis. The plot on the left hand side shows the source data. The plots on the right hand side show the regression lines generated by the command based on the selected current values (see figure 960, page 787).




Figure 960 The additional controls of the Hydrodynamic analysis command

The currents are selected using the vertical green line show in the plot on the left hand side. By default, the line is drawn at the position of the first data point (corresponding to index 1) of each curve. To reposition the line, click the line and while holding the mouse button, slide the line across the plot area (see figure 961, page 787).




Figure 961 Moving the vertical line to specifying the current

The selected current are updated as the line is moved. Releasing the mouse button validates the selection of the limiting currents (see figure 962, page 788).




Figure 962 The current values are updated

The **Hydrodynamic analysis** command automatically carries out a Levich analysis (which is normally carried out on the mass-transport limited current values) and a Koutecký-Levich analysis (which is normally carried out in the mixed kinetic - mass-transport regime). A linear regression is carried out on both these analysis methods and the results are displayed below the corresponding plots on the right-hand side.

For the Levich plot, the Slope and Intercept are provided. For the Koutecký-Levich plot, the same information is provided, as well as the extrapolated kinetic current, i_k , obtained from the intercept on the plot.

12.7 Baseline correction

The **Baseline correction** command provides additional controls that can be used when the command is used to analyze data. To use the **Baseline correction** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands and possible plots on which these command can be applied (see figure 963, page 789).




Figure 963 Adding a Baseline correction command to the Differential pulse command

The **Baseline correction** command is added to the procedure editor. Clicking the command shows the properties in the dedicated panel on the right hand side (see figure 964, page 790).




Figure 964 The Baseline correction command is added to the procedure



NOTICE

For more information on the properties of the **Baseline correction** command, please refer to *Chapter 7.8.13*.

Clicking the **More** button opens a new screen in which the additional controls of the **Baseline correction** command are shown for the scope of data analysis. The plot on the left hand side shows the source data. The properties of the **Baseline correction** command are all set to their default values (see figure 965, page 791).




Figure 965 The additional controls of the Baseline correction command

The mode and properties of the **Baseline correction** command can be adjusted in the **Properties** panel. Using the mouse, it is possible to click the plot to define a point defining the baseline. Depending on the mode selected for the **Baseline correction** command, two or more data points are necessary to define the baseline. When sufficient data points have been defined on the plot, the baseline will be drawn (see figure 966, page 791).




Figure 966 The baseline is drawn as soon as enough data points have been specified



NOTICE

The coordinates of the selected points are added to the table on the right-hand side. The table allows data points to be modified manually or removed.

Once the baseline is defined, the residual plot is automatically created in the **Plots** frame (see figure 967, page 792).




Figure 967 The residual plot is automatically created in the Plots frame when the baseline is defined

It is possible to add extra points to define the baseline by clicking additional points on the plot. Each new point added to the plot forces the baseline to be recalculated (see figure 968, page 792).




Figure 968 The baseline is update each time a point is added to the plot

Each change to the drawn baseline in turn forces the residual plot to be updated in the **Plots** frame (see figure 969, page 793).




Figure 969 Changing the baseline triggers the residual plot to be updated



NOTICE

Adding extra markers to a specific area of the plot increases the relative importance of that specific area of the plot in the baseline correction.

12.7.1 Zooming in/out

If needed, it is possible to use the controls located in the top right corner of the plot to zoom in (🔍) or out (🔍) or to rescale (🔍) the plot. It is also possible to use the controls provided in the **View** menu or the associated keyboard shortcuts (see figure 970, page 794).




Figure 970 It is possible to zoom in or out

Clicking the button or pressing the [F4] key rescales the complete plot (see figure 971, page 794).




Figure 971 Rescaling the plot



NOTICE

When working with a mouse fitted with a wheel, it is possible to zoom in or out using the wheel.

12.7.2 Fine tuning the baseline correction

If needed, it is possible to fine tune the location of the points using the table located in the **Properties** panel. To edit the location of one of the points, click the X or Y cell of the point to edit and click it again to edit the value (see figure 972, page 795).




Figure 972 Editing the location of a point

Type the new value in the selected cell (see figure 973, page 795).




Figure 973 Fine tuning the location of the selected point

Clicking away from the cell or pressing the **[Enter]** key or **[Tab]** key will validate the new location of the point.

If needed, a point marker can be deleted. To delete a point, click the cell located at the left of the X and Y cell of the point. This will select the complete row of the table. Press the **[Delete]** key to delete this point (see figure 974, page 796).




Figure 974 Selecting the point to delete

The point will be removed (see figure 975, page 796).




Figure 975 The selected point is deleted

12.8 Corrosion rate analysis

The **Corrosion rate analysis** command provides additional controls that can be used when the command is used to analyze data.



CAUTION

The **Corrosion rate analysis** command is intended to be used on current data (WE(1).Current) plotted against potential data (Potential applied).

To use the **Corrosion rate analysis** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed,

providing a list of commands and possible plots on which these command can be applied (see figure 976, page 797).




Figure 976 Adding a Corrosion rate analysis command to the linear polarization data

The **Corrosion rate analysis** command is added to the procedure editor. Clicking the command shows the properties and the results in the dedicated panel on the right hand side (see figure 977, page 798).




Figure 977 The Corrosion rate analysis command is added to the procedure



NOTICE

For more information on the properties of the **Corrosion rate analysis** command, please refer to *Chapter 7.8.14*.

Clicking the **More** button opens a new screen in which the additional controls of the **Corrosion rate analysis** command are shown for the scope of data analysis. The plot on the left hand side shows the source data, plotted on a logarithmic scale. The properties of the **Corrosion rate analysis** command are all set to their default values (see figure 978, page 799).




Figure 978 The additional controls of the Corrosion rate analysis command

Since the **Corrosion rate** command has two different modes, each mode provides dedicated controls. The *Mode* drop-down list can be used to change the mode of the command.

12.8.1 Tafel Analysis

In *Tafel Analysis* mode, it is necessary to define two points on the *anodic* part of the Tafel plot and two points on the *cathodic* part of the Tafel plot.

To define a point, click on the plot. The software will automatically select the closest point of the measured data. When two points are defined, the anodic Tafel slope will be drawn on the plot (see figure 979, page 800).




Figure 979 Defining the points for the anodic Tafel slope

The same can be done for the cathodic branch. When the two points are defined, the cathodic Tafel slope is plotted and the intercept of both lines is used to determine the corrosion potential, the exchange current and current density, the polarization resistance and the corrosion rate (see figure 980, page 800).




Figure 980 The two slopes are used to determine the corrosion rate

Depending on the position of the **Perform fit** toggle, the **Corrosion rate analysis** command will either:

1. **Perform fit off:** the command will perform the calculations of the command based on the location of the intercept. This will lead to an approximation of the corrosion data (see figure 980, page 800).
2. **Perform fit on:** the intercept will be used as a starting point for the fitting of the data using the Butler-Volmer equation. The complete curve will be fitted using this equation and the corrosion data will be determined by the results of the fit. This leads to a more accurate determination of the corrosion date, as shown in Figure 981.




Figure 981 The measured data is fitted with the Butler-Volmer equation

Once the four points required by the **Corrosion rate analysis** command have been specified, it is possible to manually adjust the location of these points by clicking a vertical line defining the location of a point and dragging the line left or right (see figure 982, page 802).




Figure 982 Adjusting the points of the Corrosion rate analysis command

As soon as the point is relocated, the calculation of the **Corrosion rate analysis** command will be updated and the new results will be displayed in the Results sub-panel (see figure 983, page 803).




Figure 983 The calculation is refreshed as soon as one of the points is modified



NOTICE

It is also possible to finetune the location of the points used in the **Corrosion rate analysis** command by using the **Selected points** table, located below the **Properties** panel.

Finally, clicking the **Reset** button resets all the properties of the **Corrosion rate analysis** command back to the default values and clears the selected points.

12.8.2 Polarization Resistance

In *Polarization Resistance* mode, no inputs are required. The calculations are automatically carried out with the specified settings. The potential range in which the analysis is carried out is shown in the plot (see figure 984, page 804).




Figure 984 The Polarization Resistance analysis is carried in out in the highlighted range

If needed, the **Range** value can be adjusted. When the value is modified, the calculation is updated and the range is adjusted on the plot (see figure 985, page 804).




Figure 985 The plot is updated when the Range value is change

Changing the value of any other property used by the command will also force the calculation to update.

12.9 Electrochemical circle fit

The **Electrochemical circle fit** command provides additional controls that can be used when the command is used to analyze data.



CAUTION

The **Electrochemical circle fit** command is intended to be used on impedance spectroscopy data.

To use the **Electrochemical circle fit** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands and possible plots on which these command can be applied (see figure 986, page 805).




Figure 986 Adding a Electrochemical circle fit command to impedance data

The **Electrochemical circle fit** command is added to the procedure editor. Clicking the command shows the properties in the dedicated panel on the right hand side (see figure 987, page 806).




Figure 987 The Electrochemical circle fit command is added to the procedure



NOTICE

For more information on the properties of the **Electrochemical circle fit** command, please refer to *Chapter 7.9.1*.

Clicking the **More** button opens a new screen in which the additional controls of the **Electrochemical circle fit** command are shown for the scope of data analysis. The plot on the left hand side shows the source data, presented in a Nyquist plot. The results of the **Electrochemical circle fit** command are shown in the panel on the right hand side (see figure 987, page 806).




Figure 988 The additional controls of the Electrochemical circle fit command



NOTICE

The coordinates of the selected points are added to the table on the right-hand side. The table allows data points to be modified manually or removed.

To use the **Electrochemical circle fit** command, it is necessary to define three or more points along a visible semi-circle in the Nyquist plot to draw a half-circle and determine the properties of the apparent time constant of this part of the plot.

To define a point, click on the plot. The software will automatically select the closest point of the measured data (if the *Snap to data* option is on) or the point will be located where the mouse is clicked (in the *Snap to data* option is off). When three points are defined, the semi-circle will be drawn on the plot and the results will be updated in the panel on the right-hand side (see figure 989, page 808).




Figure 989 The semi-circle is drawn when three or more points are selected

It is possible to add extra points to the plot. The calculation will automatically be refreshed whenever a point is added to the plot.



NOTICE

Adding extra points to a specific area of the plot increases the relative importance of that specific area of the plot in the electrochemical circle fit.



NOTICE

It is also possible to finetune the location of the points used in the **Electrochemical circle fit** command by using the **Selected points** table, located below the **Results** panel.

Finally, clicking the **Reset** button resets all the properties of the **Electrochemical circle fit** command back to the default values and clears the selected points.

12.9.1 Zooming in/out

If needed, it is possible to use the controls located in the top right corner of the plot to zoom in (🔍) or out (🔍) or to rescale (🔍) the plot. It is also possible to use the controls provided in the **View** menu or the associated keyboard shortcuts (see figure 990, page 809).




Figure 990 It is possible to zoom in or out

Clicking the 🔍 button or pressing the **[F4]** key rescales the complete plot (see figure 991, page 809).




Figure 991 Rescaling the plot



NOTICE

When working with a mouse fitted with a wheel, it is possible to zoom in or out using the wheel.

12.9.2 Fine tuning the baseline correction

If needed, it is possible to fine tune the location of the points using the table located in the **Properties** panel. To edit the location of one of the points, click the X or Y cell of the point to edit and click it again to edit the value (see figure 992, page 810).




Figure 992 Editing the location of a point

Type the new value in the selected cell (see figure 993, page 811).




Figure 993 Fine tuning the location of the selected point

Clicking away from the cell or pressing the **[Enter]** key or **[Tab]** key will validate the new location of the point.

If needed, a point marker can be deleted. To delete a point, click the cell located at the left of the X and Y cell of the point. This will select the complete row of the table. Press the **[Delete]** key to delete this point (see figure 994, page 811).




Figure 994 Selecting the point to delete

The point will be removed (see figure 995, page 812).




Figure 995 The selected point is deleted

12.9.3 Copy as equivalent circuit

It is possible to copy the results from the **Electrochemical circle fit** analysis tool to the clipboard as an equivalent circuit. This circuit can then be used in the **Fit and Simulation** command or analysis tool *Fit and simulation* (see chapter 12.10, page 813).

To copy the results, click the **Copy** button located in the **Results** panel (see figure 996, page 812).




Figure 996 Copy the results as an equivalent circuit

The results of the **Electrochemical circle fit** analysis tool will be copied to the clipboard as a R(RQ) equivalent circuit.

12.10 Fit and simulation

The **Fit and simulation** command provides additional controls that can be used when the command is used to analyze data.



CAUTION

The **Fit and simulation** command is intended to be used on impedance spectroscopy data.

To use the **Fit and simulation** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands and possible plots on which these command can be applied (*see figure 997, page 813*).




Figure 997 Adding a Fit and simulation command to impedance data

The **Fit and simulation** command is added to the procedure editor. Clicking the command shows the properties in the dedicated panel on the right hand side (*see figure 998, page 814*).




Figure 998 The Fit and simulation command is added to the procedure



NOTICE

For more information on the properties of the **Fit and simulation** command, please refer to *Chapter 7.9.2*.

For data analysis purposes, it is possible to use the **Fit and simulation** tool in two different ways:

- Direct fitting or simulation *Direct fitting or simulation* (see chapter 12.10.1, page 814)
- Fitting or simulating using the dedicated editor *Fitting or simulation using the dedicated editor* (see chapter 12.10.2, page 816)

12.10.1 Direct fitting or simulation

It is possible to use the **Fit and simulation** tool to directly fit or simulate the impedance data. Using this method, it is only necessary to specify the equivalent circuit to use, as a CDC string in the **Circuit description** field of the **Properties** panel (see figure 999, page 815).

Properties ➔

Fit and simulation

Command name

Circuit description ✕ Edit

Maximum number of iterations

Maximum change in χ^2 (scaled)

Max iterations without improvement

Fitting style ▼

Use weight factor

Fit or Simulation ▼

Measurement data format ▼

More

Results

Convergence No
Number of iterations 0
 χ^2 -

Figure 999 Typing a CDC string in Properties panel

When the string is validated, by pressing the **[Enter]** key or by unselecting the input field, the fitting or simulation will start immediately, using the default values for all the circuit elements and using all the other properties specified in the **Properties** panel (see figure 1000, page 816).




Figure 1000 The data is fitted or simulated using the specified properties and the default element values

The fitting or the simulation of the data is automatically updated whenever one of the properties provided in the **Properties** panel is modified.



NOTICE

More information on the CDC format of the equivalent circuits can be found in *Chapter 7.9.2.3*.

12.10.2 Fitting or simulation using the dedicated editor

The direct fitting method provided by the **Fit and simulation** command uses the default values for the circuit element. For a more customized analysis of the data, it is possible to use the dedicated editor instead. To use the dedicated editor, click the **Edit** button next to the **Circuit description** field of the **Properties** panel (see figure 1001, page 817).

The screenshot shows a 'Properties' dialog box with a title bar and a close button. The 'Fit and simulation' section contains several input fields and a checkbox. The 'Command name' field is set to 'Fit and simulation'. The 'Circuit description' field is empty, with an 'Edit' button to its right. The 'Maximum number of iterations' field is set to 300. The 'Maximum change in χ^2 (scaled)' field is set to 0,001. The 'Max iterations without improvement' field is set to 50. The 'Fitting style' dropdown menu is set to 'Impedance'. The 'Use weight factor' checkbox is checked. The 'Fit or Simulation' dropdown menu is set to 'Fit'. The 'Measurement data format' dropdown menu is set to 'Impedance'. A 'More' button is located at the bottom right of the 'Fit and simulation' section. The 'Results' section shows 'Convergence No', 'Number of iterations 0', and ' χ^2 -'.

Properties ✕

Fit and simulation

Command name

Circuit description Edit

Maximum number of iterations

Maximum change in χ^2 (scaled)

Max iterations without improvement

Fitting style ▼

Use weight factor

Fit or Simulation ▼

Measurement data format ▼

More

Results

Convergence No

Number of iterations 0

χ^2 -

Figure 1001 Opening the dedicated editor

The dedicated **Equivalent Circuit Editor** will be displayed. This editor provides the means to draw the equivalent circuit using the supported element and to specify the properties of each element (see figure 1002, page 818).




Figure 1002 The equivalent circuit can be specified in the dedicated editor



NOTICE

More information on the use of the **Equivalent Circuit Editor** can be found in *Chapter 7.9.2*.

When all of the properties have been defined, it is possible to run the calculation in two different ways:

- **By closing the editor:** clicking the OK button in the **Equivalent Circuit Editor** window closes the editor and triggers the calculation to run using the specified properties. The fitted or simulated data will be plotted (see figure 1003, page 819).




Figure 1003 The data is fitted or simulated

- **By using the Tools menu in the Equivalent Circuit Editor:** the Tools menu provides the possibility to *Run Fit and simulation* (**[F5]** shortcut key) or *Resume Fit and simulation* (**[F9]** shortcut key), as shown in Figure 1004.




Figure 1004 Running the calculation within the Equivalent Circuit Editor

While the calculation is running, a progress dialog will be shown (see figure 1005, page 821).




Figure 1005 A progress dialog is shown during the calculation

If needed, the calculation can be stopped by clicking the button.

At the end of the calculation, if the command is used to fit the data, the fitted values are shown for each element in the **Equivalent Circuit Editor** (see figure 1006, page 822).




Figure 1006 After the calculation, the calculated values are shown in the editor

If needed, the circuit can be modified, as shown in Figure 1006. At any time, it is possible to restart the calculation, using the **Tools** menu or the **[F5]** shortcut key or resume the calculation, using the same menu or the **[F9]** shortcut key (see figure 1007, page 823).




Figure 1007 Resuming the calculation

Using the *Run Fit and simulation* (**[F5]**) option will trigger the calculation to restart using the *Start* value of each element. Using the *Resume Fit and simulation* (**[F9]**) option will trigger the calculation to continue using the *Fitted* values of each element.



NOTICE

It is possible to assign the *Fitted* value for a circuit element as a new *Start* value by clicking the  button in the **Properties** panel .

12.10.3 Viewing the result

When the calculation is complete, it is possible to view the details by clicking the  button in the **Properties** panel (see figure 1008, page 824).




Figure 1009 Generating the report

The report will display the fitted values for each element as well as the estimated error and the total χ^2 value (see figure 1010, page 825).

Element	Parameter	Value	Estimated Error (%)
R1	R	75,878	0,265
Q1	Y0	1,5069E-05	1,063
	N	0,87619	0,187
R2	R	283	0,227
W1	Y0	0,0016259	0,176
	χ^2	0,00030595	

Figure 1010 The Circuit Report



NOTICE

The data in the **Circuit report** can be exported to ASCII using the **File** menu or can be copied to the Clipboard using the **Edit** menu.

The **Circuit Report** provides the following information:

- **Element:** this is the identification of the circuit element. If a unique name has been specified in the **Equivalent Circuit Editor**, this name will be used instead.
- **Parameter:** indicates the fitted or calculated property of the circuit element.
- **Value:** indicates the fitted or calculated property of the circuit element property.
- **Estimated error:** indicates the estimated error for the element property. This value is indicated in %.



NOTICE

The **Estimated error** is calculated by testing marginal variations of the fitted or calculated value near the convergence. For example, if the best value for a particular resistor is 100 Ohms, the value is increased/decreased until the goodness of fit starts to decrease. If 98 and 102 Ohms produces a very similar goodness of fit, but 97 and 103 Ohms produces a poorer fit, the Error is reported as $2/100 * 100 = 2\%$. Very large error estimates are typically a result of an incorrect model - often one that contains more elements than are represented by the data. If the model contains too many elements, the 'extra' element has no effect on the goodness of fit.

13 Data handling

When data has been measured, it is possible to use the data handling commands provided in NOVA to process and handle the data. To apply data handling command to data acquired in NOVA, it is necessary to add the required command to the measured procedure and apply the function of these commands on the measure data.



NOTICE

Data handling commands can be added to the initial procedure or the procedure after the measurement is finished.

To add a data handling command to a measured procedure, two methods can be used:

- Drag and drop the data handling command in the procedure
- Use the contextual shortcut  button, located in the top right corner of the procedure editor

The functionality of the data handling commands is explained in the previous chapters and will not be detailed again in this chapter. This chapter focuses on the use of these commands on **measured** data. Only the commands that provide controls that are used in a specific way on data are detailed in this chapter.

The following commands are detailed:

- Get item
- Shrink data

13.1 Get item

The **Get item** command provides additional controls that can be used when the command is used to analyze data. To use the **Get item** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands that can be applied on the selected command (*see figure 1011, page 828*).




Figure 1011 Adding a Get item command to the CV staircase command

The **Get item** command is added to the procedure editor. The **Edit links** screen will be shown immediately after the command is added (see figure 1012, page 828).




Figure 1012 The Edit links screen is automatically shown when the Get item command is added to the procedure

Using the method described in *Chapter 10.13*, the links required for the **Get item** command can be edited (see figure 1013, page 829).




Figure 1013 The links required by the **Get item** command can be edited

Clicking the  button closes the **Edit links** screen and returns to the procedure editor. The properties of the **Get item** command can now be edited in the **Properties** panel (see figure 1014, page 829).




Figure 1014 The properties of the **Get item** command can be set



If the links required by the **Get item** command are not properly, an error will be displayed in the procedure editor (see figure 1015, page 830).




Figure 1015 The procedure validation will trigger an error when the links are not set properly



NOTICE

For more information on the properties of the **Get item** command, please refer to *Chapter 7.7.4*.

13.2 Shrink data

The **Shrink data** command provides additional controls that can be used when the command is used to handle data. To use the **Shrink data** command, this command can be added to the procedure as a command, using the drag and drop method, or by using the  button. In the latter case, a pop-out menu is displayed, providing a list of commands that can be applied on the selected command (see figure 1016, page 831).




Figure 1016 Adding a Shrink data command to the CV staircase command

The **Shrink data** command is added to the procedure editor. The **Edit links** screen will be shown immediately after the command is added (see figure 1017, page 831).




Figure 1017 The Edit links screen is automatically shown when the Shrink data command is added to the procedure



Using the method described in *Chapter 10.13*, the links required for the **Shrink data** command can be edited (see *figure 1018*, page 832).




Figure 1018 Setting the links for the **Shrink data** command

Clicking the  button closes the **Edit links** screen and returns to the procedure editor. The properties of the **Shrink data** command can now be edited in the **Properties** panel (see *figure 1019*, page 832).




Figure 1019 The properties of the **Shrink data** command can be set

If the links required by the **Shrink data** command are not properly, a warning will be displayed in the procedure editor (see figure 1015, page 830).




Figure 1020 The procedure validation will trigger a warning when the links are not set properly



NOTICE

For more information on the properties of the **Shrink data** command, please refer to *Chapter 7.7.8*.




Figure 1021 Right-click a plot to create a new overlay



NOTICE

It is possible to right-click a plot from a saved file or from an ongoing measurement.

A new overlay will be created in a new tab and the data from the source dataset will be added to this overlay (see figure 1022, page 836).



Figure 1022 The data is added to the overlay



NOTICE

By default, the Overlay tab starts at number 1 and overlays will be incremented until NOVA is closed.

The information provided in the Overlay tab is distributed in three different panels:

- **Datasets panel:** this panel lists all of the datasets added to the overlay.
- **Overlay panel:** this panel provides a plot of the data from the datasets added to the overlay.
- **Properties panel:** this panel provides the possibility to change the appearance of the overlay.

14.2 Adding data to an overlay

To add a new dataset to an existing overlay, right-click a plot from the new source dataset and select the *Add to overlay X* option, where X is the number of the target overlay (see figure 1023, page 837).




Figure 1023 Adding a dataset to an existing overlay

The new dataset will be added to the target overlay. The information in the **Datasets** panel will be updated, indicating that a new dataset is available in the Overlay tab (see figure 1024, page 838).




Figure 1024 The new dataset is added to the Overlay

The data from the new dataset will be added to the Overlay panel on the right-hand side (see figure 1024, page 838).

14.3 Changing overlay plot properties

The **Properties** panel can be used to change the appearance of data series plotted in the overlay, and to give a name to the **Overlay** plot (see figure 1025, page 839).




Figure 1025 Overlay screen with plot properties panel

The **Properties** panel can be used to define the properties of the data series shown on the overlay plot (see figure 1026, page 839).

Figure 1026 The plot properties panel

The following properties can be edited in the **Properties** panel:

- **Plot title:** the title of the overlay plot.
- **Series title:** the title of the data series. Each data series may be given its own identifying title, this title will automatically appear in the overlay legend.
- **Point style:** defines the point style, color and size used by the plot, using dedicated drop-down lists.
- **Line style:** defines the line style, color and size used by the plot, using dedicated drop-down lists.



NOTICE

Changes to **Properties** such as the **Series title** or **Point style** are relevant only to the appearance of this data in the **Overlay** view. The changes are not propagated back to the source data file.

14.4 Changing overlay plot settings

It is possible to adjust the way the data is plotted in the **Overlay** panel and to change the signal used on the X, Y and Z axis at any time. To change the axes settings, right-click on one of the axes and select a new signal from the popout menu (see figure 1027, page 840).




Figure 1027 Changing the Y axis signal for the overlay



NOTICE

The popout menu shows the available common signals provided by all of the datasets in the overlay.

When a new signal is selected, all of the datasets in the Overlay panel will be replotted, using the new signal (see figure 1028, page 841).




Figure 1028 Changing the X axis signal for the overlay

It is possible to repeat this for each axis, as shown in Figure 1028. Each time a signal is modified, all the data shown in the Overlay panel will be updated (see figure 1029, page 842).




Figure 1029 The data is replotted after the new signal is selected

14.5 Hiding and showing plots

The Legend box, shown in the top right corner of the plot of the **Overlay** panel can be used to show or hide plots. For each plot, a checkbox is available (see figure 1030, page 843).




Figure 1030 Checkboxes are provided in the Legend box of the Overlay panel

Using this control, it is possible to show or hide any of the available plots (see figure 1031, page 844).




Figure 1031 The checkboxes can be used to hide or show plots in the overlay

At any time, it is possible to use this control to show hidden plots or hide shown plots (see figure 1032, page 845).




Figure 1032 The hidden and shown plots can be adjusted at any time

14.6 Remove data from overlay

It is possible to remove a dataset from an overlay by selecting the dataset in the **Datasets** panel and clicking the  button in the top right corner of the panel (see figure 1033, page 846).




Figure 1033 Select the dataset to remove

The dataset will be removed from the overlay and the plot displayed in the **Overlay** panel will be updated (see figure 1034, page 847).




Figure 1034 The selected dataset is removed from the overlay

14.7 Additional Overlay controls

Additional controls are available in the **Overlay** panel, through the dedicated buttons in the top right corner (see figure 1035, page 848).




Figure 1035 Additional controls are available in the top right corner of the Overlay frame



NOTICE

These controls are the same as the controls available normal plots.

The following controls are available:

- **3D view (3D button):** toggles the 3D plot on or off *Toggle the 3D view (see chapter 11.8.2, page 725).*
- **Zoom in (🔍 button):** zooms in on the plot *Zooming options (see chapter 11.8.5, page 728).*
- **Fit view (🔍 button):** fits all the data on the plot *Zooming options (see chapter 11.8.5, page 728).*
- **Zoom out (🔍 button):** zooms out on the plot *Zooming options (see chapter 11.8.5, page 728).*

- **Print plot** (🖨️): prints the plot *Print plot* (see chapter 11.8.6, page 729).
- **Export image** (📄 button): export the data to an image file *Export plot to image file* (see chapter 11.8.7, page 731).




Figure 1037 The procedure scheduler

The procedure scheduler provides two panels:

- **Procedures panel:** a panel that displays all open procedures, all recent procedures and a search box which can be used to search any location defined in the **Library** for a given procedure.
- **New schedule panel:** this panel list all available instrument connected to the computer and the procedure schedule for each instrument.



NOTICE

When a new procedure scheduler is started, all connected instrument are automatically listed in the **New schedule** panel. Instruments that are busy are also listed in the new schedule panel. These instruments will not be able to start a measurement until the current measurement is finished.

Using the controls provided in the procedure scheduler tab, it is possible to carry out the following tasks:

- Remove instruments from the procedure scheduler.
- Add a procedure to an instrument schedule.
- Create a synchronization point.
- Run the procedure scheduler.
- Inspect data from a running procedure.

15.2 Creating a procedure schedule

To create a procedure schedule for one or more instruments in the procedure scheduler, it is necessary to add the required procedures from the **Procedures** panel on the left hand side to the **New schedule** panel on the right hand side. Adding procedures is performed using the *drag and drop* method (see figure 1040, page 853).




Figure 1040 Creating a procedure schedule

It is possible to add procedures from three different sources to a procedure schedule:

- **Open procedures:** these are all the procedures currently open in NOVA.
- **Recent procedures:** these are the five last saved procedures.
- **Search Library:** this search option can be used to search for any procedure in the **Library**.

15.2.1 Open procedures

The procedures listed under *Open procedures* in the **Procedures** panel are all procedures that are currently open for editing in NOVA. Any of these open procedures can be dragged over to the **New schedule** panel in order to add it a procedure schedule for one of the available instrument (see figure 1041, page 854).




Figure 1041 Dragging the open procedure to the schedule

Using the drag and drop method, select an open procedure and drag it to an instrument schedule. Release the mouse button to add the procedure to the schedule. The procedure will be added to the instrument schedule, identified by a white box next to the instrument tile (see figure 1042, page 854).




Figure 1042 The procedure is added to the schedule



NOTICE

A tooltip indicates the expected duration of the procedure in the schedule.

15.2.2 Recent procedures

The procedures listed under *Recent procedures* in the **Procedures** panel are the five last saved procedure. Any of these procedures can be dragged over to the **New schedule** panel in order to add it a procedure schedule for one of the available instrument (see figure 1043, page 855).




Figure 1043 Dragging the recent procedure to the schedule

Using the drag and drop method, select a recent procedure and drag it to an instrument schedule. Release the mouse button to add the procedure to the schedule. The procedure will be added to the instrument schedule, identified by a white box next to the instrument tile (see figure 1044, page 855).




Figure 1044 The procedure is added to the schedule



NOTICE

The procedure boxes in the scheduler are scaled with respect to one another to indicate the relative difference in expected duration.

15.2.3 Search Library

The input field located under *Search Library* can be used to search any procedure in any of the locations specified in the **Library** that contains the terms specified in the field. Typing anything in this field will run a search in the background and the list of procedures will updated while typing (see figure 1045, page 856).




Figure 1045 Searching from procedures in the Library



NOTICE

NOVA searches for procedures that contain the specified search criterium in the **Name** of the procedure or in the **Remarks** field.



NOTICE

More information on the **Library** is available *Chapter 6*.

The procedures listed under *Search Library* in the **Procedures** panel are all the procedures that match the specified search criterium specified in the input field. Any of these procedure can be dragged over to the **New**

schedule panel in order to add it a procedure schedule for one of the available instrument.

Using the drag and drop method, select an searched procedure and drag it to an instrument schedule. Release the mouse button to add the procedure to the schedule. The procedure will be added to the instrument schedule, identified by a white box next to the instrument tile (see figure 1046, page 857).




Figure 1046 The procedure is added to the schedule



NOTICE

Any time a procedure is added to an instrument schedule, the size of the boxes representing these procedures will be adjusted in order to indicate their respective duration.

15.2.4 Remove procedure

To remove a procedure from a schedule, select the procedure to remove and press the **[Delete]** key (see figure 1047, page 858).




Figure 1047 Select the procedure to remove

The procedure will be removed and the schedule will be rearranged. If needed, the size of the procedure boxes will be readjusted (see figure 1048, page 858).




Figure 1048 The procedure is removed from the schedule

15.3 Using synchronization points

It is possible to force instruments involved in a procedure schedule to synchronize their measurements. This can be done by adding a synchronization point. To add a synchronization point to the procedure schedule, click the  button, located in the top right corner of the **New schedule** panel (see figure 1049, page 858).




Figure 1049 Adding a synchronization line

A vertical synchronization line will be added to the procedure schedule, after the last procedure in the schedule (see figure 1050, page 859).




Figure 1050 The synchronization line is added to the schedule

Using the drag and drop method, it is possible to relocate the procedures in the schedule on either side of the synchronization line (see figure 1051, page 859).




Figure 1051 Relocating procedures with respect to the synchronization line



NOTICE

It is also possible to add new procedures the schedule, on either side of the synchronization line.

All procedures located after a synchronization will be synchronized. This means that the all the procedures located on the left hand side of the synchronization line must be finished before starting the next procedure in the schedule. In the example shown in Figure 1052, the Chrono amperometry procedure fast and the Chrono potentiometry fast procedure will be synchronized.




Figure 1052 Creating a synchronized schedule

It is possible to relocate the synchronization line by clicking the line. The synchronization will be highlighted (see figure 1053, page 860).




Figure 1053 Selecting the synchronization line

Using the mouse, it is possible to drag the line to the left or the right to adjust its position in the schedule (see figure 1054, page 860).




Figure 1054 It is possible to relocate the synchronization line

Clicking the **[Delete]** key, when the synchronization line is selected, will delete the line (see figure 1055, page 861).




Figure 1055 Pressing the Delete key will delete the synchronization line



NOTICE

To synchronize measurements it is also possible to use the **Synchronization** command *Synchronization* (see chapter 7.2.10, page 254).

15.4 Naming and saving the schedule

For bookkeeping purposes, it is possible to provide a name to the schedule and save the schedule. To rename the schedule, click the **New schedule** name in the top left corner of the panel (see figure 1056, page 861).




Figure 1056 Renaming the schedule

An input field will be displayed (see figure 1057, page 862).




Figure 1057 A new name can be specified

A name of the procedure schedule can be specified (see figure 1058, page 862).




Figure 1058 Specifying the new name of the schedule

Press the **[Enter]** key or click away from the input field to validate the new name of the procedure schedule. The name will be updated in the top left corner of the panel (see figure 1059, page 862).




Figure 1059 The schedule name is updated

Once a name has been provided, it is possible to save the schedule by selecting the *Save My schedule* option from the **File** menu, or by using the **[CTRL] + [S]** keyboard shortcut (see figure 1060, page 863).




Figure 1060 Saving the schedule

The schedule will be saved in the default Schedules *location* defined in the **Library**. By default, this location is mapped to the **My Documents** \NOVA 2.X folder on the computer. It is also possible to specify the save location of the schedule by using the *Save My schedule As...* option from the **File** menu (see figure 1061, page 863).




Figure 1061 Saving the schedule in a specific location




Figure 1064 All the procedures are validated when the procedure schedule is started

If only warnings are displayed, it is possible to proceed with the schedule. If errors are displayed, it is not possible to proceed with the schedule. The errors will first need to be corrected.

15.5.1 Starting the complete procedure schedule

To start the complete procedure schedule, click the ▶ button in the top right corner of the panel (see figure 1065, page 865).




Figure 1065 Starting the complete schedule

After validating the procedures, the schedule will start on the all the instruments specified in the schedule (see figure 1066, page 866).




Figure 1066 The procedure schedule is started on all instruments
The procedure schedule will be executed as specified on all instruments.

15.5.2 Starting the schedule sequentially

To start the schedule sequentially, select one of the instruments in the schedule and click the  button in the top right corner of the panel (see figure 1067, page 866).




Figure 1067 Starting the procedure schedule for one instrument



NOTICE

The  button is only visible if a single instrument is selected.

The schedule of the selected instrument will start, as shown in (see figure 1068, page 867). It is possible to repeat this for the other instruments in the procedure schedule (see figure 1068, page 867).




Figure 1068 Starting the procedure schedule of the other instrument

The procedure schedule will start for the other instrument (see figure 1069, page 867).




Figure 1069 Both instrument are running

The procedure schedule will be executed as specified on all measuring instruments.

15.5.3 Procedure schedule control

At any point, it is possible to control the procedure schedule. The following actions are possible:

- **Pause one of the instrument:** select one of the measuring instruments and click the  button to pause that instrument. The schedule will be paused for that instrument (see figure 1070, page 867).




Figure 1070 Pausing one of the instruments



- **Stop one of the instruments:** select one of the measuring instruments and click the  button to stop that instrument. The schedule will be stopped for that instrument (see figure 1071, page 868).




Figure 1071 Stopping one of the instruments

- **Pause all the instruments:** click the  button to pause all instruments (see figure 1072, page 868).




Figure 1072 Pausing all of the instruments

- **Stop all the instruments:** click the  button to stop all instruments (see figure 1073, page 868).




Figure 1073 Stopping all of the instruments



NOTICE

Pause measurements can be resumed by clicking the ► or ►| buttons in the top right corner of the panel.

15.6 Inspecting procedures or data

At any time, it is possible to inspect and edit a procedure used in a procedure schedule or to inspect data measured by a procedure used in a procedure schedule.

To inspect or edit a procedure, double click on the white box for this procedure in the procedure scheduler (see figure 1074, page 869).




Figure 1074 Double click the procedure to open or edit the procedure. The procedure will be opened in a new tab (see figure 1075, page 870).




Figure 1075 The procedure is opened in a new tab

The procedure can be edited if required. Modifications that are saved will be automatically carried over to the procedure scheduler.

To inspect data recorded during the procedure scheduler, double click on the white box for a running or finished procedure in the procedure scheduler (see figure 1076, page 870).




Figure 1076 Click a running or finished procedure to inspect the measured data

The running procedure or the measured data will be opened in a new tab (see figure 1077, page 871).




Figure 1077 The data is opened in a new tab

15.7 Schedule zooming

The schedule editor frame has a limited width. If needed, the size of the items in the schedule editor frame can be adjusted with the controls located in the top right corner of the frame (see figure 1078, page 871).




Figure 1078 Zoom controls are provided in the schedule editor

Using this function will either scale the size of the items and the text up or down (between 200 % and 50 % of the original size), as shown in Figure 1079.




Figure 1079 Zooming the schedule editor out

The following zooming controls are available:

- **Zoom out:** decreases the scaling of the items and text shown on screen. The  button or **[CTRL] + [-]** keyboard shortcut can be used to do this.
- **Zoom to 100%:** resets the scaling of the items and text shown on screen to the default size. The  button or **[F4]** keyboard shortcut can be used to do this.
- **Zoom in:** increases the scaling of the items and text shown on screen. The  button or **[CTRL] + [=]** keyboard shortcut can be used to do this.

16 Hardware description

This chapter provides extended information on the Autolab hardware, extension modules and accessories.

Table 14 provides a list of the currently available instruments and modules, supported in NOVA.

Table 14 Overview of the available and supported instruments and modules

Instruments	Modules
PGSTAT100N	ADC10M
PGSTAT128N	BA
PGSTAT302N	Booster10A
PGSTAT302F	Booster20A
PGSTAT101	ECD
PGSTAT204	ECI10M
Multi Autolab M101	ECN
Multi Autolab M204	EQCM
PGSTAT302N MBA	FI20
PGSTAT128N MBA	IME303
	IME663
	MUX
	pX1000
	SCAN250

Table 15 provides a list of the phased out instruments and modules, supported in NOVA.

Table 15 Overview of the phased out and supported instruments and modules

Instruments	Modules
PGSTAT10	ADC750 (replaced by ADC10M)
PGSTAT12	ARRAY (replaced by BA)
PGSTAT20	BIPOT (replaced by BA)
PGSTAT30	FRA2 (replaced by FRA32M)



Instruments	Modules
PGSTAT100	pX (replaced by pX1000)
PGSTAT100 floating option	SCANGEN (replaced by SCAN250)
PGSTAT302	
μAutolab II	
μAutolab III	

Table 16 provides a list of the legacy instruments and modules, unsupported in NOVA.

Table 16 Overview of the legacy and unsupported instruments and modules

Instruments	Modules
PSTAT10	ADC124
μAutolab	DAC124
Multi Autolab PSTAT10	DAC168
	FRA

16.1 General considerations on the use of the Autolab potentiostat/galvanostat systems

This chapter provides general information on the use of the Autolab potentiostat/galvanostat. The information provided in this chapter applies to all instrument, unless otherwise specified. It is highly recommended to review this information before using the Autolab potentiostat/galvanostat.

16.1.1 Electrode connections

The Autolab instruments are supplied with cell cables providing connections for **three** or **four** electrodes, depending on the type of instrument. The electrode connections are provided through 4 mm male banana connectors.

These electrode connections are labeled as follows:

- **Working (or indicator electrode):** WE (red)
- **Sense electrode:** S (red)
- **Reference electrode:** RE (blue)
- **Counter electrode:** CE (black)

An additional green ground connector is provided for connections to a Faraday cage.





NOTICE

The μ Autolab type II, μ Autolab type III, PGSTAT10 and PGSTAT10 are not fitted with the Sense electrode.

16.1.1.1 Three electrode connections

Instruments with three electrode connectors can be connected to the electrochemical cell in two different ways:


- Two electrode mode:** in this mode, the counter electrode (CE) and reference electrode (RE) are connected together to one electrode while the working electrode (WE) is connected to the other electrode. The current is measured between the CE and the WE and the potential difference is measured between the RE and the WE. This mode is commonly used for the characterization of energy storage and conversion devices like batteries, fuel cells, solar cells and supercapacitors.



NOTICE

For high current applications it is highly recommended to separately connect the RE and CE to the same electrode. Furthermore, it is recommended to place the RE as close as possible to the electrodes in the cell. This will reduce ohmic losses coming from the connections.

- Three electrode mode:** in this mode, the counter electrode (CE) and reference electrode (RE) are connected to a counter and reference electrode, respectively. The working electrode (WE) is connected to the working electrode. The current is measured between the CE and the WE and the potential difference is measured between the RE and the WE. This mode is commonly used for the characterization most electrochemical cells in which a separate reference electrode is used.





NOTICE

It is common practice to place the reference electrode as close as possible to the working electrode to reduce the uncompensated resistance and reduce the ohmic losses arising from this resistance. This can be achieved by physically placing the reference electrode close to the working electrode or by using a *Luggin-Haber* capillary.

16.1.1.2 Four electrode connections

Instruments with four electrode connectors can be connected to the electrochemical cell in three different ways:

- Two electrode mode:** in this mode, the counter electrode (CE) and reference electrode (RE) are connected together to one electrode while the working electrode (WE) and sense electrode (S) are connected to the other electrode. The current is measured between the CE and the WE and the potential difference is measured between the RE and the S. This mode is commonly used for the characterization of energy storage and conversion devices like batteries, fuel cells, solar cells and supercapacitors.



NOTICE

For high current applications it is highly recommended to separately connect the RE and CE to the same electrode and to do the same with the WE and S on the other electrode. Furthermore, it is recommended to place the RE and S as close as possible to the electrodes in the cell. This will reduce ohmic losses coming from the connections.


- Three electrode mode:** in this mode, the counter electrode (CE) and reference electrode (RE) are connected to a counter and reference electrode, respectively. The working electrode (WE) and sense electrode (S) are connected to the working electrode. The current is measured between the CE and the WE and the potential difference is measured between the RE and the S. This mode is commonly used for the characterization most electrochemical cells in which a separate reference electrode is used.



NOTICE

It is common practice to place the reference electrode as close as possible to the working electrode to reduce the uncompensated resistance and reduce the ohmic losses arising from this resistance. This can be achieved by physically placing the reference electrode close to the working electrode or by using a *Luggin-Haber* capillary.

- Four electrode mode:** in this mode, the counter electrode (CE) and reference electrode (RE) are connected to a counter and reference electrode, respectively on side of the electrochemical cell. The working electrode (WE) and sense electrode (S) are connected to a second set of working electrode and reference electrode on the other side of the electrochemical cell. Both sides of the cell are separated by a membrane or by using non miscible solvent. The current is measured between the CE and the WE and the potential difference is measured between the RE and the S. This mode is commonly used for the characterization of the liquid-liquid interface.



16.1.2 Operating principles of the Autolab PGSTAT

The Autolab instrument combined with the software is a computer-controlled electrochemical measurement system. It consists of a data-acquisition system and a potentiostat/galvanostat. The basic working principle is schematically represented in *Figure 1080*.




Figure 1080 Schematic representation of the Autolab potentiostat/galvanostat

The Autolab system is fitted with the following common **digital** control components:

- **USB interface:** the interface between the Autolab and the host computer.
- **Embedded PC with real-time operating system:** a dedicated controller embedded into the instrument, which is responsible for timing and interfacing the host application and the instrument controls.
- **Decoder and DIO:** a data decoder and digital input/output interface.

The digital components are interfaced through the Autolab modules to the **analog** potentiostat/galvanostat circuit. The latter consists of the following components:

- **Summation point (Σ):** a circuit used to add the control signals required to generate the waveform used in electrochemical measurements.
- **Control amplifier (CA):** a circuit used to amplify the output of the summation point.
- **Voltage follower (VF):** a circuit used to measure the potential.
- **Current follower (CF):** a circuit used to measure the current.

The arrangement of these analog circuits with respect to the electrochemical cell are represented in *Figure 1081* for a four electrode Autolab system and in *Figure 1082* for a three electrode Autolab system.




Figure 1081 Schematic representation of the analog circuits of the Autolab in a four electrode system




Figure 1082 Schematic representation of the analog circuits of the Autolab in a three electrode system

The summation point (Σ) is an **adder** circuit that feeds the input of the control amplifier (**CA**). Each of the inputs of the summation point is divided by a hardware-defined value (see figure 1083, page 880).




Figure 1083 Schematic representation of the summation point of the Autolab

It is connected to the output of the several key modules of the Autolab:

- **DAC164 (or on-board DAC):** the digital-to-analog converters of the Autolab. Depending on the type of instrument, the following DAC inputs are available:
 - **Offset DAC:** used to generate an offset. This signal is divided by **2**.
 - **Scanning DAC:** used to generate steps and scans. This signal is divided by **2**.
 - **AC voltammetry DAC:** used for AC voltammetry only. This signal is divided by **10**.
- **FRA32M DSG or FRA2 DSG:** the digital waveform generator of the optional **FRA32M** or **FRA2** module *FRA32M module (see chapter 16.3.2.13, page 1112)*. This signal is divided by **10**.
- **SCAN250 or SCANGEN output:** the analog scan output of the optional linear **SCAN250** or **SCANGEN** module *SCAN250 module (see chapter 16.3.2.19, page 1169)*. This signal is divided by **1**.
- **E_{in}:** the external input provided through the monitor cable. This signal is divided by **1**.
- **PGSTAT feedback:** the feedback from the voltage follower (**VF**), in potentiostatic mode or the feedback from the current follower (**CF**), in galvanostatic mode.
- **iR compensation feedback:** the feedback from the iR compensation circuit, when in use in potentiostatic mode.



CAUTION

Some of the summation point inputs are not available on all Autolab instruments.

The control amplifier provides the output voltage on the counter electrode (CE) with respect to the working electrode (WE) required to keep the potential difference between the reference electrode (RE) and the sense (S) or the potential difference between the reference electrode (RE) and the working electrode (WE) at the user defined value, in potentiostatic mode, or the user required current between the counter electrode (CE) and the working electrode (WE) in galvanostatic mode.

The output of the control amplifier can be manually or remotely disconnected from the electrochemical cell through a cell ON/OFF switch. The voltage follower (**VF**) is used to measure the potential difference between the reference electrode and the sense and the current follower (**CF**). The current follower has several current ranges providing different current-to-voltage conversion factors.

The output of the current follower (**CF**) and the voltage follower (**VF**) are fed back to the analog-to-digital converter modules of the Autolab:

- **ADC164 (or on-board ADC):** general purpose analog-to-digital converter of the Autolab instrument.
- **FRA32M ADCs or FRA2 ADCs:** two synchronized analog-to-digital converters located on the *optional FRA32M* or **FRA2** module used for impedance spectroscopy measurements *FRA32M module (see chapter 16.3.2.13, page 1112)*.
- **ADC10M or ADC750:** two synchronized analog-to-digital converters located on the *optional ADC10M* or **ADC750** module for ultra-fast sampling *ADC10M module (see chapter 16.3.2.1, page 998)*.
- **FI20 (or on-board integrator):** an *optional* filter and integrator module (**FI20**) or on-board integrator which can be used to convert the current to charge and filter the current signal *FI20 module (see chapter 16.3.2.11, page 1083)*.

Furthermore, the output of the voltage follower (VF) or the current follower (**CF**) is fed back to the summation point to close the feedback loop in potentiostatic or galvanostatic mode, respectively.

The **ADC164 (or on-board ADC)** provides the possibility of measuring analog signals. The input sensitivity is software-controlled, with ranges of ± 10 V (gain 1), ± 1 V (gain 10) and ± 0.1 V (gain 100). The resolution of the measurement is 1 in 65536 (16 bits). Analog signals can be measured with a rate of up to 60 kHz. The ADC164 is used to measure the output of the



voltage follower (**VF**) and current follower (**CF**) of the potentiostat/galvanostat.

The **DAC164 (or on-board)** generates analog output signals. The output is software-controlled within a range of ± 10 V. The resolution of the DAC164 is 1 in 65536 ($300 \mu\text{V}$). In the Autolab PGSTAT two channels of the DAC (*Scanning DAC* and *Offset DAC*) are used to control the analog input signal of the potentiostat/galvanostat. The μ Autolab type II and μ Autolab type III only uses a *Scanning DAC* to control the analog input. The values of the DACs are added up in the potentiostat and divided by 2. This results in an output of ± 10 V with a resolution of $150 \mu\text{V}$.

In practice this means that the potential range available with the Autolab PGSTAT during an electrochemical experiment is ± 5 V with respect to the offset potential generated by the *Offset DAC*. The available potential range is therefore -10 V to 10 V with the Autolab PGSTAT and -5 V to 5 V with the μ Autolab type II and μ Autolab type III.

The AC voltammetry DAC, if present, is hardwired to the summation point and it is divided by 10. This input is used for measurements involving a small amplitude modulation (like AC voltammetry).

16.1.2.1 Event timing

The embedded controller of the Autolab is equipped with a 1 MHz timer that is used by the software to control the timing of events during measurements. The shortest interval time on the embedded controller is $1 \mu\text{s}$. When a procedure is started in NOVA, the procedure is first uploaded from the host PC to the embedded PC, through the USB connection. The measurement can then be started.

Depending on the type of command that NOVA encounters during the measurements two timing protocols are used:

1. **Real-time control:** all measurement commands in NOVA are timed using the embedded processor timer. Whenever NOVA encounters a measurement command, it will be executed using the timing provided by the embedded computer of the Autolab. If several measurement commands are located in sequence, the sequence is executed without interruption. This ensures that the measurement commands in the sequence are executed with the smallest possible time gap. The actual time difference between two consecutive commands depends on the hardware changes required during the transition between the two commands. Switching current ranges or using the cell switch are time consuming steps since they involve mechanical relays which require a fixed settling time. Taking into account these hardware defined interval times, the effective time gap between two consecutive measurement commands will always be ≤ 10 ms.

2. **Host computer control:** all the other commands in NOVA are timed using the timer of the host computer. Since the host computer is also involved in other Windows activity, accurate timing of events cannot be guaranteed and the effective interval time between two consecutive host commands will depend entirely on the amount of activity on the host computer. Depending on the command sequence, the time gap can be as short as ~ 1 s (transition between host command to measurement command) or several seconds (transition between measurement command and host command). Transfer of large amounts of measured data points is particularly time consuming.



NOTICE

To reduce the time gap between commands timed by the host computer it is recommended to close all unnecessary Windows applications while using NOVA.

16.1.2.2 Consequences of the digital base of the Autolab

The digital nature of the instrument control has consequences for the measurements. The consequences for the different techniques are the following:

- The minimum potential step or pulse in all techniques is 150 μV (16 Bit DAC).
- All potential steps are rounded up or down to the nearest possible multiple of 150 μV .
- In staircase potential (or current) scans, the interval time, Δt , or time between two consecutive steps is given by:

$$\Delta t = \frac{E_{\text{step}}}{\bar{v}}$$

Where E_{step} is the potential step (or current step) and \bar{v} is the scan rate.

The response of the electrochemical cell is recorded digitally. Therefore the resolution of the measurements is also limited. The actual resolution depends on the technique and on the amplitude of the signal. Since the analog-to-digital converter is equipped with a software programmable amplifier, the absolute resolution depends on the gain of the amplifier. The gains used are 1, 10 and 100 times the input signal.

NOVA automatically selects the best possible gain during a measurement. Gain 10 and 100 are used when the signal is small enough.

Depending on the gain, the resolution for potential, current and external analog signal are listed in *Table 17*.

Table 17 The resolution of the measurable signals



Signal	Potential	Current ([CR] active current range)	External
Gain 1	$\frac{20V}{2^{16} \cdot 1}$	$\frac{20[CR]}{2^{16} \cdot 1}$	$\frac{20V}{2^{16} \cdot 1}$
Gain 10	$\frac{20V}{2^{16} \cdot 10}$	$\frac{20[CR]}{2^{16} \cdot 10}$	$\frac{20V}{2^{16} \cdot 10}$
Gain 100	$\frac{20V}{2^{16} \cdot 100}$	$\frac{20[CR]}{2^{16} \cdot 100}$	$\frac{20V}{2^{16} \cdot 100}$

The effect of the limited resolution can be seen, for instance when low currents are measured at a high current range. In such cases a lower current range has to be applied, if possible. When automatic current ranging is used, the most suitable current range is selected automatically.

Care must be taken when using this option in the following situations:

- Square wave voltammetry measurements at high frequency.
- Cyclic and linear sweep voltammetry measurements at high scan rates.

Switching of the current range takes about 0.5 ms to 2 ms. Therefore an erroneous point can be measured when the current range is switched. Most of the time, this error can be corrected by smoothing the plot afterwards.

16.1.2.3 Bandwidth settings

The control amplifier of the Autolab is equipped with three different bandwidth settings:

- High stability
- High speed
- Ultra-high speed



NOTICE

The Ultra-high speed mode is not available for the PGSTAT302F, the PGSTAT10, the PGSTAT20, the μ Autolab type II and μ Autolab type III.

The bandwidth setting can be specified using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236).

This property defines the bandwidth of the instrument control amplifier. The three settings provided by the Autolab control command can be used to reach the required bandwidth while maintaining stability of the potentiostatic or galvanostatic control loop. The normal mode of operation is **High stability** (see figure 1084, page 885).




Figure 1084 The instrument bandwidth is defined in the Autolab control command

It is also possible to define the bandwidth using the **Autolab display** panel (see figure 1085, page 886).




Figure 1085 The instrument bandwidth can be specified directly from the Autolab display panel



NOTICE

The **High stability** mode is the power-up default of the instrument

The **High stability** mode is suitable for electrochemical measurements carried out at low frequencies or low scan rates. In this mode of operation, the instrument control amplifier will use the slowest possible feedback. This usually ensures a stable control loop and low noise levels on the measured potential and current. In this mode, the bandwidth of the control amplifier is limited to 10 kHz.

Faster measurements may require a higher bandwidth setting. In **High speed** mode, the control amplifier bandwidth is extended to 125 kHz

while in **Ultra-high speed** mode, the control amplifier bandwidth is extended to 1.25 MHz (for the PGSTAT302N) and to 500 kHz (for the PGSTAT128N and PGSTAT100N).

With these settings, the risk of oscillations is higher than in **High stability** mode. This is especially the case with electrochemical cells exhibiting a high capacitance. There is a significant oscillation risk in these modes of operation. Additionally, the noise in the measured potential and current signals will be higher than in **High stability** mode.



NOTICE

The specified or active bandwidth setting is not changed by any measurement command **except** the FRA measurement command which automatically selects the most suitable bandwidth setting in function of the applied frequency.



CAUTION

The higher the bandwidth, the higher the noise and probability of oscillation. When working with a high bandwidth setting (**High speed** or **Ultra high speed**), it is necessary to pay attention to adequate shielding of the cell and electrode connectors. The use of a Faraday cage is recommended in these cases.

16.1.2.3.1 Input impedance and stability

The voltage follower (**VF**) input contains a small capacitive load. If the capacitive part of the impedance between CE and RE is comparatively large, phase shifts will occur which can lead to instability problems when working in potentiostatic mode. If the impedance between the CE and the RE cannot be changed and oscillations are observed, it is recommended to select the **High stability** mode to increase the system stability.

In general, the use of **High stability** leads to a more stable control loop, compared to **High speed** or **Ultra-high speed** and a significantly lower bandwidth.

To make use of the full potentiostat bandwidth (**Ultra-high speed** mode), the impedance between CE and RE has to be lower than 35 k Ω . This value is derived by testing. In galvanostat mode, this large impedance between CE and RE, will usually not lead to stability problems, because of the current feedback regulation.



16.1.2.3.2 Galvanostat and iR compensation bandwidth limitations

For galvanostatic measurements on low current ranges, the bandwidth limiting factor becomes the current follower (**CF**) rather than the control amplifier. The same applies to potentiostatic measurements with the iR compensation circuit on. In both cases, the bandwidth of this circuit directly determines the maximum bandwidth of the control loop.

When the iR compensation circuit is not used in potentiostatic mode, the bandwidth limitation of the current follower (**CF**) does not directly influence the control loop bandwidth but it influences the measurement of the current signal.

For stability reasons, the following guidelines are provided when working in galvanostatic mode or in potentiostatic mode with the iR compensation circuit on:

- The use of the **High speed** mode is only recommended for current range of 10 μA and higher.
- The use of the **Ultra-high speed** mode is only recommended for current ranges of 1 mA and higher.

A general indication of the maximum available bandwidth for galvanostatic measurements and measurements with the iR compensation circuit on can be found in:

- *Table 18* for the N Series Autolab instruments.
- *Table 19* for the 7 Series Autolab instruments.
- *Table 20* for the PGSTAT101, PGSTAT204, M101 and M204.
- *Table 21* for the $\mu\text{Autolab}$ type II and $\mu\text{Autolab}$ type III, PGSTAT10 and PGSTAT20.

Table 18 Bandwidth overview of the N Series Autolab instruments

Instrument	PGSTAT128N, PGSTAT100N		PGSTAT302N	
	GSTAT	iR compensation	GSTAT	iR compensation
1 A - 1 mA	> 500 kHz	> 500 kHz	> 1.25 MHz	> 1.25 MHz
100 μA	125 kHz	500 kHz	125 kHz	1 MHz
10 μA	100 kHz	100 kHz	100 kHz	100 kHz
1 μA	10 kHz	10 kHz	10 kHz	10 kHz
100 nA	1 kHz	1 kHz	1 kHz	1 kHz
10 nA	100 Hz	100 Hz	100 Hz	100 Hz

Table 19 Bandwidth overview of the 7 Series Autolab instruments

Instrument	PGSTAT12, PGSTAT100		PGSTAT30, PGSTAT302	
	GSTAT	iR compensation	GSTAT	iR compensation
1 A - 1 mA	> 500 kHz	> 500 kHz	> 1.25 MHz	> 1.25 MHz
100 μ A	125 kHz	500 kHz	125 kHz	1 MHz
10 μ A	100 kHz	100 kHz	100 kHz	100 kHz
1 μ A	10 kHz	10 kHz	10 kHz	10 kHz
100 nA	1 kHz	1 kHz	1 kHz	1 kHz
10 nA	100 Hz	100 Hz	100 Hz	100 Hz

Table 20 Bandwidth overview of the Autolab PGSTAT101, M101, PGSTAT204 and M204

Instrument	PGSTAT101, M101		PGSTAT204, M204	
	GSTAT	iR compensation	GSTAT	iR compensation
100 mA			> 1 MHz	> 1 MHz
10 mA - 1 mA	> 1 MHz	> 1 MHz	> 1 MHz	> 1 MHz
100 μ A	1 MHz	1 MHz	> 1 MHz	> 1 MHz
10 μ A	10 kHz	75 kHz	10 kHz	50 kHz
1 μ A	10 kHz	20 kHz	10 kHz	50 kHz
100 nA	400 Hz	4 kHz	500 Hz	500 Hz
10 nA	400 Hz	400 Hz	500 Hz	500 Hz

Table 21 Bandwidth overview of the Autolab PGSTAT10, PGSTAT20, μ Autolab type II and μ Autolab type III

Instrument	PGSTAT10, μ Autolab type II, μ Autolab type III		PGSTAT20	
	GSTAT	iR compensation	GSTAT	iR compensation
1 A - 10 mA			> 1 MHz	> 1 MHz
10 mA - 1 mA	> 1 MHz		> 1 MHz	> 1 MHz



Instrument	PGSTAT10, μ Autolab type II, μ Autolab type III		PGSTAT20	
	GSTAT	iR compensation	GSTAT	iR compensation
100 μ A	500 kHz		500 kHz	500 kHz
10 μ A	50 kHz		50 kHz	50 kHz
1 μ A	5 kHz		5 kHz	5 kHz
100 nA	400 Hz		400 Hz	400 Hz
10 nA	20 Hz			

When a bandwidth conflict is detected in NOVA, a warning is provided in order to provide information on this conflict (see figure 1086, page 890).




Figure 1086 Warnings are provided when bandwidth limitations are encountered



NOTICE

It is possible to ignore this warning and to proceed with the measurement. This can however lead to instabilities or invalid measurements. It is therefore not recommended to adjust the procedure properties.

16.1.2.3.3 Oscillation detection and protection

The N Series Autolab instruments and the 7 Series Autolab instruments are fitted with a detector for large-amplitude oscillation. The detector will spot any signal swing that causes the control amplifier to produce both a positive and a negative voltage overload within $\sim 200 \mu\text{s}$. Thus, large oscillations at frequencies $> 2.5 \text{ kHz}$ will be detected.

Upon oscillation, the **OSC** indicator on the PGSTAT front panel will be activated (see figure 1092, page 903). The **Vovl** warning will also be shown in the Autolab display.

When an oscillation is detected, the cell will be automatically disconnected for safety reasons and the **OSC** indicator will blink on the Autolab front panel. The **Autolab display** panel will display that the cell is set to *Manually off* and the **Oscillation** warning indicator will be lit in the **Warnings** sub-panel (see figure 1087, page 891).




Figure 1087 The oscillation status is reported in the Autolab display panel

The **Cell ON** button (item 1 in Figure 1092) or the **CELL ENABLE** button (item 1 in Figure 1119), located on the right-hand side of the instrument front panel will blink.

The cell may be switched on again by pressing the **Cell ON/CELL ENABLE** button. If oscillation resumes, the cell will be switched off as soon as



the button is released. Holding the button pressed in, provides an opportunity to observe the system during oscillation. Some cells that cause ringing when switching the cell on or changing the current range can falsely trigger the oscillation detector. If this happens, the oscillation protection may be switched off in the software in order to prevent an accidental disconnection of the cell.

The oscillation protection feature can be enabled or disabled in the software, using the **Autolab control** command (see figure 1088, page 892).




Figure 1088 The oscillation protection can be enabled or disabled in the Autolab control command



CAUTION

It is **not** recommended to switch off the oscillation protection circuit.

16.1.2.4 Current range linearity

Each current range on the instrument is characterized by a specific linearity limit and this specification determines the maximum current that can be applied in galvanostatic mode. This limit also determines the maximum current that can be measured in potentiostatic mode in a given current range.

The procedure validation provides an **error** message when the specified current exceeds the linearity limit in a galvanostatic experiment (see figure 1089, page 893).




Figure 1089 An error is displayed when the applied current exceeds the linearity limit of the active current range



NOTICE

It is **not** possible to start the measurement when an error is shown.

Whenever this limit is exceeded during a potentiostat measurement, a current overload **warning** message is shown after the measurement finishes (see figure 1090, page 893).




Figure 1090 A warning is displayed if the linearity limit is exceeded during a measurement

An overview of the current range linearity can be found in:

- Table 22 for the N Series Autolab instrument.
- Table 23 for the 7 Series Autolab instruments.
- Table 24 for the PGSTAT101, PGSTAT204, M101 and M204.
- Table 25 for the μAutolab type II and μAutolab type III, PGSTAT10 and PGSTAT20.

Table 22 Linearity limit for the N Series Autolab instruments

Current range	PGSTAT128N	PGSTAT302N	PGSTAT100N
1 A	0.8	2	



Current range	PGSTAT128N	PGSTAT302N	PGSTAT100N
100 mA	3	3	2.5
10 mA - 1 mA	3	3	3
100 μ A - 1 μ A	3	3	3
100 nA - 10 nA	3	3	3

Table 23 Linearity limit for the 7 Series Autolab instruments

Current range	PGSTAT12	PGSTAT30/30 2	PGSTAT100
1 A		1/2	
100 mA	2.5	3	2.5
100 mA - 1 mA	3	3	3
100 μ A - 1 μ A	3	3	3
100 nA - 10 nA	3	3	3

Table 24 Linearity limit for the Autolab PGSTAT101, M101, PGSTAT204 and M204

Current range	PGSTAT101, M101	PGSTAT204, M204
100 mA		4
10 mA	10	7
1 mA	7	7
100 μ A - 1 μ A	7	7
100 nA - 10 nA	7	7

Table 25 Linearity limit for the Autolab PGSTAT10, PGSTAT20, μ Autolab type II and μ Autolab type III

Current range	PGSTAT10	PGSTAT20	μ Autolab type II, μ Autolab type III
1 A		1	
100 mA		4	
10 mA	5	4	5
1 mA	4	4	4
100 μ A - 1 μ A	4	4	4

Current range	PGSTAT10	PGSTAT20	μ Autolab type II, μ Autolab type III
100 nA	4	4	4
10 nA			4

The values reported in the tables indicate how many times the current range value can be applied or measured for each current range. For example, for in the 1 mA current range, the linearity limit is ± 3 mA.

16.1.2.5 Maximum input voltage

The differential electrometer input contains an input protection circuitry that becomes active after crossing the ± 10 V limit. This is implemented to avoid electrometer damage. Please note that the **Vovl** indicator, on the front panel of the instrument, will not light up for this type of voltage overload.

The measured voltage will be cutoff at an absolute value of ± 10.00 V.

Depending on the cell properties, galvanostatic control of the cell could lead to a potential difference between the reference electrode (RE) and the sense electrode (S) larger than 10 V. This situation will trigger the cut-off of the measured voltage to prevent overloading the differential amplifier.

In this case it is possible to connect the Autolab **Voltage Multiplier** to extend the measurable range of the differential amplifier to ± 100 V.

16.1.2.6 Active cells

Energy storage and conversion devices like batteries and fuel cells are capable of delivering power to the Autolab potentiostat/galvanostat. This is allowed only to a maximum active cell power, P_{MAX} . The values for P_{MAX} depend on the instrument type and are reported in *Table 26*.

Table 26 Maximum power rating for the different Autolab instruments

Instrument	Maximum power, P_{MAX} (W)
PGSTAT128N	8
PGSTAT302N	20
PGSTAT100N	2.5
PGSTAT12	2.5
PGSTAT302	20
PGSTAT30	10



Instrument	Maximum power, P _{MAX} (W)
PGSTAT100	2.5
PGSTAT101/M101	1
PGSTAT204/M204	4
PGSTAT10	1
PGSTAT20	10
μAutolab type II, μAutolab III	0.5
Booster10A	100
Booster20A	200

This means that cells showing an absolute voltage, $|V_{\text{Cell}}|$, of less than 10 V between the working electrode (WE) and counter electrode (CE) are intrinsically safe. They may drive the PGSTAT control amplifier into current limit but will not overload the amplifier. On the other hand, cells that have an absolute voltage higher than 10 V between WE and CE may only deliver a maximum current, i_{MAX} given by:

$$i_{\text{MAX}} = \frac{P_{\text{MAX}}}{|V_{\text{MAX}}|}$$



NOTICE

Instruments that can be connected to the optional **Booster10A** or **Booster20A** can work with active cell power values of 100 W and 200 W, respectively. More information on the Booster10A and Booster 20A can be found in *Chapter 16.3.2.5* and *Chapter 16.3.2.6*.

16.1.2.7 Grounded cells

The measurement circuitry of the Autolab is internally connected to protective earth (P.E.). This can be an obstacle when measurement is desired of a cell that is itself in contact with P.E.. In such a case, undefined currents will flow through the loop that is formed when the electrode connections from the PGSTAT are linked to the cell and measurements will not be possible.

Please note that not only a short circuit or a resistance can make a connection to earth, but also a capacitance is capable of providing a conductive path (for AC signals). The earth connection between the cell and P.E. should always be broken.

If there is no possibility of doing this, please contact Metrohm Autolab for a custom solution, if available.

16.1.3 Environmental conditions

The PGSTAT may be used at temperatures of 0 to 40 degrees Celsius. The instrument is calibrated at 25 degrees Celsius and will show minimum errors at that temperature. The ventilation holes on the bottom plate and on the rear panel may never be obstructed, nor should the instrument be placed in direct sunlight or near other sources of heat.

16.1.3.1 Temperature overload

As a safety precaution, the PGSTAT is equipped with a circuit that monitors the temperature of the internal power electronics. A temperature overload will be displayed as a blinking indicator in the manual cell switch, with the cell automatically turned off. You will not be able to turn the cell back on until the temperature inside the instrument has fallen to an acceptable level. It can then be switched on again by pressing the manual cell switch button on the front panel.

During normal operation the temperature should never become extremely high and no temperature overload will occur. If this does happen, the origin of the temperature overload should be identified:

1. Is the room temperature unusually high?
2. Was the PGSTAT oscillating?
3. Is the voltage selector for mains power set to the right value?
4. Is the fan turning and are all the ventilation holes unobstructed?
5. Was the cell delivering a considerable amount of power to the PGSTAT?
6. Are the WE and CE cables shorted in PSTAT mode?



NOTICE

If a temperature overload takes place repeatedly, for no obvious reason, Metrohm Autolab recommends having the instrument checked by their service department.

16.1.4 Noise considerations

When measuring low level currents, some precautions should be taken in order to minimize noise. The personal computer must be placed as far away as possible from the electrochemical cell and the cell cables. The cell cables should not cross other electrical cables. Other equipment with power supplies can also cause noise. For instance, the interface for mercury electrodes IME should also be placed with some care. If possible place the computer between the PGSTAT and other equipments. Avoid using unshielded extension cables to the electrodes. The use of a Faraday cage is also advised.



If the cell system has a ground connector, it can be connected to the analog ground connector at the front of the PGSTAT. If a Faraday cage is used, it should be connected to this ground connector. Some experiments concerning optimization of the signal-to-noise ratio can readily indicate whether or not a configuration is satisfactory.

When investigating the sources of noise, it is recommended to consider following items:

1. Problems with the reference electrode
2. Problems with unshielded cables
3. Faraday cage
4. Grounding of the instrument
5. Magnetic stirrer
6. Position of the cell with respect to the instrument and accessories
7. Measurements in a glove box



NOTICE

The **Check cell** tool, provided in the **Instrument control panel**, can be used to evaluate the noise levels. More information can be found in *Chapter 5.2.2.4*.

16.1.4.1 Problems with reference electrodes

If the reference electrode is not filled properly with electrolyte solution or when it has, for other reasons, a very high impedance, it may introduce noise in electrochemical measurement. In most cases the applied potential is not the same as the measured potential. Refer to the user manual provided by the reference electrode supplier for more information on the proper care of your reference electrode.

16.1.4.2 Problems with unshielded cables

It is not advisable to use unshielded electrode cables. Make the connections to the electrodes as close as possible to the electrode itself. Avoid the use of unshielded extension cables to the electrodes.

16.1.4.3 Faraday cage

The use of a Faraday cage is always recommended. It protects the cell from external noise interference. Connect the cage to the green ground connector embedded in the cell of the Autolab.

16.1.4.4 Grounding of the instrument

Not properly grounding of the Autolab and computer will decrease the signal-to-noise ratio. Always use a grounded power outlet and grounded power cables. Be sure to connect the Autolab and computer to the same power ground. This means they should be connected to the same power outlet.

16.1.4.5 Magnetic stirrer

In some cases a magnetic stirrer can cause noise problems. Try the measurements with the stirrer on and off and monitor the current. If the stirrer causes a lot of noise please try to find another way of stirring.

16.1.4.6 Optimizing the position of the instrument

The signal-to-noise ratio can often be improved by changing the positions of the cell, computer and ancillary equipment relative to the Autolab. In general, the electrochemical cell should be placed as far as possible from the computer and other devices, without extending the cell cables with unshielded cables. If the noise level remains too high, a Faraday cage may be necessary.

16.1.4.7 Measurements in a glove box

When the cell needs to be placed into a glove box, it is highly recommended to use **isolated** feedthrough that allows the Autolab cell cables to be connected to the cell inside the glove box. If necessary, the cell cables of the Autolab can be fitted with male BNC connectors rather than 4 mm banana connectors. This allows using BNC feedthroughs. Contact your Autolab distributor for more information about this modification.



CAUTION

The shielding of the reference electrode (RE) and sense electrode (S) cable on the Autolab is driven (or guarded). Use **isolated** cable feedthroughs for these cables in order to extend the driven shield inside the glove box. The shield of these cables must not be connected to the ground of the glove box.

16.1.5 Cleaning and inspection

It is recommended to clean the Autolab instrument and the accessories on a regular basis. This can be done with a damp cloth, optionally using a mild detergent. Never use an excessive amount of water; it may never enter into the instrument. As a precaution, disconnect Autolab from the mains when cleaning it. Also perform an inspection of the instrument and all of the connecting cables. If you find any cables with damaged insula-

16.2.1.1 Scope of delivery

The N Series Autolab systems are supplied with the following items:

- Autolab potentiostat/galvanostat
- ADC164 (installed)
- DAC164 (installed)
- Cell cable (WE/CE/GND)
- Differential amplifier (RE/S)
- Monitor cable
- Power cable
- BNC cable (50 cm)
- USB cable
- Set of four alligator clips
- Autolab dummy cell

16.2.1.2 Instrument power-up state

The power-up state of the instrument is hardware defined. The following settings are automatically selected whenever the instrument is powered on or whenever the connection to the instrument is reset by the software.

- Cell: off
- Mode: potentiostatic
- Control bandwidth: high stability
- iR compensation: off
- Current range: 10 mA
- Optional modules: off
- DIO ports: write mode, low state
- Summation point inputs: off
- Oscillation protection: on

16.2.1.3 N Series Autolab front panel

The front panel of the N Series Autolab provides a number of connections, controls and indicators (*see figure 1091, page 902*).




Figure 1091 Overview of the front panel of the N Series Autolab

- | | |
|--|--|
| <p>1 On/Off button
For switching the Autolab on or off.</p> | <p>2 ADC164 →1
Analog input for recording external signals (ADC164 →1).</p> |
| <p>3 ADC164 →2
Analog input for recording external signals (ADC164 →2).</p> | <p>4 DAC164 ←2
Analog output for controlling external signals (DAC164 ←2).</p> |
| <p>5 Monitor cable connector ⇄
For connecting the monitor cable.</p> | <p>6 CE/WE connector
For connecting the Autolab cell cable, providing connections to the counter electrode (CE), working electrode (WE) and ground.</p> |
| <p>7 RE/S connector
For connecting the Autolab differential amplifier, providing connections to the reference electrode (RE) and sense electrode (S).</p> | <p>8 Ground connector
Additional ground connector for connecting external devices to the Autolab ground.</p> |
| <p>9 Cell ON button
For enabling and disabling the cell.</p> | <p>10 Display
Display indicating real-time information on the measured current and potential and instrumental settings.</p> |
| <p>11 DAC164 ←1
Analog output for controlling external signals (DAC164 ←1).</p> | |

The display (item 10 in *Figure 1091*) is used to provide information about the Autolab to the user. *Figure 1092* shows a detail of this display.




Figure 1092 Overview of the display of the Autolab

1 Voltage indicator

Displays the measured voltage.

2 I ovl indicator

Indicates that a current overload is detected when lit.

3 T ovl

Indicates that a temperature overload is detected when lit.

4 V ovl

Indicates that a voltage overload is detected when lit.

5 BOOST mode

Indicates that a connected Booster is active when lit.

6 CELL ON

Indicates that the cell is on when lit.

7 Operation mode indicators

Indicate the operation settings of the Autolab. From top to bottom:

PSTAT indicates that the Autolab is operating in potentiostatic mode when lit.

GSTAT indicates that the Autolab is operating in galvanostatic mode when lit.

iR-C indicates that the ohmic drop compensation is on when lit

HSTAB indicates that the Autolab is operating in high stability mode when lit.

ECD indicated that the ECD module is on when lit.

OSC indicates that oscillations are detected when lit.

8 Booster current range

Indicate that a current range provided by a Booster is active when lit.

9 Autolab current ranges

The current range indicator which is lit corresponds to the active current of the Autolab.

11 Current indicator

Displays the measured current.

10 ECD current ranges

Additional current ranges provided by the ECD module. These current ranges extend the ranges of the Autolab.

**NOTICE**

The Voltage and Current values shown in the display are provided with an accuracy of 0.5 %. These values are provided for information only.

16.2.1.4 Autolab N Series back plane

The back plane of the Autolab N Series provides a number of connections, shown in *Figure 1093*.




Figure 1093 Overview of the back plane of the Autolab N Series

1 DIO P1 connector

Digital input/output connector **P1** for sending and receiving external TTL triggers.

2 DIO P2 connector

Digital input/output connector **P2** for sending and receiving external TTL triggers.

3 USB connector

Type B USB plug for connecting the USB cable to the host computer.

4 Fan

Required for cooling the Autolab during operation.

5 Mains connection socket

For connecting the Autolab to the mains supply.

6 Mains voltage indicator

Indicates the mains voltage settings of the Autolab.

7 Earth plug

For connections to the protective earth

8 GND plug

For connections to the Autolab ground.



CAUTION

Make sure that the mains voltage indicator is set properly before switching the Autolab on.

16.2.1.5 Connections for analog signals

The N Series Autolab instruments provide connections for analog signals through two different types of connectors:

- BNC connectors directly located on the front panel of the instrument *Front panel connections for analog signals (see chapter 16.2.1.5.1, page 905)*.
- BNC connectors located on the monitor cable *Monitor cable connections for analog signals (see chapter 16.2.1.5.2, page 905)*.



CAUTION

Avoid creating ground loops when connecting the Autolab to external signals as this will degrade the performance of the instrument.

16.2.1.5.1 Front panel connections for analog signals

The **ADC164** module and the **DAC164** module, installed in all the N Series Autolab instruments, are fitted with two analog inputs and two analog outputs, respectively (*see figure 1091, page 902*).

- **ADC164:** the ADC164 inputs, labeled →1 and →2 on the front panel, can be used to record any analog signal with a ± 10 V value range. The input impedance of the two analog inputs is ≥ 1 G Ω . More information on the ADC164 is provided in *Chapter 16.3.1.1*.
- **DAC164** the DAC164 outputs, labeled ←1 and ←2 on the front panel, can be used to generate any analog signal with a ± 10 V value range. The output impedance of these two inputs is 50 Ω . Corrections should be made with loads smaller than 100 k Ω . Because of dissipation, the minimum load impedance should be 200 Ω . More information on the DAC164 is provided in *Chapter 16.3.1.2*.

16.2.1.5.2 Monitor cable connections for analog signals

The **monitor cable**, supplied with the instrument, provides additional connections for analog signals, through BNC connectors. All the connections are with respect to the Autolab ground directly and indirectly with respect to the protective earth (*see figure 1094, page 906*).



Figure 1094 The monitor cable provided with the N Series Autolab instruments

To use the monitor cable, connect the cable to the matching connector located on the front panel of the Autolab. This connector, labeled \Rightarrow , is located below the front panel display (item 5 in Figure 1091).

The following connections are provided through the monitor cable:

- **Eout:** this output corresponds to the differential potential of the reference electrode (RE) with respect to the sense electrode (S). The output voltage will vary between ± 10 V. The output impedance is 50Ω , so a correction should be made if a load smaller than $100 \text{ k}\Omega$ is connected to this output. The minimum load is 200Ω .
- **iout:** this output corresponds to the output of the current-to-voltage converter circuit of the Autolab. The output corresponds to the measured current divided by the current range. The output voltage will vary between ± 10 V. The output impedance is 50Ω , so a correction should be made if a load smaller than $100 \text{ k}\Omega$ is connected to this output. The minimum load is 200Ω .
- **Ein:** this input corresponds to an analog voltage input, directly connected to the summation point of the Autolab. This input is disabled by default. When it is enabled, it can be used to control the Autolab through an external waveform generator. In potentiostatic mode, 1 V provided on this input will add 1 V to the applied potential. In galvanostat mode, 1 V provided on this input will add an extra current equal to 1 multiplied by the current range. In both cases, the converted signal is added to the value already applied by the potentiostat or the galvanostat circuit. The input range is ± 10 V and the input impedance is $1 \text{ k}\Omega$ when the connection is enabled, so a correction should be made when the source impedance is larger than 1Ω .

The **Ein** input is enabled and disabled using the **Autolab control** command (see figure 1095, page 907).




Figure 1095 The **Ein** connection is enable or disabled in the Autolab control command



CAUTION

Do not leave the **Ein** connection enabled unnecessarily to prevent noise pickup by the Autolab.

16.2.1.6 N Series Autolab restrictions

Restrictions apply when using the N Series Autolab potentiostat/galvanostat:

- **Intended use:** the Autolab potentiostat/galvanostat is intended to be used for electrochemical research only.
- **Service:** there are **no** serviceable parts inside. Servicing of the instrument can only be carried out by qualified personnel.



CAUTION

All attempts to service the instrument will lead to the immediate voiding of any warranty.



WARNING

The PGSTAT100N is fitted with a control amplifier capable of generating up to 100 V potential difference between the counter electrode (CE) and the working electrode (WE). Take all necessary precautions when working with this instrument and use the supplied warning laminated sheet to warn others.

16.2.1.7 N Series Autolab testing

NOVA is shipped with a procedure which can be used, alongside the **Diagnostics** application, to verify that the instrument is working as expected.



NOTICE

For more information on the **Diagnostics** application, please refer to *Chapter 17*.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestCV** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestCV.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1096*.





Figure 1096 The data measured by the TestCV procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestCV automatic evaluation of the data requires the following tests to succeed:

1. The residual current, determined by subtracting the expected current from the measured current, must be smaller or equal than ± 20 nA.
2. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$.
3. The intercept of the measured current versus the applied potential must be equal to ± 4 mV divided by $1000100 \Omega \pm 5 \%$.

All three conditions must be valid for the test to succeed.



16.2.1.8 Autolab N Series specifications

The specifications of the Autolab N Series are provided in *Table 27*.

Table 27 Specifications of the Autolab 7 Series instruments

Instrument	PGSTAT128N	PGSTAT302N	PGSTAT100N
Maximum current	± 800 mA	± 2 A	± 250 mA
Compliance voltage	± 12 V	± 30 V	± 100 V
Potential range	± 10 V		
Applied potential accuracy	± 0.2 % ± 2 mV		
Applied potential resolution	150 μ V		
Measured potential resolution	300 nV (gain 1000)		
Current ranges	10 nA to 1 A, 9 decades		10 nA to 100 mA, 8 decades
Current accuracy	± 0.2 % of current range		
Applied current resolution	0.015 % of current range		
Measured current resolution	0.00003 % of current range (gain 1000)		
Potentiostat bandwidth	500 kHz	1 MHz	400 kHz
Potentiostat rise/fall time	< 250 ns		< 500 ns
Input impedance of electrometer	> 1 T Ω , 8 pF		> 100 G Ω , 8 pF
Input bias current	< 1 pA		
Electrometer bandwidth	> 4 MHz		

Instrument	PGSTAT128N	PGSTAT302N	PGSTAT100N
iR compensation	2 Ω - 200 M Ω		200 m Ω - 200 M Ω
iR compensation resolution	0.025 %		
Analog output	Potential and current		
Analog voltage input	Yes		
External inputs	2		
External outputs	2		
Digital input/output	48		
Interface	USB (internal or external)		
Warm-up time	30 minutes		
Pollution degree	2		
Installation category	II		
External dimensions (without cables and accessories)	52x42x16 cm ³		
Weight	16 kg	18 kg	21 kg
Power requirements	180 W	300 W	247 W
Power supply	100 - 240 V \pm 10% in four ranges 100 V: [90 - 121 V] 120 V: [104 - 139 V] 230 V: [198 - 242 V] 240 V: [207 - 264 V]		
Power line frequency	47-63 Hz		



Instrument	PGSTAT128N	PGSTAT302N	PGSTAT100N
Fuse	100 V, 120 V: 3.15 A (slow-slow) 230 V, 240 V: 1.6 A (slow-slow)		100 V, 120 V: 3.15 A (slow-slow) 230 V, 240 V: 1.25 A (slow-slow)
Operating environment	0 °C to 40 °C, 80 % relative humidity without derating		
Storage environment	-10 °C to 60 °C		

16.2.2 Autolab F Series (AUT8) instrument

The Autolab F Series is a special version of the Autolab N Series. A single instrument is available in this series:

- **Autolab PGSTAT302F:** modular PGSTAT with 30 V compliance and 2 A maximum with floating option.

The Autolab **PGSTAT302F** is a special version of the modular **PGSTAT302N** potentiostat/galvanostat produced by Metrohm Autolab. This instrument, identified by a serial number starting with **AUT8**, is based on a modular concept that allows the instrument to be complemented by internal or external extension modules.



NOTICE

The Autolab **PGSTAT302F** derives from the **PGSTAT302N**. This chapter will only provide details on instrumental properties that deviate from the properties of the **PGSTAT302N**. Please refer to *Chapter 16.2.1* for additional information on the common properties of the instruments.

The **PGSTAT302F** is designed to be operated in two different modes:

- **Normal mode (grounded):** in this mode, the **PGSTAT302F** operates like a normal **PGSTAT302N**. In this mode, the electrochemical and working electrodes are floating with respect to the grounded instrument.
- **Floating mode:** in this mode, the **PGSTAT302F** can be used to control the potential of grounded working electrodes or can work with electrochemical cells connected to ground. In this configuration, the Autolab is floating with respect to the working electrode sample or with respect to the cell.





NOTICE

Special precautions must be taken with the cell connections when the **PGSTAT302F** is used in floating mode. Only the working electrode can be connected to ground, all other electrodes must be isolated from ground. External equipments connected to the **PGSTAT302F** must be isolated when the instrument is used in floating mode. Keep in mind that grounding of external equipment can occur through connections to a computer, if applicable (for example through a USB or RS232 cable).



CAUTION

The floating mode of the **PGSTAT302F** must only be used on grounded working electrodes or grounded cells. The working electrode or the cell can be grounded using the green ground connector embedded in the CE/WE cable of the **PGSTAT302F**.



CAUTION

Instrument performance can be substantially degraded when the **PGSTAT302F** is operated in floating mode. The instrument specifications provided by Metrohm Autolab can only be achieved when the **PGSTAT302F** is used in normal mode.

Unlike the **PGSTAT302N** from which it is derived, the **PGSTAT302F** only accommodate the following module:

- **FRA32M**: impedance spectroscopy module *FRA32M module (see chapter 16.3.2.13, page 1112)*.

16.2.2.1 Scope of delivery

The F Series Autolab systems are supplied with the following items:

- Autolab potentiostat/galvanostat
- ADC164 (installed)
- DAC164 (installed)
- Cell cable (WE/CE/GND), only suitable for PGSTAT302F
- Differential amplifier (RE/S), only suitable for PGSTAT302F
- Monitor cable
- Power cable
- BNC cable (50 cm)



- USB cable
- Set of four alligator clips
- Autolab dummy cell



CAUTION

The cables supplied with the **PGSTAT302F** can only be used in combination with this type of instrument.

16.2.2.2 Instrument power-up state

The power-up state of the instrument is hardware defined. The following settings are automatically selected whenever the instrument is powered on or whenever the connection to the instrument is reset by the software.

- Cell: off
- Mode: potentiostatic
- Control bandwidth: high stability
- iR compensation: off
- Current range: 10 mA
- Optional modules: off
- DIO ports: write mode, low state
- Summation point inputs: off
- Oscillation protection: on



CAUTION

In floating mode, the **Current overload** warning (**iOVL**) may be lit when the cell is off. This warning can be ignored.

16.2.2.3 Connections for analog signals

The Autolab PGSTAT302F provides connections for analog signals through two different types of connectors:

- BNC connectors directly located on the front panel of the instrument *Front panel connections for analog signals (see chapter 16.2.2.3.1, page 915).*
- BNC connectors located on the monitor cable *Monitor cable connections for analog signals (see chapter 16.2.2.3.2, page 915).*



CAUTION

Avoid creating ground loops when connecting the Autolab to external signals as this will degrade the performance of the instrument.

16.2.2.3.1 Front panel connections for analog signals

The **ADC164** module and the **DAC164** module, installed in all the **PGSTAT302F**, are fitted with two analog inputs and two analog outputs, respectively (see figure 1091, page 902).

- **ADC164:** the ADC164 inputs, labeled →1 and →2 on the front panel, can be used to record any analog signal with a ± 10 V value range. The input impedance of the two analog inputs is $\geq 1\text{G}\Omega$. More information on the ADC164 is provided in *Chapter 16.3.1.1*.
- **DAC164** the DAC164 outputs, labeled ←1 and ←2 on the front panel, can be used to generate any analog signal with a ± 10 V value range. The output impedance of these two inputs is $50\ \Omega$. Corrections should be made with loads smaller than $100\ \text{k}\Omega$. Because of dissipation, the minimum load impedance should be $200\ \Omega$. More information on the DAC164 is provided in *Chapter 16.3.1.2*.



CAUTION

All the signals are with respect to Autolab ground and indirectly to protective earth when the **PGSTAT302F** is operated in *normal* mode. These connectors are **floating** when the **PGSTAT302F** is operated in *floating* mode. Connected equipment may not be connected to ground and the shield of the BNC cables may not be connected to safety ground.

16.2.2.3.2 Monitor cable connections for analog signals

The **monitor cable**, supplied with the **PGSTAT302F**, provides additional connections for analog signals, through BNC connectors (see figure 1097, page 916).



Figure 1097 The monitor cable supplied with the PGSTAT302F

To use the monitor cable, connect the cable to the matching connector located on the front panel of the Autolab. This connector, labeled \Rightarrow , is located below the front panel display (item 5 in Figure 1091).

The following connections are provided through the monitor cable:

- **Eout:** this output corresponds to the inverted differential potential of the reference electrode (RE) with respect to the sense electrode (S). The output voltage will vary between ± 10 V. The output impedance is 50Ω , so a correction should be made if a load smaller than $100 \text{ k}\Omega$ is connected to this output. The minimum load is 200Ω .
- **iout:** this output corresponds to the inverted output of the current-to-voltage converter circuit of the Autolab. The output corresponds to the measured current divided by the current range. The output voltage will vary between ± 10 V. The output impedance is 50Ω , so a correction should be made if a load smaller than $100 \text{ k}\Omega$ is connected to this output. The minimum load is 200Ω .
- **Ein:** this input corresponds to an analog voltage input, directly connected to the summation point of the Autolab. This input is disabled by default. When it is enabled, it can be used to control the Autolab through an external waveform generator. In potentiostatic mode, 1 V provided on this input will add -1 V to the applied potential. In galvanostat mode, 1 V provided on this input will add an extra current equal to -1 multiplied by the current range. In both cases, the converted signal is added to the value already applied by the potentiostat or the galvanostat circuit. The input range is ± 10 V and the input impedance is $1 \text{ k}\Omega$ when the connection is enabled, so a correction should be made when the source impedance is larger than 1Ω .



CAUTION

All the signals are with respect to Autolab ground and indirectly to protective earth when the **PGSTAT302F** is operated in *normal* mode. These connectors are **floating** when the **PGSTAT302F** is operated in *floating* mode. Connected equipment may not be connected to ground and the shield of the BNC cables may not be connected to safety ground.

The **Ein** input is enabled and disabled using the **Autolab control** command (see figure 1095, page 907).




Figure 1098 The **Ein** connection is enable or disabled in the Autolab control command



CAUTION

Do not leave the **Ein** connection enabled unnecessarily to prevent noise pickup by the Autolab.

16.2.2.4 Grounded cells and grounded electrodes

The **PGSTAT302F** can be operated in two different modes:

- **Normal mode:** this mode corresponds to the operating mode using in all the PGSTAT instruments. This mode is suitable for working on electrochemical cells or electrodes that are floating with respect to the instrument.
- **Floating mode:** this mode is only available on the **PGSTAT302F**. In this mode, measurement circuitry of the Autolab is internally disconnected to protective earth (P.E.). This allows the instrument to be used in combination with a grounded working electrode or a grounded cell.

The **PGSTAT302F** can be set to either normal mode or floating mode using a dedicated short-circuit plug on the back plane of the instrument (see figure 1099, page 918).



Figure 1099 The PGSTAT302F can be set to normal mode (left) or to floating mode (right) using the provided short-circuit plug

When the short-circuit plug is connected as shown above, the instrument operates in *normal* mode. When the short-circuit plug is disconnected from the back panel, the instrument operates in *floating* mode.

16.2.2.5 Autolab PGSTAT302F restrictions

Restrictions apply when using the Autolab **PGSTAT302F** potentiostat/galvanostat:

- **Intended use:** the Autolab potentiostat/galvanostat is intended to be used for electrochemical research only.
- **Service:** there are **no** serviceable parts inside. Servicing of the instrument can only be carried out by qualified personnel.



CAUTION

All attempts to service the instrument will lead to the immediate voiding of any warranty.

- **Compliance voltage limitation:** the control amplifier of the PGSTAT302F has an output range of ± 30 V. In combination with the default cell cables, supplied with the instrument, the output range of the instrument is reduced to ± 10 V. An optional set of cell cables can be used to increase the output range to ± 30 V. These optional cables cannot be used in *floating* mode.

16.2.2.6 Autolab PGSTAT302F testing

The Autolab PGSTAT302F can be tested using the following procedures:

1. Using the TestCV procedure for the Autolab PGSTAT302F in normal (grounded) mode). Please refer to *Chapter 16.2.2.6.1* for more information.
2. Using the TestCV PGSTAT302F procedure for the Autolab PGSTAT302F in floating mode. Please refer to *Chapter 16.2.2.6.2* for more information.

16.2.2.6.1 Autolab PGSTAT302F testing in Normal mode

NOVA is shipped with a procedure which can be used, alongside the **Diagnostics** application, to verify that the instrument is working as expected.



NOTICE

For more information on the **Diagnostics** application, please refer to *Chapter 17*.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestCV** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestCV.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1100*.





Figure 1100 The data measured by the TestCV procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestCV automatic evaluation of the data requires the following tests to succeed:

1. The residual current, determined by subtracting the expected current from the measured current, must be smaller or equal than ± 20 nA.
2. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$.
3. The intercept of the measured current versus the applied potential must be equal to ± 4 mV divided by $1000100 \Omega \pm 5 \%$.

All three conditions must be valid for the test to succeed.

16.2.2.6.2 Autolab PGSTAT302F testing in Floating mode

NOVA is shipped with a procedure which can be used, alongside the **Diagnostics** application, to verify that the instrument is working as expected.



NOTICE

For more information on the **Diagnostics** application, please refer to *Chapter 17*.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestCV PGSTAT302F** procedure, provided in the NOVA 2.X installation folder (Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestCV PGSTAT302F.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1101*.





Figure 1101 The data measured by the TestCV procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestCV automatic evaluation of the data requires the following tests to succeed:

1. The residual current, determined by subtracting the expected current from the measured current, must be smaller or equal than ± 20 nA.
2. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$.
3. The intercept of the measured current versus the applied potential must be equal to ± 4 mV divided by $1000100 \Omega \pm 5 \%$.

All three conditions must be valid for the test to succeed.

16.2.2.7 Autolab F Series specifications

The specifications of the F Series Autolab are provided in *Table 28*.

Table 28 Specifications of the F Series Autolab instruments

Instrument	PGSTAT302F
Maximum current	± 2 A
Compliance voltage	± 30 V, ± 10 V (with default cables)
Potential range	± 10 V
Applied potential accuracy	± 0.2 % ± 2 mV
Applied potential resolution	150 μ V
Measured potential resolution	300 nV (gain 1000)
Current ranges	10 nA to 1 A, 9 decades
Current accuracy	± 0.2 % of current range
Applied current resolution	0.015 % of current range
Measured current resolution	0.00003 % of current range (gain 1000)
Potentiostat bandwidth	100 kHz
Potentiostat rise/fall time	< 250 ns
Input impedance of electrometer	> 1 T Ω , 8 pF
Input bias current	< 1 pA
Electrometer bandwidth	> 4 MHz
iR compensation	2 Ω - 200 M Ω



Instrument	PGSTAT302F
iR compensation resolution	0.025 %
Analog output	Potential and current
Analog voltage input	Yes
External inputs	2
External outputs	2
Digital input/output	48
Interface	USB
Warm-up time	30 minutes
Pollution degree	2
Installation category	II
External dimensions (without cables and accessories)	52x42x16 cm ³
Weight	18 kg
Power requirements	300 W
Power supply	100 - 240 V \pm 10% in four ranges 100 V: [90 - 121 V] 120 V: [104 - 139 V] 230 V: [198 - 242 V] 240 V: [207 - 264 V]
Power line frequency	47-63 Hz
Fuse	100 V, 120 V: 3.15 A (slow-slow) 230 V, 240 V: 1.6 A (slow-slow)
Operating environment	0 °C to 40 °C, 80 % relative humidity without derating

Instrument	PGSTAT302F
Storage environment	-10 °C to 60 °C

16.2.3 Autolab MBA N Series (AUT8) instruments

The Autolab MBA N Series is a special version of the modular potentiostat/galvanostat produced by Metrohm Autolab. These instruments, identified by a serial number starting with **AUT8**, are based on a modular concept that allows the instrument to be complemented by internal or external extension modules.

The following instruments belong to the Autolab MBA N Series:

- **Autolab PGSTAT302N MBA:** modular PGSTAT with 30 V compliance and 2 A maximum current.
- **Autolab PGSTAT128N MBA:** modular PGSTAT with 12 V compliance and 800 mA maximum current.



NOTICE

The **Autolab MBA N Series** derives from the **Autolab N Series**. This chapter will only provide details on instrumental properties that deviate from the properties of the **Autolab N Series**. Please refer to *Chapter 16.2.1* for additional information on the common properties of the instruments.

Unlike the Autolab N Series instruments which can accommodate a wide range of internal extension modules, the Autolab N MBA Series can only accommodate the following modules:

- **FRA32M or FRA2 module:** impedance spectroscopy module (please refer to *Chapter 16.3.2.13* and *Chapter 16.3.2.12* for more information).
- **BA module:** dual mode bipotentiostat module. Up to five BA modules can be placed in each MBA instrument *BA module (see chapter 16.3.2.3, page 1011)*.

The BA modules installed in the MBA instrument can be used to control up to five additional working electrodes, sharing a common counter electrode and reference electrode. These instruments can therefore be used to work with sensor arrays or with electrochemical cells in which more than one working electrode is located.

Each of the five BA modules are identified by a specified MBA module label (*see figure 1102, page 926*).




Figure 1102 The module labels used to identify the BA modules installed in a MBA instrument

16.2.3.1 N MBA Series Autolab front panel



NOTICE

The front panel of the **Autolab MBA N Series** instrument is arranged differently from the **Autolab N Series** instruments.

The front panel of the N MBA Series Autolab provides a number of connections, controls and indicators (see figure 1103, page 926).




Figure 1103 Overview of the front panel of the N MBA Series Autolab

- | | |
|--|--|
| <p>1 On/Off button
For switching the Autolab on or off.</p> | <p>2 ADC164 →1
Analog input for recording external signals (ADC164 →1).</p> |
|--|--|

- | | |
|--|--|
| <p>3 ADC164 →2
Analog input for recording external signals (ADC164 →2).</p> | <p>4 DAC164 ←2
Analog output for controlling external signals (DAC164 ←2).</p> |
| <p>5 Monitor cable connector ⇌
For connecting the monitor cable.</p> | <p>6 CE/WE connector
For connecting the Autolab cell cable, providing connections to the counter electrode (CE), working electrode (WE) and ground.</p> |
| <p>7 RE/S connector
For connecting the Autolab differential amplifier, providing connections to the reference electrode (RE) and sense electrode (S).</p> | <p>8 Ground connector
Additional ground connector for connecting external devices to the Autolab ground.</p> |
| <p>9 Cell ON button
For enabling and disabling the cell.</p> | <p>10 Display
Display indicating real-time information on the measured current and potential and instrumental settings.</p> |
| <p>11 DAC164 ←1
Analog output for controlling external signals (DAC164 ←1).</p> | |

The display (item 10 in *Figure 1103*) is identical to the display of the Autolab N Series (see *figure 1092, page 903*).

16.2.4 Autolab Compact Series (AUT4/AUT5) instruments

The Autolab Compact Series provides potentiostat/galvanostat instruments with a very small footprint.

The following instruments belong to the Autolab Compact Series:

- **Autolab PGSTAT101:** PGSTAT with 10 V compliance and 100 mA maximum current, identified with a serial number starting with **AUT4**. The PGSTAT101 can be complemented by **external** extension modules only.
- **Autolab PGSTAT204:** modular PGSTAT with 20 V compliance and 400 mA maximum current, identified with a serial number starting with **AUT5**. The PGSTAT204 can be complemented by a selection of **internal** and **external** extension modules.

16.2.4.1 Compact Series Autolab scope of delivery

The Compact Series Autolab systems are supplied with the following items:

- Autolab potentiostat/galvanostat
- Cell cable (RE/S/WE/CE/GND)
- Power cable
- USB cable
- Set of four alligator clips



The PGSTAT204 is also supplied with the Autolab dummy cell.

The monitor cable, used to interface to external devices, is also available for the Compact Series Autolab, as an option.

16.2.4.2 Instrument power-up state

The power-up state of the instrument is hardware defined. The following settings are automatically selected whenever the instrument is powered on or whenever the connection to the instrument is reset by the software.

- Cell: off
- Mode: potentiostatic
- Control bandwidth: high stability
- iR compensation: off
- Current range: 1 μ A
- Optional modules: off
- DIO ports: low state
- Summation point inputs: off
- Internal dummy cell: off

16.2.4.3 Autolab Compact Series front panel

The front panel of the Autolab PGSTAT101 provides a number of connections and indicators (*see figure 1104, page 929*).




Figure 1104 The front panel of the PGSTAT101

1 Status LED

Dual color LED used to indicate the status of the PGSTAT101. When the LED is red, an overload is detected. When the LED is green, the cell is on and no overloads are detected. When the LED is switched off, the cell is off and no overloads are detected.

2 I/O connector

Used to connect the optional monitor cable providing connections for E_{out} , I_{out} , V_{out} and V_{in} .

3 CELL connector

Used to connect the cell cable providing connections for the counter (CE), reference (RE), working (WE) and sense (S) electrode as well as the ground.

The front panel of the Autolab PGSTAT204 provides a number of connections and indicators (see figure 1105, page 930).




Figure 1105 The front panel of the PGSTAT204

1 DIO connector

For connecting a DIO cable or interfacing to external devices through the digital inputs and outputs.

2 I/O connector

Used to connect the optional monitor cable providing connections for E_{out} , I_{out} , V_{out} and V_{in} .

3 Status LED

Dual color LED used to indicate the status of the PGSTAT204. When the LED is red, an overload is detected. When the LED is green, the cell is on and no overloads are detected. When the LED is switched off, the cell is off and no overloads are detected.

4 CELL connector

Used to connect the cell cable providing connections for the counter (CE), reference (RE), working (WE) and sense (S) electrode as well as the ground. The cell is represented by the symbol \oplus .

16.2.4.4 Autolab Compact Series back plane

The back plane of the Autolab PGSTAT101 provides a number of connections and controls, shown in *Figure 1106*.




Figure 1106 The back plane of the PGSTAT101

- | | |
|--|--|
| <p>1 Fan
Required for cooling the Autolab during operation.</p> | <p>2 DIO port
For connecting a DIO cable or interfacing to external devices through the digital inputs and outputs.</p> |
| <p>3 USB connector
Type B USB plug for connecting the USB cable to the host computer.</p> | <p>4 Fuse holder
Holds the mains connection socket fuse.</p> |
| <p>5 Mains connection socket
For connecting the Autolab to the mains supply.</p> | <p>6 On/Off switch
For switching the Autolab on or off.</p> |

The back plane of the Autolab PGSTAT204 provides a number of connections and controls, shown in *Figure 1107*.




Figure 1107 The back plane of the PGSTAT204

- | | |
|---|--|
| <p>1 Fan
Required for cooling the Autolab during operation.</p> | <p>2 USB connector
Type B USB plug for connecting the USB cable to the host computer.</p> |
| <p>3 Mains connection socket
For connecting the Autolab to the mains supply.</p> | <p>4 Fuse holder
Holds the mains connection socket fuse.</p> |
| <p>5 On/Off switch
For switching the Autolab on or off.</p> | |

16.2.4.5 Connections for analog signals

The Autolab Compact Series instruments provide connections for analog signals through an optional monitor cable.



NOTICE

The monitor cable is not supplied with the instrument and it need to be ordered separately. A dedicated cable is available for each instrument type.



CAUTION

Avoid creating ground loops when connecting the Autolab to external signals as this will degrade the performance of the instrument.

16.2.4.5.1 Monitor cable for Autolab PGSTAT101

The **monitor cable** for **PGSTAT101** provides additional connections for analog signals, through BNC connectors. All the connections are with respect to the Autolab ground directly and indirectly with respect to the protective earth (see figure 1108, page 933).



Figure 1108 The monitor cable for PGSTAT101

To use the monitor cable, connect the cable to the matching connector located on the front panel of the Autolab. This connector, labeled **I/O**, is located on the front panel of the instrument (item 2 in Figure 1104).

The following connections are provided through the monitor cable:

- **Eout:** this output corresponds to the differential potential of the reference electrode (RE) with respect to the sense electrode (S). The output voltage will vary between ± 10 V. The output impedance is $1\text{ k}\Omega$, so a correction should be made if a load smaller than $2\text{ M}\Omega$ is connected to this output. The minimum load is $4000\ \Omega$.
- **iout:** this output corresponds to the inverted output of the current-to-voltage converter circuit of the Autolab. The output corresponds to the inverted measured current divided by the current range. The output voltage will vary between ± 10 V. The output impedance is $50\ \Omega$, so a correction should be made if a load smaller than $100\text{ k}\Omega$ is connected to this output. The minimum load is $200\ \Omega$.



- **Vout:** this output corresponds to the output of the on-board DAC of the PGSTAT101. It can be used to generate any analog signal with a ± 10 V value range. The output impedance is 1Ω and this output is capable of supplying a current of up to 5 mA. Corrections should be made with loads smaller than $2 \text{ k}\Omega$. More information on the on-board DAC is provided in *Chapter 16.3.1.2*.
- **Vin:** this input corresponds to the input of the on-board ADC of the PGSTAT101. It can be used to record any analog signal with a ± 10 V value range. The input impedance is $\geq 1 \text{ G}\Omega$. More information on the on-board ADC is provided in *Chapter 16.3.1.1*.

16.2.4.5.2 Monitor cable for Autolab PGSTAT204

The **monitor cable** for **PGSTAT204** provides additional connections for analog signals, through BNC connectors. All the connections are with respect to the Autolab ground directly and indirectly with respect to the protective earth (*see figure 1109, page 934*).




Figure 1109 The monitor cable for PGSTAT204

To use the monitor cable, connect the cable to the matching connector located on the front panel of the Autolab. This connector, labeled **I/O**, is located on the front panel of the instrument (item 2 in *Figure 1105*).

The following connections are provided through the monitor cable:

- **Eout:** this output corresponds to the differential potential of the reference electrode (RE) with respect to the sense electrode (S). The output voltage will vary between ± 10 V. The output impedance is $1 \text{ k}\Omega$, so a correction should be made if a load smaller than $2 \text{ M}\Omega$ is connected to this output. The minimum load is 4000Ω .

- **ious:** this output corresponds to the inverted output of the current-to-voltage converter circuit of the Autolab. The output corresponds to the inverted measured current divided by the current range. The output voltage will vary between ± 10 V. The output impedance is 50Ω , so a correction should be made if a load smaller than $100 \text{ k}\Omega$ is connected to this output. The minimum load is 200Ω .
- **Vout:** this output corresponds to the output of the on-board DAC of the PGSTAT204. It can be used to generate any analog signal with a ± 10 V value range. The output impedance is 1Ω and this output is capable of supplying a current of up to 5 mA. Corrections should be made with loads smaller than $2 \text{ k}\Omega$. More information on the on-board DAC is provided in *Chapter 16.3.1.2*.
- **Vin:** this input corresponds to the input of the on-board ADC of the PGSTAT204. It can be used to record any analog signal with a ± 10 V value range. The input impedance is $\geq 1 \text{ G}\Omega$. More information on the on-board ADC is provided in *Chapter 16.3.1.1*.

16.2.4.6 Compact Autolab Series restrictions

Restrictions apply when using the Compact Autolab Series potentiostat/galvanostat:

- **Intended use:** the Compact Autolab Series potentiostat/galvanostat is intended to be used for electrochemical research only.
- **Service:** there are **no** serviceable parts inside. Servicing of the instrument can only be carried out by qualified personnel.



CAUTION

All attempts to service the instrument will lead to the immediate voiding of any warranty.

16.2.4.7 Compact Series Autolab testing

The Autolab PGSTAT101 and PGSTAT204 can be tested using the following procedures:

1. For the Autolab PGSTAT101, a dedicated test with the internal dummy cell is available. Please refer to *Chapter 16.2.4.7.1* for more information.
2. For the Autolab PGSTAT204, the standard TestCV procedure with the Autolab dummy cell. Please refer to *Chapter 16.2.4.7.2* for more information.

TestCV PGSTAT101

The following problems were encountered during validation.

⚠ The internal dummy cell is on.

OK

Cancel

4 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1110*.





Figure 1110 The data measured by the TestCV procedure

5 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestCV automatic evaluation of the data requires the following tests to succeed:

will be shown. The measured data should look as shown in *Figure 1111*.





Figure 1111 The data measured by the TestCV procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestCV automatic evaluation of the data requires the following tests to succeed:

1. The residual current, determined by subtracting the expected current from the measured current, must be smaller or equal than ± 20 nA.
2. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5\%$.
3. The intercept of the measured current versus the applied potential must be equal to ± 4 mV divided by $1000100 \Omega \pm 5\%$.

All three conditions must be valid for the test to succeed.

16.2.4.8 Autolab Compact Series specifications

The specifications of the Autolab Compact Series are provided in *Table 29*.

Table 29 Specifications of the Autolab Compact Series instruments

Instrument	PGSTAT101	PGSTAT204
Maximum current	± 100 mA	± 400 mA
Compliance voltage	± 10 V	± 20 V
Potential range	± 10 V	



Instrument	PGSTAT101	PGSTAT204
Applied potential accuracy	$\pm 0.2 \% \pm 2 \text{ mV}$	
Applied potential resolution	150 μV	
Measured potential resolution	3 μV (gain 100)	
Current ranges	10 nA to 10 mA, 7 decades	10 nA to 100 mA, 8 decades
Current accuracy	$\pm 0.2 \%$ of current range	
Applied current resolution	0.015 % of current range	
Measured current resolution	0.00003 % of current range (gain 1000)	
Potentiostat bandwidth	1 MHz	
Potentiostat rise/fall time	< 300 ns	
Input impedance of electrometer	> 100 G Ω , 8 pF	
Input bias current	< 1 pA	
Electrometer bandwidth	> 4 MHz	
iR compensation	20 m Ω - 200 M Ω	200 m Ω - 200 M Ω
iR compensation resolution	0.025 %	
Analog output	Potential and current	
Analog voltage input	Yes	
External inputs	2	
External outputs	2	
Digital input/output	12	
Interface	USB	
Warm-up time	30 minutes	
Pollution degree	2	
Installation category	II	

Instrument	PGSTAT101	PGSTAT204
External dimensions (without cables and accessories)	9x21x15 cm ³	15x26x20 cm ³
Weight	2.1 kg	4.1 kg
Power requirements	40 W	75 W
Power supply	100 - 240 V \pm 10% in four ranges (auto select)	
Power line frequency	47-63 Hz	
Fuse	2 A (slow-slow)	3.5 A (slow-slow)
Operating environment	0 °C to 40 °C, 80 % relative humidity without derating	
Storage environment	-10 °C to 60 °C	

16.2.5 Multi Autolab Series (MAC8/MAC9) instruments

The Multi Autolab Series provides cabinets that can accommodate up to 12 potentiostat/galvanostat channels or a combination of potentiostat/galvanostat modules with expansion modules.

The following instruments belong to the Multi Autolab Series:

- **Multi Autolab M101:** Multi Autolab cabinet designed to accommodate up to 12 M101 potentiostat/galvanostat modules with 10 V compliance and 100 mA maximum current. The cabinet is identified with a serial number starting with **MAC8**. The M101 can be complemented by **internal** or **external** extension modules only.
- **Multi Autolab M204:** Multi Autolab cabinet designed to accommodate up to 12 M204 potentiostat/galvanostat modules with 20 V compliance and 400 mA maximum current. The cabinet is identified with a serial number starting with **MAC9**. The M204 can be complemented by **internal** or **external** extension modules only.



CAUTION

M101 modules can only be installed in a M101 Multi Autolab cabinet.

M204 modules can only be installed in a M204 Multi Autolab cabinet.

The Multi Autolab cabinet is fitted with twelve module bays, labeled 1 to 6 (from left to right) and A to F (from left to right). Potentiostat/galvanostat modules can be installed in any available module bays. Internal expansion




Figure 1113 The front panel of the Multi Autolab M101

1 On/Off button

For switching the Multi Autolab on or off.

2 DIO connector

For connecting a DIO cable or interfacing to external devices through the digital inputs and outputs.

3 I/O connector

Used to connect the optional monitor cable providing connections for E_{out} , I_{out} , V_{out} and V_{in} .

4 Status LED

Dual color LED used to indicate the status of the M101. When the LED is red, an overload is detected. When the LED is green, the cell is on and no overloads are detected. When the LED is switched off, the cell is off and no overloads are detected.

5 CELL connector

Used to connect the cell cable providing connections for the counter (CE), reference (RE), working (WE) and sense (S) electrode as well as the ground. The cell is represented by the symbol Φ .

The front panel of the Multi Autolab M204 provides a number of connections and indicators (see figure 1114, page 944).




Figure 1114 The front panel of the Multi Autolab M204

1 On/Off button

For switching the Multi Autolab on or off.

2 DIO connector

For connecting a DIO cable or interfacing to external devices through the digital inputs and outputs.

3 I/O connector

Used to connect the optional monitor cable providing connections for E_{out} , I_{out} , V_{out} and V_{in} .

4 Status LED

Dual color LED used to indicate the status of the M204. When the LED is red, an overload is detected. When the LED is green, the cell is on and no overloads are detected. When the LED is switched off, the cell is off and no overloads are detected.

5 CELL connector

Used to connect the cell cable providing connections for the counter (CE), reference (RE), working (WE) and sense (S) electrode as well as the ground. The cell is represented by the symbol Φ .

16.2.5.4 Multi Autolab Series back plane

The back plane of the Multi Autolab provides a number of connections and controls, shown in *Figure 1115*.




Figure 1115 The back plane of the Multi Autolab

- | | |
|--|--|
| <p>1 USB hub main connector
Type B USB plug for connecting the USB cable to the host computer providing connections to all modules.</p> | <p>2 USB hub sub-1 connector
Type B USB plug for connecting the USB cable to the host computer providing connections module bays 3, C, 4 and D.</p> |
| <p>3 USB hub sub-2 connector
Type B USB plug for connecting the USB cable to the host computer providing connections module bays 5, E, 6 and F.</p> | <p>4 Fan
Required for cooling the Multi Autolab during operation.</p> |
| <p>5 Ground plug
For grounding the Multi Autolab cabinet</p> | <p>6 Fuse holders
Holds the mains connection socket fuses.</p> |
| <p>7 Mains connection socket
For connecting the Autolab to the mains supply.</p> | |

16.2.5.5 Connections for analog signals

The Multi Autolab Series instruments provide connections for analog signals through an optional monitor cable.



NOTICE

The monitor cable is not supplied with the instrument and it need to be ordered separately.



CAUTION

Avoid creating ground loops when connecting the Autolab to external signals as this will degrade the performance of the instrument.

The **monitor cable** for **Multi Autolab M101** and **Multi Autolab M204** provides additional connections for analog signals, through BNC connectors. All the connections are with respect to the Autolab ground directly and indirectly with respect to the protective earth (see figure 1109, page 934).




Figure 1116 The monitor cable for Multi Autolab M101 and Multi Autolab M204

To use the monitor cable, connect the cable to the matching connector located on the front panel of the M101 or M204 module in the Multi Autolab instrument. This connector, labeled **I/O**, is located on the front panel of the potentiostat/galvanostat module (item 3 in Figure 1113 and item 3 in Figure 1114).

The following connections are provided through the monitor cable:

- **Eout:** this output corresponds to the differential potential of the reference electrode (RE) with respect to the sense electrode (S). The output voltage will vary between ± 10 V. The output impedance is $1\text{ k}\Omega$, so a correction should be made if a load smaller than $2\text{ M}\Omega$ is connected to this output. The minimum load is $4000\ \Omega$.
- **Iout:** this output corresponds to the inverted output of the current-to-voltage converter circuit of the Autolab. The output corresponds to the inverted measured current (for the M101) or the measured current (for the M204) divided by the current range. The output voltage will vary between ± 10 V. The output impedance is $50\ \Omega$, so a correction should be made if a load smaller than $100\text{ k}\Omega$ is connected to this output. The minimum load is $200\ \Omega$.
- **Vout:** this output corresponds to the output of the on-board DAC of the M101 or M204. It can be used to generate any analog signal with a ± 10 V value range. The output impedance is $1\ \Omega$ and this output is capable of supplying a current of up to 5 mA . Corrections should be made with loads smaller than $2\text{ k}\Omega$. More information on the on-board DAC is provided in Chapter 16.3.1.2.

- **Vin:** this input corresponds to the input of the on-board ADC of the M101 or M204. It can be used to record any analog signal with a ± 10 V value range. The input impedance is ≥ 1 G Ω . More information on the ADC164 is provided in *Chapter 16.3.1.1*.

16.2.5.6 Multi Autolab Series connection hub

The Multi Autolab cabinet is fitted with an internal USB hub that allows the Multi Autolab to be shared in between a maximum of three computers running NOVA.

When the Multi Autolab is connected to the **main** USB connector, the host computer controls all twelve module bays. If an additional host computer is connected to either one of the sub USB connectors (**sub-1** or **sub-2**), that computer will take over the control of the module bays indicated in the labels below the sub USB connectors (see *figure 1117, page 947*).




Figure 1117 The Multi Autolab connection hub




CAUTION

Never change the USB connections while an experiment is running on any of the Multi Autolab channels since this could lead to a loss of control of an experiment and a loss of data.

Depending on the connection to the USB hub, the host computer can control the following module bays:

- **main USB:** all the module bays in the Multi Autolab cabinet are controlled.
- **sub-1 USB:** module bays 3, C, 4 and D are controlled.
- **sub-2 USB:** module bays 5, E, 6 and F are controlled.

For example, connecting a computer to the sub-2 USB connector on the back plane will provide control over the modules installed in bays 5, E, 6 and F, exclusively. Connecting an additional computer to the main USB connector will provide control over all the remaining modules (1, A, 2, B, 3, C, 4, D). If a third PC is connected to the sub-1 USB connector, control over the modules in bays 3, C, 4 and D will be transferred from the com-



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1118*.





Figure 1118 The data measured by the TestCV procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestCV automatic evaluation of the data requires the following tests to succeed:



1. The residual current, determined by subtracting the expected current from the measured current, must be smaller or equal than ± 20 nA.
2. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5\%$.
3. The intercept of the measured current versus the applied potential must be equal to ± 4 mV divided by $1000100 \Omega \pm 5\%$.

All three conditions must be valid for the test to succeed.

16.2.5.9 Multi Autolab Series specifications

The specifications of the Multi Autolab Series are provided in *Table 30*.

Table 30 Specifications of the Multi Autolab Series instruments

Instrument	M101	M204
Maximum current	± 100 mA	± 400 mA
Compliance voltage	± 10 V	± 20 V
Potential range	± 10 V	
Applied potential accuracy	$\pm 0.2\% \pm 2$ mV	
Applied potential resolution	150 μ V	
Measured potential resolution	300 nV (gain 1000)	
Current ranges	10 nA to 10 mA, 7 decades	10 nA to 100 mA, 8 decades
Current accuracy	$\pm 0.2\%$ of current range	
Applied current resolution	0.015% of current range	
Measured current resolution	0.00003% of current range (gain 1000)	
Potentiostat bandwidth	1 MHz	
Potentiostat rise/fall time	< 300 ns	
Input impedance of electrometer	> 100 G Ω , 8 pF	
Input bias current	< 1 pA	
Electrometer bandwidth	> 4 MHz	
iR compensation	20 m Ω - 200 M Ω	200 m Ω - 200 M Ω

Instrument	M101	M204
iR compensation resolution	0.025 %	
Analog output	Potential and current	
Analog voltage input	Yes	
External inputs	2	
External outputs	2	
Digital input/output	12	
Interface	USB	
Warm-up time	30 minutes	
Pollution degree	2	
Installation category	II	
External dimensions (without cables and accessories)	52x42x17 cm ³	
Weight	13 kg	14 kg
Power requirements	200 W	700 W
Power supply	100 - 240 V \pm 10% in four ranges (auto select)	
Power line frequency	47-63 Hz	
Fuse	8 A (slow-slow)	8 A (slow-slow)
Operating environment	0 °C to 40 °C, 80 % relative humidity without derating	
Storage environment	-10 °C to 60 °C	

16.2.6 Autolab 7 Series (AUT7) instruments

The Autolab 7 Series is the predecessor version of the Autolab N Series modular potentiostat/galvanostat produced by Metrohm Autolab *Autolab N Series (AUT8) instruments* (see chapter 16.2.1, page 900). These instruments, identified by a serial number starting with **AUT7**, are based on modular concept that allows the instrument to be complemented by internal or external extension modules.

- DIO ports: write mode, low state
- Summation point inputs: off
- Oscillation protection: on

16.2.6.3 7 Series Autolab front panel

The front panel of the 7 Series Autolab provides a number of connections, controls and indicators (see figure 1119, page 953).




Figure 1119 Overview of the front panel of the 7 Series Autolab

- | | |
|--|--|
| <p>1 On/Off button
For switching the Autolab on or off.</p> | <p>2 ADC164 →2
Analog input for recording external signals (ADC164 →2).</p> |
| <p>3 DAC164 ←2
Analog output for controlling external signals (DAC164 ←2).</p> | <p>4 Monitor cable connector ⇌
For connecting the monitor cable.</p> |
| <p>5 CE/WE connector
For connecting the Autolab cell cable, providing connections to the counter electrode (CE) and working electrode (WE).</p> | <p>6 RE/S connector
For connecting the Autolab differential amplifier, providing connections to the reference electrode (RE) and sense electrode (S).</p> |
| <p>7 CELL ENABLE button
For enabling and disabling the cell.</p> | <p>8 Display switch mode button
For switching between the voltage and the current on the display.</p> |
| <p>9 Display
Display indicating real-time information on the measured current and potential and instrumental settings.</p> | <p>10 Ground connector
Additional ground connector for connecting external devices to the Autolab ground.</p> |
| <p>11 DAC164 ←1
Analog output for controlling external signals (DAC164 ←1).</p> | <p>12 ADC164 →1
Analog input for recording external signals (ADC164 →1).</p> |

The display (item 9 in Figure 1119) is used to provide information about the Autolab to the user. Figure 1120 shows a detail of this display.



Figure 1120 Overview of the display of the Autolab

1 Booster current range

Indicate that a current range provided by a Booster is active when lit.

3 ECD current range

Indicates that the ECD module is used when lit. The actual current range is corresponds to the active Autolab current range divided by 1000.

5 CELL ON

Indicates that the cell is on when lit.

7 OSC indicator

Indicates that oscillations are detected when lit.

9 I ovl indicator

Indicates that a current overload is detected when lit.

11 Voltage/Current indicator

Displays the measured voltage or current.

2 Autolab current ranges

The current range indicator which is lit corresponds to the active current of the Autolab.

4 Display switch mode button

For switching between the voltage and the current on the display.

6 Operation mode indicators

Indicate the operation settings of the Autolab. From top to bottom:

PSTAT indicates that the Autolab is operating in potentiostatic mode when lit.

GSTAT indicates that the Autolab is operating in galvanostatic mode when lit.

iR-C indicates that the ohmic drop compensation is on when lit

HSTAB indicates that the Autolab is operating in high stability mode when lit.

8 V ovl indicator

Indicates that a voltage overload is detected when lit.

10 Unit indicator

Indicates the units used for the value shown on the display.

16.2.6.4 Autolab 7 Series back plane

The back plane of the Autolab 7 Series provides a number of connections, shown in *Figure 1121*.




Figure 1121 Overview of the back plane of the Autolab 7 Series

- | | |
|--|--|
| <p>1 DIO P1 connector
Digital input/output connector P1 for sending and receiving external TTL triggers.</p> | <p>2 DIO P2 connector
Digital input/output connector P2 for sending and receiving external TTL triggers.</p> |
| <p>3 USB connector
Type B USB plug for connecting the USB cable to the host computer.</p> | <p>4 USB Hub
For connecting additional USB devices.</p> |
| <p>5 Fan
Required for cooling the Autolab during operation.</p> | <p>6 Mains connection socket
For connecting the Autolab to the mains supply.</p> |
| <p>7 Mains voltage indicator
Indicates the mains voltage settings of the Autolab.</p> | <p>8 Earth plug
For connections to the protective earth</p> |
| <p>9 GND plug
For connections to the Autolab ground.</p> | |



CAUTION

Make sure that the mains voltage indicator is set properly before switching the Autolab on.

Some of the first Autolab 7 Series instruments are not fitted with an internal USB interface. These instruments are controlled through an external USB interface adapter. The back plane of these instruments is different (see *figure 1122, page 956*).




Figure 1122 Overview of the back plane of the Autolab 7 Series (without USB)

1 PC INTERFACE connector

For connecting the external USB interface adapter.

2 DIO P1 connector

Digital input/output connector **P1** for sending and receiving external TTL triggers.

3 DIO P2 connector

Digital input/output connector **P2** for sending and receiving external TTL triggers.



NOTICE

All other items present on the back plane of *Figure 1122* are the same as in *Figure 1121*.

16.2.6.5 Connections for analog signals

The Autolab 7 Series instruments provide connections for analog signals through two different types of connectors:

- BNC connectors directly located on the front panel of the instrument *Front panel connections for analog signals (see chapter 16.2.6.5.1, page 957)*.
- BNC connectors located on the monitor cable *Monitor cable connections for analog signals (see chapter 16.2.6.5.2, page 957)*.



CAUTION

Avoid creating ground loops when connecting the Autolab to external signals as this will degrade the performance of the instrument.

16.2.6.5.1 Front panel connections for analog signals

The **ADC164** module and the **DAC164** module, installed in all the 7 Series Autolab instruments, are fitted with two analog inputs and two analog outputs, respectively (see *figure 1091*, page 902).

- **ADC164:** the ADC164 inputs, labeled $\rightarrow 1$ and $\rightarrow 2$ on the front panel, can be used to record any analog signal with a ± 10 V value range. The input impedance of the two analog inputs is ≥ 1 G Ω . More information on the ADC164 is provided in *Chapter 16.3.1.1*.
- **DAC164** the DAC164 outputs, labeled $\leftarrow 1$ and $\leftarrow 2$ on the front panel, can be used to generate any analog signal with a ± 10 V value range. The output impedance of these two inputs is 50 Ω . Corrections should be made with loads smaller than 100 k Ω . Because of dissipation, the minimum load impedance should be 200 Ω . More information on the DAC164 is provided in *Chapter 16.3.1.2*.

16.2.6.5.2 Monitor cable connections for analog signals

The **monitor cable**, supplied with the instrument, provides additional connections for analog signals, through BNC connectors. All the connections are with respect to the Autolab ground directly and indirectly with respect to the protective earth (see *figure 1123*, page 957).



Figure 1123 The monitor cable provided with the 7 Series Autolab instruments

To use the monitor cable, connect the cable to the matching connector located on the front panel of the Autolab. This connector, labeled \rightleftharpoons , is located below the front panel display (item 4 in *Figure 1119*).

The following connections are provided through the monitor cable:



CAUTION

Do not leave the **Ein** connection enabled unnecessarily to prevent noise pickup by the Autolab.

16.2.6.6 7 Series Autolab restrictions

Restrictions apply when using the 7 Series Autolab potentiostat/galvanostat:

- **Intended use:** the Autolab potentiostat/galvanostat is intended to be used for electrochemical research only.
- **Service:** there are **no** serviceable parts inside. Servicing of the instrument can only be carried out by qualified personnel.



CAUTION

All attempts to service the instrument will lead to the immediate voiding of any warranty.



WARNING

The PGSTAT100 is fitted with a control amplifier capable of generating up to 100 V potential difference between the counter electrode (CE) and the working electrode (WE). Take all necessary precautions when working with this instrument and use the supplied warning laminated sheet to warn others.

16.2.6.7 7 Series Autolab testing

NOVA is shipped with a procedure which can be used, alongside the **Diagnostics** application, to verify that the instrument is working as expected.



NOTICE

For more information on the **Diagnostics** application, please refer to *Chapter 17*.

Follow the steps described below to run the test procedure.




1 Load the procedure

Load the **TestCV** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestCV.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1125*.





Figure 1125 The data measured by the TestCV procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestCV automatic evaluation of the data requires the following tests to succeed:

1. The residual current, determined by subtracting the expected current from the measured current, must be smaller or equal than ± 20 nA.
2. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5\%$.
3. The intercept of the measured current versus the applied potential must be equal to ± 4 mV divided by $1000100 \Omega \pm 5\%$.

All three conditions must be valid for the test to succeed.

16.2.6.8 7 Series Autolab specifications

The specifications of the 7 Series Autolab are provided in *Table 31*.

Table 31 Specifications of the 7 Series Autolab instruments

Instrument	PGSTAT12	PGSTAT302/30	PGSTAT100
Maximum current	± 250 mA	± 2 A/ ± 1 A	± 250 mA
Compliance voltage	± 12 V	± 30 V	± 100 V
Potential range	± 10 V		
Applied potential accuracy	$\pm 0.2\% \pm 2$ mV		
Applied potential resolution	150 μ V		
Measured potential resolution	300 nV (gain 1000)		
Current ranges	10 nA to 100 mA, 8 decades	10 nA to 1 A, 9 decades	10 nA to 100 mA, 8 decades
Current accuracy	$\pm 0.2\%$ of current range		
Applied current resolution	0.015% of current range		



Instrument	PGSTAT12	PGSTAT302/30	PGSTAT100
Measured current resolution	0.00003 % of current range (gain 1000)		
Potentiostat bandwidth	500 kHz	1 MHz	400 kHz
Potentiostat rise/fall time	< 500 ns	< 250 ns	< 500 ns
Input impedance of electrometer	> 100 G Ω , 8 pF	> 1 T Ω , 8 pF	> 100 G Ω , 8 pF
Input bias current	< 1 pA		
Electrometer bandwidth	> 4 MHz		
iR compensation	2 Ω - 200 M Ω		200 m Ω - 200 M Ω
iR compensation resolution	0.025 %		
Analog output	Potential and current		
Analog voltage input	Yes		
External inputs	2		
External outputs	2		
Digital input/output	48		
Interface	USB		
Warm-up time	30 minutes		
Pollution degree	2		
Installation category	II		
External dimensions (without cables and accessories)	52x42x17 cm ³		

Instrument	PGSTAT12	PGSTAT302/30	PGSTAT100
Weight	22 kg	25 kg	25 kg
Power requirements	247 W	247 W	247 W
Power supply	100 - 240 V \pm 10% in four ranges 100 V: [90 - 121 V] 120 V: [104 - 139 V] 230 V: [198 - 242 V] 240 V: [207 - 264 V]		
Power line frequency	47-63 Hz		
Fuse	100 V, 120 V: 3.15 A (slow-slow) 230 V, 240 V: 1.6 A (slow-slow)	100 V, 120 V: 3.15 A (slow-slow) 230 V, 240 V: 1.25 A (slow-slow)	
Operating environment	0 °C to 40 °C, 80 % relative humidity without derating		
Storage environment	-10 °C to 60 °C		

16.2.7 μ Autolab Series instruments

The μ Autolab Series is the predecessor version of the Autolab Compact Series potentiostat/galvanostat produced by Metrohm Autolab *Autolab Compact Series (AUT4/AUT5) instruments* (see chapter 16.2.4, page 927). These instruments are identified by a serial number starting with **μ 2AUT7** (for the μ Autolab type II) or **μ 3AUT7** (for the μ Autolab type III).



NOTICE

The μ Autolab Series instruments are no longer available.

The following instruments belong to the μ Autolab Series:

- **μ Autolab type II:** compact potentiostat/galvanostat with 12 V compliance and 80 mA current. It is now replaced by the **PGSTAT101** or **PGSTAT204**.




Figure 1126 Overview of the front panel of the μ Autolab type II and type III

1 On/Off button

For switching the μ Autolab on or off.

2 Cell cable connector

For connecting the cell cable providing connections to the counter electrode (CE), reference electrode (RE), working electrode (WE) and ground.

3 CELL ENABLE button

For enabling and disabling the cell.

4 Display

Display indicating real-time information on the measured current and potential and instrumental settings.

The display (item 4 in Figure 1126) is used to provide information about the Autolab to the user. Figure 1127 shows a detail of this display.




Figure 1127 Overview of the display of the μ Autolab type II and type III

1 μ Autolab current ranges

The current range indicator which is lit corresponds to the active current of the μ Autolab.

2 PSTAT indicator

Indicates that the μ Autolab is operating in potentiostatic mode when lit.

3 FRA indicator

Indicates that the FRA2 module is in use when lit. This indicator is only available for the μ Autolab type III fitted with the **FRA2** module.

5 V ovl indicator

Indicates that a voltage overload is detected when lit.

7 HSTAB indicator

Indicates that the μ Autolab is operating in high stability mode when lit.

9 CELL ON indicator

Indicates that the cell is on when lit.

4 GSTAT indicator

Indicates that the μ Autolab is operating in galvanostatic mode when lit.

6 I ovl indicator

Indicates that a current overload is detected when lit.

8 CELL ENABLE indicator

Indicate the cell is enabled when lit.

16.2.7.4 μ Autolab Series back plane

The back plane of the μ Autolab type III provides a number of connections, shown in *Figure 1128*.




Figure 1128 Overview of the back plane of the μ Autolab type III

1 USB hub

For connecting additional USB devices.

2 USB connector

Type B USB plug for connecting the USB cable to the host computer.

3 DIO P1 connector

Digital input/output connector **P1** for sending and receiving external TTL triggers.

4 DIO P2 connector

Digital input/output connector **P2** for sending and receiving external TTL triggers.

5 Mains connection socket

For connecting the μ Autolab to the mains supply.

6 Fuse holder

Holds the mains connection socket fuse.

7 Fan

Required for cooling the μ Autolab during operation.

8 I out BNC connector

Connector providing the output of the current to voltage converter of the μ Autolab.

9 E out BNC connector

Connector providing the output of the voltage follower of the μ Autolab.

11 V in BNC connector

Connector providing the input of the on-board ADC of the μ Autolab.

10 V out BNC connector

Connector providing the output of the on-board DAC of the μ Autolab.

The back plane of the μ Autolab type II provides a number of connections, shown in (see figure 1129, page 967)




Figure 1129 Overview of the back plane of the μ Autolab type II

1 PC INTERFACE connector

For connecting the external USB interface adapter.

2 DIO P1 connector

Digital input/output connector **P1** for sending and receiving external TTL triggers.

3 DIO P2 connector

Digital input/output connector **P2** for sending and receiving external TTL triggers.

4 Mains connection socket

For connecting the μ Autolab to the mains supply.

5 Fuse holder

Holds the mains connection socket fuse.

6 Fan

Required for cooling the μ Autolab during operation.

7 V in BNC connector

Connector providing the input of the on-board ADC of the μ Autolab.

8 V out BNC connector

Connector providing the output of the on-board DAC of the μ Autolab.

9 -E out BNC connector

Connector providing the inverted output of the voltage follower of the μ Autolab.

10 I out BNC connector

Connector providing the output of the current to voltage converter of the μ Autolab.

16.2.7.5 Connections for analog signals

Four connectors, located on the back plane of the μ Autolab type II and μ Autolab type III, can be used as connections for analog signals *μ Autolab Series back plane* (see chapter 16.2.7.4, page 966). All the connections are provided through BNC connectors. All the connections are with respect to the μ Autolab ground directly and indirectly with respect to the protective earth.

16.2.7.7 μ Autolab Series testing

NOVA is shipped with a procedure which can be used, alongside the **Diagnostics** application, to verify that the instrument is working as expected.



NOTICE

For more information on the **Diagnostics** application, please refer to *Chapter 17*.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestCV** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestCV.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1130*.





Figure 1130 The data measured by the TestCV procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestCV automatic evaluation of the data requires the following tests to succeed:

1. The residual current, determined by subtracting the expected current from the measured current, must be smaller or equal than ± 20 nA.
2. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$.
3. The intercept of the measured current versus the applied potential must be equal to ± 4 mV divided by $1000100 \Omega \pm 5 \%$.

All three conditions must be valid for the test to succeed.

16.2.7.8 μ Autolab Series specifications

The specifications of the μ Autolab Series are provided in Table 32.

Table 32 Specifications of the μ Autolab Series instruments

Instrument	μ Autolab type II	μ Autolab type III
Maximum current	± 80 mA	
Compliance voltage	± 12 V	

Instrument	μAutolab type II	μAutolab type III
Potential range	± 5 V	
Applied potential accuracy	± 0.2 % ± 2 mV	
Applied potential resolution	150 μV	
Measured potential resolution	300 nV (gain 1000)	
Current ranges	10 nA to 10 mA, 7 decades	
Current accuracy	± 0.2 % of current range	
Applied current resolution	0.015 % of current range	
Measured current resolution	0.00003 % of current range (gain 1000)	
Potentiostat bandwidth	500 kHz	
Potentiostat rise/fall time	< 500 ns	
Input impedance of electrometer	> 100 GΩ, 8 pF	
Input bias current	< 1 pA	
Electrometer bandwidth	> 4 MHz	
Analog output	Potential and current	
Analog voltage input	No	
External inputs	1	
External outputs	1	
Digital input/output	48	



Instrument	μ Autolab type II	μ Autolab type III
Interface	External USB	Internal USB
Warm-up time	30 minutes	
Pollution degree	2	
Installation category	II	
External dimensions (without cables and accessories)	27x27x9 cm ³	
Weight	3.6 kg	3.6 kg or 4.4 kg with FRA2
Power requirements	75 W	144 W
Power supply	100 - 240 V \pm 10% in four ranges (auto select)	
Power line frequency	47-63 Hz	
Fuse	1.6 A (slow-slow)	
Operating environment	0 °C to 40 °C, 80 % relative humidity without derating	
Storage environment	-10 °C to 60 °C	

16.3 Module description

This chapter describes the extension modules available for the Autolab potentiostat/galvanostat instruments. The modules are grouped into two groups:

- **Common modules:** these modules are included standard in all Autolab systems.
- **Optional modules:** these internal or external optional modules can be installed in the Autolab or connected to the Autolab to extend the functionality of the instrument.



NOTICE

Some of the modules described in this chapter are no longer available but are still supported in the software. Whenever applicable, the successor module is specified.

16.3.1 Common modules

Common modules are always present in all Autolab instruments. These modules are either present as exchangeable modules or built-in modules. Depending on the type of Autolab system, these modules can be slightly different. However, the functionality these modules provide is common to all instruments.

The following common modules are available:

- **ADC164 or on-board ADC:** the *analog-to-digital* converter module of the Autolab. It is used to convert measured values into digital words that can be recorded by the host computer *ADC164 or on-board ADC* (see chapter 16.3.1.1, page 973).
- **DAC164 or on-board DAC:** the *digital-to-analog* converter module of the Autolab. It is used to convert digital words generated by the host computer into analog values that can be used to control the Autolab or external devices connected to the Autolab *DAC164 or on-board DAC* (see chapter 16.3.1.2, page 980).
- **DIO48 or DIO12:** the digital input/output module of the Autolab. This module can be used to send or receive TTL (*Transistor-Transistor Logic*) triggers in order to synchronize the Autolab control with external devices *TTL Triggers* (see chapter 16.3.1.3, page 986).

16.3.1.1 ADC164 or on-board ADC

The ADC164 or on-board ADC is the *analog-to-digital* converter used by the Autolab instrument to perform all analog control actions during measurements. The ADC164 or on-board ADC used by the Autolab is a multi-channel *analog-to-digital* converter. Each channel is fitted with a 16 bit converter, with an input range of ± 10 V.

The resolution of the ADC164 or the on-board ADC is given by:

$$\frac{20\text{ V}}{2^{16}} = 305.175\mu\text{V}$$

The ADC164 or on-board ADC is also fitted with a gain circuit with three settings: gain 1, gain 10 and gain 100. Each of these circuits divides the resolution by a factor equal to the gain. This means that in gain 100, the resolution is 3.05 μV .



The ADC164 provides two inputs for external signals, while the on-board ADC provides an input for a single external signal. One or two signals are provided in the Sampler (see figure 1131, page 974).




Figure 1131 The ADC164 module provides the External(1).External 1 and External(1).External 2 signals

In the case of the ADC164 the following signals are provided:

- External(1).External 1
- External(1).External 2

In the case of the on-board ADC, the following signal is provided:

- External(1).External 1



NOTICE

The names of the signals can be modified in the hardware setup *ADC164 and on-board ADC hardware setup* (see chapter 16.3.1.1.3, page 975).

16.3.1.1.1 ADC164 module front panel connections

The ADC164 module is fitted with two female BNC connectors, labeled →1 and →2 (see figure 1132, page 975).




Figure 1132 The front panel label of the ADC164 module

These two connectors provide inputs that can be used to record external signals. They have an input range of ± 10 V and an input impedance of 50Ω .

16.3.1.1.2 On-board ADC connections

Instruments that are fitted with an on-board ADC provide connections either through a dedicated connector located on the back plane of the instrument or through an optional monitor cable.

- **For the μ Autolab type II and μ Autolab type III:** the on-board ADC provides a connection, labeled **V in**, on the back plane of the instrument. Please refer to *Chapter 16.2.7.4* for more information.
- **For the PGSTAT101, PGSTAT204, M101 and M204:** the on-board ADC provides a connection, labeled **Vin**, through the optional monitor cable. Please refer to *Chapter 16.2.4.5.1*, *Chapter 16.2.4.5.2* and *Chapter 16.2.5.5* for more information.

16.3.1.1.3 ADC164 and on-board ADC hardware setup

To use the ADC164 or the on-board ADC for the measurement of external signals, the hardware setup needs to be adjusted. The External devices checkbox, provided in the Additional modules panel, adds the ADC164 or on-board ADC to the hardware setup for the purpose of recording external signals (*see figure 1133, page 976*).




Figure 1133 The External Devices module adds the ADC164 or on-board ADC to the hardware setup

For each available input, the following properties can be defined:

- **Signal name:** the name of the signal to record.
- **Signal unit:** the units of the signal to record.
- **Conversion slope:** the slope of the conversion function used to convert the signal.
- **Conversion offset:** the offset of the conversion function used to convert the signal.

Predefined settings are available using the drop-down list provided for the **Signal name** property (see figure 1134, page 977).

Properties

Vin

Signal name: Autolab Default (dropdown menu)
Signal unit: Autolab Default (dropdown menu)
Conversion slope: (input field)
Conversion offset: 0 V (input field)

Vout

Signal name: Autolab Default (dropdown menu)
Signal unit: V (input field)
Conversion slope: 1 V/V (input field)
Conversion offset: 0 V (input field)
Enable limits:
Upper limit: 10 V (input field)
Lower limit: -10 V (input field)

Figure 1134 Predefined settings are available

All the properties are automatically adjusted when one of the predefined setting is selected (see figure 1135, page 978).

Properties

Vin

Signal name TDI electronic load

Signal unit A

Conversion slope 30 A/V

Conversion offset 0 A

Vout

Signal name Autolab Default

Signal unit V

Conversion slope 1 V/V

Conversion offset 0 V

Enable limits

Upper limit 10 V

Lower limit -10 V

Figure 1136 Saving a new preset in the hardware setup



NOTICE

Once a new preset is saved, it can be reused with other instruments connected to the computer.



NOTICE

Clicking the button deletes the preset from the computer. It is not possible to delete predefined presets.

16.3.1.1.4 ADC164 and on-board ADC settings

The ADC164 and on-board ADC have no user-definable settings, except the Sampler settings *Sampler* (see chapter 9.1, page 616).

The Sampler settings define which signals are sampled during an electrochemical measurement.

The on-board ADC located in the Autolab PGSTAT101, PGSTAT204 and in the M101 and M204 modules of the Multi Autolab systems have an addi-

16.3.1.2.1 DAC164 module front panel connections

The DAC164 module is fitted with two female BNC connectors, labeled ←1 and ←2 (see figure 1138, page 981).




Figure 1138 The front panel label of the DAC164 module

These two connectors provide outputs that can be used to generate signals suitable for controlling external devices. They have an output range of ± 10 V and an output impedance of 50Ω .



CAUTION

DAC164 ←2 is reserved for AC voltammetry measurements. This means that this output cannot be used to control an external device, unless the instrument is modified.

16.3.1.2.2 On-board DAC connections

Instruments that are fitted with an on-board DAC provide connections either through a dedicated connector located on the back plane of the instrument or through an optional monitor cable.

- **For the μ Autolab type II and μ Autolab type III:** the on-board DAC provides a connection, labeled **V out**, on the back plane of the instrument. Please refer to *Chapter 16.2.7.4* for more information.
- **For the PGSTAT101, PGSTAT204, M101 and M204:** the on-board DAC provides a connection, labeled **Vout**, through the optional monitor cable. Please refer to *Chapter 16.2.4.5.1*, *Chapter 16.2.4.5.2* and *Chapter 16.2.5.5* for more information.

16.3.1.2.3 DAC164 and on-board DAC hardware setup

To use the DAC164 or the on-board DAC for the generating signals to control external devices, the hardware setup needs to be adjusted. The External devices checkbox, provided in the Additional modules panel, adds the DAC164 or on-board DAC to the hardware setup for the purpose of generating external signals (see figure 1139, page 982).





Figure 1139 The External Devices module adds the DAC164 or on-board DAC to the hardware setup

For each available output, the following properties can be defined:

- **Signal name:** the name of the signal to generate.
- **Signal unit:** the units of the signal to generate.
- **Conversion slope:** the slope of the conversion function used to generate the signal.
- **Conversion offset:** the offset of the conversion function used to generate the signal.
- **Enable limits:** a toggle that is provided to enable or disable limits for the generated signal.
- **Upper limit:** the upper limit for the generated signal. This limit is only used if the *Enable limits* property is set to on.
- **Lower limit:** the lower limit for the generated signal. This limit is only used if the *Enable limits* property is set to on.

Predefined settings are available using the drop-down list provided for the **Signal name** property (see figure 1140, page 983).



Properties

Vin ⏏ ✕

Signal name Autolab Default ▾

Signal unit V

Conversion slope 1 V/V

Conversion offset 0 V

Vout ⏏ ✕

Signal name Autolab Default ▾

Signal unit Autolab Default

Conversion slope Autolab LED Driver

Conversion offset Autolab R(R)DE

Conversion offset 0 V

Enable limits

Upper limit 10 V

Lower limit -10 V

Figure 1140 Predefined settings are available

All the properties are automatically adjusted when one of the predefined setting is selected (see figure 1141, page 984).

Properties

Vin

Signal name Autolab Default

Signal unit V

Conversion slope 1 V/V

Conversion offset 0 V

Vout

Signal name TDI electronic load

Signal unit A

Conversion slope 30 A/V

Conversion offset 0 A

Enable limits

Upper limit 140 A

Lower limit 0 A

Save preset

Figure 1142 Saving a new preset in the hardware setup



NOTICE

Once a new preset is saved, it can be reused with other instruments connected to the computer.



NOTICE

Clicking the button deletes the preset from the computer. It is not possible to delete predefined presets.

16.3.1.2.4 DAC164 and on-board DAC settings

The DAC164 or on-board DAC module settings are completely defined in the NOVA software. The following user-definable settings are available, through the **Autolab control** command (see figure 1143, page 986):

- **DAC164 ←1/Vout:** this setting defines the output of the DAC164 output 1 (DAC164 ←1) or on-board DAC (Vout), unconverted, as a voltage in the ± 10 V range.

The Autolab is able to send and receive triggers, using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236) and the **Wait** command *Wait for DIO* (see chapter 7.2.4.2, page 242), respectively.

Every Autolab instrument is equipped with one or two digital input/output connectors (DIO) that can be used to receive or send a digital TTL trigger. Depending on the instrument type, two different connector layouts are available:

- **For all other Autolab instruments:** two programmable, 25 pin SUB-D connectors located on the front panel or the back plane of the instrument are available for TTL triggering. Both connectors are identified as a **DIO48**.
- **For the PGSTAT101 or M101 module and the PGSTAT204 and M204 module:** a single, female, 15 pin SUB-D connector located on the front panel or the back plane of the instrument or module is available for TTL triggering. This connector is identified as **DIO12**.



CAUTION

There is a chance of introducing a ground loop when connecting external devices to the Autolab DIO. This can result in higher than expected noise levels during measurements. It is recommended to disconnect external devices from the DIO connector(s) of the Autolab when TTL triggering is not required.



CAUTION

Although the Autolab **PGSTAT302F** is fitted with two DIO ports on the back plane, these ports **cannot** be used for TTL triggering.

16.3.1.3.1 DIO48 type connectors

The DIO48 connectors for TTL triggering consist of two, 25 pin, female SUB-D connectors. Each connector has a total of 24 user-addressable input/output pins, grouped in three sections:

- **Section A:** pins 1 to 8.
- **Section B:** pins 17 to 14.
- **Section C:** pins 9 to 16.




Figure 1145 The DIO12 connector layout

The pin layout is detailed in Table 33.

Table 33 Inputs and outputs of the DIO12 connector

Assignment	Pin number	Section
Input 1	1	B
Input 2	9	B
Input 3	2	B
Input 4	10	B
Output 1	12	A
Output 2	5	A
Output 3	13	A
Output 4	6	A
Output 5	14	A
Output 6	7	A
Output 7	15	A
Output 8	8	A
Digital ground	4	
Digital ground	11	
Isolated ground	3	



CAUTION

The write lines of the DIO12 connector are capable of supplying a maximum current of **200 mA**. Suitable pull-down resistors should be placed in the write lines of the DIO12. A typical value for the pull-down resistance is about 1 k Ω . Please refer to the user manual of the external device connected to the instrument for more information.

16.3.1.4 Dummy cell

For testing purposes and for illustration purposes, the Autolab potentiostat/galvanostat systems are supplied with a **dummy cell**. These cells can be used to carry out specially designed tests or diagnostics of the instrument, as explained in *Chapter 17* and can be used to perform measurements on a known circuit.

Depending on the type of instrument, the following dummy cells are available:

- **Dummy cell 2:** all Autolab potentiostat/galvanostat systems are supplied with the standard Autolab dummy cell *Autolab Dummy cell 2* (see chapter 16.3.1.4.1, page 991).
- **Internal dummy cell:** all the PGSTAT101, M101, PGSTAT204 and M204 Autolab systems are supplied with an internal dummy cell *Internal dummy cell* (see chapter 16.3.1.4.2, page 993).
- **Option ECI10M dummy cell:** this dummy cell is an optional item that can be used in combination with the ECI10M module *ECI10M optional dummy cell* (see chapter 16.3.1.4.3, page 995).
- **Booster10A test cell:** all BOOSTER10A systems are supplied with a dedicated high power test cell *Booster10A test cell* (see chapter 16.3.1.4.4, page 996).
- **Booster20A test cell:** all BOOSTER20A systems are supplied with a dedicated high power test cell *Booster20A test cell* (see chapter 16.3.1.4.5, page 997).



CAUTION

All the **dummy cells** supplied with or available for the Autolab instruments are uncalibrated. These cells **cannot** be used to verify that the instrument is reaching all of the specifications. These cells should only be used to carry out qualitative measurements unless otherwise specified in this manual.

16.3.1.4.1 Autolab Dummy cell 2

The Autolab dummy cell 2 is the standard dummy cell, supplied with all instruments except the PGSTAT101. The Dummy cell is shown in *Figure 1148*, schematically.




Figure 1148 The Autolab dummy cell 2

This dummy cell is fitted with five circuits, consisting of resistors and capacitors. The actual values and the tolerances of these components are shown in Figure 1149.




Figure 1149 Component values and tolerances used in the Autolab dummy cell 2

All resistors have a tolerance of 5% and all capacitors have a tolerance of 10%.



CAUTION

The **Dummy cell 2** is not calibrated. This cell **cannot** be used to verify that the instrument is reaching all of the specifications. This cell should only be used to carry out qualitative measurements unless otherwise specified in this manual.

16.3.1.4.2 Internal dummy cell

The Autolab PGSTAT101, M101, PGSTAT204 and M204 are all fitted with an internal dummy cell. This cell is built inside the instrument and cannot be removed from the instrument. To use the internal dummy cell, it is necessary to connect the cell cable to the instrument and to shorten the electrode connectors as shown in (see figure 1150, page 993).




Figure 1150 The electrode connections used in combination with the internal dummy cell

The internal dummy cell can be activated using the dedicated switch available through the **Autolab control** command (see figure 1151, page 994).



Figure 1151 The internal dummy cell can be activated using the dedicated toggle in the Autolab control command

This dummy cell is fitted with a single circuit, consisting of two resistors and one capacitor. The actual values and the tolerances of these components are shown in Figure 1152.




Figure 1152 Component values and tolerances used in the internal dummy cell

Both resistors have a tolerance of 5% and the capacitor have a tolerance of 10%.



CAUTION

The **Internal dummy cell** is not calibrated. This cell **cannot** be used to verify that the instrument is reaching all of the specifications. This cell should only be used to carry out qualitative measurements unless otherwise specified in this manual.



16.3.1.4.3 ECI10M optional dummy cell

The ECI10M test cell is an optional dummy cell, designed for measurements in combination with the ECI10M module. The Dummy cell is shown in *Figure 1153*, schematically.




Figure 1153 The ECI10M test cell



NOTICE

This dummy cell must be directly connected to the front panel of the ECI10M external interface!

This dummy cell is fitted with a single circuit, consisting of resistors and capacitors. The actual values and the tolerances of these components are shown in *Figure 1154*.




Figure 1154 Component values and tolerances used in the ECI10M test cell

All resistors have a tolerance of 5% and all capacitors have a tolerance of 10%.



CAUTION

The **ECI10M test cell** is not calibrated. This cell **cannot** be used to verify that the instrument is reaching all of the specifications. This cell should only be used to carry out qualitative measurements unless otherwise specified in this manual.

16.3.1.4.5 Booster20A test cell

The BOOSTER20A systems are supplied with a high power test cell. This cell is mounted on a heat sink to dissipate heat while the cell is used (see figure 1156, page 997).




Figure 1156 The Booster20A test cell



NOTICE

More information on the BOOSTER20A is available in *Chapter 16.3.2.6*.

The dummy cell supplied with the BOOSTER20A is fitted with a single resistor of 50 mΩ. This resistor has a tolerance of 5% and can dissipate up to 50 W of power.



CAUTION

The **Booster20A test cell** is not calibrated. This cell **cannot** be used to verify that the instrument is reaching all of the specifications. This cell should only be used to carry out qualitative measurements unless otherwise specified in this manual.

16.3.2.1.2 ADC10M scope of delivery

The ADC10M module is supplied with the following items:

- ADC10M module
- ADC10M module label

16.3.2.1.3 ADC10M hardware setup

To use the **ADC10M** module, the hardware setup needs to be adjusted. The checkbox for the module needs to be ticked (*see figure 1157, page 999*).




Figure 1157 The ADC10M module is selected in the hardware setup

16.3.2.1.4 ADC10M module settings

The ADC10M module can be used in combination with the Cyclic voltammetry linear scan and the Chrono methods command.

In the case of the Cyclic voltammetry linear scan command, all the module settings are automatically controlled by the measurement command.

- **Channel 2**
 - **Measure external:** specifies the input signal for Channel 2 of the ADC10M module, using the provided  toggle. When this toggle is off, the WE(1).Current signal is measured through Channel 2 of the ADC10M. When this toggle is on, the signal provided on the →2 BNC input on the front panel of the ADC10M module is sampled *ADC10M module front panel connections (see chapter 16.3.2.1.6, page 1002)*.
 - **Gain:** specifies the amplification gain for the signal measured by Channel 2 of the ADC10M module, using the drop-down list. The default gain is 1 and optional gains 5, 10 and 20 are available.
 - **Filter:** specifies if a filter should be applied on the signal measured on Channel 2 of the ADC10M, using the provided  toggle. When this filter is on, the bandwidth of Channel 2 of the ADC10M is reduced to 200 kHz.
- **Other**
 - **High bandwidth:** specifies if the high bandwidth mode of the ADC10M should be used, using the provided  toggle. When this setting is off, the bandwidth of both ADC10M channels is set to 600 kHz. When this setting is on, the bandwidth of both ADC10M channels is increased to 1.2 MHz.



NOTICE

Select the gain carefully to avoid exceeding the measurable range of the ADC10M.




NOTICE

The Filters provided for Channel 1 and Channel 2 can be used to overrule the settings defined by the High bandwidth toggle.

16.3.2.1.5 ADC10M module restrictions

Restrictions apply when using the ADC10M module:

- **No real-time data display:** the ADC10M is fitted with an on-board memory that can be used to store up to 1,024,000 data points. When the ADC10M module is used in an experiment, each new data point is stored in the on-board memory of the module until the experiment is finished. At the end of the measurement, all the stored data points are transferred to the computer for data analysis.



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test uses a high-speed chrono methods measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1160*.





Figure 1160 The results of the TestADC procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestADC automatic evaluation of the data requires the following tests to succeed:

Two versions of the ADC750 are available: the ADC750 (revision 5) and the ADC750 (revision 4). The former is identified as ADC750 in the hardware setup while the latter is identified as ADC750r4.

The module revision is not indicated on the front panel of the instrument. It is therefore recommended to declare the ADC750 module in the hardware setup (see figure 1161, page 1005).




Figure 1161 Declaring the ADC750 in the hardware setup

If an error message is shown after adjusting the hardware setup, then the ADC750 is a revision 4 module. The hardware setup should then be adjusted accordingly (see figure 1162, page 1006).




Figure 1163 Additional settings of the ADC750 are available in the Chrono methods command

The following settings are available:

- **Channel 1**

- **Measure external:** specifies the input signal for Channel 1 of the ADC750 module, using the provided toggle. When this toggle is off, the WE(1).Potential signal is measured through Channel 1 of the ADC750. When this toggle is on, the signal provided on the →1 BNC input on the front panel of the ADC750 module is sampled *ADC750 module front panel connections (see chapter 16.3.2.2.5, page 1008)*.
- **Gain:** specifies the amplification gain for the signal measured by Channel 1 of the ADC750 module, using the drop-down list. The default gain is 1 and optional gains 10 and 100 are available.

- **Channel 2**

- **Measure external:** specifies the input signal for Channel 2 of the ADC750 module, using the provided toggle. When this toggle is off, the WE(1).Current signal is measured through Channel 2 of the ADC750. When this toggle is on, the signal provided on the →2 BNC input on the front panel of the ADC750 module is sampled *ADC750 module front panel connections (see chapter 16.3.2.2.5, page 1008)*.
- **Gain:** specifies the amplification gain for the signal measured by Channel 2 of the ADC750 module, using the drop-down list. The default gain is 1 and optional gains 10 and 100 are available.

16.3.2.2.6 ADC750 module testing

NOVA is shipped with a procedure which can be used to verify that the **ADC750** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestADC** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestADC.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (c).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test uses a high-speed chrono methods measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1165*.





Figure 1165 The results of the TestADC procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestADC automatic evaluation of the data requires the following tests to succeed:

1. The applied potential in Step 1 of the measurement must be $0\text{ V} \pm 10\text{ mV}$.
2. The applied potential in Step 2 of the measurement must be $0.3\text{ V} \pm 10\text{ mV}$.
3. The applied potential in Step 3 of the measurement must be $-0.3\text{ V} \pm 10\text{ mV}$.
4. The applied potential in Step 4 of the measurement must be $0\text{ V} \pm 10\text{ mV}$.

All four conditions must be valid for the test to succeed.

16.3.2.2.7 ADC750 module specifications

The specifications of the ADC750 module are provided in *Table 35*.

Table 35 Specifications of the ADC750 module

Specification	Value
Number of channels	2
Maximum sampling rate	750,000 samples/second
Shortest interval time	1.33 μs

Specification	Value
ADC resolution	12 bit
Maximum resolution, potential	500 μ V (gain 10)
Maximum resolution, current	0.0025 % of current range (gain 10)
Maximum number of points	512,000
Input range	± 10 V
Input impedance	≥ 5 G Ω

16.3.2.3 BA module

The BA module is an extension module for the Autolab PGSTAT and the Multi Autolab. This module provides a second working electrode, **WE(2)**. The BA module has two different operation modes:

- **BIPOT mode:** in this mode, the potential of the second working electrode, WE(2), is defined with respect to the common reference electrode.
- **Scanning BIPOT:** in this mode, the potential of the second working electrode, WE(2), is defined with respect to the main working electrode, WE(1).



NOTICE

The BA module only works in potentiostatic mode. The main potentiostat can be set to galvanostatic mode.

The BA module adds the following signal to the Sampler (*see figure 1166, page 1012*):

- **WE(2).Current (A):** this signal corresponds to the current flowing through the second working electrode, WE(2).
- **WE(2).Charge (C):** this signal corresponds to the charge obtained by numerical integration of the measured WE(2).Current.




Figure 1167 The BA module is selected in the hardware setup

The toggles provided in the **Properties** panel can be used to specify the number of **BA** modules installed in the instrument.



NOTICE

It is only possible to have more than one **BA** module in combination with **MBA** instruments *Autolab MBA N Series (AUT8) instruments* (see chapter 16.2.3, page 925).

16.3.2.3.4 BA module settings

The BA module settings are completely defined in the NOVA software. The following user-definable settings are available, through the **Autolab control** command (see figure 1168, page 1014):

- **Cell:** a  toggle control that can be used to switch the WE(2) on or off.



CAUTION

When the electrode control is set to independent, a specific order must be respected to avoid current leakage between WE(1) and WE(2). Always set the cell of WE(2) to ON after the cell has been set to ON for WE(1). Always set the cell of WE(2) to OFF before the cell has been set to OFF for WE(1).

16.3.2.3.5 BA module restrictions

No restrictions apply when using the BA module.

16.3.2.3.6 BA module front panel connections

The BA module is fitted with a single female BNC connector, labeled ←I (in the case of the Autolab potentiostat/galvanostat instrument) or with two female BNC connectors, labeled ←I and WE, respectively (*see figure 1169, page 1015*).




Figure 1169 The front panel labels of the BA module (left: BA module in PGSTAT, right: BA module in Multi Autolab)



NOTICE

The WE connector provided in the case of the Multi Autolab instruments is used to connect the WE2 cable.

The signal provided through the ←I connector on the front panel corresponds to the output of the current-to-voltage converter located on the



BA module. The output signal is a voltage, referred to the instrument ground, corresponding to the converted current according to:

$$E_{out}(\leftarrow I) = \frac{i_{(WE2)}}{[CR]}$$

Where $E_{out}(\leftarrow I)$ corresponds to the output voltage signal of the module, in V, $i_{(WE2)}$ corresponds to the current measured by the BA module, in BA and [CR] is the active current range of the BA module.



NOTICE

The front panel $\leftarrow I$ BNC output is provided for information purposes only.

16.3.2.3.7 BA module testing

NOVA is shipped with a procedure which can be used to verify that the **BA** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestBA** procedure, provided in the NOVA 2.X installation folder (Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestBA.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a) and the second working electrode (BA) to the Autolab dummy cell (b).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. Both modes of the BA module are tested. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1170*.





Figure 1170 The results of the TestBA procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestBA automatic evaluation of the data requires the following tests to succeed:

1. The average WE(2).Current measured in Bipotentiostat mode must be equal to $1 \text{ V} \pm 5 \text{ mV}/1000100 \text{ } \Omega \pm 5 \%$.
2. The intercept of the WE(2).Current measured in Bipotentiostat mode must be equal to $1 \text{ V} \pm 5 \text{ mV}/1000100 \text{ } \Omega \pm 5 \%$.
3. The inverted slope of the measured WE(2).Current versus the applied potential must be equal to $1000100 \text{ } \Omega \pm 5 \%$.

- **WE(2).Current (A):** this signal corresponds to the current flowing through the second working electrode, WE(2).
- **WE(2).Charge (C):** this signal corresponds to the charge obtained by numerical integration of the measured WE(2).Current.




Figure 1171 The BIPOT/ARRAY module provides the WE(2).Charge and WE(2).Current signals

16.3.2.4.1 BIPOT/ARRAY module compatibility

The BIPOT/ARRAY module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- PGSTAT100N and PGSTAT100
- PGSTAT20
- PGSTAT10



NOTICE

The BIPOT/ARRAY module is **not** compatible with the Autolab instruments not listed above.

16.3.2.4.2 BIPOT/ARRAY module scope of delivery

The BIPOT/ARRAY module is supplied with the following items:

- BIPOT or ARRAY module
- BIPOT or ARRAY module label
- WE2 connection cable

- **Electrode control:** defines the relationship between the cell switch of WE(2) and WE(1). When this setting is set to Linked to WE(1), the cell switch of WE(2) will automatically be set to the same status of the cell switch of WE(1). Using the Cell command, both cell switches will be toggled at the same time. When this setting is set to Independent, the cell switch of WE(2) is decoupled from the cell switch of WE(1). In that case, the cell switch of the WE(2) must be set manually, using the control Cell control provided by the Autolab control command. This setting only affects the transition from cell off to cell on.
- **Potential (V):** defines the potential difference between WE(2) and the reference electrode for the BIPOT module or between WE(2) and WE(1) in for the ARRAY mode.




Figure 1173 The BIPOT/ARRAY module settings are defined in the Autolab control command



NOTICE

The default startup Electrode control setting is linked to WE(1).



NOTICE

The Cell command switches the main potentiostat and the second working electrode off automatically, regardless of the Electrode control setting.



NOTICE

The front panel ←I BNC output is provided for information purposes only.

16.3.2.4.7 BIPOT/ARRAY module testing

Two test procedures are provided for testing the BIPOT/ARRAY module:

- For the BIPOT module, please refer to *Chapter 16.3.2.4.7.1*.
- For the ARRAY module, please refer to *Chapter 16.3.2.4.7.2*.

16.3.2.4.7.1 BIPOT module testing

NOVA is shipped with a procedure which can be used to verify that the **BIPOT** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestBIPOT** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestBIPOT.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a) and the second working electrode (BIPOT) to the Autolab dummy cell (b).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1175*.





Figure 1175 The results of the TestBIPOT procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestBIPOT automatic evaluation of the data requires the following tests to succeed:

1. The average WE(2).Current measured in Bipotentiostat mode must be equal to $1 \text{ V} \pm 5 \text{ mV}/1000100 \Omega \pm 5 \%$.
2. The intercept of the WE(2).Current measured in Bipotentiostat mode must be equal to $1 \text{ V} \pm 5 \text{ mV}/1000100 \Omega \pm 5 \%$.

Both conditions must be valid for the test to succeed.

16.3.2.4.7.2 ARRAY module testing

NOVA is shipped with a procedure which can be used to verify that the **ARRAY** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestARRAY** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestARRAY.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a) and the second working electrode (ARRAY) to the Autolab dummy cell (b).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1176*.

16.3.2.4.8 BIPOT/ARRAY module specifications

The specifications of the BIPOT module are provided in *Table 37*.

Table 37 Specifications of the BIPOT module

Specification	Value
Operation mode	BIPOT (hardware defined)
Control DAC	DAC164
Maximum current	± 35 mA
Current ranges	100 nA to 10 mA (6 ranges)
Current accuracy	± 0.2 % of current range
Potential range	± 5 V
Potential accuracy	± 2 mV

The specifications of the ARRAY module are provided in *Table 38*.

Table 38 Specifications of the ARRAY module

Specification	Value
Operation mode	Scanning BIPOT (hardware defined)
Control DAC	DAC164
Maximum current	± 35 mA
Current ranges	100 nA to 10 mA (6 ranges)
Current accuracy	± 0.2 % of current range
Potential range	± 5 V
Potential accuracy	± 2 mV

16.3.2.5 Booster10A

The Booster10A is an **external** extension module for the Autolab potentiostat/galvanostat. This module extends the maximum current of the Autolab system to which they are connected to 10 A. The Booster10A can be used with all the measurement commands provided in NOVA.

16.3.2.5.1 Booster10A module compatibility

The Booster10A is compatible with the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N
- PGSTAT100N and PGSTAT100
- PGSTAT204 and M204
- PGSTAT20

16.3.2.5.3 Booster10A module settings

The Booster10A extends the available current ranges of the controlling instrument by adding a single 10 A current range. The module does not provide specific settings. The Booster10A can be used in two different modes, when connected to the Autolab:

- **Bypass mode:** in this mode, the Booster10A is connected to the Autolab but the extra current range provided by the Booster10A is not used. The Booster10A is bypassed and only provides connections to the electrochemical cell.
- **Operation mode:** in this mode, the extra current range provided by the Booster10A is used.

The mode of operation is controlled by the active current range. The user can specify which current range is used during a measurement or at any time by using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236) or the **Autolab display** panel *Autolab display panel* (see chapter 5.2.3, page 127).

Figure 1178 shows the additional current range provided in the Autolab control command.




Figure 1178 The Booster10A is controlled by the current range

7 Cell on indicator LED
Indicates that the cell is on when lit.

8 Manual cell On/Off indicator LED
The LED is lit when the cell is enabled.

9 Cell On/Off switch
For manually enabling or disabling the cell.

16.3.2.5.6 Booster10A module back plane connections

The back plane of the Booster10A provides a number of connections, shown in *Figure 1180*.




Figure 1180 The back plane of the Booster10A

1 Air flow holes
Required for cooling the Booster10A during operation.

2 Mains voltage indicator
Indicates the mains voltage settings of the Booster10A.

3 Mains connection socket
For connecting the Booster10A to the mains supply.

4 GND plug
For grounding purposes.

5 DIO connector (digital control)
For connecting the digital control cable to interface with the Autolab PGSTAT.

6 Fan
Required for cooling the Booster10A during operation.



CAUTION

Make sure that the mains voltage indicator is set properly before switching the Booster10A on.

16.3.2.5.7 Booster10A installation and configuration

The Booster10A can be used in combination with any compatible instrument. The installation and configuration can be carried out by the end-user at any time.

Depending on the type of instrument it is connected to, the Booster10A has to be installed and configured according to a specific procedure:

1. For all the compatible Autolab instruments, except the PGSTAT204 and the M204, please refer to *Chapter 16.3.2.5.7.1*.

Hardware setup

Autolab module

Main module: PGSTAT302N

Power supply frequency: 50 Hz

C1: 1.6E-11

C2: 3E-13

Automatic configuration

Additional modules

- FRA32M
- FRA2
- ADC10M
- ADC750
- ADC750r4
- SCAN250
- SCANGEN
- BA
- BIPO/ARRAY
- ECD
- FI20 - Filter
- FI20 - Integrator
- Booster20A
- Booster10A**
- EQCM
- pX1000
- pX
- ECN
- External Devices
- IME303
- IME663
- MUX

Properties

DIO connector: P1

5 Connect the cell

Use the CE and WE connectors provided by the Booster10A and the RE and S connectors provided by the Autolab to the electrochemical cell.

16.3.2.5.7.2 **Booster10A installation and configuration (PGSTAT204 and M204 only)**

For the Autolab PGSTAT204 and the M204 module, a special set of cables is required for connecting the Booster10A. The set of cables includes two cables:

1. **DIO adapter cable:** a female 25 pin SUB-D to male 15 pin SUB-D adapter cable.
2. **Cell adapter cable:** a dedicated cell cable assembly providing an interface between the PGSTAT204 or M204 module and the Booster10A (see figure 1181, page 1034).




Figure 1181 The cell adapter cable assembly

1 Booster10A PGSTAT cable connector

For connecting the PGSTAT cable from the Booster10A (item 6 in *Figure 1179, page 1030*).

2 PGSTAT204/M204 cell cable connector

For connecting to the PGSTAT204 or M204 module cell connector (item 4 in *Figure 1105, page 930* or item 5 in *Figure 1114, page 944*).

3 RE/S cable

Cable providing reference electrode (RE) and sense electrode (S) to connect to the electro-chemical cell.

The following steps describe how to install and configure the Booster10A for the PGSTAT204 and the M204 module.

1 Remove the cell cable

Unscrew and remove the cell cable from the PGSTAT204 front panel or the M204 module front panel.



NOTICE

It is recommended to store this cable carefully for future use.

2 Connect the Booster10A PGSTAT cable to the Booster10A PGSTAT cable connector

Connect the PGSTAT cable, located on the front panel of the Booster10A (item 6 in *Figure 1179, page 1030*) to the Booster10A PGSTAT cable connector (item 1 in *Figure 1181*).

3 Connect the DIO adapter cable

Connect the DIO adapter cable, supplied with the Booster10A, to the digital control connector located on the backplane of the Booster10A (item 5 in *Figure 1180, page 1031*). Connect the other end of the cable to the DIO connector located on the front panel of the Autolab PGSTAT204 or the M204 module (item 1 in *Figure 1105, page 930* or item 2 in *Figure 1114, page 944*).

4 Connect the PGSTAT204/M204 cell cable connector to the PGSTAT204 or M204

Connect the PGSTAT204/M204 cell cable connector (item 2 in *Figure 1181*) to the cell cable connector located on the front panel of the PGSTAT204 or M204 (item 4 in *Figure 1105, page 930* or item 5 in *Figure 1114, page 944*).

5 Specify the hardware setup

Adjust the hardware setup.

The screenshot displays the 'Hardware setup' window. It is divided into three main sections: 'Autolab module', 'Additional modules', and 'Properties'.

- Autolab module:** Contains two dropdown menus. The first is labeled 'Main module' and is set to 'PGSTAT204'. The second is labeled 'Power supply frequency' and is set to '50 Hz'. Below these is a button labeled 'Automatic configuration'.
- Additional modules:** A list of modules with checkboxes. The 'Booster10A' checkbox is checked and highlighted in green. Other modules listed include FRA32M, BA, FI20 - Integrator, EQCM, pX1000, External Devices, IME303, IME663, and MUX.
- Properties:** A section on the right side of the window, currently empty.

will be shown. The measured data should look as shown in *Figure 1182*.





Figure 1182 The results of the TestBooster10A procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestBooster10A automatic evaluation of the data requires the following tests to succeed:

1. The inverted slope of the measured current versus the applied potential must be equal to $100 \text{ m}\Omega \pm 5 \%$.
2. The intercept of the measured current versus the applied potential must be equal to $\pm 4 \text{ mV}$ divided by $100 \text{ m}\Omega \pm 5 \%$.

Both conditions must be valid for the test to succeed.




Figure 1183 The 50 mΩ dummy cell supplied with the Booster20A

16.3.2.6.3 Booster20A module settings

The Booster20A extends the available current ranges of the controlling instrument by adding a single 20 A current range. The module does not provide specific settings. The Booster20A can be used in two different modes, when connected to the Autolab:

- **Bypass mode:** in this mode, the Booster20A is connected to the Autolab but the extra current range provided by the Booster20A is not used. The Booster20A is bypassed and only provides connections to the electrochemical cell.
- **Operation mode:** in this mode, the extra current range provided by the Booster20A is used.

The mode of operation is controlled by the active current range. The user can specify which current range is used during a measurement or at any time by using the **Autolab control** command *Autolab control* (see chapter 7.2.1, page 236) or the **Autolab display** panel *Autolab display panel* (see chapter 5.2.3, page 127).

Figure 1184 shows the additional current range provided in the Autolab control command.




Figure 1185 Overview of the front panel of the Booster20A

- | | |
|---|---|
| <p>1 On/Off button
For switching the Booster on or off.</p> | <p>2 CE (Counter electrode) cable connector
For connecting the CE cable.</p> |
| <p>3 WE (Working electrode) cable connector
For connecting the WE cable.</p> | <p>4 To PGSTAT connector
For connecting the analog control cable between the Booster and the Autolab PGSTAT.</p> |
| <p>5 Emergency stop connector
For connecting the emergency stop button.</p> | <p>6 Cell enable button
For enabling or disabling the cell.</p> |
| <p>7 Display
For indications and warnings.</p> | |

The display (item 7 in *Figure 1185*) is used to provide information about the Booster20A to the user. *Figure 1186* shows a detail of this display.




Figure 1186 Overview of the display of the Booster20A

- | | |
|---|---|
| <p>1 I ovl LED
Indicates that a current overload is detected when lit.</p> | <p>2 V ovl LED
Indicates that a voltage overload is detected when lit.</p> |
|---|---|



3 T ovl LED

Indicates that a temperature overload is detected when lit.

5 PGSTAT MODE LED

Indicates that the Booster20A is in **Bypass** mode when lit.

7 CELL ON

Indicates that the cell is on when lit.

4 SYS FAIL and UNIT FAIL LEDs

Indicates that the Booster20A is malfunctioning when lit.

6 UNIT ENABLE

Indicates that the Booster20A is in operation when lit.

16.3.2.6.6 Emergency stop button

The Booster20A is supplied with an emergency stop button (see figure 1187, page 1042).



Figure 1187 The emergency switch of the Booster20A

The emergency stop button is connected to the front panel of the Booster20A (item 5 in Figure 1185).

The emergency stop button can be pressed at any time to immediately disconnect the Booster20A from the electrochemical cell. The stop button remains engaged until it is disengaged by the user. To disengage the emergency stop button, rotate the red knob counter-clockwise until it releases.



NOTICE

Pressing the emergency stop button does not stop the NOVA measurement.

16.3.2.6.7 Booster20A module back plane connections

The back plane of the Booster20A provides a number of connections, shown in *Figure 1188*.




Figure 1188 The back plane of the Booster20A

1 DIO connector (TO AUTOLAB)

For connecting the digital control cable to interface with the Autolab PGSTAT.

2 Fans

Required for cooling the Booster20A during operation.

3 Fuse holders

Fuse holders containing the fuses protecting the Booster20A.

4 Mains connection socket

For connecting the Booster10A to the mains supply.

5 Earth plug

For connections to the protective earth.

6 GND plug

For connections to the ground.

16.3.2.6.8 Booster20A installation and configuration

The Booster20A can be used in combination with any compatible instrument. The installation and configuration can be carried out by the end-user at any time. The following steps describe how to install and configure the Booster20A.

1 Remove the CE/WE cable

Unscrew and remove the CE/WE cable from the PGSTAT front panel.



NOTICE

It is recommended to store this cable carefully for future use.

6 Connect the WE cable

Connect the WE cable to the dedicated connector on the front panel of the Booster20A (item 3 in *Figure 1185*).

7 Connect the CE cable

Connect the CE cable to the dedicated connector on the front panel of the Booster20A (item 2 in *Figure 1185*).

8 Connect the cell

Use the CE and WE connectors provided by the Booster20A and the RE and S connectors provided by the Autolab to the electrochemical cell.

16.3.2.6.9 Booster20A testing

NOVA is shipped with a procedure which can be used to verify that the **Booster20A** is working as expected.

Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestBooster20A** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestBoosert20A.nox)

2 Connect the Autolab dummy cell

Connect the Autolab and Booster20A to the Booster20A test cell.





Figure 1189 The results of the TestBooster20A procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestBooster20A automatic evaluation of the data requires the following tests to succeed:

1. The inverted slope of the measured current versus the applied potential must be equal to $50 \text{ m}\Omega \pm 5 \%$.
2. The intercept of the measured current versus the applied potential must be equal to $\pm 4 \text{ mV}$ divided by $50 \text{ m}\Omega \pm 5 \%$.

Both conditions must be valid for the test to succeed.

16.3.2.7.2 ECD module scope of delivery

The ECD module is supplied with the following items:

- ECD module
- ECD module label

16.3.2.7.3 ECD hardware setup

To use the **ECD** module, the hardware setup needs to be adjusted. The checkbox for the module needs to be ticked (*see figure 1190, page 1049*).




Figure 1190 The ECD module is selected in the hardware setup

16.3.2.7.4 ECD module settings

The ECD module settings are completely defined in the NOVA software. The following user-definable settings are available, through the **Autolab control** command (*see figure 1191, page 1050*):

- **ECD:** a toggle which can be used to switch the ECD module on or off.




Figure 1192 The current ranges provided by the ECD module can be selected in the Autolab control command

16.3.2.7.5 ECD module restrictions

Restrictions apply when using the ECD module. The offset compensation of the ECD module is controlled directly from the DAC164 located in the instrument. This forces the following restrictions:

- **BIPOT/ARRAY module:** the BIPOT or ARRAY module cannot be used at the same time as the ECD module.
- **RDE or RRDE:** the remote control option of the rotating disc electrode (RDE) or rotating ring-disc electrode (RRDE) is not possible when the ECD module is used.

Additionally, the following restrictions apply to the available current ranges and the automatic current ranging option (see figure 1193, page 1052):

- **Current range restrictions:** when the ECD module is used, the 1 mA and higher current ranges are no longer usable by the Autolab. An **error** message is provided by NOVA when this situation is encountered.
- **Automatic current ranging restriction:** when the ECD on-board filter is used, the automatic current ranging option is not available. An **error** message is provided by NOVA when this situation is encountered.




Figure 1193 The ECD current ranges are available in the Automatic current ranging option

The ECD module is also not available for use in combination with galvanostatic measurement and measurements involving the FRA module.

Finally, the current ranges provided by the ECD module have a limited bandwidth. When the interval time specified in the procedure is too small with respect to the bandwidth of the ECD current ranges, a warning is provided by the procedure validation. The bandwidth values of the ECD module are reported in Table 41.

Table 41 Overview of the bandwidth of the ECD module current ranges

Current range	Bandwidth
100 µ and 10 µA	2000 Hz
1 µA	1000 Hz
100 nA	250 Hz
10 nA	100 Hz
1 nA	50 Hz
100 pA	10 Hz



CAUTION

It is possible to force the procedure to continue despite the bandwidth warning. This is not recommended since the measured data could be affected by this limitation and be invalid.

16.3.2.7.6 ECD module front panel connections

The ECD module is fitted with a single female BNC connector, labeled ←I (see figure 1194, page 1053).




Figure 1194 The front panel label of the ECD module

The signal provided through the ←I connector on the front panel corresponds to the output of the current-to-voltage converter located on the ECD module. The output signal is a voltage, referred to the instrument ground, corresponding to the converted current according to:

$$E_{\text{out}}(\leftarrow I) = \frac{i_{(\text{ECD})}}{[\text{CR}]}$$

Where $E_{\text{out}}(\leftarrow I)$ corresponds to the output voltage signal of the module, in V, $i_{(\text{ECD})}$ corresponds to the current measured by the ECD module, in A and [CR] is the active current range of the ECD module.



NOTICE

The front panel ←I BNC output is provided for information purposes only.

16.3.2.7.7 ECD module testing

NOVA is shipped with a procedure which can be used to verify that the **ECD** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestECD** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestECD.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1195*.





Figure 1195 The results of the TestECD procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestECD automatic evaluation of the data requires the following tests to succeed:

1. The residual current, determined by subtracting the expected current from the measured current, must be smaller or equal than ± 20 nA.
2. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$.
3. The intercept of the measured current versus the applied potential must be equal to ± 4 mV divided by $1000100 \Omega \pm 5 \%$.

All three conditions must be valid for the test to succeed.

16.3.2.7.8 ECD module specifications

The specifications of the ECD module are provided in *Table 42*.

Table 42 Specifications of the ECD module

Specification	Value
Current ranges	100 μ A to 100 pA, in 7 ranges
Current accuracy	$\pm 0.5 \%$ of current range
Current offset	± 2 pA
On-board filter	3 rd order Sallen-Key
Filter time constants	10 ms, 100 ms and 500 ms
Maximum offset compensation	± 1 μ A

16.3.2.8 ECI10M module

The ECI10M is an extension module for the Autolab PGSTAT. With the ECI10M module, it is possible to perform electrochemical impedance spectroscopy measurements up to a frequency of 10 MHz.

The ECI10M module is an auxiliary external high bandwidth potentiostat/galvanostat that is used instead of the main Autolab PGSTAT during high frequency impedance measurements.

The ECI10M consist of a module installed in the Autolab and an external interface to be placed close to the cell.




Figure 1196 The ECI10M module is selected in the hardware setup

16.3.2.8.4 ECI10M module settings



CAUTION

The ECI10M cannot be used at the same time as the Autolab PGSTAT. The ECI10M and the Autolab PGSTAT are both electrochemical interfaces to the electrochemical cell. To use the ECI10M module, it is necessary to switch the control of the electrochemical cell from the Autolab PGSTAT to the ECI10M.




Figure 1198 The electrochemical interface can be directly selected through the Dashboard



CAUTION

Switching the **Electrochemical interface** from PGSTAT to ECI10M or the other way around using the provided drop-down list takes a couple of seconds. During this time, the Autolab system cannot be used.

At any time, a tooltip shows the active electrochemical interface, in bold (see figure 1199, page 1059).




Figure 1199 A tooltip shows the active electrochemical interface in bold



Follow these steps to switch the control of the cell from the Autolab PGSTAT to the ECI10M module:

1. Make sure that the cell switch of the Autolab PGSTAT is in the **off** position.
2. Disconnect the Autolab PGSTAT cell connections (CE, RE, WE and S) from the electrochemical cell.
3. Switch the **Electrochemical interface** to **ECI10M**.
4. Connect the cables supplied with the ECI10M to the electrochemical cell.

Follow these steps to switch the control of the cell from the ECI10M module to the Autolab PGSTAT:

1. Make sure that the cell switch of the ECI10M module is in the **off** position.
2. Disconnect the cables supplied with the ECI10M from the electrochemical cell.
3. Switch the **Electrochemical interface** to **PGSTAT**.
4. Connect the Autolab PGSTAT cell connections (CE, RE, WE and S) to the electrochemical cell.

When the ECI10M module is used as the **Electrochemical interface**, the following user-definable settings are available, through the **Autolab control** command (*see figure 1200, page 1061*):

- **Cell:** a  toggle that can be used to switch the cell of the ECI10M module on or off.
- **Mode:** specifies the mode of operation of the ECI10M module (potentiostatic, galvanostatic), using the provided drop-down list.
- **Current range:** specifies the active current range of the ECI10M module, using the provided drop-down list.
- **Bandwidth:** specifies the bandwidth of the ECI10M module (high stability, high speed, ultra high speed), using the provided drop-down list.
- **FRA32M input:** a  toggle that can be used to switch the FRA32M input for the summation point of the ECI10M on or off.
- **Reference potential:** an input slider that can be used to set the reference potential of the ECI10M module. This value can be set in the range of ± 10 V.
- **Offset potential/current:** an input slider that can be used to set the offset potential or current of the ECI10M module. This value can be set in the range of ± 5 V (in potentiostatic mode) and in the range of ± 5 times the active current range (in galvanostatic mode).




Figure 1200 The ECI10M module settings are defined in the Autolab control command

16.3.2.8.5 ECI10M bandwidth and contour map

The bandwidth for each current range of the **ECI10M** module are reported in *Table 43*.

Table 43 Bandwidth overview of the ECI10M module

Current range	Bandwidth
100 mA - 1 mA	10 MHz
100 μ A	4 MHz
10 μ A	150 kHz
1 μ A	15 kHz
100 nA	1.5 kHz
10 nA	150 Hz

A typical contour map for the **ECI10M** module in combination with the **PGSTAT302N** potentiostat/galvanostat with **FRA32M** module is shown in *Figure 1201*.




Figure 1201 Contour map of the ECI10M/PGSTAT302N/FRA32M combination

The map reported in Figure 1201 shows a dark green area, which corresponds to the area of the map where an error of $\pm 0.3^\circ$ on the measured phase angle and $\pm 0.3\%$ on the measured impedance value is expected. The light green area corresponds to the area of the map where an error of $\pm 10^\circ$ on the measured phase angle and $\pm 10\%$ on the measured impedance value is expected.



NOTICE

The contour map is determined empirically with an amplitude of 10 mV, in potentiostatic mode.

16.3.2.8.6 ECI10M module restrictions and precautions

Restrictions apply when using the ECI10M module:

- **FRA32M required:** in order to perform electrochemical impedance measurements up to 10 MHz using the ECI10M the **FRA32M** module must be installed in the instrument *FRA32M module* (see chapter 16.3.2.13, page 1112).
- **Concurrent use:** when the ECI10M is in use, the Autolab potentiostat/galvanostat is bypassed and cannot be used. The reverse applies to the ECI10M when the Autolab potentiostat/galvanostat is in use.
- **Module incompatibility:** while the ECI10M is in use, the additional internal **optional modules** installed in the Autolab cannot be used, with the exception of the FRA32M module.

Precautions apply when using the ECI10M module at high frequency:

- **Capacitive leakage:** for high frequency measurements with the ECI10M is it recommended to use unshielded cables to connect the ECI10M external interface to the electrochemical cell. Shielded cables have a capacitance of about 50 pF per meter. At very high frequency, the coaxial assembly of these cables will add a significant parallel capacitance to the impedance of the cell.
- **Cell proximity:** if possible, it is recommended to reduce the length of the cables between the ECI10M external interface and the electrochemical cell as much as possible. The default length of the cables supplied with the ECI10M is 25 cm and the device is optimized for this length of cable.

16.3.2.8.7 ECI10M front panel connection

The ECI10M module is fitted with a single female, DVI connector. A matching cable is used to connect to the external ECI10M external interface (see figure 1202, page 1063).




Figure 1202 The front panel labels of the ECI10M

The external interface of the ECI10M is fitted with a cable that directly connects to the front panel of the ECI10M module.

The front panel of the ECI10M external interface provides a number of connections and indicators, shown in Figure 1203.




Figure 1203 Overview of the front panel of the ECI10M

- | | |
|--|---|
| <p>1 RE connection (core)
4 mm banana connection for the Reference electrode (RE).</p> | <p>2 Air flow holes
For cooling the ECI10M during operation.</p> |
| <p>3 Ground plug
4 mm banana connector for grounding purposes.</p> | <p>4 S connection (core)
4 mm banana connection for the Sense electrode (S).</p> |
| <p>5 WE connection (core)
4 mm banana connection for the Working electrode (WE).</p> | <p>6 WE connection (shield)
4 mm banana connection for shield of the Working electrode (WE), if applicable.</p> |
| <p>7 S connection (shield)
4 mm banana connection for shield of the Sense electrode (S), if applicable.</p> | <p>8 Status LED
Used to indicate the status of the ECI10M (off, green or red).</p> |
| <p>9 RE connection (shield)
4 mm banana connection for shield of the Reference electrode (RE), if applicable.</p> | <p>10 CE connection (shield)
4 mm banana connection for shield of the Counter electrode (CE), if applicable.</p> |
| <p>11 CE connection (core)
4 mm banana connection for the Counter electrode (CE).</p> | |



CAUTION

Never connect the electrode connectors from the PGSTAT instrument to the front panel of the ECI10M!

16.3.2.8.8 ECI10M module testing

NOVA is shipped with a procedure which can be used to verify that the **ECI10M** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestECI10M** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestECI10M.nox)

2 Connect the Autolab dummy cell

Connect the **ECI10M** to the Autolab dummy cell circuit (e).



3 Switch on the ECI10M

Using the **Electrochemical interface** drop-down list provided in the **Instrument control** panel, switch the control of the electrochemical cell to the ECI10M (see figure 1204, page 1066).




Figure 1204 Switch on the ECI10M module in the Instrument control panel



NOTICE

Please allow the **ECI10M** module to warm up to 30 minutes.

4 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out an impedance spectroscopy measurement. During the measurement, the data is fitted using a (RC) equivalent circuit. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in Figure 1205.





Figure 1205 The data measured by the TestECI10M procedure

5 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestECI10M automatic evaluation of the data requires the following tests to succeed:

1. The fitted resistance must be equal to $10000 \Omega \pm 5\%$.

This condition must be valid for the test to succeed.

16.3.2.8.9 ECI10M module specifications

The specifications of the ECI10M module are provided in *Table 44*.

Table 44 Specifications of the ECI10M

Specification	Value
Compliance voltage	± 10 V
Applied voltage range	± 10 V
Maximum current	± 100 mA
Current ranges	100 mA to 10 nA, in 8 decades
Electrode connections	2, 3 and 4
Frequency range	10 μ Hz - 10 MHz
Frequency resolution	0.003 %
Maximum AC amplitude	700 mV (RMS)
Bandwidth	15 MHz
Potential accuracy	± 0.2 %
Potential resolution	0.3 μ V
Current accuracy	± 0.2 %
Current resolution	0.0003 % (of current range)
Input impedance	> 100 G Ω

16.3.2.9 ECN module

The ECN is an extension module for the Autolab PGSTAT. The ECN module provides the means to perform Electrochemical Noise (ECN) measurements. Electrochemical noise corresponds to seemingly random fluctuations in current and potential generated by stochastic phenomena occurring at the electrochemical interface.



The fluctuations of potential and current signals that arise directly from the electrochemical reactions, taking place on the electrode surface, can be measured using the optional ECN module. The ECN measurement is non-invasive because no external perturbation is applied and the current is monitored through a so-called Zero Resistance Ammeter (ZRA).

ECN can be used to monitor localized corrosion (pitting), uniform corrosion through measurement of the Noise Resistance, and the deterioration of paints on metal substrates. The same technique can also be used to monitor galvanic coupling in the presence of an electrolyte.



NOTICE

Electrochemical noise measurements are also possible without the use of the dedicated ECN module. However, the resolution that can be achieved with the ECN is significantly better than without the use of this module.

The ECN module provides a dedicated differential amplifier, which can be used instead of the default differential amplifier provided by the Autolab. The ECN offers the possibility to carry out DC potential compensation and a four times additional amplification of the measured potential, leading to a maximum resolution of 760 nV.

The ECN module adds the following signal to the Sampler (*see figure 1206, page 1068*):

- **ECN(1).Potential (V):** this signal corresponds to the potential measured by the ECN module.




Figure 1206 The ECN module provides the ECN(1).Potential signal

16.3.2.9.1 **ECN module compatibility**

The ECN module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- PGSTAT10 and PGSTAT20



NOTICE

The ECN module is **not** compatible with the Autolab instruments not listed above.

16.3.2.9.2 **ECN module scope of delivery**

The ECN module is supplied with the following items:

- ECN module
- ECN module label
- BNC to banana connection cable

16.3.2.9.3 **ECN hardware setup**

To use the **ECN** module, the hardware setup needs to be adjusted. The checkbox for the module needs to be ticked (*see figure 1207, page 1070*).




Figure 1208 The ECN module setting is defined in the Autolab control command

16.3.2.9.5 ECN module restrictions

Restrictions apply when using the ECN module during electrochemical noise measurements:

- **Current range restrictions:** during electrochemical noise measurements in combination with the ECN module, the current ranges of 1 mA and higher cannot be used. An **error** message is shown by NOVA when this situation is encountered.
- **Cell switch:** electrochemical noise measurements are normally carried out at open circuit potential. Whenever the ECN(1).Potential signal is sampled, the cell should be switched off. NOVA provides a **warning** message when ECN measurements are carried out with the cell switched on.

16.3.2.9.6 ECN module front panel connections

The ECN module is fitted with two female BNC connector, labeled ←V and →E, respectively (see figure 1209, page 1071).




Figure 1209 The front panel label of the ECN module

The signal provided through the ←V connector on the front panel corresponds to the output of the voltage follower located on the ECN module.

- Connect the working electrode (WE) provided by the Autolab cell cable to electrode #1.
- Connect the green ground connector provided by the Autolab cell cable to electrode #2.
- Connect the red connector of the ECN cable to electrode #1.
- Connect the black connector of the ECN cable to the reference electrode (if present), or to electrode #2.

16.3.2.9.8 ECN module testing

NOVA is shipped with a procedure which can be used to verify that the **ECN** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestECN** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestECN.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a) and the ECN module to the Autolab dummy cell (a).



NOTICE

Connect the cables from the ECN module to the dummy cell first.

Both conditions must be valid for the test to succeed.

16.3.2.9.9 ECN module specifications

The specifications of the ECN module are provided in *Table 45*.

Table 45 Specifications of the ECN module

Specification	Value
Input range	± 2.5 V
Maximum potential resolution	760 nV (Gain 100)
Potential offset compensation range	± 2.5 V
Potential accuracy	300 μ V
Input impedance	100 G Ω
Input bias current	< 25 fA

16.3.2.10 EQCM module

The EQCM module is an extension module for the Autolab PGSTAT and the Multi Autolab. The EQCM module provides the means to perform Electrochemical Quartz Crystal Microbalance measurements.

The EQCM module measures a mass change per unit area by measuring the change in resonant frequency of a quartz crystal. Quartz crystals belong to a group of materials displaying the so-called piezoelectric effect. When a properly cut crystal (AT-cut) is exposed to an AC current, the crystal starts to oscillate at its resonant frequency and a standing shear wave is generated.

In first approximation, the resonant frequency depends on the thickness of the crystal. As mass is deposited on the surface of the crystal, the thickness increases; consequently the frequency of oscillation decreases from the initial value. With some simplifying assumptions, this frequency change can be quantified and correlated precisely to the mass change using Sauerbrey's equation:

$$\Delta f = -\frac{2f_0^2}{A\sqrt{\rho_q\mu_q}} \cdot \Delta m$$

Where Δf is the change in oscillating frequency, f_0 is the nominal resonant frequency of the crystal (6 MHz), Δm is the change in mass, in g/cm², A is the area of the crystal in cm², ρ_q is the density of quartz, in g/cm³ and μ_q is the shear modulus of quartz, in g/cm·s².

For a 6 MHz crystal, the same equation can be reduced to:

$$-\Delta f = \Delta m \cdot C_f$$

Where C_f is 0.0815 Hz/ng/cm².



NOTICE

For more information on the EQCM module please refer to the EQCM User Manual. More information on the validity and the application of the Sauerbrey equation can be found in the peer-reviewed literature.

The EQCM module adds the following signal to the Sampler (*see figure 1212, page 1076*):

- **EQCM(1).Temperature (°C):** this signal corresponds to the temperature recorded by the sensor located at the bottom of the EQCM cell.
- **EQCM(1).Driving force (V):** this value represents the amount of energy required to sustain the oscillation of the crystal. When the loading of the crystal increases, the driving force also increases. In air, the typical driving force is close to 0 V. In water, the driving force is about 0.85 V.
- **EQCM(1).ΔFrequency (Hz):** this signal corresponds to the relative change in oscillation frequency of the quartz crystal. This variation is expressed with respect to an arbitrary, user-defined reference frequency (zero Hz).




Figure 1212 The EQCM module provides the EQCM(1).ΔFrequency, EQCM(1).Driving force and EQCM(1).Temperature signals

16.3.2.10.1 EQCM module compatibility

The EQCM module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- PGSTAT100N and PGSTAT100
- M101
- PGSTAT204/M204



NOTICE

The EQCM module is **not** compatible with the Autolab instruments not listed above.

16.3.2.10.2 EQCM module scope of delivery

The EQCM module is supplied with the following items:

- EQCM module
- EQCM module label
- EQCM oscillator
- EQCM connection cable
- EQCM cell
- Ag/AgCl reference electrode
- Au counter electrode
- 6 MHz, double-sided, Au coated quartz crystal (2)
- EQCM trimmer
- EQCM User Manual

16.3.2.10.3 EQCM hardware setup

To use the **EQCM** module, the hardware setup needs to be adjusted. The checkbox for the module needs to be ticked (*see figure 1213, page 1078*).




Figure 1214 The EQCM module setting is defined in the Autolab control command



NOTICE

The EQCM oscillator will be powered as long as the EQCM module is on. It is recommended to allow the EQCM oscillator to warm up before starting electrochemical measurements and to keep it on when needed to avoid temperature changes.

16.3.2.10.5 EQCM module restrictions

Restrictions apply when using the EQCM module:

- **Sampling rate:** The EQCM module is capable of providing one new set of values for the measured signals (Δ Frequency, Driving force, Temperature) with an interval time of 20 ms (50 samples/s). When the sampling rate specified in the procedure is smaller than 20 ms, the EQCM module is not able to provide new data points quickly enough. In practice this means that last available data point provided by the EQCM module will be measured several times, until a new data point is available.
- **Dynamic frequency range:** The EQCM module uses a dynamic window of frequency of 1000 Hz. This ensures that the frequency is measured with the highest possible resolution. When changes in frequency larger than 1000 Hz are measured, the measurement window is adjusted downwards or upwards depending on the direction of the frequency change. This software rewindowing takes 100 ms. During this adjustment, the EQCM is not able to supply new data points and this can be seen in the measured data

1 Load the procedure

Load the **TestEQCM** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestEQCM.nox)

2 Insert an EQCM crystal in the EQCM cell

insert a 6 MHz EQCM crystal in the EQCM cell. Fill the cell with 2 ml of water and check for leakage. Connect the cell to the EQCM oscillator and the oscillator to the Autolab PGSTAT using the provided cable. Leave the cell connectors from the PGSTAT disconnected.



NOTICE

For more information on the EQCM hardware, please refer to the EQCM User Manual.

3 Start the procedure

Start the procedure and follow the instructions on-screen.



NOTICE

Please allow the **EQCM** module to warm up to 15 minutes.

After 15 minutes, the test can be continued. The Determine EQCM zero frequency window will be displayed. Using the provided adjustment tool, rotate the trimmer on the EQCM oscillator in order to minimize the driving force and zero the Δ Frequency signal (as explained in the EQCM User Manual).

The test carries out a time resolved measurement measurement. The measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1216*.





Figure 1216 The results of the TestEQCM procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestEQCM automatic evaluation of the data requires the following tests to succeed:

1. The EQCM(1).Driving force must stable within ± 0.025 V.
2. The average value of the EQCM(1). Δ Frequency must be larger or equal to -5 Hz and must be smaller or equal to 5 Hz.

Both condition must be valid for the test to succeed.

16.3.2.10.8 EQCM module specifications

The specifications of the EQCM module are provided in *Table 46*.

Table 46 Specifications of the EQCM module

Specification	Value
Oscillation frequency	6 MHz
Frequency resolution	0.07 Hz
Relative accuracy	1 Hz
Absolute accuracy	10 Hz
Frequency range	80 kHz
Temperature sensor accuracy	1 °C
Temperature sensor resolution	0.1 °C

16.3.2.11 FI20 module

The FI20 module is an extension module for the Autolab potentiostat/galvanostat. This module provides two electronic circuits that can be used to process the current measured by the instrument. Each of these circuits fulfills a specific role:

- **Filter circuit:** the filter circuit is designed to filter the current during electrochemical measurements. The filter can be used to remove noise on the measurements in cases where it is impossible to remove the noise by the use of proper shielding of the cell and electrodes or by using a Faraday cage. The module uses a third order Sallen-Key filter, with three different filter time constants (0.1 s, 1 s and 5 s).
- **Integrator circuit:** the integrator circuit provides the means to integrate the measured current. The integrator can be used to perform chronocoulometric experiments and the so-called cyclic or linear sweep voltammetry current integration. The module consists of an analog integrator fitted with four different integration time constants (0.01 s, 0.1 s, 1 s, 10 s).



NOTICE

The **integrator** circuit is present by default (as on-board integrator) on the following instruments: PGSTAT10, PGSTAT20, μ Autolab II, μ Autolab III, PGSTAT101, M101, PGSTAT204 and M204.



CAUTION

The FI20 and the on-board integrator can only be used to process the current measured by the Autolab potentiostat/galvanostat (WE(1).Current).

The FI20 module and the on-board integrator add the following signal to the Sampler (see figure 1217, page 1084):

- **Integrator(1).Charge (C):** this signal corresponds to the measured charge.
- **Integrator(1).Integrated current (A):** this signal corresponds to the converted equivalent current obtained by deriving the measured charge over the interval time used in the measurement. This signal can be used in order to perform so-called current integration cyclic and linear sweep voltammetry measurements. In first approximation, at low scan rates, the results obtained with the current integration method can be compared to the results obtained with a linear scan generator.



NOTICE

The FI20 module is **not** compatible with the Autolab instruments not listed above.

The on-board integrator is fitted in the following instruments:

- μ Autolab II and μ Autolab III
- PGSTAT101 and M101
- PGSTAT204 and M204

16.3.2.11.2 FI20 module scope of delivery

- FI20 module
- FI20 module label

16.3.2.11.3 FI20 hardware setup

To use the **FI20** module, the hardware setup needs to be adjusted. The checkbox for the module needs to be ticked (*see figure 1218, page 1085*).

The screenshot shows a software interface for hardware configuration. It is divided into two main sections: 'Additional modules' and 'Properties'.

Additional modules: A list of modules with checkboxes. The 'FI20 - Integrator' checkbox is checked and highlighted with a green background. A mouse cursor is pointing at this checkbox. Other modules listed include FRA32M, ECI10M, FRA2, ADC10M, ADC750, ADC750r4, SCAN250, SCANGEN, BA, BIPOT/ARRAY, ECD, FI20 - Filter, Booster20A, Booster10A, EQCM, pX1000, pX, ECN, External Devices, IME303, IME663, and MUX.

Properties: A section containing a 'Calibration factor' input field with the value '1'.

Figure 1218 The FI20 module is selected in the hardware setup

16.3.2.11.4 FI20 module and on-board integrator settings

The settings of the FI20 module and the on-board integrator are completely defined in the NOVA software.

The following user-definable settings are available for the **filter** circuit, the **integrator** circuit of the FI20 module and the on-board integrator, through the **Autolab control** command (see figure 1220, page 1087):

- In the **Filter** sub-panel:
 - **Filter time:** a drop-down control that can be used to select the time constant for the filter circuit. The filter time can be set to off, 100 ms, 1 s and 5 s.
- In the **Integrator** sub-panel:
 - **Integration time:** a drop-down control that can be used to select the time constant for the integrator circuit. The integration time can be set to 10 ms, 100 ms, 1 s, 10 s and **hold**.
 - **Discharge integrator:** a that can be used to discharge the integrator.
 - **Integrator drift:** a value field which can be used to manually set the integrator drift, in C/s.




Figure 1220 The settings of the filter circuit and the integrator circuit of the FI20 module and the on-board integrator are defined in the Autolab control command

The **Filter time** setting defines the time constant of the filter. The effect of a filter constant of x s is that x seconds after a potential perturbation has been applied, the current response can be measured correctly. A high time constant results in heavy filtering.



NOTICE

It is highly recommended to set the Discharge integrator property to On at the end of each measurement involving the **integrator** circuit of the FI20 module or the on-board integrator. Leaving the integrator charged could introduce extra noise during consequent measurements.

The **Integrator drift** can be manually specified. The drift is defined as the charge accumulation due to leakage current, in C/s. Setting the drift allows to compensate for the non-ideality of the Autolab. This value can be set manually or it can be determined automatically using the provided drift determination tool, available in the instrument control panel *Reset integrator drift* (see chapter 5.2.2.6, page 126).



NOTICE

The drift depends on the active current range. It is recommended to determine the drift whenever the current range is modified.

16.3.2.11.5 FI20 module restrictions

Restrictions apply when using the FI20 module or the on-board integrator.

The following restrictions apply when using the **filter** circuit of the FI20 module:

- **Automatic current ranging:** the Automatic Current Ranging option cannot be used with the **filter** circuit of the FI20 module is used. An **error** message is provided by NOVA when this situation is encountered.
- **Bandwidth limit:** the **filter** circuit of the FI20 module has a limited bandwidth. Depending on the interval time, a **warning** can be provided by NOVA indicating that the selected filter time constant is too slow to measure properly.

The following restrictions apply when using the **integrator** circuit of the FI20 module or the on-board integrator:

- **Automatic current ranging:** the Automatic Current Ranging option cannot be used with the **integrator** circuit of the FI20 module or the on-board integrator is used. An **error** message is provided by NOVA when this situation is encountered.

16.3.2.11.7 FI20 module and on-board integrator module testing

Three test procedures are provided for testing the FI20 module and the on-board integrator:

- For the Filter circuit of the FI20 module, please refer to *Chapter 16.3.2.11.7.1*.
- For the Integrator circuit of the FI20 module and the on-board integrator of all instrument equipped with an on-board integrator except the PGSTAT101, please refer to *Chapter 16.3.2.11.7.2*.
- For the on-board integrator of the PGSTAT101, please refer to *Chapter 16.3.2.11.7.3*.

16.3.2.11.7.1 FI20 module filter test

NOVA is shipped with a procedure which can be used to verify that the filter circuit of the **FI20** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestFI20-Filter** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestFI20-Filter.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1222*.

16.3.2.11.7.2 FI20 module integrator test

NOVA is shipped with a procedure which can be used to verify that the integrator circuit of the **FI20** module and the on-board integrator of all equipped instruments, except the Autolab PGSTAT101, is working as expected.



NOTICE

To test the on-board integrator of the PGSTAT101, a dedicated test is provided. Please refer to *Chapter 16.3.2.11.7.3* for more information.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestFI20-Integrator** procedure, provided in the NOVA 2.X installation folder (Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestFI20-Integrator.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Reset the integrator drift

Open the instrument control panel, set the current range to 10 μA and reset the integrator drift using the provided button (see figure 1223, page 1094).




Figure 1223 Resetting the integrator drift

4 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in (see figure 1224, page 1094).





Figure 1224 The results of the TestFl20-Integrator procedure

5 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestFI20-Integrator automatic evaluation of the data requires the following tests to succeed:

1. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$.
2. The calculated capacitance, determined from the integrated current signal, must be equal to $1 \mu\text{F} \pm 10 \%$.

Both conditions must be valid for the test to succeed.

16.3.2.11.7.3

On-board integrator test (PGSTAT101 only)

NOVA is shipped with a procedure which can be used to verify that the on-board integrator of the Autolab PGSTAT101 is working as expected.



CAUTION

This test is only suitable for the Autolab PGSTAT101! For the other instruments, please refer to *Chapter 16.3.2.11.7.2*.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestFI20-Integrator-PGSTAT101** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestFI20-Integrator-PGSTAT101.nox)

2 Connect the electrode cables

Connect the counter electrode (CE) and reference electrode (RE) together and the working electrode (WE) and sense electrode (S) together.



3 Reset the integrator drift

Open the instrument control panel, set the current range to 10 μA and reset the integrator drift using the provided button (see figure 1225, page 1096).




Figure 1225 Resetting the integrator drift

4 Ignore the warning

Ignore the warning message shown when the procedure is loaded and when the procedure is started.

TestFI20 integrator-PGSTAT101

The following problems were encountered during validation.

 The internal dummy cell is on.

OK

Cancel

5 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1226*.





Figure 1226 The results of the TestFI20-Integrator-PGSTAT101 procedure

6 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestFI20-Integrator-PGSTAT101 automatic evaluation of the data requires the following tests to succeed:

1. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$.
2. The calculated capacitance, determined from the integrated current signal, must be equal to $1 \mu\text{F} \pm 10 \%$.

- μ Autolab III



NOTICE

The FRA2 module is **not** compatible with the Autolab instruments not listed above.

16.3.2.12.2 FRA2 module scope of delivery

The FRA2 module is supplied with the following items:

- FRA2 module
- FRA2 module labels

16.3.2.12.3 FRA2 hardware setup

To use the **FRA2** module, the hardware setup needs to be adjusted. The checkbox for the module needs to be ticked. The **FRA2** module also has a number of additional properties that need to be specified correctly in the **Properties** panel (*see figure 1227, page 1100*).




Figure 1227 The FRA2 module is selected in the hardware setup

The following additional properties must be defined:

- **FRA2 calibration:** the calibration file (*fra2cal.ini*) for the FRA2 module. The **Import** button is provided to specify the location of this file.
- **FRA2 offset DAC range:** specifies the offset DAC range of the FRA2 module using the provided drop-down list. This property can be set to either 5 V or 10 V.



CAUTION

The serial number specified in the **FRA2 calibration file** must match the serial number of the connected instrument.



CAUTION

The **FRA2 offset DAC range** property must be specified carefully. Failure to set this value properly may result in faulty data at frequencies of 25 Hz and lower (refer to front panel labels of the FRA2 module on the instrument *FRA2 module front panel connections* (see chapter 16.3.2.12.9, page 1109))

16.3.2.12.4 C1 and C2 calibration factors

When the **FRA2** module is used in combination with the Autolab, the **C1** and **C2** calibration factors need to be determined.



NOTICE

The **C1** and **C2** calibration factors are predetermined when the **FRA2** module is preinstalled. These factors must be determined experimentally when a **FRA2** module is installed into an existing instrument. This determination must only be carried out upon installation of the module.

Two procedures are supplied with NOVA to determine these calibration factors:

- PGSTAT C1 calibration
- PGSTAT C2 calibration


The determination of **C1** and **C2** requires the following items:

- Autolab Dummy cell
- Faraday cage


Typical values are indicated in *Table 49*.

Table 49 Typical values for C1 and C2

Instrument type	C1	C2
PGSTAT302N	1.6 E-11	3.0 E-13
PGSTAT302F	1.6 E-11	1.0 E-12
PGSTAT128N (serial number \leq AUT84179)	2.6 E-11	1.0 E-12
PGSTAT128N (serial number $>$ AUT84179)	1.6 E-11	1.0 E-12
PGSTAT100N	1.6 E-11	5.0 E-13



- 4 Start the measurement and wait until it finishes. Ignore the warning message displayed at the beginning of the measurement.
- 5 During the measurement, the data will be plotted as a Bode plot and should be similar to the example shown.
- 6 The measured data is automatically fitted and a message is shown at the end, displaying the measured **C1** value.



- 7 Open the instrument hardware setup and type the measured value in the **C1** field.


Autolab module

Main module


Power supply frequency

C1

C2



- 7 Open the instrument hardware setup and type the measured value in the **C2** field.



- 8 Close the hardware setup and wait for the Autolab to be reinitialized using the updated Hardware setup.

16.3.2.12.5 FRA2 module settings

The FRA2 module settings are completely defined in the NOVA software. The following user-definable settings are available, through the **Autolab control** command (see figure 1228, page 1106):

- **FRA2:** a toggle that can be used to switch the output of the FRA2 on or off.
- **Frequency:** the output frequency of the FRA2 module, in Hz.
- **Amplitude:** the amplitude of the FRA2 module output, in V, specified as a TOP value.

FRA2

Properties

Frequency Hz

Amplitude mV

Input connection

Wave type

Integration time ms

Minimum cycles to integrate

FRA2

Results

Elapsed time -

E(DC) (V)	E(AC) (V)	i(DC) (A)	i(AC) (A)	% E	% i
999,7 m	10,91 m	943,2 μ	10,30 μ	35	41

Freq. (Hz)	Z (Ω)	-Phase ($^{\circ}$)	Z' (Ω)	-Z'' (Ω)
1,000	1,058 k	0,3161	1,058 k	5,838

Figure 1229 The FRA2 manual control panel

As soon as the FRA2 is switched on using the provided toggle, the module will start a manual measurement using the properties specified in the panel. The measured values will be displayed after the measurement in the Results sub-panel.



NOTICE

The measurement will continue until the module is switched off using the provided toggle.

16.3.2.12.8 FRA2 module restrictions

No restrictions apply when using the FRA2 module.

16.3.2.12.9 FRA2 module front panel connections

The FRA2 module has twice the size of a normal Autolab module. The module consists of a function generator module and transfer function analyzer. The FRA2 module is fitted with three female BNC connectors, labeled $\rightarrow X$, $\rightarrow Y$ (on the transfer function analyzer front panel) and $\leftarrow V$ (on the function generator front panel).

Two versions of the FRA2 module exist:

- **FRA2 module:** the default FRA2 module with a 5 V input range. This module is identified with the module labels shown in *Figure 1231*.
- **FRA2V10 module:** a modified version of the FRA2 module with a 10 V input range. This module is identified with the module labels shown in *Figure 1232*.




Figure 1231 The front panel labels of the FRA2 module (5 V input range version)




Figure 1232 The front panel labels of the FRA2 module (10 V input range version)

The two connectors, labeled $\rightarrow X$ and $\rightarrow Y$ are input connectors that can be used to analyze external transfer functions. They have an input range of ± 5 V or ± 10 V and an input impedance of 50Ω .

The signal provided through the $\leftarrow V$ connector on the front panel corresponds to the output of the sinewave generator of the FRA2. Whenever the FRA2 module is used, the voltage provided at this output corresponds to the signal generated by the module (either single sine or multi sine).



The output signal is a voltage, referred to the instrument ground, corresponding to the applied amplitude, multiplied by 10 (when the instrument is working in potentiostatic mode), or the converted amplitude, multiplied by 10 (when the instrument is working in galvanostatic mode):

$$E_{\text{out}}(\leftarrow V) = E_{(\text{FRA2})} \cdot 10$$

$$E_{\text{out}}(\leftarrow V) = \frac{i_{(\text{FRA2})}}{[\text{CR}]} \cdot 10$$

Where $E_{\text{out}}(\leftarrow V)$ corresponds to the output voltage signal of the module, in V, $E_{(\text{FRA2})}$ and $i_{(\text{FRA2})}$ corresponds to the specified amplitude, in V or A, respectively and $[\text{CR}]$ is the active current range of the FRA2 module.



NOTICE

The front panel $\leftarrow V$ BNC output is provided for information purposes only except for impedance measurement involving external transfer functions.

16.3.2.12.10

FRA2 module testing

NOVA is shipped with a procedure which can be used to verify that the **FRA2** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestFRA** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestFRA.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (c).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out an impedance spectroscopy measurement. During the measurement, the data is fitted using a R(RC) equivalent circuit. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1233*.





Figure 1233 The data measured by the TestFRA procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestFRA automatic evaluation of the data requires the following tests to succeed:

1. The fitted series resistance must be equal to $100 \Omega \pm 5\%$.
2. The fitted parallel resistance must be equal to $1000 \Omega \pm 5\%$.
3. The fitted parallel capacitance must be equal to $1 \mu\text{F} \pm 10\%$.
4. The calculated χ^2 must be smaller or equal to 0.01.

16.3.2.13.2 FRA32M module scope of delivery

The FRA32M module is supplied with the following items:

- FRA32M module
- FRA32M module label

16.3.2.13.3 C1 and C2 calibration factors

When the **FRA32M** module is used in combination with the Autolab, the **C1** and **C2** calibration factors need to be determined.



NOTICE

The **C1** and **C2** calibration factors are predetermined when the **FRA32M** module is preinstalled. These factors must be determined experimentally when a **FRA32M** module is installed into an existing instrument. This determination must only be carried out upon installation of the module.



NOTICE

On some instruments, the value of **C1** and **C2** is already determined and stored in the on-board processor of the instrument. In this case, the values reported in the Hardware setup are not 0. For these instruments, it is not necessary to determine **C1** and **C2**.

Two procedures are supplied with NOVA to determine these calibration factors:

- PGSTAT C1 calibration
- PGSTAT C2 calibration

The determination of **C1** and **C2** requires the following items:

- Autolab Dummy cell
- Faraday cage



CAUTION

The determination of the **C1** and **C2** calibration factors is not required for the PGSTAT204 and for the M101 and M204 modules used in combination with the **FRA32M** module in the Multi Autolab instrument.


Typical values are indicated in *Table 51*.




CAUTION


Do **not** connect the ground connector from the PGSTAT to the Autolab Dummy cell. Place the dummy cell in the Faraday cage.

- 1 Start the NOVA software and allow the instrument to warm up for at least 30 minutes.
- 2 Open the **PGSTAT C1 calibration** procedure.
- 3 Connect the Autolab Dummy cell as shown. Connect the ground lead from the PGSTAT to the Faraday cage.




- 4 Start the measurement and wait until it finishes. Ignore the warning message displayed at the beginning of the measurement.
- 5 During the measurement, the data will be plotted as a Bode plot and should be similar to the example shown.
- 6 The measured data is automatically fitted and a message is shown at the end, displaying the measured **C1** value.





- 4 Start the measurement and wait until it finishes. Ignore the warning message displayed at the beginning of the measurement.
- 5 During the measurement, the data will be plotted as a Bode plot and should be similar to the example shown.
- 6 The measured data is automatically fitted and a message is shown at the end, displaying the measured **C2** value.



- 7 Open the instrument hardware setup and type the measured value in the **C2** field.

Autolab module

Main module PGSTAT302N

Power supply frequency 50 Hz

C1 1,61E-11

C2 5,43E-13

Automatic configuration

- 8 Close the hardware setup and wait for the Autolab to be reinitialized using the updated Hardware setup.

- **Input connection:** specifies if the measurement should be carried out internally (through the PGSTAT) or externally, using the external inputs provided on the front panel of the FRA32M module.
- **Wave type (Single sine, 5 sines or 15 sines):** specifies the type of signal used during the measurement. The choice is provided between the default single sine or the multi sine wave types.
- **Integration time:** specifies the time during which the signal is measured, in s.
- **Minimum number of cycles to integrate:** specifies the minimum number of cycles to integrate during the measurement.
- **FRA32M:** a toggle that can be used to switch the FRA32M module on or off.

FRA32M

Properties

Frequency Hz

Amplitude mV

Input connection ▼

Wave type ▼

Integration time ms

Minimum cycles to integrate

FRA32M

Results

Elapsed integration time 0.8 s (1 s)

E(DC) (V)	E(AC) (V)	i(DC) (A)	i(AC) (A)	% E	% i
999.7 m	10.78 m	-300.6 n	1.002 μ	64	10

Freq. (Hz)	Z (Ω)	-Phase ($^{\circ}$)	Z' (Ω)	-Z'' (Ω)
1.000	163.1 k	81.56	23.94 k	161.3 k

Figure 1235 The FRA32M manual control panel

As soon as the FRA32M is switched on using the provided toggle, the module will start a manual measurement using the properties specified in the panel. The measured values will be displayed after the measurement in the **Results** sub-panel.



NOTICE

The measurement will continue until the module is switched off using the provided  toggle.



NOTICE

For impedance measurements using the potentiostat/galvanostat, it is necessary to set the potential or current, specify the DC potential or current and select the appropriate current range using the instrument manual control panel *Autolab display panel* (see chapter 5.2.3, page 127).

16.3.2.13.6 FRA32M contour map

A typical contour map for the **FRA32M** module in combination with the **PGSTAT302N** potentiostat/galvanostat is shown in *Figure 1236*.




Figure 1236 Contour map of the PGSTAT302N/FRA32M combination

The map reported in *Figure 1236* shows a dark green area, which corresponds to the area of the map where an error of $\pm 0.3^\circ$ on the measured phase angle and $\pm 0.3\%$ on the measured impedance value is expected. The light green area corresponds to the area of the map where an error of $\pm 5^\circ$ on the measured phase angle and $\pm 2\%$ on the measured impedance value is expected.



NOTICE

The contour map is determined empirically with the maximum possible amplitude, in potentiostatic mode.

16.3.2.13.7 FRA32M module restrictions

No restrictions apply when using the FRA32M module.

16.3.2.13.8 FRA32M module front panel connections

The FRA32M module is fitted with three female SMB connectors, labeled $\rightarrow X$, $\rightarrow Y$ and $\leftarrow V$, from top to bottom (see figure 1237, page 1121).




Figure 1237 The front panel labels of the FRA32M module (left: FRA32M module in PGSTAT, right: FRA32M module in Multi Autolab)

The two connectors, labeled $\rightarrow X$ and $\rightarrow Y$ are input connectors that can be used to analyze external transfer functions. They have an input range of ± 10 V and an input impedance of 50Ω .

The signal provided through the $\leftarrow V$ connector on the front panel corresponds to the output of the sinewave generator of the FRA32M. Whenever the FRA32M module is used, the voltage provided at this output corresponds to the signal generated by the module (either single sine or multi sine).

The output signal is a voltage, referred to the instrument ground, corresponding to the applied amplitude, multiplied by 10 (when the instrument is working in potentiostatic mode), or the converted amplitude, multiplied by 10 (when the instrument is working in galvanostatic mode):



$$E_{out}(\leftarrow V) = E_{(FRA32M)} \cdot 10$$

$$E_{out}(\leftarrow V) = \frac{i_{(FRA32M)}}{[CR]} \cdot 10$$

Where $E_{out}(\leftarrow V)$ corresponds to the output voltage signal of the module, in V, $E_{(FRA32M)}$ and $i_{(FRA32M)}$ corresponds to the specified amplitude, in V or A, respectively and $[CR]$ is the active current range of the FRA32M module.



NOTICE

The front panel $\leftarrow V$ SMB output is provided for information purposes only except for impedance measurement involving external transfer functions.

16.3.2.13.9 FRA32M module testing

NOVA is shipped with a procedure which can be used to verify that the **FRA32M** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestFRA** procedure, provided in the NOVA 2.X installation folder (Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestFRA.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (c).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out an impedance spectroscopy measurement. During the measurement, the data is fitted using a R(RC) equivalent circuit. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1238*.





Figure 1238 The data measured by the TestFRA procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestFRA automatic evaluation of the data requires the following tests to succeed:

1. The fitted series resistance must be equal to $100 \Omega \pm 5\%$.
2. The fitted parallel resistance must be equal to $1000 \Omega \pm 5\%$.
3. The fitted parallel capacitance must be equal to $1 \mu\text{F} \pm 10\%$.
4. The calculated χ^2 must be smaller or equal to 0.01.

All four conditions must be valid for the test to succeed.

16.3.2.13.10 FRA32M module specifications

The specifications of the FRA32M module are provided in *Table 52*.

Table 52 Specifications of the FRA32M module

Specification	Value
Frequency range	10 μHz - 32 MHz
Frequency range in combination with Autolab PGSTAT	10 μHz - 1 MHz



Specification	Value
Frequency range in combination with ECI10M module	40 Hz - 10 MHz
Frequency resolution	0.003 %
Input range	± 10 V
Output amplitude, potentiostatic mode	0.2 mV to 350 mV (RMS)
Output amplitude, galvanostatic mode	0.0002 - 0.35 times current range, A (RMS)
Output amplitude, external mode	2 mV - 3.5 V (RMS)
Input resolution	14 bit

16.3.2.14 IME303 module

The IME303 is an external interface to the Princeton Applied Research PAR303(A) Stand.

The Princeton Applied Research PAR303(A) Stand is a polarographic stand which provides the means to perform electrochemical measurements using a mercury drop electrode. The mercury drops are formed at the very end of a narrow glass capillary. The Princeton Applied Research PAR303(A) Stand can operate in three different modes:

- **Dropping Mercury Electrode (DME):** in this mode mercury drops form at the end of the capillary. The drop grows until the weight of the drop exceeds the surface tension and the drop falls into the solution, leading to a new drop at end of the capillary.
- **Static Drop Mercury Electrode (SDME):** in this mode mercury drops are formed at the end of the capillary. The drop grows until a hardware-controlled size and it remains at the end of the capillary until the tapper, built into the Princeton Applied Research PAR303(A) Stand, is activated. This dislodges the mercury drop, which falls into the solution, leading to a new, identical drop at the end of the capillary.
- **Hanging Drop Mercury Electrode (HDME):** in this mode a single mercury drop is formed at the end of the capillary. The drop grows until a hardware-controlled size and it remains at the end of the capillary. The tapper, can be activated if needed to dislodge this drop and create a new drop at the end of the capillary.

The IME303 provides remote controls for the Princeton Applied Research PAR303(A) Stand. The following actions can be controlled through the IME303:

- **Stirrer on/off:** the optional stirrer connected to the Princeton Applied Research PAR303(A) Stand can be remotely switched on or off.

- **Purge on/off:** the purge function of the Princeton Applied Research PAR303(A) Stand can be remotely switched on or off.
- **Create new drop:** the tapper of the Princeton Applied Research PAR303(A) Stand can be remotely activated. This knocks the current mercury drop from the electrode and created a new drop.



NOTICE

For more information on the Princeton Applied Research PAR303(A) Stand, please consult the corresponding User Manual.



CAUTION

Take all necessary precautions when working with mercury. It is highly recommended to consult the Material Safety Data Sheet (MSDS) before operating the Princeton Applied Research PAR303(A) Stand. It is also recommended to dispose of the mercury waste properly.

16.3.2.14.1 IME303 module compatibility

The IME303 module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- PGSTAT100N and PGSTAT100
- PGSTAT101, M101, PGSTAT204 and M204
- PGSTAT20 and PGSTAT10
- μ Autolab II and μ Autolab III



NOTICE

The IME303 module is **not** compatible with the Autolab instruments not listed above.

16.3.2.14.2 IME303 module scope of delivery

Depending on the type of instrument it is connected to, the IME303 is the following items:




Figure 1239 The IME303 module settings



NOTICE

A 500 ms settling time is used each time the tapper is activated. This settling time can be adjusted in the hardware setup.

16.3.2.14.4 IME303 module manual control

The IME303 can be manually controlled, using the **Autolab display** provided in the instrument control panel (*see figure 1240, page 1128*). The dedicated manual control panel can be used to perform the following tasks:

- **Purge (on/off toggle):** this control can be used to switch the nitrogen purge of the Princeton Applied Research PAR303(A) Stand on or off.
- **Stirrer (on/off toggle):** this control can be used to switch the stirrer of the Princeton Applied Research PAR303(A) Stand on or off if the stirrer is installed.
- **Create new drop (button):** this control can be used to create a new drop by activating the tapper of the Princeton Applied Research PAR303(A) Stand.




Figure 1240 Manual control of the IME303 provided in the Autolab display



NOTICE

The manual control panel provided for the IME303 can be used to set or display the current settings of the IME303.

16.3.2.14.5 IME303 module restrictions

No restrictions apply when using the IME303 module.

16.3.2.14.6 IME303 module front panel controls

The front panel of the IME303 provides a number of controls and indicators, shown in Figure 1241.




Figure 1241 Overview of the front panel of the IME303

1 Power On/Off button

For switching the IME303 on or off.

2 New drop indicator LED

Flashes when the tapper of the Princeton Applied Research PAR303(A) Stand is activated to create a new drop.

3 Purge LED

Indicates that the nitrogen purge is on when lit.

4 Stirrer LED

Indicates that the stirrer is on when lit.

16.3.2.14.7 IME303 module back plane connections

The back plane of the IME663 provides a number of connections, shown in *Figure 1242*.




Figure 1242 Overview of the back plane of the IME303

- | | |
|--|--|
| <p>1 PAR303(A) connector
For connecting the IME303 to the J1 or Remote connector located on the back plane of the Princeton Applied Research PAR303(A) Stand.</p> | <p>2 DIO connector
For connecting the IME303 to the DIO connector of the Autolab.</p> |
| <p>3 SOLENOID - connector
For connecting the negative pole of a third party tapper.</p> | <p>4 SOLENOID + connector
For connecting the positive pole of a third party tapper.</p> |
| <p>5 Mains voltage indicator
Indicates the mains voltage settings of the IME303.</p> | <p>5 Mains connection socket
For connecting the IME303 to the mains supply.</p> |



NOTICE

A third party drop tapper can be controlled with the IME303. Use the SOLENOID marked banana sockets located on the back plane of the IME303. Every time a new drop is created (through the a command or manual control of the IME303) a 15 V pulse will be generated between the + and - connectors.



CAUTION

Make sure that the mains voltage indicator is set properly before switching the IME303 on.

The IME663 provides manual and remote controls for the Metrohm 663 VA Stand. The following actions can be controlled through the IME663:

- **Stirrer on/off:** the stirrer of the Metrohm 663 VA Stand can be manually or remotely switched on or off.
- **Purge on/off:** the purge function of the Metrohm 663 VA Stand can be remotely switched on or off.
- **Create new drop:** the tapper of the Metrohm 663 VA Stand can be remotely activated. This knocks the current mercury drop from the electrode and created a new drop.



NOTICE

Remote control of the Metrohm 663 VA Stand only works when the stand is set to **SDME** mode.



NOTICE

For more information on the Metrohm 663 VA Stand, please consult the corresponding User Manual.



CAUTION

Take all necessary precautions when working with mercury. It is highly recommended to consult the Material Safety Data Sheet (MSDS) before operating the Metrohm 663 VA Stand. It is also recommended to dispose of the mercury waste properly.

16.3.2.15.1 IME663 module compatibility

The IME663 module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- PGSTAT100N and PGSTAT100
- PGSTAT101, M101, PGSTAT204 and M204
- PGSTAT20 and PGSTAT10
- μ Autolab II and μ Autolab III

16.3.2.15.3 IME663 module settings

The IME663 module settings are defined in the NOVA software. The following user-definable setting is available, through the **Autolab control** command (see figure 1244, page 1133):

- **Purge:** this control can be used to switch the nitrogen purge of the Metrohm 663 VA Stand on or off.
- **Stirrer:** this control can be used to switch the stirrer of the Metrohm 663 VA Stand on or off if the stirrer switch on the IME663 is set to remote control.
- **Number of new drops:** this control can be used to create the specified number of new drops by activating the tapper of the Metrohm 663 VA Stand as many times.




Figure 1244 The IME663 module settings



NOTICE

A 500 ms settling time is used each time the tapper is activated. This settling time can be adjusted in the hardware setup.

16.3.2.15.4 IME663 module manual control

The IME663 can be manually controlled, using the **Autolab display** provided in the instrument control panel (see figure 1245, page 1134). The dedicated manual control panel can be used to perform the following tasks:

- **Purge (on/off toggle):** this control can be used to switch the nitrogen purge of the Metrohm 663 VA Stand on or off.




Figure 1246 Overview of the front panel of the IME663

1 Power On/Off button

For switching the IME663 on or off.

2 New drop indicator LED

Flashes when the tapper of the Metrohm 663 VA Stand is activated to create a new drop.

3 Purge LED

Indicates that the nitrogen purge is on when lit.

4 Stirrer On/Remote button

For switching the stirrer on or to software control. When the button is engaged, the stirrer is on. When the button is disengaged, the stirrer is controlled by the software.

16.3.2.15.7 IME663 module back plane connections

The back plane of the IME663 provides a number of connections, shown in Figure 1247.




Figure 1247 Overview of the back plane of the IME663

1 VA STAND connector


For connecting the IME663 to the cartridge connector of the Metrohm 663 VA Stand (labeled H) and the ground and instrument front panel (in the case of PGSTAT10, PGSTAT20, μ Autolab II and μ Autolab III).

2 DIO connector

For connecting the IME663 to the DIO connector of the Autolab.

led in the Metrohm 663 VA Stand. The control provides the choice between four positions:

- **HDME:** sets the MME/MME PRO to *hanging mercury drop electrode* mode. When this mode is selected, a drop can be manually created by using the **man** switch located next to the rotary knob. Activating this switch triggers the built-in tapper once.
- **0:** the MME/MME PRO is switched off.
- **SDME:** sets the MME/MME PRO to *static mercury drop electrode* mode. When this mode is selected, the MME/MME PRO can be remotely controlled by NOVA using the IME663.
- **DME:** sets the MME/MME PRO to *dropping mercury electrode* mode.



4 deaeration

This switch can be used to manually switch the N₂ purge on or off. This control can be used to overrule the purge control provided through the NOVA software controls.




deaeration

16.3.2.15.9 IME663 and Metrohm 663 VA Stand installation

The IME663 and Metrohm 663 VA Stand can be used in combination with any compatible instrument. The installation and configuration can be carried out by the end-user at any time.

Depending on the type of instrument it is connected to, the IME663 and Metrohm 663 VA Stand have to be installed according to a specific procedure:

1. For the PGSTAT10, PGSTAT20, μ Autolab II and μ Autolab III, please refer to *Chapter 16.3.2.15.9.1*.
2. For the PGSTAT101, M101, PGSTAT204 and M204, please refer to *Chapter 16.3.2.15.9.2*.
3. For all other instruments, please refer to *Chapter 16.3.2.15.9.3*.



16.3.2.15.9.2 IME663 and Metrohm 663 VA Stand installation (PGSTAT101, M101, PGSTAT204, M204)

The following steps describe how to install the IME663 and Metrohm 663 VA Stand in combination with an Autolab PGSTAT. These steps apply to the PGSTAT101, M101, PGSTAT204 and M204.

1 Connect the DIO cable

Connect the DIO cable, supplied with the IME663, to the DIO connector, located on the back plane of the IME663 (item 2 in *Figure 1247*). Connect the other end of the cable to the DIO connector located on the front panel of the Autolab PGSTAT.

2 Connect the Stirrer cable

Connect the Stirrer cable, supplied with the IME663, to the STIRRER connector, located on the back plane of the IME663 (item 3 in *Figure 1247*). Connect the other end of the cable to the E connector located on the back plane of the Metrohm 663 VA Stand.

3 Connect the VA Stand cable

Connect the VA Stand cable, supplied with the IME663, to the VA STAND connector, located on the back plane of the IME663 (item 1 in *Figure 1247*). Connect the other end to the H cartridge connector located on the back plane of the Metrohm 663 VA Stand.

2 Connect the Stirrer cable

Connect the Stirrer cable, supplied with the IME663, to the STIRRER connector, located on the back plane of the IME663 (item 3 in *Figure 1247*). Connect the other end of the cable to the E connector located on the back plane of the Metrohm 663 VA Stand.

3 Connect the VA Stand cable to the Metrohm 663 VA Stand

Connect the VA Stand cable, supplied with the IME663, to the VA STAND connector, located on the back plane of the IME663 (item 1 in *Figure 1247*). Connect the other end to the H cartridge connector located on the back plane of the Metrohm 663 VA Stand.

4 Remove the cell cable

Remove the cell cable from the front panel of the Autolab instrument.



NOTICE


It is recommended to store this cable carefully for future use.

5 Connect the VA Stand cable to the Autolab

Connect the cell cable connector, embedded in the VA Stand cable to the front panel of the Autolab.

6 Specify the hardware setup

Adjust the hardware setup and specify the DIO connector (P1 or P2) in the **Properties** panel.



16.3.2.16 MUX module

The MUX is an optional module for the Autolab PGSTAT and the Multi Autolab. With the MUX module, it is possible to multiplex the PGSTAT electrode connections. Depending on the type of MUX used, measurements on several electrochemical cells or working electrodes can be done sequentially

The MUX module is available in three standard configurations:

- MUX-MULTI4:** in this configuration, shown in *Figure 1248*, the WE, S, RE and CE leads from the PGSTAT are multiplexed to 4 (or more) sets of connections (see Figure 1). This means that sequential measurements can be performed on multiple independent electrochemical cells. It is possible to add up to 16 MUX-MULTI4 boxes in series in order to multiplex up to 64 independent electrochemical cells.




Figure 1248 The MUX-MULTI4 switch box



NOTICE

The MUX-MULTI4 is typically used for sequential measurements on individual electrochemical cells or for measurements requiring different electrode arrangements in the same cell (like Van der Pauw measurements).

- MUX-SCNR8:** in this configuration, shown in *Figure 1249*, the RE and S leads from the PGSTAT are multiplexed to 8 (or more) sets of connections (see *Figure 2*). This means that the differential amplifier of the PGSTAT can be multiplexed across as many individual cells. It is possible to add up to 16 MUX-SCNR8 boxes in series in order to multiplex up to 128 sets of RE and S electrodes.




Figure 1249 The MUX-SCNR8 switch box



NOTICE

The MUX-SCNR8 is intended to be used for measurements on individual cells in a series of cells. The most common example is the measurement of cell voltages of individual cells or sections of cells in a fuel cell stack.

- MUX-SCNR16:** in this configuration, shown in *Figure 1250*, the WE lead from the PGSTAT is multiplexed to 16 (or more) WE connections (see *Figure 3*). This means that sequential measurements can be performed on multiple working electrodes sharing a common reference electrode and counter electrode. It is possible to add up to 16 MUX-SCNR16 boxes in series in order to multiplex up to 255 working electrodes.




Figure 1250 The MUX-SCNR16 switch box



NOTICE

The MUX-SCNR16 is intended to be used for measurements on individual working electrodes contained in a single electrochemical cell. Corrosion and sensor measurements usually benefit from this hardware extension.



CAUTION

Whenever the active channel of the MUX is changed, a 8 ms settling time is added before the next command is executed.

16.3.2.16.1 MUX module compatibility

The MUX module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- M101
- PGSTAT204/M204
- PGSTAT20
- PGSTAT10



NOTICE

The MUX module is **not** compatible with the Autolab instruments not listed above.

16.3.2.16.2 MUX module scope of delivery

The MUX module is supplied with the following **common** items:

- MUX module.
- MUX module label.
- Control cable.
- Connection plugs.

Depending on the type of MUX, the following **specific** items are provided:

- For the **MUX-MULTI4**
 - Electrode connection cable for WE, S, CE and RE.
 - MULTI4 switch box.
 - 16 SMB to banana cables (4 labeled WE1 to WE4, 4 labeled S1 to S4, 4 labeled CE1 to CE4 and 4 labeled RE1 to RE4).
 - 16 alligator clips (8 red and 8 black).
- For the **MUX-SCNR8**
 - Electrode connection cable for S and RE.
 - SCNR8 switch box.
 - 16 SMB to banana cables (8 labeled RE1 to RE8 and 8 labeled S1 to S8).
 - 16 alligator clips (8 red and 8 black).
- For the **MUX-SCNR16**
 - Electrode connection cable for WE.
 - SCNR16 switch box.
 - 16 SMB to banana cables (labeled WE1 to WE16).
 - 16 red alligator clips.

16.3.2.16.3 MUX hardware setup

To use the **MUX** module, the hardware setup needs to be adjusted. The checkbox for the module needs to be ticked (*see figure 1251, page 1146*).




Figure 1252 The MUX module settings



NOTICE

When the Channel is set to 0, no MUX channel is selected and the MUX is bypassed.

16.3.2.16.5 MUX module manual control

The MUX can be manually controlled, using the **Autolab display** provided in the instrument control panel (see figure 1253, page 1147). The dedicated manual control panel can be used to perform the following tasks:

- **Active channel:** this control can be used to define the active channel of the MUX. A value between 0 and the number of channels defined in the hardware setup can be specified.




Figure 1253 Manual control of the MUX provided in the Autolab display

The value can be specified directly as an integer or through the provided slider control (see figure 1254, page 1148).




Figure 1255 The front panel labels of the MUX module (left: MUX module in PGSTAT, right: MUX module in Multi Autolab)

The MUX-MULTI4 switch box has the following connections (see figure 1256, page 1149):




Figure 1256 The connections provided by the MUX-MULTI4 switch box

1 Digital control connector

A 15 pin male SUB-D connector used to connect the digital control cable from the MUX module.

2 Electrode connection cable connector

A 9 pin male SUB-D connector used to connect the electrode connection cable providing inputs for the WE, S, CE and RE from the Autolab PGSTAT.

3 Digital control extension connector

A 15 pin female SUB-D connector used to extend the digital control of the MUX to the adjacent switch box.

4 Electrode connection cable extension connector

A 9 pin female SUB-D connector used to extend the electrode connections from the Autolab PGSTAT to the adjacent switch box.



5 WE1 to WE4 connections

Output connections for attaching the WE1 to WE4 connections cables provided with the MULTI4 switch box.

6 S1 to S4 connections

Output connections for attaching the S1 to S4 connections cables provided with the MULTI4 switch box.

7 CE1 to CE4 connections

Output connections for attaching the CE1 to CE4 connections cables provided with the MULTI4 switch box.

8 RE1 to RE4 connections

Output connections for attaching the RE1 to RE4 connections cables provided with the MULTI4 switch box.

The MUX-SCNR8 switch box has the following connections (see figure 1257, page 1150):




Figure 1257 The connections provided by the MUX-SCNR8 switch box

1 Digital control connector

A 15 pin male SUB-D connector used to connect the digital control cable from the MUX module.

2 Electrode connection cable connector

A 9 pin male SUB-D connector used to connect the electrode connection cable providing inputs for the RE and S from the Autolab PGSTAT.

3 Digital control extension connector

A 15 pin female SUB-D connector used to extend the digital control of the MUX to the adjacent switch box.

4 Electrode connection cable extension connector

A 9 pin female SUB-D connector used to extend the electrode connections from the Autolab PGSTAT to the adjacent switch box.

5 RE1 to RE8 connections

Output connections for attaching the RE1 to RE8 connections cables provided with the SCNR8 switch box.

6 S1 to S8 connections

Output connections for attaching the S1 to S8 connections cables provided with the SCNR8 switch box.

The MUX-SCNR16 switch box has the following connections (see figure 1258, page 1151):






Figure 1258 The connections provided by the MUX-SCNR16 switch box

1 Digital control connector

A 15 pin male SUB-D connector used to connect the digital control cable from the MUX module.

2 Electrode connection cable connector

A 9 pin male SUB-D connector used to connect the electrode connection cable providing inputs for the WE from the Autolab PGSTAT.

3 Digital control extension connector

A 15 pin female SUB-D connector used to extend the digital control of the MUX to the adjacent switch box.

4 Electrode connection cable extension connector

A 9 pin female SUB-D connector used to extend the electrode connections from the Autolab PGSTAT to the adjacent switch box.

5 WE1 to WE16 connections

Output connections for attaching the WE1 to WE16 connections cables provided with the SCNR16 switch box.

16.3.2.16.8 MUX module testing

Two test procedures are provided for testing the MUX module:

- For the MULTI4 and the SCNR16 module option, please refer to *Chapter 16.3.2.16.8.1*.
- For the SCNR8 module option, please refer to *Chapter 16.3.2.16.8.2*.

16.3.2.16.8.1 MUX-SCNR16 and MUX-MULTI4 module testing

NOVA is shipped with a procedure which can be used to verify that the MUX module in combination with the **SCNR16** and the **MULTI4** option is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure


Load the **TestMUX-SCNR16-MULTI4** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestMUX-SCNR16-MULTI4.nox)

2 Connect the Autolab dummy cell

For the SCNR16 option, connect the MUX channel 1 and 2 to the Autolab dummy cell circuit (a) and (c).



For the MULTI4 option, connect the MUX channel 1 and 2 to the Autolab dummy cell circuit (a) and (c).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1259*.





Figure 1259 The results of the TestMUX-SCNR16-MULTI4 procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestMUX-SCNR16-MULTI4 automatic evaluation of the data requires the following tests to succeed:

1. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$ for channel 1 and must be equal to $1100 \Omega \pm 5 \%$ for channel 2.
2. The intercept of the measured current versus the applied potential must be equal to $\pm 4 \text{ mV}$ divided by $1000100 \Omega \pm 5 \%$ for channel 1 and must be equal to $\pm 4 \text{ mV}$ divided by $1100 \Omega \pm 5 \%$ for channel 2.

All four conditions must be valid for the test to succeed.





Figure 1260 The results of the TestMUX-SCNR8 procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestMUX-SCNR8 automatic evaluation of the data requires the following tests to succeed:

1. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5\%$ for both channels.
2. The intercept of the measured current versus the applied potential must be equal to $\pm 4 \text{ mV}$ divided by $1000100 \Omega \pm 5\%$ for both channels.

All four conditions must be valid for the test to succeed.

16.3.2.16.9 MUX module specifications

The MUX module is available in three versions:

- **MUX-MULTI4:** can be used to multiplex all four electrodes provided by the Autolab: CE, RE, S and WE. Using this module, it is possible to work on different electrochemical cells, sequentially.
- **MUX-SCNR8:** can be used to multiplex two independent electrodes provided by the Autolab. The most common use of this special module is to multiplex the S and RE connection from the differential amplifier in order to measure the potential difference across several electrochemical interfaces, sequentially.



NOTICE

The pX module is no longer available and it is now replaced by its successor module, the pX1000.

The pX module adds the following signals to the signal sampler (see figure 1261, page 1157):

- **pX.Voltage (V):** this signal corresponds to the potential difference measured by the pX module.
- **pX.pH:** this signal corresponds to the converted pH value, determined based on the voltage value measured by the pX.

The screenshot shows the 'Sampler' interface with a table of signal options. The 'pX(1).Voltage' and 'pX(1).pH' signals are highlighted in green, indicating they are selected. A mouse cursor is pointing at the 'pX(1).pH' signal.

Signal	Sample	Average	d/dt
WE(1).Current	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
WE(1).Potential	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Power	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Resistance	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
WE(1).Charge	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
pX(1).Voltage	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
pX(1).pH	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Time	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Sample alternating

Figure 1261 The pX module provides the pX(1).Voltage and pX(1).pH signals

16.3.2.17.1 pX module compatibility

The pX module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- PGSTAT100N and PGSTAT100
- PGSTAT20
- PGSTAT10



NOTICE

The pX module is **not** compatible with the Autolab instruments not listed above.

16.3.2.17.4 pX module settings

The pX module settings are controlled using the parts supplied with the module. When the sensor connected to the pX module is used outside of the cell containing the working electrode, a 50 Ohm BNC resistor must be connected to the \ominus G BNC connector located on the pX module, above the \rightarrow E input plug. This resistor is used to ground the pX module.

When the sensor connected to the pX module is located inside the electrochemical cell, the 50 Ohm BNC resistor must be removed.



NOTICE

If the values measured through the pX module are unusually noisy, please verify that you are using the correct configuration.



CAUTION

When the sensor connected to the pX module is located in the electrochemical cell, the potential difference measured by the module during a measurement will be affected by the changes in electrical field in the solution generated by the polarization of the electrodes in the cell. Depending on several experimental factors like temperature, conductivity, ionic strength, cell and electrode geometry, etc., the effect on the potential difference measured by the pX module can be more or less pronounced. Trial and error may be required to determine the experimental conditions required to minimize these effects.

16.3.2.17.5 pX module restrictions

No restrictions apply when using the pX module.

16.3.2.17.6 pX module front panel connections

The pX module is fitted with two female BNC connectors, labeled \ominus G and \rightarrow E, respectively (see figure 1263, page 1159).




Figure 1263 The front panel label of the pX module

Both connectors located on the front panel of the pX module are input connectors. The \odot G connector is used ground the input circuit of the pX module, using the provided 50 Ω plug. The \rightarrow E is used to connect the pH sensor electrode. A BNC to LEMO converter is supplied with the module.



NOTICE

The pX module must be grounded when the pH sensor is not located in the same cell as the working electrode.

16.3.2.17.7 pX module testing

NOVA is shipped with a procedure which can be used to verify that the pX module is working as expected.



CAUTION

The pX module must be grounded using the provided 50 Ω resistor plug.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestpX** procedure, provided in the NOVA 2.X installation folder (Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestpX.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a) and the pX module to the Autolab dummy cell (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1264*.





Figure 1264 The results of the TestpX procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestpX automatic evaluation of the data requires the following tests to succeed:

1. The inverted slope of the measured $pX(1).Potential$ versus the applied potential must be equal to $1 \pm 5\%$.
2. The intercept of the measured $pX(1).Potential$ versus the applied potential must be equal to ± 0.004 .




Figure 1265 The pX1000 module provides the pX1000(1).Temperature, pX1000(1).Voltage and pX1000(1).pH signals

16.3.2.18.1 pX1000 module compatibility

The pX1000 module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- PGSTAT100N and PGSTAT100
- M101
- PGSTAT204/M204



NOTICE

The pX1000 module is **not** compatible with the Autolab instruments not listed above.

16.3.2.18.2 pX1000 module scope of delivery

The pX1000 module is supplied with the following items:

- pX1000 module
- pX1000 module label
- V+ connection cable
- V- connection cable
- Lemo connector cover

- **Cell setup:** this drop-down control can be used to define the location of the sensor connected to the pX1000 module. In the Autolab potentiostat, the working electrode is connected to virtual ground. When the sensor connected to the pX1000 module is located in the same cell as the working electrode, the inputs of the pX1000 must be floating, to avoid leakage or ground loops.




Figure 1267 The pX1000 module settings



NOTICE

If the values measured through the pX1000 module are unusually noisy, please verify that you are using the correct configuration.



CAUTION

When the sensor connected to the pX1000 module is located in the electrochemical cell, the potential difference measured by the module during a measurement will be affected by the changes in electrical field in the solution generated by the polarization of the electrodes in the cell. Depending on several experimental factors like temperature, conductivity, ionic strength, cell and electrode geometry, etc., the effect on the potential difference measured by the pX1000 module can be more or less pronounced. Trial and error may be required to determine the experimental conditions required to minimize these effects.

16.3.2.18.5 pX1000 module restrictions


No restrictions apply when using the pX1000 module.

1 Load the procedure

Load the **TestpX1000** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestpX1000.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a) and the pX1000 module to the Autolab dummy cell (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1269*.

16.3.2.18.8 pX1000 module specifications

The specifications of the pX1000 module are provided in *Table 55*.

Table 55 Specifications of the pX1000 module

Specification	Value
Input impedance	> 1000 G Ω
Input range	± 10 V
Module configuration	Software control
Pt1000/NTC temperature sensor	Yes
Potential resolution	30 μ V (gain 10)
Potential accuracy	± 2 mV
Temperature resolution	0.015 $^{\circ}$ C
Temperature accuracy	± 0.5 $^{\circ}$ C

16.3.2.19 SCAN250 module

The SCAN250 module is an extension module for the Autolab PGSTAT. This module provides a true linear scan generator which can be used to perform linear scan cyclic voltammetry measurements. Linear scan cyclic voltammetry measurement use a perfectly smooth potential scan, defined only by a scan range, in mV and a scan rate, in V/s. Unlike the cyclic voltammetry staircase method, the sampling of the electrochemical signals is decoupled from the scan generation.

Using the SCAN250 module it is possible to measure both the faradic and the capacitive currents during a cyclic voltammetry. Therefore, this particular form of cyclic voltammetry is well suited for monitoring changes in the double layer structure or electrochemical reactions characterized by fast electron transfer kinetics.



NOTICE

The SCAN250 can generate a potential scan between 10 mV/s and 250 kV/s.

Depending on the hardware configuration, the SCAN250 can be used in two different scan rate ranges:

- **Without the ADC10M or the ADC750 module:** if the instrument is **not** fitted with the ADC10M and ADC750 module it is possible to use the SCAN250 module up to 250 V/s.




Figure 1270 The SCAN250 module is selected in the hardware setup

16.3.2.19.4 SCAN250 module settings

The SCAN250 module does not have any user-definable settings. All the module settings are directly controlled through the **CV linear scan** command.


16.3.2.19.5 SCAN250 module restrictions

Restrictions apply when using the SCAN250 module:

- **Potentiostatic control only:** the SCAN250 module can only be used in Potentiostatic mode.
- **Vertex overshoot or undershoot:** The SCAN250 module is fitted with a vertex detection circuit which is used to synchronize the reversal of the scan direction with the detection of the vertex. At very high scan rates, this circuit can cause the linear scan generator module to stop the scan at a potential value that no longer matches the experimental parameters. This potential difference depends on the scan rate.

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1272*.




Figure 1272 The results of the TestSCAN procedure



NOTICE

The SCANGEN module is no longer available and it is now replaced by its successor module, the SCAN250.

Using the SCANGEN module it is possible to measure both the faradic and the capacitive currents during a cyclic voltammetry. Therefore, this particular form of cyclic voltammetry is well suited for monitoring changes in the double layer structure or electrochemical reactions characterized by fast electron transfer kinetics.



NOTICE

The SCANGEN can generate a potential scan between 10 mV/s and 10 kV/s.

Depending on the hardware configuration, the SCANGEN can be used in two different scan rate ranges:

- **Without the ADC10M or the ADC750 module:** if the instrument is **not** fitted with the ADC10M and ADC750 module it is possible to use the SCANGEN module up to 250 V/s.
- **With the ADC10M or ADC750 module:** if the instrument is fitted with the ADC10M *ADC10M module (see chapter 16.3.2.1, page 998)* or the ADC750 *ADC750 module (see chapter 16.3.2.2, page 1004)* module it is possible to use the SCANGEN module up to 10 kV/s.



NOTICE

The SCANGEN module only works in potentiostatic mode.

16.3.2.20.1 SCANGEN module compatibility

The SCANGEN module is available for the following instruments:

- PGSTAT302N, PGSTAT302 and PGSTAT30
- PGSTAT128N and PGSTAT12
- PGSTAT100N and PGSTAT100
- PGSTAT20
- PGSTAT10

16.3.2.20.4 SCANGEN module settings

The SCANGEN module does not have any user-definable settings. All the module settings are directly controlled through the **CV linear scan** command.

16.3.2.20.5 SCANGEN module restrictions

Restrictions apply when using the SCANGEN module:

- **Potentiostatic control only:** the SCANGEN module can only be used in Potentiostatic mode.
- **Vertex overshoot or undershoot:** The SCANGEN module is fitted with a vertex detection circuit which is used to synchronize the reversal of the scan direction when the vertex is reached. At very high scan rates, this circuit can cause the linear scan generator module to stop the scan at a potential value that no longer matches the experimental parameters. This potential difference depends on the scan rate.

16.3.2.20.6 SCANGEN module front panel connections

The SCANGEN module is fitted with a single female BNC connector, labeled $\leftarrow V$ (see figure 1274, page 1177).




Figure 1274 The front panel label of the SCANGEN module

The signal provided through the $\leftarrow V$ connector on the front panel corresponds to the output of scan generator signal of the SCANGEN module. The output signal is a voltage, referred to the instrument ground, corresponding to the specified potential scan.



NOTICE

The output of the SCANGEN module does not include the offset potential which is created by the DAC164 module of the Autolab potentiostat/galvanostat.



NOTICE

The front panel ←V BNC output is provided for information purposes only.

16.3.2.20.7 SCANGEN module test

NOVA is shipped with a procedure which can be used to verify that the integrator circuit of the **SCANGEN** module is working as expected.


Follow the steps described below to run the test procedure.

1 Load the procedure

Load the **TestSCAN** procedure, provided in the NOVA 2.X installation folder (\Metrohm Autolab\NOVA 2.X\SharedDatabases\Module test\TestSCAN.nox)

2 Connect the Autolab dummy cell

Connect the PGSTAT to the Autolab dummy cell circuit (a).



3 Start the procedure

Start the procedure and follow the instructions on-screen. The test carries out a cyclic voltammetry measurement. At the end of the measurement, the measured data will be processed and a message will be shown. The measured data should look as shown in *Figure 1275*.





Figure 1275 The results of the TestSCAN procedure

4 Test evaluation

If the test is successful, a message will be shown at the end of the measurement.



The TestSCAN automatic evaluation of the data requires the following tests to succeed:

1. The inverted slope of the measured current versus the applied potential must be equal to $1000100 \Omega \pm 5 \%$.
2. The calculated capacitance, determined from the integrated current signal, must be equal to $1 \mu\text{F} \pm 10 \%$.

Both conditions must be valid for the test to succeed.

**16.3.2.20.8 SCANGEN module specifications**

The specifications of the SCANGEN module are provided in *Table 57*.

Table 57 Specifications of the SCANGEN module

Specification	Value
Scan range	± 5 V
Vertex resolution	2.5 mV
Vertex accuracy	± 5 mV
Output offset	± 1 mV
Minimum scan rate	10 mV/s
Maximum scan rate	10 kV/s

17 Diagnostics

NOVA is supplied with a **diagnostics** tool that can be used to test the Autolab instrument. This tool is provided as a standalone application and can be accessed from the start menu, in the Autolab group (Start menu – All programs – Autolab – Tools).

The diagnostics tool can be used to troubleshoot an instrument or perform individual tests to verify the correct operation of the instrument. Depending on the instrument type, the following items are required:

- **μAutolab type II, μAutolab type III and μAutolab type III/FRA2:** the standard Autolab dummy cell. For the diagnostics test, the circuit (a) is used.
- **PGSTAT101 and M101 module:** the internal dummy cell is used during the test, no additional items are required.
- **PGSTAT204 and M204 module:** the standard Autolab dummy cell. For the diagnostics test, the circuit (a) is used.
- **Other Autolab PGSTAT instruments:** the standard Autolab dummy cell and a 50 cm BNC cable. For the diagnostics test, the circuit (a) is used. The BNC cable must be connected between the ADC164 channel 2 and the DAC164 channel 2 on the front panel of the instrument .



NOTICE

The **PGSTAT302F** must be tested in normal (grounded) mode.

17.1 Connecting the instrument

The Diagnostics application supports multiple Autolab instruments. When the application starts it detects all available instruments connected to the computer.



NOTICE

If only one instrument is detected, the Diagnostics application will start immediately.

If more than one instrument is detected, a selection menu is displayed (see *figure 1276, page 1182*).



NOTICE

Instruments with serial number beginning with AUT9 or with μ 2AUT7, connected through an external USB interface, are identified by the serial number of the interface, USB7XXXX (see figure 1276, page 1182).

Select one of the available instruments and click the OK button to continue.



NOTICE

The Diagnostics application can only test one instrument or one potentiostat/galvanostat module at a time.

17.2 Running the Diagnostics

When the application is ready, a series of tests can be performed on the connected instrument. Pressing the start button will initiate all the selected tests. A visual reminder will be displayed at the beginning of the test, illustrating the connections required for the test (see figure 1278, page 1184).




Figure 1278 A visual reminder is shown at the beginning of the Diagnostics test

During the test, the progress will be displayed and a successful test will be indicated by a green symbol and a progress bar (see figure 1279, page 1185).




Figure 1279 The diagnostics report after all the tests have been performed successfully

If one or more of the tests fails, a red symbol will be used to indicate which test failed and what the problem is (see figure 1280, page 1185).




Figure 1280 A failed test will be indicated in the diagnostics tool

The following options are available through the Menu in the Diagnostics application:

File	
Save Report As...	Used to save the Diagnostics report as a TXT file.
Print Report	Used to print the Diagnostics report.
Exit	Closes the Diagnostics application.
Edit	
Select All Tests	Activates all available tests provided by the Diagnostics application.
Deselect Optional Tests	Deactivates all optional tests provided by the Diagnostics application.
Hardware Setup	Opens the Hardware Setup window.
Select Instrument	Opens the instrument selection dialog window.
Tools	
Checksum error wizard	Opens the checksum error wizard. This option is only used for service purposes.

17.5 Firmware update

For some instruments, a firmware update may be required. If this is the case an update message will be displayed (see figure 1283, page 1187).




Figure 1283 An update message is displayed when the outdated firmware is detected



Clicking the Yes button when prompted will update the firmware. The firmware update is permanent. The Firmware Update window will close automatically at the end of the update and the diagnostics test will continue.



CAUTION

Do **not** switch off the instrument or disconnect the instrument during the firmware update since this will damage the instrument.

18 Warranty and conformity

This chapter provides information about warranty, safety and conformity of the Autolab product range.

18.1 Warranty

The warranty on Autolab products is limited to defects or failures that are traceable to material, construction or manufacturing errors, which occur within 36 months from the day of delivery. In this case, the defects or failures will be rectified by Metrohm Autolab free of charge. Transport costs are to be paid by the customer, if applicable.



Glass breakage in the case of electrodes, cells or other parts is not covered by the warranty. Consumables (electrodes, QCM crystals, etc.) are not covered by the warranty.

If damage of the packaging is evident on receipt of the goods or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report is demanded. Lack of an official damage report releases Metrohm Autolab from any liability to pay compensation.

If any instruments or parts have to be returned, the original packaging should be used. This applies to all instruments, electrodes, cells and other parts. If the original packaging is not available it can be ordered at Metrohm Autolab or at your local distributor. For damage that arises as a result of non-compliance with these instructions, no warranty responsibility whatsoever will be accepted by Metrohm Autolab.

Do not modify the cell cable or the differential amplifier cable connectors. These cables are designed for the best possible operation. Modifications of these connections, i.e. with other connectors, will lead to the loss of any warranty.



CAUTION

There are no user-serviceable parts inside. Servicing must only be done by qualified personnel.



WARNING

Removal of the protective external panels results in exposure to potentially dangerous voltages; this must only be done by qualified personnel. The instrument must be disconnected from all power sources before removing the protective panels.

- **Fuses:** Replace blown fuses only with size and rating stipulated on or near the fuse panel holder and in the manual.
- **Cords:** Replace or repair faulty or frayed insulation on power cords and control cables. Replace control cables only with original spare parts. When replacing power cord, use only an approved type that conforms to local regulations. Do not move the instrument with power cords connected.
- **Grounding:** Ensure that power cords are plugged into the correct mains voltage source and always use a wall outlet with protective earth. Check all connected equipment and accessories for proper grounding.



WARNING

The PGSTAT100N is fitted with a control amplifier capable of generating up to 100 V potential difference between the counter electrode (CE) and the working electrode (WE).

- Take all necessary precautions when working with this instrument and use the supplied warning laminated sheet to warn others.

Means of carrying the instrument

This section describes the safe method for carrying and lifting the Autolab instruments. These guidelines must be followed when moving the N and F series potentiostats, as well as the Booster 10A and Booster20A.



CAUTION

The potentiostats can weigh up to 20 kg; the boosters can weigh up to 25 kg. If you cannot move the instrument alone, find someone who can assist you.

- Stand with your feet positioned one shoulder width apart. Be sure the ground is level and free of all obstacles, then put one foot next to the instrument and the other in front of it..
- Bend with your knees down to the instrument. Keep your back straight by tucking your chin into your chest.
- Grab the instrument with the handles on the side.
- Shift your weight into your feet and lift using your leg muscles; push up from your legs as well. Make sure your heels are firmly planted into the ground.
- Keep your arms and elbows close to your body and the box/object close to your chest.
- Choose the shortest, most direct route to your destination and avoid twisting your body as you walk.
- When you reach your destination, lower the instrument by bending at the knees. Do not bend from the waist, as this can strain your lower back muscles.

18.4 EU Declaration of conformity

This chapter provides the following certificates of conformity:

- For the Autolab N and F series instruments, namely: PGSTAT128N, PGSTAT302N, PGSTAT100N, and PGSTAT302F *EU Declaration of Conformity (see chapter 18.4.1, page 1192)*.
- For the Autolab compact series instruments, namely: PGSTAT101 and PGSTAT204 *EU Declaration of Conformity (see chapter 18.4.2, page 1195)*.
- For the Autolab multichannel instruments, namely: MAC101 and MAC204 *Declaration of Conformity (see chapter 18.4.3, page 1197)*.
- For Autolab accesories, namely: Motorcontroller, IME663, and Optical bench *Declaration of Conformity (see chapter 18.4.4, page 1199)*.

18.4.1 EU Declaration of Conformity

This declaration attests the compliance of the instrument with the standard specifications for electrical instruments and accessories and with the standard specifications for safety and system validation of the manufacturing company.

Name of commodity

Autolab PGSTAT128N, PGSTAT302N, PGSTAT100N,

PGSTAT302F used or not used in conjunction with one or more of the modules listed in Appendix A: Mutations

Modular Potentiostat/galvanostat laboratory measurement equipment used for measuring electrochemical cells.

Directives

These instruments have the CE marking and comply with the following EU directives:



- 2014/35/EU (Low Voltage Directive, LVD)
- 2014/30/EU (EMC Directive, EMC)
- 2011/65/EU (Directive for certain hazardous substances, RoHS)

Safety specifications

These instruments comply with safety specifications according to the following standards:

- | | |
|-------------------------|---|
| Design and type testing | <ul style="list-style-type: none"> ▪ EN 61010-1:2016
<i>Safety requirements for electrical equipment for measurement, control, and laboratory use.</i> |
| Testing in production | <p>Every instrument is routine-tested in the production division according to:</p> <ul style="list-style-type: none"> ▪ EN/IEC 61010-1 Appendix F
<i>Check of the protective conductor connection and the insulation from power circuit.</i> |

Electromagnetic compatibility (EMC)

These instruments have been built and have undergone final type testing according to the standards:

- | | |
|--------------|--|
| Requirements | <ul style="list-style-type: none"> ▪ EN 61326-1:2013
<i>Electrical equipment for measurement, control, and laboratory use - general EMC requirements.</i> |
| Emission: | <p>Standards fulfilled</p> <ul style="list-style-type: none"> ▪ EN 61000-3-2:2014 ▪ EN 61000-3-3:2013 ▪ EN 61326-1:2013 |
| Immunity: | <p>Standards fulfilled</p> <ul style="list-style-type: none"> ▪ EN 61326-1:2013 |

Manufacturer

Metrohm Autolab B.V., Kanaalweg 29-G, 3526 KM Utrecht

Metrohm Autolab B.V. is holder of the TÜV-certificate of the quality system ISO 9001:2015 for quality assurance in development, production,



sales, and repair of measuring instruments and accessories in electrochemistry including technical support (registration number 7528).

Utrecht, November 30th, 2017

J. J. M. Coenen
Head of R&D

A. Idzerda
Head of Production

Appendix A: mutations

Appendix A: mutations

The possible combinations with modules and instruments are shown below in a matrix.

Modules	PGSTAT 128N	PGSTA T302N	PGSTA T100N	PGSTA T302F
BOOSTER10A*	•	•	•	n.a.
BOOSTER20A*	n.a.	•	n.a.	n.a.
FRA32M	•	•	•	•
ECI10M	•	•	n.a.	n.a.
ADC10M	•	•	•	n.a.
SCAN250	•	•	•	n.a.
MUX*	•	•	n.a.	n.a.
BA	•	•	•	n.a.
ECN⁺	•	•	n.a.	n.a.
ECD	•	•	•	n.a.
FI20	•	•	•	n.a.
PX1000⁺	•	•	n.a.	n.a.
EQCM	•	•	n.a.	n.a.
ADC164	•	•	•	•
DAC164	•	•	•	•

(*) (+) These modules are mutually exclusive.

ADC164 and DAC164 modules are not optional.

18.4.2 EU Declaration of Conformity

This declaration attests the compliance of the instrument with the standard specifications for electrical instruments and accessories and with the standard specifications for safety and system validation of the manufacturing company.

Name of commodity

Autolab PGSTAT101, PGSTAT204 used or not used in conjunction with one or more of the modules listed in Appendix A: Mutations
Modular Potentiostat/galvanostat laboratory measurement equipment used for measuring electrochemical cells.

Directives

These instruments have the CE marking and comply with the following EU directives:



- 2014/35/EU (Low Voltage Directive, LVD)
- 2014/30/EU (EMC Directive, EMC)
- 2011/65/EU (Directive for certain hazardous substances, RoHS)

Safety specifications

These instruments comply with safety specifications according to the following standards:

- | | |
|-------------------------|---|
| Design and type testing | <ul style="list-style-type: none"> ▪ EN 61010-1:2016
<i>Safety requirements for electrical equipment for measurement, control, and laboratory use.</i> |
| Testing in production | <p>Every instrument is routine-tested in the production division according to:</p> <ul style="list-style-type: none"> ▪ EN/IEC 61010-1 Appendix F
<i>Check of the protective conductor connection and the insulation from power circuit.</i> |

Electromagnetic compatibility (EMC)

These instruments have been built and have undergone final type testing according to the standards:

- | | |
|--------------|--|
| Requirements | <ul style="list-style-type: none"> ▪ EN 61326-1:2013
<i>Electrical equipment for measurement, control, and laboratory use - general EMC requirements.</i> |
| Emission: | <p>Standards fulfilled</p> <ul style="list-style-type: none"> ▪ EN 61000-3-2:2014 ▪ EN 61000-3-3:2013 ▪ EN 61326-1:2013 |



Immunity: Standards fulfilled
 ■ EN 61326-1:2013

Manufacturer

Metrohm Autolab B.V., Kanaalweg 29-G, 3526 KM Utrecht

Metrohm Autolab B.V. is holder of the TÜV-certificate of the quality system ISO 9001:2015 for quality assurance in development, production, sales, and repair of measuring instruments and accessories in electrochemistry including technical support (registration number 7528).

Utrecht, November 30th, 2017

J. J. M. Coenen
 Head of R&D

A. Idzerda
 Head of Production

Appendix A: mutations

Appendix A: mutations

The possible combinations with modules and instruments are shown below in a matrix.

Modules	PGSTAT101	PGSTAT204
M101 / SYSTAT	•	
M204	n.a.	•
BOOSTER10A	n.a.	•
BOOSTER20A	n.a.	n.a.
FRA32M	n.a.	•
ECI10M	n.a.	n.a.
ADC10M	n.a.	n.a.
SCAN250	n.a.	n.a.
MUX	n.a.	•
BA	n.a.	•
ECN	n.a.	n.a.
ECD	n.a.	n.a.
FI20	n.a.	n.a.
PX1000	n.a.	•

EQCM	n.a.	•
-------------	------	---

M101 or M204 module is not optional.

Only one other module, beside the M101/M204, can exist in the PGSTAT204 at one time

18.4.3 Declaration of Conformity

This declaration attests the compliance of the instrument with the standard specifications for electrical instruments and accessories and with the standard specifications for safety and system validation of the manufacturing company.

Name of commodity

Autolab MAC101, MAC204 used or not used in conjunction with one or more of the modules listed in Appendix A: Mutations

Multi-channel Potentiostat/galvanostat laboratory measurement equipment used for measuring electrochemical cells.

Directives

These instruments have the CE marking and comply with the following EU directives:



- 2014/35/EU (Low Voltage Directive, LVD)
- 2014/30/EU (EMC Directive, EMC)
- 2011/65/EU (Directive for certain hazardous substances, RoHS)

Safety specifications

These instruments comply with safety specifications according to the following standards:

Design and type testing

- EN 61010-1:2016
Safety requirements for electrical equipment for measurement, control, and laboratory use.

Testing in production

- Every instrument is routine-tested in the production division according to:
- EN/IEC 61010-1 Appendix F
Check of the protective conductor connection and the insulation from power circuit.

Electromagnetic compatibility (EMC)

These instruments have been built and have undergone final type testing according to the standards:

Requirements

- EN 61326-1:2013
Electrical equipment for measurement, control, and laboratory use - general EMC requirements.



- Emission: Standards fulfilled
- EN 61000-3-2:2014
 - EN 61000-3-3:2013
 - EN 61326-1:2013
- Immunity: Standards fulfilled
- EN 61326-1:2013

Manufacturer

Metrohm Autolab B.V., Kanaalweg 29-G, 3526 KM Utrecht

Metrohm Autolab B.V. is holder of the TÜV-certificate of the quality system ISO 9001:2015 for quality assurance in development, production, sales, and repair of measuring instruments and accessories in electrochemistry including technical support (registration number 7528).

Utrecht, November 30th, 2017

J. J. M. Coenen
Head of R&D

A. Idzerda
Head of Production

Appendix A: mutations

Appendix A: mutations

The possible combinations with modules and instruments are shown below in a matrix.

Modules	MAC101	MAC204
M101 / SYSTAT	•	n.a.
M204	n.a.	•
BOOSTER10A	n.a.	•
BOOSTER20A	n.a.	n.a.
FRA32M	•	•
ECI10M	n.a.	n.a.
ADC10M	n.a.	n.a.
SCAN250	n.a.	n.a.
MUX	•	•
BA	•	•
ECN	n.a.	n.a.

ECD	n.a.	n.a.
FI20	n.a.	n.a.
PX1000	•	•
EQCM	•	•

Only one other module per M101/M204 module can exist in the multi-channel instrument with a max to 12 modules in total.

18.4.4 Declaration of Conformity

This declaration attests the compliance of the instrument with the standard specifications for electrical instruments and accessories and with the standard specifications for safety and system validation of the manufacturing company.

Name of commodity

Autolab Motorcontroller Controller for high-end rotating disk electrodes for measuring electrochemical cells.

IME663 Controller for VA-stand.

Optical bench Highly focused LED setup for measuring the behaviour of photovoltaic cells.

Directives

The Autolab Motorcontroller and IME663 have the CE marking and comply with the following EU directives:



- 2014/35/EU (Low Voltage Directive, LVD)
- 2014/30/EU (EMC Directive, EMC)
- 2011/65/EU (Directive for certain hazardous substances, RoHS)

The optical bench has the CE marking and complies with the following EU directives:



- 2014/30/EU (EMC Directive, EMC)
- 2011/65/EU (Directive for certain hazardous substances, RoHS)

Safety specifications

These instruments comply with safety specifications according to the following standards:

Design and type testing

- EN 61010-1:2016
Safety requirements for electrical equipment for measurement, control, and laboratory use.
- EN 62471:2006
Photobiological safety of lamps and lamp systems. Applicable to the Optical bench.



Testing in production
 Every instrument is routine-tested in the production division according to:

- EN/IEC 61010-1 Appendix F
Check of the protective conductor connection and the insulation from power circuit.

Electromagnetic compatibility (EMC)

These instruments have been built and have undergone final type testing according to the standards:

- Requirements
- EN 61326-1:2013
Electrical equipment for measurement, control, and laboratory use - general EMC requirements.
- Emission:
- Standards fulfilled
- EN 61000-3-2:2014
 - EN 61000-3-3:2013
 - EN 61326-1:2013
- Immunity:
- Standards fulfilled
- EN 61326-1:2013

Manufacturer

Metrohm Autolab B.V., Kanaalweg 29-G, 3526 KM Utrecht

Metrohm Autolab B.V. is holder of the TÜV-certificate of the quality system ISO 9001:2015 for quality assurance in development, production, sales, and repair of measuring instruments and accessories in electrochemistry including technical support (registration number 7528).

Utrecht, November 30th, 2017

J. J. M. Coenen
 Head of R&D

A. Idzerda
 Head of Production

Appendix A: mutations

Appendix A: mutations

The possible combinations with modules and instruments are shown below in a matrix.

Modules	IME663	Motorcontroller	LED Bank
VA stand	•	n.a.	n.a.

RDE rotator	n.a.	•	n.a.
RRDE rotator	n.a.	•	n.a.

18.5 Environmental protection

The pictogram shown in *Figure 1284* located on the product(s) and / or accompanying documents means that used electrical and electronic equipment (WEEE) should not be mixed with general household waste. For proper treatment, recovery and recycling, please take this product(s) to designated collection points where it will be accepted free of charge.




Figure 1284 The WEEE logo shown on Autolab products

Alternatively, in some countries, you may be able to return your products to your local retailer upon purchase of an equivalent new product. Disposing of this product correctly will help save valuable resources and prevent any potential negative effects on human health and the environment, which could otherwise arise from inappropriate waste handling.

Please contact your local authority for further details of your nearest designated collection point. Penalties may be applicable for incorrect disposal of this waste, in accordance with you national legislation.



Index

Numbers/Symbols

3D plots	725
μAutolab Series	
Back plane	966
Connections for analog signals	967
Front panel	964
Instrument power-up state	964
Instrument test	969
Restrictions	968
Scope of delivery	964
Specifications	970
μAutolab type II	963
μAutolab type III	963
μAutolab type III/FRA2	963

A

AC voltammetry	551
Active cells	895
ADC10M	
Front panel connections ..	1002
Hardware setup	999
Module compatibility	998
Module restrictions	1001
Module settings	999
Module specifications	1004
Module test	1002
Scope of delivery	999
ADC164	
Front panel connections ...	974
Hardware setup	975
Module	973
Module restrictions	980
Module settings	979
ADC750	
Front panel connections ..	1008
Module compatibility	1006
Module restrictions	1008
Module settings	1006
Module specifications	1010
Module test	1009
Scope of delivery	1006
Add analysis	727
Adding commands	
Double click method	666
Drag and drop method	660
Additional components	
Installation	5
Autolab	
Hardware installation	10
Autolab 7 Series	
Back plane	955
Connections for analog signals	956, 957
Front panel	953
Instrument power-up state	952
Instrument test	959
Monitor cable	957
Restrictions	959
Scope of delivery	952
Specifications	961

Autolab Compact Series	
Back plane	930
Connections for analog signals	932
Front panel	928
Instrument power-up state	928
Instrument test	935
Restrictions	935
Scope of delivery	927
Specifications	939
Autolab control panel	96
Autolab display	
FRA2 manual control	1106
FRA32M manual control ..	1118
IME303 manual control ..	1127
IME663 manual control ..	1133
MUX manual control	1147
Autolab display panel	
Docking	131
Instrument Properties	128
Instrument signals	129
Instrument warnings	130
Undocking	131
Autolab F Series	
Connection for analog signals	914
Connections for analog signals	915
Floating mode	918
Grounded mode	918
Instrument test	919
Monitor cable	915
Restrictions	918
Scope of delivery	913
Specifications	923
Autolab N MBA Series	
Front panel	926

- Autolab N Series
 Back plane 904
 Connections for analog signals
 905
 Front panel 901
 Instrument power-up state
 901
 Instrument test 908
 Monitor cable 905
 Restrictions 907
 Scope of delivery 901
 Specifications 910
 Autolab PGSTAT302F
 Instrument power-up state
 914
 Autolab RHD Microcell HC
 Installation 9
 Support 8
 Autolab RHD Microcell HC control
 panel 134
 Autolab spectrophotometer
 Installation 8
 Support 7
 Autolab Spectrophotometer con-
 trol panel
 Absorbance 144
 Data display 144
 Data export 146
 Hardware setup 137
 Manual control 140
 Step through data 145
 Transmittance 144
 Avantes devices
 Installation 8
 Support 7
- B**
- BA
 Front panel connections . 1015
 Hardware setup 1012
 Module compatibility 1012
 Module restrictions 1015
 Module settings 1013
 Module specifications 1018
 Module test 1016
 Scope of delivery 1012
 BIPOT/ARRAY
 Front panel connections . 1022
 Hardware setup 1020
 Module compatibility 1019
 Module restrictions 1022
 Module settings 1020
 Module specifications 1027
 Module test 1023
 Scope of delivery 1019
 Booster10A
 Back plane connections .. 1031
 Front panel controls 1030
 Installation 1031
 Module 1027
 Module compatibility 1027
 Module restrictions 1030
 Module settings 1029
 Module specifications 1038
 Module test 1036
 Scope of delivery 1028
 Booster20A
 Back plane connections .. 1043
 Emergency stop button .. 1042
 Front panel controls 1040
 Installation 1043
 Module 1038
 Module compatibility 1038
 Module restrictions 1040
 Module settings 1039
 Module specifications 1048
 Module test 1045
 Scope of delivery 1038
 Build signal
 Remove filter 344
 Select signals 341
 Using filters 339
- C**
- C1
 Calibration 1102
 C2
 Calibration 1104
 Calculate signal
 Additional functions 354
 Link using drop down list .. 349
 Linking arguments 350
 Mathematical operators 351
 Calibration
 C1 1102
 C2 1104
 Check cell 110
 Chrono amperometry ($\Delta t > 1$ ms)
 556
 Chrono amperometry fast 564
 Chrono amperometry high speed
 573
 Chrono charge discharge galvano-
 static 583
 Chrono charge discharge potention-
 static 580
 Chrono coulometry ($\Delta t > 1$ ms)
 559
 Chrono coulometry fast 567
 Chrono methods
 Level 298
 Repeat 300
 Repeat (unsampled) 301
 Sequence editor 295
 Step 296
 Chrono potentiometry ($\Delta t > 1$ ms)
 561
 Chrono potentiometry fast 570
 Chrono potentiometry high speed
 576
 Cleaning and inspection 899
 Command
 Send email 205
 Stirrer 468
 Command groups
 Grouping commands 655
 Renaming groups 657
 Ungrouping commands 656
 Command options
 Automatic current ranging 618
 Automatic integration time
 641
 Counters 624
 Cutoffs 620
 Plots 633
 Sampler 616
 Value of Alpha 642

- Command stack
 Create stack 675
 Remove command from stack
 677
- Commands
 .NET 229
 .NET - Call a Method 233
 .NET - Call a static method 230
 .NET - Cast object 232
 .NET - Create object 231
 .NET - Get or set a field 233
 .NET - Get or set a static field
 231
 Analysis - general 367
 Analysis - impedance 403
 Apply 239
 Autolab control 236
 Autolab RHD control 483
 Autolab RHD control - Get
 temperature 484
 Autolab RHD control - Set tem-
 perature 485
 Avantes - DIO trigger 475
 Baseline correction 392
 Baseline correction - Exponen-
 tial 396
 Baseline correction - Linear
 395
 Baseline correction - Moving
 average 397
 Baseline correction - Polynomial
 395
 Build signal 337
 Build text 228
 Calculate charge 385
 Calculate signal 345
 Cell 240
 Chrono methods 286, 290
 Chrono methods - High speed
 293
 Chrono methods - Normal 291
 Control 204
 Convolution 381
 Convolution - FRLT differente-
 gration 384
 Convolution - GO differintegration
 383
 Convolution - Kinetic 385
 Convolution - Spherical 384
 Convolution - Time semi-deriv-
 ative 382
 Convolution - Time semi-inte-
 gral 382
 Corrosion rate analysis 398
 Cyclic and linear sweep voltam-
 metry 257
 Cyclic and linear sweep voltam-
 metry - CV linear scan 260
 Cyclic and linear sweep voltam-
- metry - CV staircase 257
 Cyclic and linear sweep voltam-
 metry - LSV staircase 262
 Data handling 332
 Derivative 376
 Dosino 455
 Dosino - Dose 456
 Dosino - Empty 458
 Dosino - Exchange 460
 Dosino - Fill 459
 Dosino - Prepare 457
 Dosino - To end 459
 ECN spectral noise analysis
 387
 ECN spectral noise analysis -
 FFT 389
 ECN spectral noise analysis -
 MEM 390
 EFM 327
 Electrochemical circle fit ... 403
 Electrochemical Frequency
 Modulation 327
 Export data 362
 External devices 472
 FFT analysis 379
 Fit and simulation 405
 FRA measurement 303
 FRA single frequency 305
 General 235
 Generate index 364
 Get item 358
 Hydrodynamic analysis 386
 Import data 358
 Include all FRA data 452
 Increment 224
 Increment - with Signal 226
 Increment - with Value 225
 Integrate 378
 Interpolate 379
 iR drop correction 391
 Kramers-Kronig test 449
 MDE control 251
 MDE control - New drop .. 253
 MDE control - Purge 252
 MDE control - Set stirrer ... 253
 Measurement - impedance
 302
 Message 204
 Metrohm devices 454
 OCP 246
 Peak search 372
 Play sound 227
 Potential scan FRA data 453
 Procedure Information 233
 R(R)DE 250
 Record signals 286
 Regression 373
 Regression - Exponential ... 376
 Regression - Linear 374
 Regression - Polynomial 375
 Remote 469
 Remote - Inputs 470
 Remote - Outputs 471
 Repeat 206

- values 208
- Repeat - Repeat n times 207
- Repeat - Timed repeat 224
- Reset EQCM delta frequency
..... 248
- RS232 - Initialize 480
- RS232 - Send 482
- RS232 control 479
- RS232 control - Close 483
- RS232 control - Receive 482
- Sample processor 460
- Sample processor - Inject .. 466
- Sample processor - Lift 463
- Sample processor - Move . 462
- Sample processor - Peristaltic
pump 467
- Sample processor - Pump . 464
- Sample processor - Stir 466
- Sample processor - Swing . 465
- Sample processor - Valve .. 463
- Set pH measurement tempera-
ture 247
- Shrink data 365
- Smooth 368
- Smooth - FFT smooth 371
- Smooth - SG smooth 369
- Spectroscopy 472
- Spectroscopy - Software trigger
..... 473
- Synchronization 254
- Synchronization - Hardware
..... 256
- Synchronization - Software
..... 255
- Voltammetric analysis 264
- Voltammetric analysis - AC vol-
tammetry 283
- Voltammetric analysis - Differ-
ential normal pulse 273
- Voltammetric analysis - Differ-
ential pulse 270
- Voltammetric analysis - Normal
pulse 267
- Voltammetric analysis - PSA
..... 279, 280
- Voltammetric analysis - PSACC
..... 282
- Voltammetric analysis - Sam-
pled DC 265
- Voltammetric analysis - Square
wave 276
- Wait 240
- Wait - For DIO 242
- Wait - For Metrohm device
..... 244
- Wait - For remote inputs ... 243
- Wait - For seconds 241
- Windower 332
- Common modules 973
- Conformity 1189
- Control amplifier bandwidth
 High stability, High speed,
 Ultra-high speed 884
- Conventions
 Scientific conventions 16
- Corrosion rate analysis
 Polarization Resistance 401
- Tafel Analysis 399
- Counters
 Autolab control 629
- Combining counters 632
- Configuration 625
- Get spectrum 631
- Pulse 627
- Shutter control 630
- Current interrupt 101
- Current range logging 739
- Cutoffs
 Combining cutoffs 623
- Configuration 621
- Cyclic voltammetry 487
- Cyclic voltammetry galvanostatic
..... 490
- Cyclic voltammetry potentiostatic
..... 487
- Cyclic voltammetry potentiostatic
current integration 493
- Cyclic voltammetry potentiostatic
linear scan 496
- Cyclic voltammetry potentiostatic
linear scan high speed 499
- D**
- DAC164
 Front panel connections ... 981
- Hardware setup 981
- Module 980
- Module restrictions 986
- Module settings 985
- Dashboard
 Actions 85
- Instrument panel 92
- Instruments 90
- Recent items 86
- What's going on 88
- Data analysis
 Baseline correction 788
- Corrosion rate analysis 796
- Electrochemical circle fit ... 805
- Fit and simulation 813
- Hydrodynamic analysis 784
- Integrate 775
- Interpolate 779
- Peak search 754
- Regression 771
- Smoothing 748
- Data grid
 Column order 742
- Current range logging 739
- Exporting 743
- Formatting 740
- Sorting 741
- Viewing 737
- Data handling
 Get item 827
- Shrink data 830
- Data overlays 834
- Database 173

Default procedures

AC voltammetry	551
Chrono amperometry ($\Delta t > 1$ ms)	556
Chrono amperometry fast	564
Chrono amperometry high speed	573
Chrono charge discharge galvanostatic	583
Chrono charge discharge potentiostatic	580
Chrono coulometry ($\Delta t > 1$ ms)	559
Chrono coulometry fast ...	566
Chrono methods	555
Chrono potentiometry ($\Delta t > 1$ ms)	561
Chrono potentiometry fast	570
Chrono potentiometry high speed	576
Cyclic voltammetry	487
Cyclic voltammetry galvanostatic	490
Cyclic voltammetry potentiostatic	487
Cyclic voltammetry potentiostatic current integration ...	493
Cyclic voltammetry potentiostatic linear scan	496
Cyclic voltammetry potentiostatic linear scan high speed	499
Differential normal pulse voltammetry	543
Differential pulse voltammetry	539
Electrochemical Frequency Modulation	612
FRA current scan	601
FRA impedance galvanostatic	594
FRA impedance potentiostatic	591
FRA potential scan	597
FRA time scan galvanostatic	608
FRA time scan potentiostatic	605
Hydrodynamic linear sweep with RDE	512
Hydrodynamic linear sweep with RRDE	518
Impedance spectroscopy ..	590
Linear polarization	509
Linear sweep voltammetry	503
Linear sweep voltammetry gal-	

vanostatic	506
Linear sweep voltammetry potentiostatic	503
Normal pulse voltammetry	535
Potentiometric stripping analysis	587
Potentiometric stripping analysis constant current	588
Sampled DC polarography	530
Spectroelectrochemical linear sweep	525
Square wave voltammetry	547
Voltammetric analysis	529
Diagnostics	
Firmware update	1187
Instrument connection ...	1181
Integrator calibration	1186
Options	1186
Running	1183
Differential normal pulse voltammetry	543
Differential pulse voltammetry	539
DIO	
Sending triggers	990
DIO12	988
DIO48	987
Disable commands	658
Dummy cells	
Booster10A test cell	996
Booster20A test cell	997
Dummy cell 2	991
EC110M dummy cell	995
Internal dummy cell	993

E

ECD	
Front panel connections .	1053
Hardware setup	1049
Module compatibility	1048
Module restrictions	1051
Module settings	1049
Module specifications	1055
Module test	1053
Scope of delivery	1049
EC110M	
Bandwidth	1061
Contour map	1061
Front panel connections .	1063
Hardware setup	1056
Module compatibility	1056
Module restrictions	1062
Module settings	1057
Module specifications	1067
Module test	1064
Precautions	1062
Scope of delivery	1056
ECN	
Front panel connections .	1071
Hardware setup	1069
Module compatibility	1069
Module connections	1072
Module restrictions	1071
Module settings	1070
Module test	1073
Scope of delivery	1069
Electrochemical Frequency Modulation	612
Electrode connections	
Four electrode connections	876
Three electrode connections	875
Enable commands	659
End of measurement	719
Environmental conditions	
Temperature overload	897
EQCM	
Front panel connections .	1080
Hardware setup	1077
Module compatibility	1077
Module restrictions	1079
Module settings	1078
Module specifications	1082
Module test	1080
Scope of delivery	1077
Event timing	882
Export plot	731

F

F120	
Front panel connections .	1090
Hardware setup	1085
Module compatibility	1084
Module restrictions	1089
Module settings	1087
Module specifications	1098
Module test	1091
Scope of delivery	1085
Fit and simulation	
Advanced editing	438
Build equivalent circuit	429
Circuit elements	410
Direct fit	814
Edit element properties	442
Equivalent Circuit Editor	816
Generate circuit from CDC	433
Importing and exporting equiv-	
alent circuits	437
Linkable element properties	
.....	444
Load pre-defined circuit	435
Save to library	447
Viewing results	823
Four electrode connections	876
FRA	
Additional properties	307
Additional properties - Adding	
frequencies	318
Additional properties - External	
.....	324
Additional properties - Manual	
modification	317
Additional properties - Multi	
sine measurements	322
Additional properties - Options	
.....	313
Additional properties - Plots	
.....	314
Additional properties - Sampler	
.....	309
Additional properties - Sorting	
table	321
Additional properties - Sum-	
mary	315
FRA current scan	601
FRA impedance galvanostatic .	594
FRA impedance potentiostatic	591
FRA potential scan	597
FRA time scan galvanostatic ...	608
FRA time scan potentiostatic ..	605

FRA2

C1 and C2 determination	
.....	1101
Contour map	1108
Front panel connections .	1109
Hardware setup	1099
Manual control	1106
Module compatibility	1098
Module restrictions	1109
Module settings	1105
Module specifications	1112
Module test	1110
Scope of delivery	1099

FRA32M

C1 and C2 determination	
.....	1113
Contour map	1120
Front panel connections .	1121
Manual control	1118
Module compatibility	1112
Module restrictions	1121
Module settings	1118
Module specifications	1123
Module test	1122
Scope of delivery	1113

G

Grounded cells	896
----------------------	-----

H

Hardware compatibility	2
Hardware description	873
Hardware setup	
Autolab module	99
Optional modules	100
Hydrodynamic linear sweep	512,
518	

I

i-Interrupt	101
IME303	
Back plane connections ..	1129
Front panel controls	1128
Manual control	1127
Module compatibility	1125
Module restrictions	1128
Module settings	1126
Scope of delivery	1125
IME663	
Back plane connections ..	1135
Front panel controls	1134
Installation	1137
Manual control	1133
Module compatibility	1131
Module restrictions	1134
Module settings	1133
Scope of delivery	1132
Impedance spectroscopy	590
Input impedance and stability .	887
Instrument	
Hardware setup	99
Instrument description	900
Instrument information panel ...	97
Instrument panel	92
Instrument Properties	
Sub-panel	128
Instrument Signals	
Sub-panel	129
Instrument Warnings	
Sub-panel	130

K

Koutecký-Levich analysis	386
--------------------------------	-----

L

Legacy hardware support	2
Levich analysis	386
Library	
Add location	177
Arranging columns order ..	192
Column visibility	187
Data repository	195
Default procedures	175
Delete files	194
Edit name	183
Edit remarks	183
Filtering files	188
Hide columns	187
Load file	182
Locating files	193
Merge data	197
Open data	182
Open procedure	182
Preview plot	186
Rating	184
Remove locationl	180
Search	200
Show columns	187
Show in Windows Explorer	
.....	193
Sorting files	191
Tagging	184
Linear polarization	509
Linear sweep voltammetry	503
Linear sweep voltammetry galvano-	
static	506
Linear sweep voltammetry poten-	
tiostatic	503
Links	
Creating	681
Editing	689
Link between more than two	
commands	683
Link between two commands	
.....	682
Linking order	686
Viewing	679
Load from Library	182

M

M101	941
M204	941
Mains connection	9
Manual control	
FRA2	1106
FRA32M	1118
IME303	1127
IME663	1133
MUX	1147
Maximum input voltage	895
Measurement	
Convert data to procedure	745
Post validation	720
Time stamp	719
Measurements	
Manual control	713
Plots frame	704
Procedure cloning	702
Procedure validation	701
Real time modifications	708
Metrohm device installation	7
Metrohm devices	
Support	6
Metrohm devices control panel	
.....	149
Module	
ADC10M	998
ADC164	973
ADC750	1004
BA	1011
BIPOT/ARRAY	1018
Booster10A	1027
Booster20A	1038
DAC164	980
ECD	1048
ECI10M	1055
ECN	1067
EQCM	1075
FI20	1083
FRA2	1098
FRA32M	1112
IME303	1124
IME663	1130
MUX	1142
On-board integrator	1083
pX	1156
pX1000	1162
SCAN250	1169
SCANGEN	1174
Module description	972
Moving commands	
Drag and drop method	669
Moving group	670

Multi Autolab measurements ..	850
Multi Autolab Series	
Back plane	944
Connection hub	947
Connections for analog signals	
.....	945
Front panel	942
Instrument power-up state	
.....	942
Instrument test	948
Restrictions	948
Scope of delivery	942
Specifications	950
MUX	
Front panel connections .	1148
Hardware setup	1145
Manual control	1147
Module compatibility	1144
Module restrictions	1148
Module settings	1146
Module specifications	1155
Module test	1151
MULTI4	1142
SCNR8	1142
SCNR16	1142
Scope of delivery	1145
My commands	
Edit command	695
Saving command	693

N

Normal pulse voltammetry	535
NOVA	
Installation	2
License	12
Requirements	1
NOVA installation	1
Number format	17
Numbering conventions	17

O

On-board ADC	
Connections	975
On-board DAC	
Connections	981
On-board integrator	
Module restrictions	1089
Module settings	1087
Module specifications	1098
Module test	1091
Operating principles	
Active cells	895
Consequences of the digital	
base of the Autolab	883
Current range linearity	892
Event timing	882
Grounded cells	896
Input impedance and stability	
.....	887
Maximum input voltage	895
Oscillation protection	890
Operating principles of the Autolab	
PGSTAT	877
Optional modules	
Properties	100
Oscillation protection	890
Overlays	
Adding data	836
Changing plot properties ..	838
Controls	847
Create overlay	834
Hiding plots	842
Plot settings	840
Removing data	845
Showing plot	842

P

Pause measurement	711
Peak search	
Automatic peak search	757
Exponential	759
Fine tuning results	768
Linear curve cursor	762
Linear free cursos	763
Linear front	764
Linear front tangent	766
Linear rear	765
Linear rear tangent	767
Manual peak search	757
Polynomial	761
Results	770
Zero base	760
PGSTAT12	951
PGSTAT30	951
PGSTAT100	951
PGSTAT100N	900
PGSTAT101	927
PGSTAT128N	900
PGSTAT128N MBA	925
PGSTAT204	927
PGSTAT302	951
PGSTAT302F	912
PGSTAT302N	900
PGSTAT302N MBA	925
pH calibration	
Adding calibration points ..	119
Editing calibration points ..	122
Printing calibration report ..	123
Removing all points	117
Removing points	122
Saving the data	125
Plot	
Export	731
Print	729
Step through data	725
Toggle 3D view	725
Zooming	728
Plot options	
Axes	639
Chart	640
Data	639
Plot preview	721
Plot properties	723
Plots	
Custom plots	634
Default plots	634
Plot properties	636

Plots frame	
Arranging plots	706
Enable and disable plots ...	715
Plot properties	722
Relocating plots	733
Positive feedback	107
Potentiometric stripping analysis	
.....	587
Potentiometric stripping analysis	
constant current	588
Powering the instrument	9
Preview plot	
Library	186
Print plot	729
Procedure editor	
Adding commands	660
Command groups	655
Command links	678
Command stacks	674
Command wrapping	653
Disable commands	658
Enable commands	658
End status Autolab	651
Global options	647
Global sampler	647
Links	678
Moving commands	668
My commands	692
New procedure	645
Procedure tracks	652
Removing commands	660
Stacking commands	674
Zooming	654
Procedure scheduler	
Add open procedures to sched-	
ule	853
Add procedures from the	
Library	856
Add recent procedures to	
schedule	855
Create schedule	853
Inspect procedure	869
Inspecting data	869
Remove instrument	852
Remove procedure from sched-	
ule	857
Running the schedule	864
Runtime control	867
Saving the schedule	861
Start schedule sequentially	
.....	866
Starting complete schedule	
.....	865
Synchronization point	858
Zooming	871



- PSA 587
- pX
- Front panel connections . 1159
 - Hardware setup 1158
 - Module compatibility 1157
 - Module restrictions 1159
 - Module settings 1159
 - Module specifications 1162
 - Module test 1160
 - Scope of delivery 1158
- pX module 114
- pX1000
- Front panel connections . 1166
 - Hardware setup 1164
 - Module compatibility 1163
 - Module restrictions 1165
 - Module settings 1164
 - Module specifications 1169
 - Module test 1166
 - Scope of delivery 1163
- pX1000 module 114
- Q**
-
- Q- 718
- Q+ 718
- Q+ and Q- determination 718
- R**
-
- Real time modifications
- Enable and disable plots ... 715
 - Measurement parameters . 708
 - Procedure control 711
 - Properties modifications ... 708
 - Q+ and Q- determination . 718
 - Reverse scan direction 712
- Record signals
- Fast options 289
 - Fast options and time-derivative threshold 289
 - Normal 288
- Release notes
- NOVA 2.0 80
 - NOVA 2.0.1 66
 - NOVA 2.0.2 58
 - NOVA 2.1 44
 - NOVA 2.1.1 34
 - NOVA 2.1.2 34
 - NOVA 2.1.3 28
 - NOVA 2.1.4 23
- Removing commands 667
- Repeat for multiple values
- Add additional columns 217
 - Add additional values 211
 - Add values using the Add range option 212
 - Clear column 222
 - Delete values 220
 - Edit column name 210
 - Import a table using the Import table option 214
 - Moving columns 219
 - Remove column 223
 - Sorting values 222
- Reverse scan direction 712
- Run
- Measurement 698
 - Procedure 698
- S**
-
- Safety information 1190
- Sampled DC polarography 530
- SCAN250
- Front panel connections . 1172
 - Hardware setup 1170
 - Module compatibility 1170
 - Module restrictions 1171
 - Module settings 1171
 - Module specifications 1174
 - Module test 1172
 - Scope of delivery 1170
- SCANGEN
- Front panel connections . 1177
 - Hardware setup 1176
 - Module compatibility 1175
 - Module restrictions 1177
 - Module settings 1177
 - Module specifications 1180
 - Module test 1178
 - Scope of delivery 1176
- Scientific conventions 16
- Search 200
- Send email
- Command 205
- Skip command 711
- Smooth
- FFT 752
 - Savitzky-Golay 750
- Software
- Compatibility 1
 - Installation 2
 - License 12
- Spare part availability 1190
- Spectroelectrochemical linear sweep 525
- Square wave voltammetry 547
- Step through data 725
- Stop measurement 711
- Supported hardware 2
- T**
-
- Three electrode connections ... 875
- Tools panel 98
- Triggers
- Sending triggers 990
- TTL Triggers
- Sending triggers 990
- U**
-
- Uncompensated resistance
- Determination 101, 107
- V**
-
- Voltammetric analysis 529



W

Warranty 1189