



All-in-One SPE platform  
for thick-film electrodes


# All-in-One SPE Platform



The **All-in-One SPE cell** (Ref. **ED-AIO-SPE-CELL**) provides an **unique multipurpose** interface with **movable add-ons** that can be easily **interchanged** for using the standard *thick-film electrodes* supplied by MicruX Technologies.

## » AIO-SPE-cell general features

The **AIO-SPE-cell** enables the use of the thick-film electrodes in **static** (*Batch-cell*) or **dynamic** (*Flow-cell*) conditions, fulfilling the requirements of **multiple** electroanalytical **applications**.



# All-in-One SPE Platform



The **All-in-One (SPE) cell** is provided in two basic versions, the *supporting base + one add-on* (Ref. **ED-AIO-SPE-CELL-1x**) or the *supporting base + two add-ons* (Ref. **ED-AIO-SPE-CELL-2x**).

## » AIO-cell specific Add-ons features



The **AIO-SPE cell** consists of a supporting base containing the accommodation for placing the standard *thick-film electrodes* (27.5 x 10.1 mm). The supporting base also includes the basic assembly parts for the different add-ons.



## » Batch-cell Add-on



The **batch-cell** add-on enables the use of the standard (27.5 x 10.1 mm) thick-film electrodes in applications (batch analysis, standard additions, etc...) in which are required larger **sample volume** up to 1 mL.



## » Flow-cell Add-on



The **flow-cell** add-on enables the use of the standard (27.5 x 10.1 mm) thick-film electrodes as EC detection system in **flowing liquids** such as FIA, HPLC, CE, etc...



- » **Wall-jet** based flow-cell. The inlet flow is perpendicular to the working electrode surface.
- » Standard **fluidic ports** (¼" - 28 UNF) with inlet channel of 0.8 mm I.D.
- » **Low dead-volume** (internal volume <5 µL). The cell volume is limited by an O-ring (8 mm I.D.).
- » **High sensitivity** electrochemical measurements.
- » **Low sample** requirements (microvolume < 100 µL).



## » AIO-SPE-cell instrument-interfacing

### » AIO-SPE-cell interfacing




A universal **CABLE connector** (ref. **ED-SPE-CABLE**) is supplied with the AIO-SPE cell. The cable enables the interface between the electrodes (up to four contact pads) with the potentiostat.

### » CABLE connector



The cable ends are available with 2 mm female or male bananas.



*Plug to instrumentation may be available in other format under previous request.*



Mora-Garay Industrial Park  
Juan de la Cierva, 2C, Bldg. # 6  
33211 · Gijón (Asturias) · SPAIN

Phone/FAX: +34 984151019

E-mail: [info@micruxfluidic.com](mailto:info@micruxfluidic.com)

Web: [www.micruxfluidic.com](http://www.micruxfluidic.com)

