




MicruX[®] HVStat

MicruX[®] HVStat

HVStat (ref. HVSTAT2010) combines in a compact and portable equipment a *high voltage power supply* and a *bipotentiostat* for dual amperometric detection.



- Miniaturization*
- Easy-handle*
- Portable*
- Adaptable*
- Low-cost*
- Low power requirements*
- Integrated lab-on-a-chip solution*

Technical specifications

- » Dimensions: 165 x 150 x 85 mm (L x W x H).
- » Battery-powered (LiPo – 3300 mAh).
- » Control PC software.
- » Interfacing: Serial RS232/ USB Adapter/ Bluetooth®.
- » LED indicators: power, Bluetooth®, cable.

TECHNICAL FEATURES: High Voltage Power Supply

» Power:	1 W
» Channels/ Outputs:	1/4
» Outputs polarity:	Positive/negative
» Output voltage:	±3000 V
» Max. output current:	0,34 mA
» Output voltage tolerance:	±3 % typical
» Ripple:	< 1%
» Operating temperature:	-20°C to +70°C
» Storage temperature:	-20°C to +105°C
» Humidity:	20% to 85% RH

TECHNICAL FEATURES: Bipotentiostat



» DC-potential range:	± 2,00 V
» Voltage reference:	4,00 V
» DC-potential resolution:	1 mV
» DC-offset error:	± 1 mV
» Accuracy:	≤ 0,1 %
» Current ranges:	1 nA to 1 µA (4 ranges)
» Maximum current:	± 2 µA
» Current resolution:	0,1 % of current range / 1 pA on lowest current range
» Electrochemical techniques:	DC amperometric detection (AD) Pulsed amperometric detection (PAD)
» Run time:	1 s - ∞ (Experiments 1 s – 1000 h)
» Interval time:	10 ms - 1000 s (RS232/ USB)
» Pulse time:	5 ms - 1000 s (RS232/ USB)
» Maximum number of points:	No limited (depending of computer memory)

Specifications are subject to change without previous notice



MicruX[®] Manager

MicruX Manager is a graphical user interface (GUI) to control the high voltage power supply and bipotentiostat of HVStat.


MicruX Manager main interface for single- and dual-mode detection



MicruX Manager experiment window: basic and advanced options




MicruX Manager results window




Separation of neurotransmitters and phenolic compounds performed using MicruX[®] HVStat instrument and different microchips electrophoresis.

:: MCE-SU8-Pt001T



Electropherogram for the separation of 100 μM DA and 500 μM DOPA using a SU-8/pyrex single-channel microchip. Conditions: Running buffer: 25 mM MES-His pH = 6.0; $V_{inj} = +750$ V for 5s, $V_{sep} = +1000$ V, $E_d = +0.75$ V (vs. Pt).

:: MCE-SU8-IDA-Pt005T



Electropherograms for the separation of 100 μM DA, 100 μM pAP and 250 μM APAP using a SU-8/pyrex single-channel microchip with an interdigitated array microelectrode. Conditions: Running buffer: 20 mM MES pH = 6.0; $V_{inj} = +750$ V for 3s, $V_{sep} = +1000$ V.

Severo Ochoa Building · Floor -1 – Room 4 & 6
Julián Clavería s/n · Oviedo (Asturias) · SPAIN

Phone/FAX: +34 984151019

E-mail: info@micruxfluidic.com

Web: www.micruxfluidic.com